Product Instructions

FTBFD9060/FTBFD9960/FTBPW9060/FTPHB9060 FTPHB9960/FTPW9060/FTPW9960

-  (EN) Flip-Top Dilution Bottle
-  (FR) Flip-Top
-  (DE) Flip-Top Verdünnungsröhrchen
-  (IT) Flip-Top flacone per diluizione
-  (ES) Frasco de Dilución Flip-Top
-  (NL) Flip-Top Fles voor verdunning
-  (SV) Flip-Top Dilution Bottle
-  (DA) Flip-Top fortyndingsflaske
-  (NO) Fortynningsflaske med vippelukk
-  (FI) Flip-Top Laimennuspullo
-  (PT) Flip Top Frasco de Diluição
-  (EL) Flip-Top Φιάλη Αραίωσης
-  (PL) Butelka do rozcieńczeń typu Flip-Top
-  (HU) Felpattintható-kupakú hígító palack
-  (CS) Ředící nádobka s uzávěrem Flip-Top
-  (RO) Sticlă de diluție Flip-Top
-  (RU) Флакон Флип-Топ для разведений
-  (TR) Flip-Top Seyreltme Şişesi
-  (JA) フリップトップ 希釈水
-  (ZH) Flip-Top 稀释瓶
-  (TH) ขวดบรรจุสารละลายสำเร็จรูปเพื่อการเจือจาง
-  (KO) Flip-Top 희석액

Product Instructions

Flip-Top Dilution Bottles

Product Description and Intended Use

Neogen® Flip-Top Dilution Bottles come pre-filled with a variety of buffers that may be used as diluents in the enumeration of microorganisms from food and food process environmental samples. The Neogen Flip-Top Dilution Bottles offer simplicity in preparing dilutions.

Table 1: Product Description and Storage Conditions

Catalog Number	Product Description	Volume	Storage
FTBFD9060	Neogen® Flip-Top Dilution Bottle with Butterfield's Buffer	90 mL	20-25°C
FTBFD9960	Neogen® Flip-Top Dilution Bottle with Butterfield's Buffer	99 mL	20-25°C
FTBPW9060	Neogen® Flip-Top Dilution Bottle with Buffered Peptone Water	90 mL	2-30°C
FTPHB9060	Neogen® Flip-Top Dilution Bottle with Phosphate Buffer	90 mL	20-25°C
FTPHB9960	Neogen® Flip-Top Dilution Bottle with Phosphate Buffer	99 mL	20-25°C
FTPW9060	Neogen® Flip-Top Dilution Bottle with Peptone Water	90 mL	20-25°C
FTPW9960	Neogen® Flip-Top Dilution Bottle with Peptone Water	99 mL	20-25°C

Safety

The user should read, understand, and follow all safety information in the instructions for the Neogen Flip-Top Dilution Bottles. Retain the safety instructions for future reference.

NOTICE: Indicates a hazardous situation, which, if not avoided, could result in property damage.

NOTICE

- Neogen has not documented Neogen Flip-Top Dilution Bottle diluents for use in industries other than food. For example, Neogen has not documented Neogen Flip-Top Dilution Bottles for testing water, pharmaceuticals, or cosmetics.
- Neogen Flip-Top Dilution Bottles have not been tested with all possible food products, food processes, testing protocols or with all possible strains of bacteria.
- Do not use Neogen Flip-Top Dilution Bottles in the diagnosis of conditions in humans or animals.
- The user must train its personnel in proper testing techniques, for example: Good Laboratory Practices⁽¹⁾, ISO 7218⁽²⁾, or ISO 17025⁽³⁾.
- For information on documentation of product performance, visit our website at www.neogen.com or contact your local Neogen representative or distributor.

User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at www.neogen.com, or contact your local Neogen representative or authorized distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results. The food sample itself may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any Neogen Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

Limitation of Warranties / Limited Remedy

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, NEOGEN DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any Neogen Food Safety Product is defective, Neogen or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. Please contact your Neogen representative or authorized Neogen distributor for any further questions.

Limitation of Neogen Liability

NEOGEN WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall Neogen's liability under any legal theory exceed the purchase price of the product alleged to be defective.

Storage

Upon receipt, store unopened Neogen Flip-Top Dilution Bottles at the indicated temperatures (See Table 1).

△ Disposal

After use, Neogen Flip-Top Dilution Bottles may contain microorganisms that may be a potential biohazard. Follow current industry standards or local regulations for disposal.

Instructions for Use

Follow all instructions carefully. Failure to do so may lead to inaccurate results.

Wear appropriate protective apparel and follow standard good laboratory practices GLP¹.

1. Remove tamper evident seal. Do not use if seal is broken. Open cap on Neogen Flip-Top Dilution Bottle.
2. Aseptically add sample to the bottle.
3. Close lid and shake to complete the dilution.

If used with Neogen® Petrifilm® Plates:

Refer to the appropriate "Neogen® Petrifilm® Plate Product Instructions".

If you have questions about specific applications or procedures, please visit our website at www.neogen.com or contact your local Neogen representative or distributor.

References:

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Practice Studies.
2. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
3. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.

Refer to the current versions of the standard methods listed above.

Explanation of Symbols

info.neogen.com/symbols

Neogen Food Safety

Neogen Corporation

620 Leshar Place
Lansing, MI 48912 USA
neogen.com

Neogen Europe Ltd.

The Dairy School
Auchincruive
Ayr, KA6 5HU
Scotland, UK

Neogen Ireland, Ltd.

Bray Business Park, Bray
Co. Wicklow
A98YV29, Ireland

Neogen Corporation

620 Leshar Place Lansing, MI 48912 USA
www.neogen.com

Neogen and Petrifilm are trademarks of Neogen Corporation.

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Referenties:

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Practice Studies.
2. ISO 7218. Microbiology of food and animal feeding stuffs - General requirements and guidance for microbiological examinations.
3. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.

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1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Practice Studies.
2. ISO 7218 Microbiology of food and animal feeding stuffs - General requirements and guidance for microbiological examinations.
3. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.

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All rights reserved.
© Neogen Corporation 2024.