700002758  
9548

# Product Instructions

Reveal® for Histamine

Reveal® para Histamina

## Product Instructions

# Reveal for Histamine

### Introduction and Intended Use

High levels of histamine may develop in a variety of fish species as they decompose. These species include tuna, mahi-mahi, marlin, bluefish, sardines, anchovy, bonito, herring, and mackerel. Ingestion of histamine may cause scombroid poisoning in humans, which may lead to a variety of symptoms, including rash, nausea, vomiting, diarrhea, hypotension, heart palpitations, and muscle weakness. Paralysis and death have also been reported in cases of scombroid poisoning.

Because of its potential for human illness, the U.S. Food and Drug Administration (FDA) has ruled that extensive refrigeration records and/or histamine testing must be included in a hazard analysis and critical control point (HACCP) program for relevant fish species. The FDA has set an action level of 35 parts per million (ppm) for histamine in domestic and imported fish.

### Intended Use

Reveal® for Histamine is intended for the qualitative analysis of histamine in scombroid species of fish, such as tuna and mahi-mahi.

### Intended User

The test kit is designed for use by quality control personnel and others familiar with histamine analysis in fish. Since technique is very important, operators should be trained by a Neogen representative or someone who has completed Neogen training.

### Assay Principles

Reveal for Histamine is a single-step lateral flow assay based on a competitive immunoassay format. The extract is wicked through a reagent zone, which contains histamine-specific antibodies conjugated to colored particles. If histamine is present, it will be captured by the particle-antibody complex. The histamine-antibody-particle complex then is wicked onto a membrane, which contains a zone of histamine conjugated to a protein carrier. This zone captures only uncomplexed antibody-particle conjugate, allowing the colored particles to concentrate and form a visible line. The membrane also contains a control zone where the antibody-particle complex present in the reagent zone is captured by an antibody, forming a visible line. The control line should always form regardless of the presence of histamine, indicating the strip is functioning properly. If histamine is present above the pre-determined limit, no test line will form in the test zone. Therefore, a positive sample will result in an absence of a line in the test zone. Alternatively, a negative sample will form a line in both the control zone and the test zone.

### Storage Requirements

Store kit components at room temperature, 18–30°C (64–86°F) to facilitate full shelf life. Test strips should remain capped in their original sample tubes until used to facilitate optimal performance.

## Materials Provided

1. 25 histamine test strips
2. 25 sample cups
3. 1 bottle sample diluent
4. 25 red dilution cups

## Materials Recommended But Not Provided

1. Graduated cylinder (700006494 | 9368)
2. Sample extraction bottles, 250 mL (700004005 | 9399)
3. Sample collection tubes (700002726 | 9421)
4. Blender (700004031 | 9493, 700002745 | 9477)
5. Scale capable of weighing 10–50 g (700004010 | 9427)
6. Pipettor, 100  $\mu$ L (700003935 | 9278)
7. Pipette tips (700004009 | 9417)
8. Pipettor, 1,000  $\mu$ L, (700003969 | 9337)
9. Pipette tips (100001373 | 9464)
10. Paper towels or equivalent absorbent material
11. Rack for Reveal sample cups (100001376 | 9475)
12. Timer (100001358 | 9426)
13. Waterproof marker
14. Distilled or deionized water

## Precautions

1. The test strips must remain inside the stay-dry tube before use.
2. Store test kit at room temperature 18–30°C (64–86°F) when not in use.
3. Do not use kit components beyond expiration date.
4. Treat all used liquids, including sample extract, and labware as if contaminated with histamine. Gloves and other protective apparel should be worn at all times.
5. To avoid cross-contamination, use clean equipment for each sample, and thoroughly wash all equipment between samples.
6. Glassware should not be used for extraction purposes. Histamine may adhere to glass, therefore using glassware may affect test results.

## Sample Preparation

1. Canned tuna: AOAC 937.07b  
Place entire contents of can/pouch (meat and liquid) into a blender. Blend until homogenous. Store samples at 2–8°C (35–46°F) until analyzed.
2. Fresh or thawed frozen raw fish: AOAC 937.07a  
Clean and eviscerate three fish. Cut three cross-sectional pieces 2.5 cm (1 inch) thick, from back of the pectoral fin, halfway to the vent, and one posterior to the vent. Debone slices and blend or grind combined samples until homogenous. Store samples at 2–8°C (35–46°F) until analyzed.

## Sample Extraction

1. Add 10 g of the homogenous mixture to a clean disposable extraction bottle containing 190 mL of distilled or deionized water.
2. Tightly cap and vigorously shake the bottle for 15–20 seconds to suspend the fish tissue in the water.
3. Wait approximately 5 minutes, and then shake the bottle for 15–20 seconds to resuspend the fish tissue.
4. Wait an additional 5 minutes, and again shake the bottle for 15–20 seconds to resuspend the fish tissue. Allow the tissue to settle to the bottom of the bottle for about 30 seconds. The sample is now ready for extract dilution.

## Test Procedure

1. Remove the appropriate number of clear sample cups and red dilution cups and place in a sample cup rack.
2. Add 100  $\mu$ L of sample extract to each red dilution cup.
3. Add 1,000  $\mu$ L of sample diluent to the red dilution cup. Mix well by pipetting up and down 5 times.
4. Transfer 100  $\mu$ L of diluted sample extract to a clear sample cup.
5. Place a new Reveal for Histamine test strip with the sample end down into the sample cup and set timer for 5 minutes. Ensure the test strip comes into contact with liquid and begins to wick.
6. Allow the test strip to develop in the sample cup for 5 minutes.
7. Remove the strip and read results visually by holding the strip up to line of sight. Strip must be read within 30 seconds after the 5 minute incubation.

## Test Interpretation

### Negative Result

If a line forms in the test zone and another line forms in the control zone within 5 minutes (2 lines total), the sample contains less than 35 ppm histamine (negative sample).

**Note:** The presence of the control line, regardless of intensity, indicates a valid test.

### Positive Results

If after the full 5 minutes there is no visible line in the test zone, but a visible line in the control zone, the sample contains at least 35 ppm histamine (positive sample).

### Invalid Results







If the control line is absent, the test is invalid.

## Notes

1. Observations should be made immediately following the 5 minute development. Observations after 5 minutes may be inaccurate due to over-development of the device.
2. Positives can be confirmed by Neogen's Reveal for Histamine, Veratox for Histamine, or HPLC.

## Validated Matrices

Tuna and mahi-mahi.

## Customer Service

Neogen Customer and Technical Services can be contacted through [neogen.com](https://neogen.com) and product training is available by request.

## SDS Information Available

Safety data sheets are available for all test kits at [neogen.com](https://neogen.com) or by calling 800.234.5333 or 517.372.9200.

## Terms and conditions

Neogen's full terms and conditions are available [online](#).

## Warranty

Neogen makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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