**NEOGEN<sup>®</sup>**  
**ANSR<sup>®</sup>**

***USER***  
***MANUAL***

**ANSR reader and software**

1. Introduction .....	2
2. Warnings and limitations.....	2
3. Materials provided .....	2
4. How it works.....	2
5. Using ANSR with a computer.....	3
6. ANSR stand-alone operation .....	10
7. Cleaning and maintenance.....	12
8. Troubleshooting .....	12
9. Limited warranty.....	12
10. Terms and conditions.....	12
11. ANSR reader specifications.....	13
12. Neogen company profile .....	14

## 1. Introduction

This manual describes setup, configuration and operation of a reader and software for use with the Neogen ANSR® pathogen detection system.

The incubation function of the reader has been set to maintain a constant temperature of  $56 \pm 1^\circ\text{C}$  to meet the specific requirements of the Neogen ANSR pathogen detection system.

This reader has been designed to operate in one of two methods:


1. With a dedicated personal computer and ANSR software package to track, interpret and analyze assay results.
2. As a stand-alone unit that displays assay results on the reader's screen.

The flexibility of the Neogen ANSR pathogen detection system allows users to evaluate up to 48 samples an hour depending on an individual laboratory's workflow needs.

For more information on the use of the reader, or the Neogen ANSR pathogen detection system, please contact Neogen Technical Services or a Neogen sales representative at 800/234-5333 or 517/372-9200.

## 2. Warnings and limitations

1. Always operate the instrument on a dry, level surface out of direct sunlight.
2. Never move the instrument while a test is in progress.
3. To avoid damaging the reader, never place objects on top of it.
4. The heater block in the reader operates at up to  $56^\circ\text{C}$  and may cause contact burns or damage to materials in contact with these areas.
5. This instrument is designed only to operate with the provided power supply plug pack. This module forms part of the system. Do not operate the system with a different power supply module. The correct power supply is required to maintain the safety and electromagnetic compatibility of the system.
6. To avoid electrical shock, do not operate the instrument or the power supply plug pack if it has been opened, damaged or exposed to moisture, condensation or rain. The external power supply plug pack is sealed with no user serviceable parts. Do not operate this module with any damaged or exposed parts.
7. Do not open or attempt to repair the reader or other accessories. This instrument does not contain serviceable parts and should be returned to the supplier for repair. Attempting to service the instrument also will void the warranty.
8. The real-time clock coin battery included in the equipment is not a user replaceable item.
9. Only use the instrument for its intended purpose and in accordance with this manual and warnings. Protection provided by the equipment may be impaired if the equipment is operated in a manner contradictory to the instructions and warnings. This instrument (including power supply) is designed to operate within the manufacturer's specifications. Do not exceed the manufacturer's specifications when in use.

**Simpler. Better. Faster.**

10. Keep connected cables clear of work areas. The main socket outlet is intended to be used with the reader's external power pack and should be readily accessible. It is recommended that the user unplug the reader when not in use.
11. If intended for connection to external equipment, ensure the interfaces of such equipment are separated from the main socket by double or reinforced insulation and present no risk of electrical shock.

## 3. Materials provided

- Neogen ANSR reader
- Computer with pre-installed software (Windows 7/10)
- Case
- 4 power cord plug adapters
- Ethernet cable (if purchased with computer)
- USB cable
- Power cord

## 4. How it works

The ANSR reader uses two channels of light, referred to as **Control** and **Test**, to analyze samples. Fluorescence is generated by periodically exposing the reagents in the tubes to the **Control** and **Test** wavelengths. The resulting fluorescence is captured by photodiodes, where it is converted to millivolts as a function of time. Mathematical algorithms are applied through the software to signal the detection of target pathogens in the sample.

## Turning the reader on and off

To turn the reader on, press the rectangular power button located on the lower left corner of the reader's front panel for **2–3 seconds**. The reader will take about **2 minutes** to boot.

To turn the reader off, press the power button for **2–3 seconds**.

## 5. Using the ANSR reader with a computer (remote mode)

### 5.1 Set up

Setting up the reader depends on the desired usage: with a computer or stand-alone. If using a computer, ensure the reader is connected to the laptop before turning on the reader. Once the reader is turned on, it will heat to its resting temperature. **NOTE:** The reader's default setting is remote mode. To switch a reader previously used for stand-alone mode to remote mode, see Section 6. The reader is in stand-alone mode if the touch screen displays three icons upon start up: a test tube, a set of gears and a magnifying glass.

### Remote mode (connected to a computer)

This mode allows the reader to be used with a computer and the Neogen ANSR software. Results will be displayed on the computer screen. **NOTE:** Prior to starting the reader, ensure connected computer is on and the ANSR software is running.

1. Plug the reader into a power source, using one of the provided outlet adapters if necessary. **NOTE:** If the computer is to be connected to a network, do so prior to connecting the reader to the computer.
2. Using either the provided USB cable or the Ethernet cable, connect the reader to the computer by plugging the cable into the appropriate port on the back of the reader. Connect the other end of the cable to the provided laptop. **NOTE:** If using an Ethernet cord to connect the reader to the computer, the computer will not be able to connect to a network via Ethernet.
3. Press the rectangular power button located in the lower left corner on the reader's front panel for **2–3 seconds**. The reader will take about **2 minutes** to boot.

*ANSR reader with case, power cord, USB cable, Ethernet cable, and alternate plug adapters*


## 5.2 Installing the software

The ANSR reader software comes pre-installed on the computer provided in the start up system.

## 5.3 Using the software

### 5.3.1 Launching the software

1. Open the reader program by double clicking on the **T16 ISO** desktop icon. 
2. Before running tests, the computer must be connected to the reader. To do this, click the **Connect** icon in the upper right corner of the software window. Select either **USB** or **Networking**, depending on the cable used to connect the reader to the computer. Select the appropriate choice and hit **Connect**. The circle in the upper right corner will turn green and indicate the reader is connected. **NOTE:** Users also may use the **Locate** button, which is in the **Connect** menu to ensure the reader is connected. Clicking **Locate** will cause the reader to emit a beep. This feature only is available when the reader is connected via a networking cable.


### 5.3.2 Software utilities

**Open:** To open a previously saved file, click **Open** and select the appropriate .json file. **NOTE:** The reader does not need to be connected to the computer in order to use the software.

**Print:** Allows users to print a summary report and/or tube graphs.

**Export:** Allows users to export all files associated with the test or only desired files. Click **Export** and select which files to export by checking or unchecking boxes located next to the options. Click **Export** once all selections are made. See example below.


**Save as:** Allows the user to select the location where they would like the data saved.

**Support:** This button displays a copy of the user manual.

**Instrument:** This button displays software and firmware version information.

### 5.4 Running a test (with computer)


1. Once the reader is set up and connected to the software, tests can be run. **IMPORTANT:** Do not open the reader's lid when the test is running.
2. To begin a test, click **Run** and select which test the reader will perform (e.g., ANSR *Salmonella*) from the drop-down menu.

3. The screen will read, "**Please wait. Setting temperature.**" It also will display the target temperature of 56°C and the instrument's current temperature. During this time, users may enter a user name, lot number and sample identification numbers. Users may manually type in sample identification numbers or may use the **Autofill** option, which is activated by entering a sample ID in the first field and clicking the drop-down arrow. **Autofill** then will fill all 16 ID fields with sequential IDs. Numbers also can be added by copying and pasting fields from Microsoft Excel. **NOTE:** If no sample ID is entered, the result will not be displayed. Also note that user name, lot number and sample IDs cannot be edited after the test run has started. The sample ID must have a number at the end to use **Autofill**.

4. Once the instrument is ready, the window will display, "**Instrument ready**". After placing the samples in the reader, press the **Start workflow** button near the bottom of the screen to begin the test. **NOTE:** For instructions on preparing samples and reagents, and on loading them into the reader, see the instructions that came with the ANSR test kit. Wait to add samples until the reader reaches 56°C. Ensure the permanent caps are placed onto the tubes to avoid contamination.

5. During the test, the user can select from a range of real-time graphs on the left side of the program window. Individual tubes

can be selected along with a total view, which will show the formation of the curve that denotes a negative or a positive result.

- Once the assay is complete, the software will display an overall graphical view of the results. See example above.
- Click **Results** on the left side of the window to view a graphic depicting whether the tubes were positive (red), negative (green) or invalid (yellow). On the left side of the window, users have the

option to click on individual tube tabs for graphs related only to one sample, all of the graphs together, or temperature data. To view the **Summary report**, which contains a summary of the test results, click **Print** and check the box for **Summary report**. Users also can print the **Tube graphs** from the **Print** menu. Once the desired selections are made, click **Print preview** to review the selections, and click **Print**.

- The data can be saved to the computer or exported. To export data, click the **Export** icon in the toolbar. Select which data files to export (e.g., the data file, a .csv file for use in Microsoft Excel, or a graph). **NOTE:** The data file is the most important and contains results and traceability information. It is saved as a .json file and only can be opened using the ANSR software. The data automatically is saved in the **Records** folder on the computer's C: drive (**Program files(x86)/Axxin/AxxinT16ISODesktop/TestResults/Records**). The default file name includes the instrument ID, date and time. The file extension is .json. This file type can only be opened using the ANSR software. The user can use the **Save As** function to save the file with a specific name and in a user-specified location.
- To begin a new test, click **New** in the toolbar. The software will autofill the user name and lot number. These may be edited. **NOTE:** Ensure data is saved before clicking **New**.

- To shut down the reader, click **Disconnect** in the upper right corner of the software window. Once the reader is disconnected from the software, the computer and reader may be shut down.

### Running high throughput samples – Two step lysis (ANSR for *Salmonella* and ANSR for *Listeria*)

**NOTE:** This procedure begins after the 37°C incubation step. See the ANSR for *Salmonella* or ANSR for *Listeria* kit insert for further instructions.

The ANSR reader can process 32 samples (*Listeria*) or 48 samples (*Salmonella*) an hour by adhering to the following procedure:

- Follow the procedures for preparing and enriching the samples in the ANSR kit insert. Ensure the reader is on and heated to 56°C.
- Place up to four rows of lysis tubes containing samples in the 80°C heater block for **20 minutes** per the kit's instructions.
- Transfer the lysed samples from the first row of tubes to reaction tubes 1–8 (from the left) in the reader. Transfer the lysed samples from the second row of tubes to reaction tubes 9–16 (from middle) in the reader. Follow pipetting, mixing and capping procedures as listed in the ANSR kit insert.
- Close the reader's lid and run the test by pressing **Start** on the software screen.
- Once the 10 or 18 minute assay is complete, save the results to a desired location. Discard the used reaction and lysis tubes once results have been obtained. **NOTE:** After each run of 16 samples, select **New** for the next set of 16. The remaining lysis tubes containing samples may remain in the heater block for up to **60 minutes**.
- Complete the assay with the next two rows of lysis tubes containing samples.
- Repeat steps 3–6 until all samples have been run.


## 5.5 Open Remote List NEW

Open Remote List allows a user to access and download test files saved on a remote server.

Users can navigate SFTP or FTP servers, and download relevant results for review. The T8 desktop software supports viewing (and listing) of T8 results only (.csv files).

### 5.5.1 Add and Remove 'Open Remote List' Profiles

To set up a remote connection, click on the Open Remote List button. From the drop-down list, select New Connection.

When prompted, fill in the remote server's details in the Add Connection window.

The following fields must be filled in:

- **Address** – IP address or domain name of the remote server. eg: 118.168.4.75 or sftp.example.com
- **Type** – The type of server the user is connecting to. This can be either SFTP or FTP.
- **Port** – The port number through which the server can be reached.
- **Destination Directory** – The user can explicitly specify which subdirectory (on the remote server) should be viewed. If this field is left blank, the root directory of the server will be displayed.
- **User Name/ Password** – The user's login credentials for the remote server.

It's useful to click on the Test Connection button to confirm the entered details are valid. An error icon will be displayed if the connection was unsuccessful, along with further details for troubleshooting.

Hovering over the error icon will provide helpful information for troubleshooting

A few common connection errors are:


- Incorrect IP address
- Incorrect server type – the desktop software supports connecting to SFTP/FTP servers
- Incorrect or invalid user login credentials – The user must have access to the server for results to be listed.

If all fields are valid and correct, Test Connection will indicate that a connection could be successfully made.

Test Connection was successful. The Open Remote List profile can now be added.

The user can now add the profile, by clicking the Add button.

Once a profile has been created, the user can add other profiles too, to browse other desktop servers.


A list of other valid Open Remote List profiles. Each will pull results from a different server.

To delete a profile, click the red 'X' icon as shown below:

### 5.5.2 Browsing, Viewing and Downloading Results

Once a valid connection has been established, the user can now browse results on the remote server.

A collection of test result files on the remote server.


The user can navigate up a folder by clicking “Open” on the Parent Folder. Similarly, clicking “Open” on a subfolder will display that subfolder’s contents.

A user can also search for test results using the search field.


The search field will look for matching filenames, and will return a list of partial matches as well.

Once a relevant result has been located, click on the Load button:


View a result by clicking Load.

The result is loaded into the desktop view:

Viewing a result remotely does not automatically download it. To do this, click on the “Export” button, and provide a local location for the file to be saved to.


## 5.6 Sync Page NEW

The T16-ISO desktop now supports synchronizing test results saved on the laptop to a remote server. This is known as the “Sync” function. This Sync is unidirectional; results can only be transmitted to the server, and not retrieved. Up to two remote Sync profiles can be configured for use.

A sync profile will look for test records to transfer to the server. It’s important to note that only results that reside in the Records folder will be synchronized. By default the records folder is created in *C:\Program Files (x86)\Axxin\Axxin 16 ISO Desktop\TestResults\Records*

In addition, the software can also be configured to act as a server, and receive test results when transmitted by another desktop sync profile. This is convenient when several instruments are in use, and it’s necessary to store a copy of all results in one centralized location.

### 5.6.1 Add an Outgoing Connection

Adding an ‘outgoing connection’ can be thought of as creating an automatic backup. When an outgoing connection is configured, any new results will be transmitted to the designated server.

To create a sync profile, click the ‘Add Outgoing Connection’ button. When prompted, enter the required server details:

**Add New Sync Destination**

Address:

Type:

Port:

Destination Directory:

User Name:

Password:

Periodic Sync:  Enabled  Disabled

The following fields need to be set:


- **Address** – The address of the remote server. (E.g.: 192.168.10.10 or *server.example.com*)
- **Type** – Selected between SFTP and FTP
- **Port** – The server access port number
- **Destination Directory** – The location where files will be stored on the server. (E.g.: *MyT16Results* will store the results in *server.example.com/MyT16Results*)
- **User Name/ Password** – Authentication credentials for the server
- **Periodic Sync** – If enabled, the software will transfer any unsynced results every 10 minutes to the server. If disabled, the user must manually initiate the Sync function.

Click Test Connection to confirm that with a valid username and password, a connection can be made to the server. If successful, a green 'success' icon will be displayed:


If unsuccessful, the software will indicate which field was invalid.

In the above screenshot, the login credentials were incorrect.

Hovering over the red exclamation icon will provide more information about the error.

## 5.6.2 Sync Operation (Outgoing Connections)


Two outgoing connection profiles can be configured for use. It's important to note that while a test is running, no Sync operations will be performed.

The `\TestResults` folder is the 'root' folder which the Sync function will search through to find files. Additionally, the Sync function will search 3 levels deep through subfolders for other files to transmit. It's important to note that test results are saved in `\TestResults\Records` by default.


As an example, any files that have been saved in `\TestResults\Records\SubFolder` will be transmitted, however, data in `\TestResults\Records\SubFolder\AnotherFolder` will be ignored.

NOTE: Sync is restricted to transmitting **only** .csv files, the main file format exported/saved by the T8 desktop software. Other files in the `\Records` folder (if any) are not transferred.

NOTE: Result files larger than 5Mb will not be transferred. During normal operation, it is unlikely that a result file larger than 5Mb will be generated.

Clicking on the 'Sync Now' button will force the desktop software to connect and transmit data to the server.


Clicking on the 'More' button will display additional profile options:

- **Remove** – Will remove the sync profile.
- **Edit** – This will allow the user to edit the sync parameters for the profile, as in Figure 61.
- **Test Connection** – Test the current connection parameters. If the connection was successful, a green tick icon will be displayed next to the test connection button. If unsuccessful, a red error icon is displayed.
- **Enable/Disable Periodic Sync** – Toggle between Periodic Sync being enabled or disabled. Once enabled, the Sync process will begin immediately to transmit any unsent results.

While files are being transmitted, the sync profile displays the text 'Running' and begins counting down the number of files remaining to be pushed to the server. The user cannot trigger a manual sync or modify the contents of the sync profile while a Sync is in progress.

Once all files have been successfully pushed to the server, the User Interface will display 'Synchronized'. If periodic send is enabled, the time the last transmission was successful will be displayed along with the time the next Periodic Sync will happen.

## 5.6.3 Add Incoming Connection (Server)


The software can be configured as a server, to receive data transmitted by a sync profile (outgoing connection) from another instance of desktop software on a separate computer.

Click Add Incoming Connection:

The following fields need to be filled out by the user:

- **File Storage Path** – The location where received results will be saved.

- **Start Server on Instrument Startup** – If this is checked, the server will be online and ready to receive files, as soon as the desktop software is powered up.
- **Connection Port** – Specify the port number over which results can be received.
- **User management** – At least 1 user needs to be created before a server can be set up. (see section below)

### 5.6.3.1 Configuring Users


To create a user, enter a user name between 5 and 20 characters.

Enter a password for this user. The password must:

- be between 6 and 20 characters
- contain at least 1 capital letter
- contain at least 1 number
- contain at least 1 non-alphanumeric character (eg: #, \$, @)
- not match the username.


Click on “Show Passwords” at the bottom of the panel to simplify user creation.

To delete a user, click on the red “X” next to their entry.


### 5.6.4 Sync Server Operation (Incoming Connection)

Once a server is added, the desktop software is now ready to receive test results.


Only 1 Sync server (Incoming Connection) configuration can be set up and active at any time. It's important to note that the sync server will only run as long as the desktop software is open; if it is closed, any profiles which are configured to transfer results to it, will be unable to do so.

A server is set up and running on these IP addresses.

The following pieces of information are displayed to the user:

- The local IP address(es) on which the server can be connected to. (Note: if attempting to connect from a computer that is not on a local network, use your global IP address instead.)
- **The Connection Type** – the server type is SFTP
- **Auto-start** – Whether the server will automatically start once the desktop software starts
- **Port** – the port on which data can be transmitted
- **Storage Location** – The folder in which received data will be stored.
- **Server state** – Online or Offline
- **Current Connections** – The number of users currently connected to the server
- **Files Received** – The total number of files received for the duration of the server session (i.e.: from the last time the desktop software was opened).
- **Last File Received** – The filename of the last file transmitted.

## 6. ANSR stand-alone operation

### 6.1 Set up

To operate the reader without the use of a computer or a network connection, plug the reader into an appropriate power outlet. Turn the reader on by pressing the rectangular power button located on the lower left corner of the reader's front panel for **2–3 seconds**. The reader will take about **2 minutes** to boot. Once ready to use, the screen will display three icons: a test tube (**Run test**), a magnifying glass (**Search**), and gears (**Settings**). The reader is now ready for use.

### 6.2 Printer set up

1. Touch the gears (**Settings**) icon on the **Home screen**.
2. Touch the down arrow followed by the **Printer Set up** icon.
3. Using the touch screen key pad, enter the IP address of the desired printer and press the green check mark to confirm.

### 6.3 Remote mode

To switch from stand-alone operation to operation with a computer (remote mode), users must ensure remote mode is activated prior to attaching the reader to a computer. For instructions on how to set up the reader for use in remote mode, see Section 5.

1. Touch the magnifying glass (**Settings**) icon on the **Home screen**.
2. Touch **Remote mode**. Touch the box in the lower right corner next to **Bootup into remote mode**. The next time the reader is turned on, it will be in remote mode.

### 6.4 Settings

#### About this device

This option provides information about the instrument name, ID, application version, IP address, MAC address, OS version, kernel version and light channel information.

1. To access **About this device**, touch the gears (**Settings**) icon on the **Home screen**.
2. Select **About this device**. The information will be displayed.

#### Help info

1. To access **Help info**, touch the gears (**Settings**) icon on the **Home screen**.
2. Select **Help info**. The information will be displayed.

#### Audio

This option allows the user to change the volume of the reader's tones.

1. To access **Audio** settings, touch the gears (**Settings**) icon on the **Home screen**.
2. Touch the **Audio** icon.
3. Use the up and down arrows to adjust the reader's tones for the touch screen and for alerts (0–low, 5–high). Press **Test** to hear the sound. To confirm the selection, press the green check mark.

#### Brightness

This option allows the user to adjust the brightness of the touch screen.

1. To access **Brightness** settings, touch the gears (**Settings**) icon on the **Home screen**.
2. Touch the **Brightness** icon.
3. Using the white up and down arrows, adjust the brightness of the screen to the desired level (30–low, 100–high). The screen brightness will change as the user touches the arrows. Touch the green check mark to confirm the selection.

### 6.5 Administrator settings

Administrator settings are password protected. For the password, please contact a Neogen Technical Services representative.

#### Import and export

1. Touch the gears (**Settings**) icon on the **Home screen**.
2. Touch the **Administrator settings** icon. Enter the password and touch the green check mark to confirm.
3. Select **Import and export functions**. Plug in a USB drive to the port located on the back of the reader.
4. Touch the **Export test results** icon. This will export results as a .json file that can only be read using the ANSR software.
5. The reader will ask if a copy of the exported results should remain on the reader or be deleted. Make the appropriate choice and touch the green check mark to confirm.

#### Export data log

This option exports a folder named with .json files (which only can be opened using ANSR software) and .csv files, which can be opened in Microsoft Excel.

1. Touch the gears (**Settings**) icon on the **Home screen**.
2. Touch the **Administrator settings** icon. Enter the password and press the green check mark to confirm.
3. Select **Import and export functions**. Plug in a USB drive to the port located on the back of the reader.
4. Select the **Export data log** icon. The reader will export a folder containing a .json file and a .csv file.

#### Import test types

To import a new test type, please contact Neogen Technical Services.

#### Time and date

This setting allows the user to adjust the date and time on the reader.

1. Touch the gears (**Settings**) icon on the **Home screen**.
2. Touch the **Administrator settings** icon. Enter the password and press the green check mark to confirm.
3. Select **Time and date**.
4. Use the touch screen to select the item to be changed (e.g., a.m. or p.m., the month, the hour) and use the white up and down arrows to scroll through the options until the desired information is reached.
5. Select the green checkmark to save the time and date settings.

## Park optics

Selecting **Park optics** will move the optical head to a position where it can be removed for service purposes. This option should only be used by Neogen Technical Services.

## Network settings

Network settings allow the user to configure the reader's Ethernet connection. Users can choose between DHCP or a static network, enter an IP address and enter a subnet mask.

The default setting is DHCP. To switch to a static IP address:

1. Touch the gears (**Settings**) icon on the **Home screen**.
2. Touch the **Administrator settings** icon. Enter the password and press the green check mark to confirm.
3. Select **Network**.
4. Select **Static** by touching the box next to the word. It will turn white once selected.
5. Using the touch screen keypad, input the IP address and the subnet mask.
6. Press the green check mark to save the settings.

## Reset to factory default

For Neogen Technical Services use only.

## Factory settings


For Neogen Technical Services use only.

## 6.6 Running a test

### Prior to placing the sample in the reader:


1. Touch the test tube (**Run test**) icon.
2. Touch the **Test type** icon and select the appropriate test (e.g., ANSR for *Salmonella*). Once the desired choice has been selected, touch the green check mark to confirm or the red X to choose a different test.
3. The screen will display the selected test type. Touch the green check mark to confirm or the red X to go back and select a different test type.
4. Enter the test name by using the touch screen. Once the name is input, press the green check mark to confirm. **NOTE:** The bar code scanner function is not yet validated.

5. The screen will display the individual test tubes. Touch the desired tubes and input the sample ID using the touch screen. Press the green check mark to confirm. **NOTE:** Users also may use the **Autofill** option to fill in sample IDs. To do this, enter the sample ID of the first test tube by using the touch screen. Once the sample ID is confirmed, touch the blue **Autofill** button at the bottom center of the screen. This will autofill the remaining test tube IDs with


sequential numbers from the first tube. This option only is available when beginning from the first tube. The sample ID must have a number at the end to use **Autofill**.

6. Once the reader reaches 56°C and the start button appears, place the samples into the reader. **NOTE:** For instructions on preparing samples and reagents, and on loading them into the reader, see the instructions that came with the ANSR test kit. Ensure the permanent caps are placed onto the tubes to avoid contamination.

7. To begin running the test, press the green **Start test** arrow on the touch screen.

8. Once the test has been completed, the touch screen will display the results. To view data and a graph for a specific tube, touch the tube on the screen and then press the graph icon in the lower left corner. To return to the overall results, press the white back arrow in the lower left corner of the screen.

9. To print results, press the **Print** icon in the lower right corner of the touch screen. **NOTE:** Printing must be set up prior to using this function. For instructions, see Section 6.2.

## 6.7 Search


The **Search** function is marked by a magnifying glass icon on the **Home screen**.

### View all test results

This function shows the test results of the most recent test run. Users can use the arrows to see other tests in order of newest to oldest.

1. Touch the magnifying glass (**Search**) icon on the **Home screen**.
2. Select **View all test results**.
3. Use the arrows to scroll through the tests and select the desired choice.

### Search filter

This option allows the user to search by test name, test type, test number and by sample ID.

1. Touch the magnifying glass (**Search**) icon on the **Home screen**.
2. Select **Search filter**.

3. Select the desired search criteria and enter the search information.

### Go to test number

This function allows the user to search by test number.

1. Touch the magnifying glass (**Search**) icon on the **Home screen**.
2. Select **Search by test number**.
3. Input the desired test number and press the green check mark to search.

## 7. Cleaning and maintenance

To clean spills or stains on the reader, dampen a lint-free wipe with isopropyl alcohol or 10% bleach solution. Gently wipe the down reader. Do not use any free liquid.

If liquid is spilled into the reader, contact Neogen Technical Services for assistance by calling 800/234-5333 or 517/372-9200.

## 8. Troubleshooting

### Q: Why won't the reader and the computer connect?

- A: Restart the reader while it is connected to the computer. Ensure the touch screen on the reader says "remote mode." Also, check to make sure the first three numbers of the version software on the computer match the first three numbers of the firmware version on the reader, which is located under **About this device**. See Section 6.3 for more information.

### Q: Why did the reader display an invalid reading?

- A: In the reader displays an invalid reading, re-run the sample. This may occur if the reader's lid is not properly closed during a test run or if the run isn't started quickly after samples have been added.

### Q: Why do error messages appear when I attempt to load test type files or update the reader firmware using a USB drive?

- A: Ensure the USB drive is formatted for a PC and not Mac.

For other questions, please contact Neogen Technical Services.

## 9. Limited warranty

Neogen Corporation warrants the ANSR product(s) to the original Customer for a period of one year after date of installation against defects in material and workmanship and defect arising from failure to conform to specifications applicable on the date of installation, and Neogen agrees to correct, either by repair or, at its election, by replacement, any such defect found on examination to have occurred, under normal use and service, during such one-year period provided Neogen Corporation is promptly notified in writing upon discovery of such defect. All of the product elements in the ANSR product and the total system are warranted to be new or equivalent to new for the full product warranty period of one year.

Neogen Corporation shall have no liability for any instrument or part failure which Neogen determines, at its sole discretion, resulted from accident, improper installation or site preparation, unauthorized modification, customer-supplied software or interfacing, inadequate maintenance, misuse, failure to maintain proper environment, misapplication of the product or failure to follow instructions set out in ANSR product literature.

## 10. Terms and conditions

For Neogen's full terms and conditions, please visit [www.neogen.com/Corporate/termsconditions.html](http://www.neogen.com/Corporate/termsconditions.html).

## 11. ANSR reader specifications

<b>Number of tubes</b>	Configured for up to 16 tubes
<b>Assays</b>	Neogen ANSR assays
<b>Measurement technology</b>	Fluorescent 3 channel: FAM, HEX and ROX
<b>Color touch screen (3.5")</b>	Simple workflow management Ease of use: Color icon driver menus and prompts
<b>Communications</b>	Cable Ethernet for data transfer USB port for data export, barcode reader and printer Mini USB port for laptop connection
<b>Data storage</b>	Up to 99 results (configurable)
<b>Power</b>	12 V DC from external AC/DC supplied plug pack Battery power (optional) DC voltage fluctuation $\pm$ 10% DC current consumption: 12 V DC, 3.0 amps
<b>Dimensions</b>	380 mm (L) x 122 mm (W) x 139 mm (H) 14.0" (L) x 4.8" (W) x 5.4" (H)
<b>Weight</b>	1.95 kg 4.1 lbs
<b>Printers</b>	Network printer
<b>Operating environment</b>	Indoor use 10–30°C 10–70% relative humidity (noncondensing) 0–2000 m altitude Pressure: 0.763 bar (2000 m)–1.01 bar (sea level) Pollution category 2 Minimum light conditions of 100 LUX Maximum light intensity of 5000 LUX
<b>Storage environment</b>	0–45°C 20–70% relative humidity (noncondensing) Pressure: 0.763 (2000 m)–1.01 bar(sea level) 0–2000 m altitude Maximum duration: 1 month (cumulative)

## Neogen company profile

Founded in 1982, Neogen Corporation (NASDAQ: NEOG) has grown to have a worldwide presence, with more than 1,600 employees in multiple U.S. and international locations. Neogen remains firmly committed to the quality standards customers expect and to consistently meeting those standards.

Neogen's Food Safety Division develops and markets one-stop solutions to simply, rapidly and accurately detect bacterial pathogens, mycotoxins, food allergens, genetic modifications, drug and pesticide residues, ruminant by-products and sanitation concerns in food and animal feed.

Neogen's testing products are simple and accurate, and supported by a highly-skilled staff. In addition to the ANSR system, Neogen also offers additional microbiological solutions:

- **NeoSeek™ food safety testing**
- **GeneQuence® DNA hybridization assays**
- **Lateral flow immunochromatonic assays**
- **Hydrophobic grid membrane filtration**
- **Neogen Culture Media**
- **Soleris® rapid microbiological testing system**

## Contact us

620 Leshler Place • Lansing, MI 48912

**800/234-5333** (USA/Canada) or 517/372-9200

E-mail: [foodsafety@neogen.com](mailto:foodsafety@neogen.com) • [foodsafety.neogen.com](http://foodsafety.neogen.com)

## Neogen ANSR Molecular Diagnostic for Foodborne Pathogen Detection – limited use label license

### SYTO® 82

SYTO® 82 contained within this product is provided under an intellectual property license from Life Technologies Corporation. The transfer of this product is contingent on the buyer using the purchased product for *in vitro* detection and analysis of (i) food, feeds and beverages, including nutraceuticals, (ii) ingredients for food, feeds and beverages, (iii) process samples from food, feed and beverage preparation, distribution and delivery, and (iv) water from any source for human consumption, all for the purpose of safety and quality assurance. The buyer must not sell or otherwise transfer this product or its components for any other use, including but not limited to, human *in vitro* diagnostics, veterinary diagnostics, human identity or paternity testing, human forensics, or *in vivo* detection of nucleic acid sequences in living persons, animals or cells. For information on purchasing a license to SYTO® 82 for purposes other than food, beverage and water safety and quality assurance, contact Life Technologies Corporation, Cell Analysis Business Unit, Business Development, 29851 Willow Creek Road, Eugene, OR 97402, Tel: (541) 465-8300. Fax: (541) 335-0354.

### Molecular beacon probes

One or more molecular beacon probes contained within this product is sold under license from PHRI Properties and may be used under PHRI Properties patent rights only for tests on food products, feeds, beverages, and water.

### NEAR™ Technology

This product utilizes the patent pending NEAR isothermal technology and is sold under license from Ionian Technologies, San Diego, CA, and may be used under Ionian Technologies patent rights only for tests on food, beverage, and water safety.

### Axxin Pty. Ltd.

© 2012 Axxin Pty. Ltd. All Rights Reserved.

This software is provided under license and may be used and/or copied only in accordance with the terms of such license. This document is provided as an operational summary to describe the use of the instrument. This document does not describe operation of the specific diagnostic test or provide any input to the requirements, safety or processing of the diagnostic test which is outside of the scope of Axxin Pty Ltd.

Axxin Pty. Ltd. and the Axxin Pty. Ltd. Logo are trademarks of Axxin Pty. Ltd. in Australia and/or other countries.

© Neogen Corporation, 2017. Neogen, ANSR, GeneQuence, Acumedia and Soleris are registered trademarks, and NeoSeek is a trademark of Neogen Corporation. All other brand and product names are trademarks or registered trademarks of their respective companies.

Lined writing area with a faint green molecular structure background.

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**Neogen Food and Animal Security**

+91 484 2306598, 2301582  
info@neogenindia.com  
www.neogenindia.com

**Natural Toxins • Foodborne Bacteria • Food Allergen and Adulteration • Sanitation  
Drug Residues • Additional Rapid Tests • Dehydrated Culture Media**

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