**U.S. English**

SKU #: 700002564

REF #: 8420

 **Veratox**<sup>®</sup>  
*for Hazelnut Allergen*  
*Quantitative Test*

Refrigerate at 2–8°C (35–46°F). Do not freeze.

# Veratox<sup>®</sup> for Hazelnut Allergen

SKU #: 700002564 | REF #: 8420

## Hazelnut Allergen

Food allergens are proteins in food that can create an immune response in sensitive individuals. Once ingested, food allergens can cause a number of reactions, ranging in severity from hives and itching to anaphylaxis. Anaphylaxis is a severe allergic reaction, involving vomiting, diarrhea, difficulty breathing, swelling of the mouth and tongue, and a rapid drop in blood pressure.

An estimated 3.5% to 4% of adults, and 6% to 8% of children are sensitive in some degree to food allergens. More than 12 million people in the United States alone are known to have a food allergy.

Food manufacturers protect those with food allergies by clearly labeling their products with a list of ingredients. Testing for the presence of hazelnut components can play an important role in allergen control plans and can help detect presence of unintended allergens in food products.

## Intended Use

Veratox for Hazelnut Allergen is intended for the quantitative analysis of hazelnut protein residues in food products, ingredients, and clean-in-place (CIP) rinses. Qualitative analysis of environmental swabs is also available.

The Neogen<sup>®</sup> Veratox Kit for Hazelnut Allergen is intended for use in the food and beverage industry by trained personnel. Neogen has not documented the use of this product in industries other than food or beverage. For example, Neogen has not documented this kit for testing pharmaceutical, cosmetic, clinical, or veterinary samples. The Neogen Veratox Kit for Hazelnut Allergen has not been evaluated with all possible food products, food processes, and testing protocols.

## Intended User

This test kit is designed for use by quality control personnel and others familiar with foods possibly contaminated by hazelnuts or hazelnut products. Since technique is very important, operators should be trained by a Neogen representative or someone who has completed the Neogen training.

## Assay Principles

The Veratox for Hazelnut Allergen test is a sandwich enzyme-linked immunosorbent assay (S-ELISA). Hazelnut protein is extracted from samples with a phosphate buffered salt (PBS) solution by shaking in a heated water bath, followed by centrifugation or filtration. Extracted hazelnut protein is sampled and added to antibody-coated wells (capture antibody) where it binds to the antibody during an incubation. Any unbound hazelnut protein is washed away and a second antibody (detector antibody), which is enzyme labeled, is added. The detector antibody binds to the already bound hazelnut protein. After a second wash, substrate is added. Color develops as a result of the presence of bound detector antibody. Red Stop reagent is added and the color of the resulting solution is observed. The test is read in a microwell reader to yield optical densities. The optical densities of the controls form a standard curve, and the sample optical densities are plotted against the curve to calculate the concentration of hazelnut expressed as parts per million (ppm) of hazelnut.

## Storage Requirements

The kit can be used until the expiration date on the label when stored refrigerated at 2–8°C (35–46°F). Do not freeze.

## Materials Provided

1. 48 antibody-coated microwells
2. 48 red-marked transfer wells
3. 5 yellow-labeled bottles of 0, 2.5, 5, 10, and 25 ppm hazelnut controls
4. 2 blue-labeled bottles of enzyme-labeled antibody conjugate
5. 1 green-labeled bottle of K-Blue<sup>®</sup> Substrate
6. 1 red-labeled bottle of Red Stop solution

- 5 foil pouches of 10 mM PBS dry powder extraction solvent. Each pouch is enough to prepare 1 L in distilled or deionized water (pH 7.4)
- 40 mL of 10 mM PBS-Tween washing reagent in a wide mouth bottle. Each bottle is enough to prepare 1 L in distilled or deionized water (pH 7.4)
- 50 g of extraction additive in a specimen cup
- Plastic scoop to measure extraction additive

### **Materials Required, Not Provided**

- Extraction kit (Neogen item 8429)
  - 20 disposable plastic extraction bottles
  - 20 sample collection tubes (12 x 75 mm) with caps
- Allergen environmental swabbing kit (Neogen item 8432S)
  - 100 sterile swabs
  - 100 dropper tips
- Shaker water bath capable of maintaining  $60 \pm 1^\circ\text{C}$  with holders for 250 mL disposable plastic bottles (Neogen items 9298, 9299)
- Whatman No. 4 filters or equivalent (Neogen item 9429)
- Centrifuge (optional)
- Pipettor, 50–200  $\mu\text{L}$ , adjustable (Neogen item 9276)
- Pipettor, 12-channel (Neogen item 9273)
- Pipette tips (Neogen item 9410, 9407, 9417)
- Timer (Neogen items 9426, 9452)
- Microwell reader with a 650 nm filter (Neogen item 9303)
- 1 L bottle to prepare washing solution (Neogen item 9472)
- 1 L heat safe bottle to prepare extract solution (Neogen item 9472)
- Paper towels or equivalent absorbent material
- Microwell holder (Neogen item 9402)
- Waterproof marker
- Wash bottle (Neogen item 9400)
- Distilled or deionized water
- 3 reagent boats for 12-channel pipettor (Neogen item 9435)
- Graduated cylinder capable of measuring 125 mL (Neogen item 9368)
- Scale capable of weighing  $5 \pm 0.1$  g (Neogen item 9427)

### **Precautions**

- Components of Veratox for Hazelnut Allergen, such as controls and extraction reagents, may contain one or more of the following potentially allergenic materials: casein, egg protein, peanut protein, soy protein, or tree nut protein. If allergic to any of these compounds, please use caution when handling this product.
- Concentrated food additives, colors, and flavors may cause interferences on ELISA test methods. Contact Neogen's Technical Services for validation information.
- Hydrolyzed and fermented proteins may not be detected using ELISA methods for allergen testing. Due to the breakdown of proteins to small peptides or amino acids, they may become undetectable by this assay, but still could be allergenic and cause an allergic reaction.
- Sponges should not be used to collect samples for allergen testing. Sample collection swabs other than Neogen swabs should be validated prior to use. General sponges and swabs may contain solutions or materials that may interfere with the test kit.
- Store test kit between  $2\text{--}8^\circ\text{C}$  ( $35\text{--}46^\circ\text{F}$ ) when not in use. Do not freeze test kits.
- Do not use kit components beyond expiration date.
- Do not mix reagents from one kit serial with reagents from a different kit serial.
- Do not run more than 24 wells per test.

9. Follow proper pipetting techniques (e.g., prime tips and use clean tips).
10. Use only incubation times specified. Others may give inaccurate results.
11. Bring kits to room temperature 18–30°C (64–86°F) prior to use.
12. Avoid prolonged storage of kits at ambient temperatures.
13. Use clean pipette tips and glassware for each sample to avoid cross-contamination. Thoroughly wash all glassware between samples.

## User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at [neogen.com](http://neogen.com) or contact your local Neogen representative or distributor for more information. As with all test methods used for food analysis, the test matrix can influence the results. When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results. The food sample itself may influence results. It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples to satisfy the user that the chosen test method meets the user's criteria. It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements. As with any test method, results obtained from use of any Neogen product do not constitute a guarantee of the quality of the matrices or processes tested.

## Procedural Notes

1. **Substrate**  
Substrate is ready for use. The substrate should be clear to light blue — discard if it has turned dark blue. Only pour the needed volume of substrate into a reagent boat. Do not return unused substrate to the bottle. Cover the reagent boat to keep the substrate protected from light until needed.
2. **Conjugate**  
The conjugate supplied with this kit is ready to use. One bottle is enough for 24 wells. Cover the reagent boat to keep the conjugate protected from direct light and contaminants.
3. **Antibody Wells**  
Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after samples are extracted, and the test procedure is set to begin.
4. **Extraction Solution**  
Prepare extraction solution by adding a foil pouch of extraction solvent, 10 mM PBS, to 1 L distilled or deionized water. Swirl to mix thoroughly. Cover and store any unused portions refrigerated at 2–8°C (35–46°F).
5. **Wash Buffer**  
Prepare the wash buffer solution by pouring all the wash buffer concentrate into an empty 1 L container. Rinse the wash buffer concentrate bottle with distilled or deionized water and pour into the 1 L container to ensure all the concentrate is used. Fill the 1 L container with additional distilled or deionized water, and swirl to assure thorough mixing. Cover and store any unused portions refrigerated at 2–8°C (35–46°F).  
**Note:** Discard unused portions of extraction solution and wash buffer when the test kit has been used completely.

## Sample Preparation and Extraction

The sample to be tested should be collected according to accepted sampling techniques (see Neogen's Food Allergen Handbook). The sample should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

1. Prepare the extraction solution as described in the procedural notes.
2. Preheat extraction solution to 60°C (140°F) by immersing the bottle containing the solution into the water bath and allowing it to reach 60°C.
3. Using your sampling and collection procedure, obtain a representative sample and grind it to a very fine particle size.
4. Transfer 5 g of sample or 5 mL of liquid sample into a 250 mL disposable plastic bottle.
5. Add one level scoop of the extraction additive to the sample bottle.
6. Pour 125 mL of the 60°C (140°F) extraction solution to the sample plastic bottle.
7. Cap the bottle to prevent contents from splashing during the extraction.
8. Extract by shaking (150 rpm) in a water bath at 60°C (140°F) for 15 minutes. Remove the bottle from the bath.

9. If a solid food product was extracted, allow 5 minutes for the sample to settle prior to filtration in step 10.
10. Filter the extract by pouring at least 5 mL through a Whatman No. 4 filter and collecting the filtrate as a sample.  
**Alternative:** Centrifuge at 14,000 rpm for 5 minutes (20 minutes for lower speeds). Use the clear supernatant as a sample.
11. Allow extracts to cool to room temperature before beginning analysis.
12. Discard extracts after completion of analysis.

### Environmental Swabbing Preparation and Extraction

1. Prepare the extraction solution as described in the procedural notes.
2. Prepare a working solution by combining 125 mL of the extraction solution to a separate container, and add one level scoop of extraction additive to a container. Prepare a fresh working solution daily.  
**Note:** Prepare 125 mL of working solution for every 25 environmental swabs to be tested.
3. Preheat the working solution to 60°C (140°F). Shake container to dissolve additive.
4. Gather the sample with a swab, using one of the following methods:  
**For dry surfaces:** Open a new swab and wet with extraction solution. Swab a 10 x 10 cm area by using a crosshatch technique. Do not use working solution to moisten swabs.  
**For wet surfaces:** Open a new swab and swab a 10 x 10 cm area by using a crosshatch technique. Do not moisten swab prior to use.
5. Return the swab to its original tube once sampling is complete. Remember to label each tube.
6. Remove the swab from its tube, and add 5 mL of working solution at 60°C to the tube. Mix by placing the swab back into the tube and shaking for 2 minutes by hand (inverting tube), or for 30 seconds with a vortex mixer.
7. Remove the swab from its tube.
8. Place a new sample dropper tip onto the tube. The solution in the tube now serves as the sample.  
**Note:** Use caution when inverting the tube, as some liquid may drip.

### Test Procedure for Quantitation

Allow the test kit and all reagents to warm to room temperature 18–30°C (64–86°F) before using.

1. Remove 1 red-marked mixing well for each sample to be tested plus 5 red-marked wells for controls, and place in the well holder.
2. Remove an equal number of antibody-coated wells. Return antibody wells which will not be used immediately to the foil pack with desiccant. Reseal the foil pack to protect the antibody. Mark one end of the strip with a "1", and place strip in the well holder with the marked end on the left.
3. Mix each reagent by swirling the reagent bottle prior to use.
4. Using a new pipette tip for each, transfer 150 µL of controls and sample extracts to the redmarked transfer wells as shown in the template below. Only run up to two 12-well strips at a time.

|    |     |     |     |     |     |     |     |     |     |     |     |         |
|----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|---------|
| 0  | 2.5 | 5   | 10  | 25  | S1  | S2  | S3  | S4  | S5  | S6  | S7  | Strip 1 |
| S8 | S9  | S10 | S11 | S12 | S13 | S14 | S15 | S16 | S17 | S18 | S19 | Strip 2 |

5. Place tips on the 12-channel pipettor and transfer 100 µL of the controls and sample extracts to the antibody-coated wells. Mix for 20 seconds by sliding the well holder back and forth on a flat surface.
6. Incubate microwells 10 minutes at room temperature 18–30°C (64–86°F). Discard the redmarked transfer wells.
7. Empty the contents of the wells into a sink. With a wash bottle fill each antibody well with the wash buffer solution and dump out. Repeat the washing 5 times, then turn the wells upside down and tap out on a paper towel until the remaining washing solution is removed.
8. Pour the needed volume of conjugate from the blue-labeled bottle into a clean reagent boat.
9. Using the 12-channel pipettor and new tips, transfer 100 µL of the conjugate into all the wells and mix for 20 seconds by sliding the well holder back and forth on a flat surface.
10. Incubate for 10 minutes at room temperature 18–30°C (64–86°F).
11. Wash all wells with the wash buffer solution as described in step 7.

12. Pour the needed volume of substrate solution from the green-labeled bottle into a clean reagent boat.
13. Place new tips on the 12-channel pipettor and transfer 100  $\mu$ L of substrate into each well and mix for 20 seconds. Do not eject tips.
14. Incubate for 10 minutes at room temperature 18–30°C (64–86°F).
15. Pour the needed volume of Red Stop solution from the red-labeled bottle into a clean reagent boat.
16. With the same tips used to dispense the substrate, transfer 100  $\mu$ L of Red Stop into each well and mix for 20 seconds.
17. Wipe the bottom of the microwells and read in a microwell reader with a 650 nm filter. Results should be read within 20 minutes after the addition of Red Stop solution.
18. Interpret the test's results using Neogen's Stat-Fax microwell reader, or an equivalent strip reader. If using a strip reader, calculate the results using Neogen's Veratox software.

## Test Procedure for Screening

Allow the test kit and all reagents to warm to room temperature 18–30°C (64–86°F) before using.

1. Remove 1 well for each sample to be tested plus 1 well for the control, and place into the well holder.
2. Choose one yellow-labeled control bottle to serve as your screening level for the test. 5 ppm is recommended.
3. Mix each reagent by swirling its bottle prior to use.
4. Add 100  $\mu$ L from the yellow-labeled control bottle to the first well. Add 100  $\mu$ L from each sample extract to a respective well as indicated in the template below. For environmental swabs, add 3 drops from swab tube with dropper tip. Mix for 20 seconds by sliding the microwell holder back and forth on a flat surface.

|         |    |    |    |    |    |
|---------|----|----|----|----|----|
| Control | S1 | S2 | S3 | S4 | S5 |
|---------|----|----|----|----|----|

5. Incubate microwells 10 minutes at room temperature 18–30°C (64–86°F).
6. Shake out the contents of the wells. Using a wash bottle filled with wash buffer solution, fill each well and shake out. Repeat 5 times. Remove excess wash buffer by turning wells upside down and vigorously tapping wells on absorbent towel.
7. Add 100  $\mu$ L from the blue-labeled conjugate bottle to each well. Mix for 20 seconds by sliding the microwell holder back and forth on a flat surface.
8. Incubate for 10 minutes at room temperature 18–30°C (64–86°F).
9. Shake out the contents of the wells. Using a wash bottle filled with wash buffer solution, fill each well and shake out. Repeat 5 times. Remove excess wash buffer by turning wells upside down and vigorously tapping wells on absorbent towel.
10. Add 100  $\mu$ L from the green-labeled substrate bottle to each well. Mix for 20 seconds by sliding the microwell holder back and forth on a flat surface.
11. Incubate for 10 minutes at room temperature 18–30°C (64–86°F).
12. Add 100  $\mu$ L from the red-labeled Red Stop bottle to each well. Mix for 20 seconds by sliding the microwell holder back and forth on a flat surface. The results are now ready to be interpreted.
13. Visually compare the color of a sample well to the color of the control well. If the sample well has more blue color than the control well, the sample tests positive for hazelnut contamination of more than the control used. If the sample well has less blue color, or more red color, than the control well, the sample contains less than the control used of hazelnut contamination.

**Alternative:** Read wells (wipe bottom of wells first) in a microwell reader with a 650 nm filter. If the sample well has an optical density (OD) higher than the control well, the sample is positive for hazelnut contamination of more than the control used. If the sample well has an OD lower than the control well, the sample contains less than the control used of hazelnut contamination.

## Performance Characteristics

Limit of quantitation: 2.5 ppm total hazelnut (described as the lowest concentration point on the calibration curve that this test can reliably detect hazelnut allergen).

Range of quantitation: 2.5–25 ppm total hazelnut (for quantitating samples above 25 ppm, contact a Neogen representative for dilution instructions).

Allergen detection: This test detects hazelnut protein, and the results are expressed as ppm of total hazelnut.

Cross reactivity: Cross-reacts with pecan and walnut. Please see Veratox for Hazelnut validation report for additional information.

Protein conversion: Test kit results are expressed as total hazelnut. To express results as protein, multiply the hazelnut result by 0.15 (e.g., 2.5 ppm total hazelnut x 0.15 = 0.375 ppm hazelnut protein). \*On average, hazelnut contains 15% protein.

\*Source: USDA National Nutrient Database for Standard Reference Release 28, #12120 — Nuts, hazelnut, or filberts.

## Customer Service

Neogen Customer and Technical Services can be contacted through [neogen.com](http://neogen.com) and product training is available by request.

## Safety Data Sheets (SDS) Information Available

SDS are available for all test kits at [neogen.com](http://neogen.com) or by calling 800.234.5333 or 517.372.9200.

## Terms and Conditions

Neogen's full terms and conditions are available online.

## Warranty

Neogen makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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