Read instructions carefully before starting test

ALERT[®]

for Total Milk Allergen

Screening Test

REFRIGERATE AT 2–8°C (35–46°F) • DO NOT FREEZE

MILK ALLERGEN

Food allergens are proteins in food that can create an immune response in sensitive individuals. Once ingested, food allergens can cause a number of reactions, ranging in severity from hives and itching to anaphylaxis. Anaphylaxis is a severe allergic reaction, involving vomiting, diarrhea, difficulty breathing, swelling of the mouth and tongue, and a rapid drop in blood pressure.

An estimated 3.5 to 4 percent of adults, and 6 to 8 percent of children, are sensitive in some degree to food allergens. More than 12 million people in the United States alone are known to have a food allergy.

Food manufacturers protect those with food allergies by clearly labeling their products with a list of ingredients. Testing for the presence of milk residue ensures food manufacturers that an unlabeled — and potentially dangerous — ingredient did not make its way into a food product.

INTENDED USE

Alert for Total Milk Allergen is intended to be used to screen food, such as juices, sauces and sorbets, or environmental surfaces for the presence of the milk components, casein and whey.

INTENDED USER

The test kit is designed for use by quality control personnel and others familiar with foods possibly contaminated by casein or whey. Since technique is very important, operators should be trained by a Neogen representative or someone who has successfully completed the Neogen training.

ASSAY PRINCIPLES

Alert for Total Milk Allergen is a sandwich enzyme-linked immunosorbent assay (S-ELISA). Milk protein residue is extracted from samples by blending with a hot buffered salt solution (PBS). Extracted milk protein is sampled and added to antibody-coated wells (capture antibody) where it binds to the antibody during an incubation. Any unbound milk protein is washed away and a second antibody (detector antibody), which is enzyme labeled, is added. The detector antibody binds to the already bound milk protein. After a second wash, substrate is added. Color develops as a result of the presence of bound detector antibody. Red Stop reagent is added and the color of the resulting solution is observed.

STORAGE REQUIREMENTS

The kit can be used until the expiration date on the label when stored refrigerated at 2–8°C (35–46°F). Do not freeze.

MATERIALS PROVIDED

1. 24 antibody-coated microwells
2. 2 yellow-labeled bottles of controls made from nonfat dried milk; 1 bottle at a concentration of 5 ppm and 1 bottle at 10 ppm (see procedural note No. 1)
3. 1 blue-labeled bottle of enzyme-labeled antibody conjugate
4. 1 green-labeled dropper bottle of K-Blue® Substrate solution
5. 1 red-labeled dropper bottle of Red Stop solution
6. 20 unlabeled empty sample dropper bottles with tips
7. 3 foil pouches of 10 mM PBS dry powder extraction solvent; each pouch contains enough powder to prepare 1 L of extraction solution
8. 1 wide-mouth bottle of 40 mL PBS-Tween wash buffer concentrate; a bottle contains enough concentrate to prepare 1 L of wash buffer
9. 50 g of extraction additive in a specimen cup
10. Plastic scoop to measure extraction additive

MATERIALS RECOMMENDED BUT NOT PROVIDED

1. Allergen Extraction Kit (Neogen item 8429)
 - a. 20 disposable extraction bottles
 - b. 20 disposable transfer pipettes
2. Allergen Environmental Swabbing Kit (Neogen item 8432S)
 - a. 100 sterile swabs
 - b. 100 dropper tips
3. Scale capable of weighing 5 ± 0.1 g (Neogen item 9427)
4. Shaker water bath, hot plate or equivalent heat source capable of maintaining $60 \pm 1^\circ\text{C}$
5. Thermometer
6. Two 1 L bottles to prepare washing and extraction solutions (Neogen item 9472)
7. Paper towels or equivalent absorbent material
8. Microwell holder (Neogen item 9402)
9. Waterproof marker
10. Wash bottle (Neogen item 9400)
11. Distilled or deionized water
12. Timer (Neogen item 9426)
13. Graduated cylinder capable of measuring 125 mL (Neogen item 9368)
14. Pipettor, 100 μL (Neogen item 9272, 9278, 9276)
15. Pipette tips (Neogen item 9410, 9407, 9417)
16. Vortex (Neogen item 9494)
17. 50 mL sample tubes (Neogen item 9381)

PRECAUTIONS

1. Samples intended to be tested for **milk must be extracted separately** from samples intended to be tested for other food allergens. The extraction additives for each type of test are designed specifically for the target food allergen.
2. Components of Alert for Total Milk Allergen, such as controls and extraction reagents, may contain one or more of the following potentially allergic materials: milk, egg protein, peanut protein or soy protein. If allergic to any of these compounds, please use caution when using this product.
3. Concentrated food additives, colors and flavors may cause interferences on ELISA test methods. Contact Neogen's technical services for validation information.
4. Hydrolyzed and fermented proteins may not be detected using ELISA methods for allergen testing. Due to the nature of the proteins they may be undetectable in the assay, but there could still be active allergenic protein residue present.
5. Infant formula: Contact Neogen for additional information concerning the testing of infant formulas.
6. Sponges should not be used for sample collection and allergen testing. Sample collection swabs other than Neogen swabs should be validated prior to use. General sponges and swabs may contain solutions or materials that may interfere with the test kit.
7. Disposable extraction bottles and transfer pipettes must be used to extract milk samples to avoid cross-contamination.
8. The testing environment should be clean and dust-free.
9. Store test kit between 2–8°C (35–46°F) when not in use. Do not freeze.
10. Do not use kit components beyond expiration date.
11. Do not mix reagents from kit serial with those from a different serial.
12. Do not run more than 6 wells per test.
13. Use only the incubation times specified. Others may give inaccurate results.
14. Bring kits to room temperature 18–30°C (64–86°F) prior to use.
15. Avoid prolonged storage of kits at ambient temperatures.

PROCEDURAL NOTES

1. **Controls.** To increase the flexibility of this test kit, two test controls are supplied: 5 ppm and 10 ppm. **Choose one control** level to screen samples against for each test performed.
2. **Extraction solution.** Prepare extraction solution by adding a foil pouch of extraction solvent, 10 mM PBS, to 1 L distilled or deionized water at pH 7.4. Cover and store unused portion refrigerated at 2–8°C (35–46°F).
3. **Wash buffer.** Prepare the wash buffer solution by pouring all the wash buffer concentrate into an empty 1 L bottle. Rinse the wash buffer concentrate container with distilled or deionized water and pour into the 1 L bottle to ensure all the concentrate is used. Fill the 1 L bottle with additional distilled or deionized water, and swirl to assure thorough mixing. Cover and store any unused portions refrigerated at 2–8°C (35–46°F).
NOTE: Discard unused portions of extraction solution and wash buffer when the test kit has been used completely.
4. **Antibody wells.** Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after samples are extracted, and the test procedure is set to begin.

SAMPLE PREPARATION AND EXTRACTION

The sample to be tested should be collected according to accepted sampling techniques (see Neogen's Food Allergen Monitoring Handbook).

NOTE: Glassware and material used for other allergen testing cannot be used for milk residue testing due to the potential of cross-contamination. For this reason, it is highly recommended that any labware, such as wash bottles, graduated cylinders or 1 L bottles, are dedicated solely for use with the milk allergen kit.

1. Prepare the extraction solution as described in the procedural notes.
2. Preheat extraction solution to 60°C (140°F).

3. Using your sampling and collection procedure, obtain a representative sample and grind it to a very fine particle size.
4. Transfer 5 g of sample, or 5 mL of liquid sample, into a 250 mL disposable plastic extraction bottle.
5. Add one level scoop of the milk extraction additive into the sample bottle. Do not use the extraction additive from another allergen test kit.
6. Pour 125 mL of the 60°C (140°F) extraction solution into the sample bottle.
7. Cap the sample bottle to prevent contents from splashing during the extraction.
8. Extract by shaking (150 rpm) in a 60°C water bath for 15 minutes.
9. Remove the bottle from the bath and let the material settle for 5 minutes.
10. Using a disposable pipette, transfer approximately 1 mL of the supernatant (the top liquid portion of the extract) to a sample dropper bottle. Then, label and place a dropper tip on the bottle.
11. Begin the test procedure once the sample has cooled to room temperature (at least 15 minutes).

ALTERNATE EXTRACTION

1. Prepare the extraction solution as described in the procedural notes.
2. Prepare a working solution by adding 125 mL of extraction solution into container and adding one scoop of extraction additive and mix.
3. Preheat the working solution to 60°C (140°F) by using a hot plate, incubator, microwave, water bath or equivalent. Use a thermometer to assure temperature is correct.
4. Transfer 1 g of sample, or 1 mL of liquid sample, into a 50 mL tube.
5. Pour 25 mL of the 60°C (140°F) working extraction solution into the sample tube and cap.
6. Extract by vortexing for 3 minutes. Use caution, the liquid is hot.
7. Let the material settle for 5 minutes.
8. Using a disposable pipette, transfer approximately 1 mL of the supernatant (the top liquid portion of the extract) to a sample dropper bottle. Then label and place a dropper tip on the bottle.
9. Begin the test procedure once the sample has cooled to room temperature (at least 15 minutes).

ENVIRONMENTAL SWABBING PREPARATION AND EXTRACTION

1. Prepare the extraction solution as described in the procedural notes.
2. Prepare a working solution by combining 125 mL of the extraction solution to a separate container, and add one level scoop of extraction additive to a container. Prepare a fresh working solution daily.
NOTE: Prepare 125 mL of working solution for every 25 environmental swabs to be tested.
3. Preheat the working solution to 60°C (140°F). Shake container to dissolve additive.
4. Gather the sample with a swab, using one of the following methods:
For dry surfaces: Open a new swab and wet with extraction solution. Swab a 10 x 10 cm area by using a crosshatch technique. (Do not use working solution to moisten swabs.)
For wet surfaces: Open a new swab and swab a 10 x 10 cm area by using a crosshatch technique. Do not moisten swab prior to use.
5. Return the swab to its original tube once sampling is complete. Remember to label each tube.
6. Remove the swab from its tube, and add 5 mL of working solution at 60°C to the tube. Mix by placing the swab back into the tube and shaking for 2 minutes by hand (inverting tube), or for 30 seconds with a vortex mixer.
7. Remove the swab from its tube.
8. Place a new sample dropper tip onto the tube. The solution in the tube now serves as the sample.
NOTE: Use caution when inverting the tube, as some liquid may drip.

TEST PROCEDURE

Allow the test kit and all reagents to warm to room temperature (18–30°C) 64–86°F before using.

1. Remove 1 well from the foil pouch for each sample to be tested plus 1 well for the control, and place into the well holder.
2. Mix each reagent by swirling the reagent bottle prior to use.
3. Pipette 100 µL from a yellow-labeled control bottle (*choose the 5 ppm or 10 ppm control supplied—see procedural note No. 1*) to the first well. Add three drops, or pipette 100 µL from each sample extract to a respective well as indicated in the template below. Mix for 20 seconds by sliding the microwell holder back and forth on a flat surface.

Control	S1	S2	S3	S4	S5
---------	----	----	----	----	----

4. Incubate for **10 minutes** at room temperature 18–30°C (64–86°F).
5. Shake out the contents of the wells. Using a wash bottle filled with wash buffer, fill each well and shake out. **Repeat 10 times.** Remove excess wash buffer by turning wells upside down and vigorously tapping wells on an absorbent towel.
6. Pipette 100 µL from the blue-labeled conjugate bottle to each well. Mix for 20 seconds by sliding the microwell holder back and forth on a flat surface.
7. Incubate for **10 minutes** at room temperature 18–30°C (64–86°F).
8. Shake out the contents of the wells. Using a wash bottle filled with wash buffer, fill each well and shake out. **Repeat 10 times.** Remove excess wash buffer by turning wells upside down and vigorously tapping wells on an absorbent towel.
9. Add 3 drops from the green-labeled substrate dropper bottle to each well. Mix for 20 seconds by sliding the microwell holder back and forth on a flat surface.
10. Incubate for **10 minutes** at room temperature 18–30°C (64–86°F).
11. Add 3 drops from the red-labeled Red Stop dropper bottle to each well. Mix for 20 seconds by sliding the microwell holder back and forth on a flat surface. The results are now ready to be interpreted.

INTERPRETATION OF RESULTS

Visually compare the color of a sample well to the color of the control well. If the sample well has **more blue color** than the control well, the sample tests positive for milk protein contamination of **more than 5 ppm (or 10 ppm, depending on the control used)**. If the sample well has **less blue color, or more red color**, than the control well, the sample contains **less than 5 ppm (10 ppm)** of milk protein contamination.

Alternative: Read wells (wipe bottom of wells first) in a microwell reader with a 650 nm filter. If the sample well has an optical density (OD) higher than the control well, the sample is positive for milk contamination of more than 5 ppm (10 ppm). If the sample well has an OD lower than the control well, the sample contains less than 5 ppm (10 ppm) of milk protein contamination.

PERFORMANCE CHARACTERISTICS

Allergen detection: This test detects casein and whey proteins from cow, goat and sheep. The results are expressed in terms of detection of non-fat dried milk (NFDm).

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Service can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's Food Safety test kits, at www.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/Corporate/termsconditions.html.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TESTING KITS AVAILABLE FROM NEOGEN

Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, crustaceans, eggs, gliadin, hazelnut, lupine, milk, mustard, peanuts, sesame, soy, walnuts

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