- 
- 
-  (EN) Dilution Pouch
 -  (FR) Poche de dilution
 -  (DE) Verdünnungsbeutel
 -  (IT) Busta per diluizione
 -  (ES) Bolsa de Dilución
 -  (NL) Zak voor verdunning
 -  (SV) Dilution Pouch
 -  (DA) Fortyndingspose
 -  (NO) Fortynningspose
 -  (FI) Laimennuspussi
 -  (PT) Saco de Diluição
 -  (EL) Σάκος Αραίωσης
 -  (PL) Worek do rozcieńczeń
 -  (HU) Hígító tasak
 -  (CS) Sáček na ředění
 -  (RO) Săculeț pentru diluție
 -  (RU) Пакет для разведений
 -  (TR) Seyreltme Poşeti
 -  (JA) 希釈水 パウチ
 -  (ZH) 稀释袋
 -  (TH) ไโดลูชั่น เพาซ์
 -  (KO) 희석 파우치

Product Instructions

Dilution Pouch

PRODUCT DESCRIPTION AND INTENDED USE

The Neogen® Dilution Pouch is a convenient alternative to the traditional preparation and sterilization of diluent media. The Neogen Dilution Pouch prefilled with 99, 225 or 450 mL of specified dilution media consists of a standalone, heavyweight, puncture-resistant pouch which can be used for the dilution of product samples that require stomaching.

Table 1: Product Description and Storage Conditions

Catalog Number	Product Description	Volume	Storage
QDBFD99	Neogen® Dilution Pouch with Butterfield's Buffer	99 mL	20-25°C
QDBFD225	Neogen® Dilution Pouch with Butterfield's Buffer	225 mL	20-25°C
QDBF225F	Neogen® Dilution Pouch with Butterfield's Buffer & Filter	225 mL	20-25°C
QDBFD450	Neogen® Dilution Pouch with Butterfield's Buffer	450 mL	20-25°C

SAFETY

The user should read, understand, and follow all safety information in the instructions for the Neogen Dilution Pouches. Retain the safety instructions for future reference.

⚠ CAUTION Indicates a hazardous situation, which, if not avoided, could result in minor or moderate injury and/or property damage.

⚠ CAUTION

- The Neogen Dilution Pouch is intended for use in a laboratory environment by professionals trained in laboratory techniques.
- Neogen has not documented the use of this product in industries other than food. For example, Neogen has not documented this product for testing water, pharmaceutical, cosmetic, clinical or veterinary samples.
- The Neogen Dilution Pouches have not been evaluated with all possible food products, food processes and food processing environments, testing protocols or with all possible strains of bacteria.
- Do not use Neogen Dilution Pouches in the diagnosis of conditions in humans or animals.
- Neogen has not validated this product using composite samples.
- The user must train its personnel in proper testing techniques, for example: Good Laboratory Practices⁽¹⁾, ISO 7218⁽²⁾, or ISO 17025⁽³⁾.
- For information on documentation of product performance, visit our website at www.neogen.com or contact your local Neogen representative or distributor.
- Consult the Safety Data Sheet for additional information.

USER RESPONSIBILITY

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at www.neogen.com, or contact your local Neogen representative or distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results. The food sample itself may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any Neogen Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

LIMITATION OF WARRANTIES / LIMITED REMEDY

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, NEOGEN DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any Neogen Food Safety Product is defective, Neogen or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. Please contact your Neogen representative or authorized Neogen distributor for any further questions.

LIMITATION OF NEOGEN LIABILITY

NEOGEN WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall Neogen's liability under any legal theory exceed the purchase price of the product alleged to be defective.

STORAGE

Upon receipt, store unopened Neogen Dilution Pouches at the indicated temperatures (See Table 1).

⚠ DISPOSAL

After use, Neogen Dilution Pouches may contain microorganisms that could be a potential biohazard. Follow current industry standards or local regulations for disposal.

INSTRUCTIONS FOR USE

Follow all instructions carefully. Failure to do so may lead to inaccurate results.

Wear appropriate protective apparel and follow standard good laboratory practices GLP¹.

1. Stand the Neogen Dilution Pouch on flat surface.
2. Put on gloves.
3. Tear the Neogen Dilution Pouch open.
4. Pull red tabs, so the Neogen Dilution Pouch stays open.
5. Aseptically add sample to the Neogen Dilution Pouch.
6. The Neogen Dilution Pouch can be paddle blended utilizing a paddle style blender to homogenize the contents.
7. Fold down the Neogen Dilution Pouch to close.
8. Fold ends of blue wires inward.
9. Follow user established procedures.

REFERENCES

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
3. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.

Refer to the current version of the standard methods listed above.

EXPLANATION OF SYMBOLS

info.neogen.com/symbols

Neogen Food Safety

Neogen Corporation

620 Leshar Place
Lansing, MI 48912 USA
Neogen.com

Neogen Europe Ltd.

The Dairy School
Auchincruive
Ayr, KA6 5HU
Scotland, UK

Neogen Ireland, Ltd.

Bray Business Park, Bray
Co. Wicklow
A98YV29, Ireland

Neogen Corporation

620 Leshar Place Lansing, MI 48912 USA
www.neogen.com

Neogen is a trademark of Neogen Corporation.
© Neogen Corporation 2024.
All rights reserved.
FS00843B

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