Read instructions carefully before starting test

Reveal[®] MAX
for Aflatoxin
Qualitative Test
With Aqueous Extraction

Bulk Kit — 50 tests

THE TOXIN

Aflatoxin is a toxic and carcinogenic substance produced by certain strains of the molds *Aspergillus flavus* and *A. parasiticus*. There are four principle types of aflatoxin: B₁, B₂, G₁ and G₂. Aflatoxin B₁ is the most frequently encountered of the group and the most toxic. The commodities most affected by aflatoxin are corn, peanuts, cottonseed, milo and the majority of tree nuts.

The effects on animals ingesting excessive amounts of the toxin range from chronic health and performance problems to death. Aflatoxin has been shown to cause liver damage or cancer, decreased milk and egg production, immune suppression and interference with reproductive efficiency.

Many countries have set maximum allowable levels of aflatoxin in food and feed. Therefore, accurate determination of the toxin's presence is of major importance to those monitoring the quality of food and feed in which aflatoxin may occur. Testing these commodities for the toxin requires careful sampling, chemical extraction, sanitation and quantitative analysis.

Regulatory limits for aflatoxin are issued by regional authorities:

FOOD

Region	Commodity	Maximum Level		
		B1	Total	M1
US	All foods except milk		20 ppb	
	Milk			0.5 ppb

Region	Commodity	Maximum Level		
		B1	Total	M1
EU	Cereals, processed	2 ppb	4 ppb	
	Dried fruits to be processed	5 ppb	10 ppb	
	Groundnuts and nuts to be processed	8 ppb	15 ppb	
	Groundnuts, nuts and dried fruit for direct human consumption	2 ppb	4 ppb	
	Maize to be processed	5 ppb	10 ppb	
	Spices	5 ppb	10 ppb	
	Almonds, pistachios and apricots for direct human consumption	8 ppb	10 ppb	
	Almonds, pistachios and apricots to be processed	12 ppb	15 ppb	
	Dietary foods for special medical purposes intended specifically for infants	0.10 ppb		
	Dried fruit for direct human consumption	2 ppb	4 ppb	
	Hazelnuts and Brazil nuts for direct human consumption	5 ppb	10 ppb	
	Hazelnuts and Brazil nuts to be processed	8 ppb	15 ppb	
	Infant formula and follow-on formula, including infant milk and follow-on milk			0.025 ppb
	Processed cereal-based foods and baby foods for infants and young children	0.10 ppb		
	Raw milk, heat-treated milk and milk for the manufacture of milk-based products	2 ppb	4 ppb	0.05 ppb
	Tree nuts for direct human consumption	5 ppb	10 ppb	

FEED

Region	Commodity	Maximum Level	
		B1	Total
US	Corn and peanut products intended for finishing beef cattle		300 ppb
	Cottonseed meal intended for beef cattle, swine or poultry		300 ppb
	Corn and peanut products intended for finishing swine of 100 lbs or greater		200 ppb
	Corn and peanut products intended for breeding beef cattle, breeding swine or mature poultry		100 ppb
	Corn, peanut products and other animal feeds and feed ingredients, excluding cottonseed meal, intended for immature animals		20 ppb
	Corn, corn products, cottonseed meal and other animal feed and feed ingredients		20 ppb

Region	Commodity	Maximum Level	
		B1	Total
EU	All feed materials	20 ppb	
	Complementary and complete feed	10 ppb	
	Compound feed for dairy cattle and calves, dairy sheep and lambs, dairy goats and kids, and young poultry animals	5 ppb	
	Compound feed for cattle, sheep, goats, pigs and poultry not listed above	20 ppb	

Please refer to the Mycotoxin Handbook for additional regulatory levels.

INTENDED USE/USER

Reveal® MAX for Aflatoxin is intended for the qualitative analysis of corn, corn products and wheat for aflatoxin. The test uses an aqueous extraction, eliminating the need for the disposal of hazardous waste. This screening kit is designed for use by quality control personnel and others familiar with commodities possibly contaminated by aflatoxin.

ASSAY PRINCIPLES

Reveal MAX for Aflatoxin is a lateral flow immunochromatographic assay based on a competitive immunoassay format. In the test, the extract is wicked through the reagent zone, which contains antibodies specific for aflatoxin conjugated to colloidal gold particles. If aflatoxin is present, it will be captured by the particle-antibody complex. The aflatoxin antibody-particle complex is then wicked onto a membrane, which contains a zone of aflatoxin conjugated to a protein carrier. This zone captures any uncomplexed aflatoxin antibody, allowing the particles to concentrate and form a visible line. If aflatoxin is present above a pre-determined limit, no line will form. Therefore, a positive sample will result in the absence of a line in the test zone. Alternatively, a negative sample will form a visible line in the test zone. The membrane also contains a control zone where an immune complex present in the reagent zone is captured by an antibody, forming a visible line. This control line will always form regardless of the presence of aflatoxin, ensuring the strip is functioning properly.

STORAGE REQUIREMENTS

Store kit components at room temperature 18–30°C (64–86°F) to ensure full shelf life. Test strips should remain capped in their original tubes until used to ensure optimal performance.

MATERIALS PROVIDED

1. 50 Reveal MAX for Aflatoxin test strips
2. 50 red dilution cups
3. 2 bottles of sample diluent
4. 50 MAX 1 Aqueous Extraction packets
5. Kit instructions for use

MATERIALS RECOMMENDED BUT NOT PROVIDED

1. Extraction Materials
 - a. Sample collection cups with lids 125 mL (Neogen item 9428, 9428B)
 - b. Sample collection tubes with caps (Neogen item 9421, 9421B)
 - c. Filter syringes, Whatman #4 filter paper or equivalent (Neogen item 9420, 9429, 9519)
OR
 - d. Centrifuge, mini (Neogen item 9330)
 - e. Microcentrifuge tubes (Neogen item 9172)
 - f. Dispensing pump or graduated cylinder (Neogen item 9448, 9447)
2. Agri-Grind grinder or equivalent (Neogen item 9401, 9453)
3. Scale capable of weighing 5–50 grams \pm 0.1 g (Neogen item 9427)
4. Timer (Neogen item 9426, 9452)
5. Reveal sample cup rack (Neogen item 9475)
6. Pipettor, 100 μ L (Neogen item 9860, 9290)
7. Pipette tips, 1–200 μ L (Neogen item 9407, 9410, 9417)
8. Distilled or deionized water

PRECAUTIONS

1. The test strips must remain inside the stay-dry tube before use.
2. Store test kit at room temperature 18–30°C (64–86°F) when not in use. Do not freeze.
3. Do not use kit components beyond expiration date.
4. Treat all used liquids, including sample extract, and labware as if contaminated with aflatoxin.
5. Commodity extracts should have a pH of 6–8 before testing. Excessively acidic or alkaline samples should be adjusted. For instructions on adjusting pH, contact a Neogen representative or Technical Services.

SAMPLE PREPARATION

The sample to be tested should be collected according to accepted sampling techniques (see FGIS sampling protocol or contact your Neogen representative). Obtain a representative sample (minimum 100 gm). Grind the sample so at least 95% of the ground material passes through a 20 mesh sieve (about the particle size of fine espresso).

SAMPLE EXTRACTION

1. Weigh out 10 g \pm 0.1 g of sample into extraction cup.
2. Add contents of one (1) MAX 1 aqueous extraction packet to the extraction cup.
3. Add 50 mL distilled or deionized water to the extraction cup.
4. Vigorously shake, using hand or mechanical means, for 1 minute, or blend for 1 minute.
5. Allow the sample to settle, then filter with a filter syringe or Whatman #4 filter paper, or you may also pipette 1 mL of sample into a 1.5 mL micro-centrifuge tube, and centrifuge for 1 minute using a micro-centrifuge.

TEST PROCEDURE

1. Place the appropriate number of red dilution cups into a sample cup rack. Label cups if necessary.
2. Using a pipettor and a new tip, add 100 μ L of sample diluent to each red dilution cup.
3. Using a new tip, add 200 μ L of sample extract to the dilution cup. Mix by pipetting up and down 5 times.
4. Place a new Reveal MAX for Aflatoxin test strip with the sample end down into the sample cup and set timer for 3 minutes. Ensure the test strip comes into contact with liquid and begins to wick.
5. Remove the strip from the sample cup after it has developed for **3 minutes**.
6. Read visually by holding strip up to line of sight. Strip must be read within **30 seconds** after the 3 minute incubation.


INTERPRETATION OF RESULTS

Results less than 20 ppb

If a line forms in the test zone and another line forms in the control zone within 3 minutes (2 lines total), the sample contains less than 20 ppb aflatoxin (negative sample). A negative sample may be determined as soon as 2 lines are visible on the strip. If there is no line in the control zone, the test is invalid and the sample should be re-tested with another strip.

Results more than 20 ppb

If after the full 3 minutes there is no visible line in the test zone, but a visible line in the control zone, the sample contains greater than 20 ppb aflatoxin (positive sample). If there is no line in the control zone, the test is invalid and the sample should be re-tested with another strip.

NOTES

- Visible reading of the test strips should be made between 3–3.5 minutes. Observations after 3.5 minutes may be inaccurate due to over-development of the device.
- The strips can be read using a Neogen AccuScan® Gold Reader.
- Positive samples can be confirmed by Neogen's quantitative Veratox® for Aflatoxin or by HPLC.

VALIDATED COMMODITIES

Corn, wheat

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at foodsafety.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/en/terms-and-conditions.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TEST KITS AVAILABLE FROM NEOGEN

Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, crustaceans, eggs, gliadin, hazelnut, milk, mustard, peanuts, sesame, soy, walnuts, multi-treenut

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