(EN) Neogen[®] Microbial Luminescence System II

(FR) Neogen[®] Microbial Luminescence System II

Microbial Luminescence System II

Contents

1. Important Information	2	6. Cleaning the System for Initial Use	18
1.1. Product Description and Intended Use	2	7. Use of the Neogen MLS II Washing/Priming the Injectors	19
1.2. User Responsibility	2	7.1. Washing the Injectors	19
1.3. Safety Information	3	7.2. Priming the Injectors	20
1.4. Environmental Operating Conditions	5	7.3. Reagent Use Calculation for the Neogen® Microbial Luminescence System Ultra-High Temperature (UHT) Screen Kits Reagents	20
1.5. Regulatory Compliance	5	8. Ultra-High Temperature (UHT) Sterility Testing	22
1.6. Package Contents	5	9. Cleaning the Instrument After Testing	22
1.7. User-Supplied Items	6	10. Weekly Cleaning	23
1.8. Accessory Items and Kits	6	11. Routine and Periodic Maintenance Summary	24
2. Technical Data	6	12. Moving the Neogen MLS II Instrument to Another Workbench	24
3. Overview of Luminescence	7	13. Instructions for Shipping or Prolonged Storage of the Neogen MLS II Instrument	24
3.1. Bioluminescence Reaction	7	14. Instructions for Installing the Neogen MLS II Instrument After Service Repair	25
3.2. Bioluminescence and the Neogen® Microbial Luminescence System (MLS) II	7	15. Instructions for Cleaning External Surfaces of the Neogen MLS II Instrument	25
3.3. Summary of a Neogen Microbial Luminescence System (MLS) II Test Assay	7	16. Contact Information	25
4. Initial Set-up of the Neogen MLS II	8	Appendix I – Troubleshooting Guide:	26
4.1. Neogen Microbial Luminescence System (MLS) II Instrument Installation	8		
4.2. Neogen Microbial Luminescence System (MLS) II Software Installation	8		
4.2.1 Change Password	9		
4.2.2 Access Control	9		
5. Neogen MLS II Software Functions	10		
5.1. Log-on window	10		
5.2. Actual Tab	10		
5.2.1 Live Window	10		
5.2.2 Actual Tab with Opened Assay Protocol	11		
5.2.3 Actual Tab with Finished Assay	11		
5.3. Lot Specific Values For Tab	12		
5.4. Load Plate Tab	13		
5.5. Report Tab	14		
5.6. Protocol Editor Tab	15		
5.7. Tools Tab	15		
5.7.1 Plate Command	15		
5.7.2 Liquid Handling	16		
5.7.3 Plate Editor	17		
5.7.4 Change Password	17		
5.7.5 Access Control	17		
5.7.6 Preferences	17		

1. Important Information

Please read, understand, and follow all safety information contained in this User's Manual and in the associated Neogen® Microbial Luminescence System (MLS) Kits inserts prior to using the Neogen® Microbial Luminescence System (MLS) II Instrument and Neogen MLS Kits. Retain this Important Information for future reference.

1.1. Product Description and Intended Use

The Neogen® Microbial Luminescence System (MLS) II is a microplate luminometer designed for microbial ATP detection in end products. This unit has been documented for use in the food and beverage industry. Contact your official Neogen Food Safety Department representative for a current listing of Neogen® Microbial Luminescence System (MLS) Reagents intended for use with Neogen MLS II Instrument. Neogen has neither designed nor documented this instrument for use with other manufacturers' products. Use with other manufacturers' designs, or use in any other application that has not been evaluated or documented by Neogen may lead to an unsafe condition. It is expected that all users have been properly trained in the use of the Neogen MLS II Instrument and the Neogen MLS Kits. The Neogen MLS II has been designed and tested only to be used with the power cable that is supplied with the instrument. If this product is not used in a manner not specified, the protection provided by the product may be impaired.

For information on documentation of product performance, visit our website at neogen.com, or contact your local Neogen representative or distributor.

1.2 User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at neogen.com, or contact your local Neogen representative or distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory techniques may influence results.

It is the user's responsibility in selecting any test method to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet customers' and suppliers' requirements.

As with any test method, results obtained from use of this product do not constitute a guarantee of the quality of the matrices or processes tested.

Limited Warranty

Neogen Food Safety Hardware ("Hardware") includes the Neogen Microbial Luminescence System (MLS) II Instrument and any associated Software and other components supplied by Neogen which are described in the applicable Installation and Use Guides ("User Documentation").

Limited Warranty and Disclaimer

Neogen warrants that for one year from the date of shipment that the Hardware will substantially perform in accordance with the User Documentation. THIS WARRANTY IS VOID IF: (A) THE HARDWARE HAS BEEN REPAIRED BY PERSONS NOT AUTHORIZED BY Neogen; OR (B) THE HARDWARE HAS BEEN ALTERED, MODIFIED, OR MISUSED; OR (C) THE HARDWARE IS USED WITH PRODUCTS, SUPPLIES, COMPONENTS OR SOFTWARE NOT SUPPLIED BY Neogen FOR USE WITH THE HARDWARE; OR (D) THE HARDWARE OR A COMPONENT IS USED FOR OTHER USES (FOR EXAMPLE USE WITH OTHER CIRCUIT BOARDS OR SOFTWARE) OR (E) THE HARDWARE HAS NOT BEEN MAINTAINED OR USED IN ACCORDANCE WITH THE USER DOCUMENTATION. UNLESS PROHIBITED BY LAW, THIS WARRANTY IS MADE IN LIEU OF ALL OTHER WARRANTIES, EXPRESS OR IMPLIED, INCLUDING BUT NOT LIMITED TO THE IMPLIED WARRANTY OF FITNESS FOR A PARTICULAR PURPOSE, THE IMPLIED WARRANTY OF MERCHANTABILITY, OR ANY IMPLIED WARRANTY ARISING OUT OF A COURSE OF DEALING OF PERFORMANCE, CUSTOM OR USAGE OF TRADE. Neogen DOES NOT WARRANT THAT THE SOFTWARE WILL FUNCTION ERROR FREE.

If within one year after shipment, the Hardware does not conform to the express warranty set forth above, Neogen's sole obligation and User's sole remedy shall be, at Neogen's option: 1) to repair or replace the non-conforming component; or, 2) refund the purchase price.

Limitation of Neogen Liability

UNLESS PROHIBITED BY LAW, Neogen WILL NOT BE LIABLE TO USER OR OTHERS FOR ANY OTHER DIRECT, INDIRECT, CONSEQUENTIAL, INCIDENTAL OR SPECIAL DAMAGES INCLUDING, FOR EXAMPLE, LOST PROFITS, BUSINESS, INVESTMENTS, OR OPPORTUNITIES EVEN IF Neogen HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES. The parties agree that Neogen's total cumulative liability to User for direct damages for all causes shall not exceed One Hundred Dollars, or the price paid for the Hardware, whichever is higher. Some states or countries may have laws which require liability rights different from those stated above. In such state or countries the minimum required liability terms shall apply.

To request service, you must request a Returned Material Authorization (RMA) number from the service center. Please call 1-800-234-5333.

1.3. Safety Information

Explanation of Signal Word Consequences

△ WARNING: Indicates a hazardous situation which, if not avoided, could result in serious injury or death

△ CAUTION: Indicates a hazardous situation which, if not avoided, could result in minor or moderate injury and/or property damage.

NOTICE: Indicates a situation which, if not avoided, could result in property damage.

Explanation of Safety and Related Symbols	
	Warning: Hazardous Voltage
	Caution: Lifting Hazard
	Attention: Read Accompanying Documentation
	CAUTION and WARNING

Explanation of Used Warning Symbols	
	Biohazard

Other Symbols	
	CE mark
	NEMKO Certification Mark for USA and Canada (Canada Standards Association)
	Catalog number
	Serial number
	Date of manufacture
	This system is subject to the European WEEE Directive 2002/96/EC. This product contains electrical and electronic components and must not be disposed of using standard refuse collection. Please consult local directives for disposal of electrical and electronic equipment.

**WARNING****To reduce the risks associated with hazardous voltage:**

- Do not service this instrument - there are no user serviceable parts. Use only Neogen authorized service personnel.
- Do not get the power supply wet.
- Use only the power cord specified for this product and certified for the country of use.
- Plug instrument into a properly grounded electrical outlet.
- Do not use a two-pin plug adaptor to connect to power.
- Ensure the instrument is positioned so that the power supply and main switch are easily accessible.
- Always turn the instrument off and disconnect it from the mains supply before cleaning, disinfection or decontamination.
- Always turn the instrument off and disconnect from outlet before positioning the pump tower.
- Do not use a damaged power cord or power supply.
- Qualified Neogen representative shall return instrument into service.
- Use only the bottles, tubing and components supplied and recommended for containing liquids on the instrument.
- Surfaces (floors, work table) must be dry when you are working with the system.
- Do not place heavy objects on the power cord.
- Do not bend, crimp, or otherwise damage the power cord.
- Replace damaged components with only the Neogen designated replacement parts.

To reduce the risks associated with release of contaminated product:

- Do not change default settings or update or modify software and software settings. These settings should not be changed except by Neogen qualified service representative.
- Ensure only qualified Neogen representative commissions the hardware after receiving instrument from service.
- The instrument must only be serviced by qualified Neogen representative.

⚠ CAUTION**To reduce the risk of exposure to biohazards:**

- Always follow standard lab safety practices, including wearing appropriate protective apparel and eye protection while handling instrument, contaminated samples and reagents.
- Ensure only qualified Neogen representative decommissions and commissions the hardware prior to shipping and after receiving.
- Avoid contact with the contents of the Neogen® MLS Microwell Plate.

To reduce the risks associated with release of contaminated product:

- Do not operate the instrument if the probes, tubing, syringe, pump or plunger are damaged, or an error message occurs on the screen.

To reduce the risks associated with environmental contamination:

- Follow applicable regulations when disposing of this device.
- Follow current industry standards for disposal of contaminated waste.

NOTICE**To avoid instrument damage:**

- Do not block or cover ventilation slots or openings.
- Do not place the instrument on an unstable or uneven surface.
- Do not use solvents such as acetone or thinner to clean the instrument; avoid abrasive cleaners.
- Clean instrument exterior with soft cloth using plain water or a mild all-purpose, or nonabrasive cleaner.
- Liquid detergents, disinfectants or decontamination liquids may not be poured into the instrument or sprayed inside the system.

To avoid instrument and/or property damage:

- Ensure that area in front of instrument loading tray is clear of obstructions.

Note: For Neogen MLS Reagent Kits safety information, please refer to specific Neogen MLS Kits Instructions for Use and Safety Data Sheet for additional reagent chemical hazards.

14. Environmental Operating Conditions

Environmental Condition	Operating Condition	Units
Altitude	Up to 2000 m (1.24 mi) above sea level. Storage: as required for air travel.	Meters
Temperature	Operating: +15 to +30°C Storage: +5 to 45°C Transport: +5 to 45°C	Celsius
Relative Humidity	20–80, Non-condensing	%
Frequency	50/60	Hertz
Sunlight	No direct sunlight. May mislead optical sensors and affect performance. Indoor use only.	
Dust	No excessive dust.	
Installation / overvoltage category	II	
Pollution Degree	2	

15. Regulatory Compliance

ROHS Information

Neogen MLS II Instrument complies with the RoHS Directive 2011/65/EU.

FCC Class A Notice

This device complies with part 15 of the FCC Rules. Operation is subject to the following two conditions (1) this device may not cause harmful interference, and (2) this device must accept any interference received, including interference that may cause undesired operation.

Note: This equipment has been tested and found to comply with the limits for a Class A digital device, pursuant to part 15 of the FCC Rules. These limits are designed to provide reasonable protection against harmful interference when the equipment is operated in a commercial environment. This equipment generates, uses, and can radiate radio frequency energy and, if not installed and used in accordance with the instruction manual, may cause harmful interference to radio communications. Operation of this equipment in a residential area is likely to cause harmful interference in which case the user will be required to correct the interference at his own expense.

IC Information

This Class (A) digital apparatus complies with Canadian ICES-003.

CAN ICES-3(A)/NMB-3(A)

EU Regulations

Warning

This is a class A Information Technology Equipment (ITE) product. In a domestic environment this product may cause radio interference in which case the user may be required to take adequate measures.

16. Package Contents

Do not unpack! The initial set-up of the Neogen MLS II Instrument shall only be performed by a Neogen authorized service personnel.

- Pump Tower Assembly (3 bags)
- Arrangement tool
- Microplate
- Bottle Holder Waste
- Pump Cover
- Reagent Bottle Holder Accessories
- Tower distribution PCR
- Polyethylene Waste Bottle
- Power Cable
- Allen Wrenches (3)
- Power Supply
- Certificate of Conformance
- Neogen MLS II User Software CD
- Neogen MLS II User Manual

1.7. User-Supplied Items

- Computer, connected through a USB port
 - Minimum requirements: Windows-based PC, Pentium 500 MHz or higher, CD-ROM drive, ≥100 MB free space, ≥64 MB RAM (128 MB recommended), Microsoft® Windows® 2000, Windows® XP, Windows® 7 or Windows® 8
- Computer monitor (SVGA color graphics card and compatible color monitor)
- Mouse
- Keyboard
- Printer (optional; Windows 2000 or XP compatible)
- Pipettes capable of dispensing 50-1000 µL
- Pipette tips
- Incubator

1.8. Accessory Items and Kits

- Neogen® Microbial Luminescence System (MLS) Ultra-High Temperature Beverage Screen Kit
- Neogen® Microbial Luminescence System (MLS) Ultra-High Temperature Dairy Screen Kit
- Neogen® Microbial Luminescence System (MLS) Injector Cleaning Kit
- Neogen® Microbial Luminescence System (MLS) Weekly Cleaning Kit
- Neogen® Microbial Luminescence System (MLS) Maintenance Solution
- Neogen® MLS Reagent Control Kit (ATP50)
- Neogen® Microbial Luminescence System (MLS) Microwell Plate
- Neogen® Microbial Luminescence System (MLS) Microwell Strip
- Neogen® Microbial Luminescence System (MLS) Microwell Strip Holders
- Neogen® Microbial Luminescence System (MLS) Short Neck Bottle Holders
- Neogen® Microbial Luminescence System (MLS) Service Contract

2. Technical Data

General Technical Data	
Neogen MLS II Dimensions	203 x 408 x 360 mm (with pump) 7.992 x 16.063 x 7.835 inches (with pump)
Neogen MLS II Weight	10.9 kg (24 lb) with 3 pumps and waste pump
Noise	67 dB with distance 1 m (3.28 ft)
Power Requirements	
Power Supply	Input: Voltage 100 V - 240 VAC +/- 10% Current: 1.6 A Frequency: 50 - 60 Hz Output: Voltage: 24 VDC Current: 5 A
Neogen MLS II	Voltage 24 VDC Current: 3 A
Modules	
Injectors	Up to 3 syringe pumps in measuring position, 10 – 300 µl in 1 µl increments
Shaker	Linear, orbital and double orbital: speed programmable
Waste Pump	Built in, software driven pump with priming function


Detection Unit	
Detector	Low noise PMT in photon counting mode
Dynamic range	>5 orders of magnitude
Plate	96 well opaque microplates and strip plates
Spectral range	280-650 nm
Performance/Detection Limit	
Crosstalk	< 2 x 10 ⁻⁵
Sensitivity	< 0.3 amol Aequorin
Packaging	
Cardboard box dimensions	515 x 570 x 310 mm 20.3 x 22.4 x 12.2 inch
Weight	12.9 kg (28.4 lb) (with 3 pumps and 1 waste pump)

3. Overview of Luminescence

Adenosine triphosphate (ATP) is the basic energy molecule in all living organisms. It is present in all microbial, plant, and animal cells. The term 'somatic' cell is commonly used for plant and animal (non-microbial) cells. ATP can also be present in materials of organic origin as 'free' ATP (outside of cells) when cells are disrupted during natural processes or manufacturing processes, such as homogenization or heat treatment.

The term bioluminescence is used to describe a reaction that produces light through a biological process. In the presence of ATP, an enzyme called luciferase, together with its substrate D-luciferin (collectively referred to as LL1 in this document) triggers a biochemical reaction that produces light as an end product. The amount of light emitted is proportional to the amount of ATP that is present. Low levels of light may be detected and measured using a sensitive luminometer.

3.1. Bioluminescence Reaction

3.2. Bioluminescence and the Neogen® Microbial Luminescence System (MLS) II

The Neogen MLS II platform utilizes ATP bioluminescence to detect potential microbial contamination in products such as those processed using ultra-high temperature (UHT) or extended shelf-life (ESL) technologies. To obtain reliable results for microbial detection in these products, they must be incubated prior to testing, allowing microbial growth to reach a detectable level. Substantial reduction of somatic and free ATP is also critical to allow accurate detection of microbial (contaminant) ATP; this is achieved using an ATP-degrading enzyme called ATPase. If no microbes are present, little or no detectable ATP will remain after the ATPase treatment. Upon addition of LL1, the signal will be low or undetectable, and therefore within the 'Pass' limit.

3.3. Summary of a Neogen Microbial Luminescence System (MLS) II Test Assay


Up to 96 samples may be loaded into an assay plate. The timing and addition of reagents is controlled by the Neogen MLS II Instrument. **Step 1:** ATPase is added to each test well to degrade somatic and free ATP; followed by a 15-minute incubation. **Step 2:** an Extractant is added to each test well to lyse the intact microbial cells to release the microbial ATP. **Step 3:** the LL1 reagent is added to the test well, creating a light-emitting reaction between the LL1 and microbial ATP. The Neogen MLS II Instrument reads the level of light that is emitted from each test sample and expresses the results in Relative Light Units (RLU). Figure 1 illustrates the steps in an UHT assay using the Neogen MLS II.

Results are available in 27 minutes or less (depending on the number of test samples), compared to 2-5 days using traditional microbiological methods.


Default Pass/Fail criteria are established, however, these criteria may depend upon the type of product and user's requirements. Appropriate Pass/Fail criteria are initially established for each type of product by running multiple samples against existing standard methods.

Step 1.

Step 2.

Step 3.


Figure 1. Steps in the bioluminescence reaction occurring during an UHT assay using the Neogen MLS II.

4. Initial Set-up of the Neogen MLS II

4.1. Neogen Microbial Luminescence System (MLS) II Instrument Installation

The initial set-up of the Neogen MLS II Instrument shall only be performed by Neogen authorized service personnel.

4.2. Neogen Microbial Luminescence System (MLS) II Software Installation

The end user understands and agrees that the Neogen MLS II Software is provided to the end user for the end user's convenience, without charge, on an «AS IS» basis. Neogen DISCLAIMS ALL WARRANTIES, EXPRESS OR IMPLIED, INCLUDING THE WARRANTIES OF MERCHANTABILITY AND FITNESS FOR A SPECIFIC PURPOSE, RELATING TO THE Neogen MLS II SOFTWARE. Neogen WILL NOT BE LIABLE TO END USER UNDER ANY CIRCUMSTANCE FOR ANY DAMAGE OR LOSS OF ANY KIND, INCLUDING DIRECT, INDIRECT, INCIDENTAL, SPECIAL, OR CONSEQUENTIAL DAMAGES, BASED UPON BREACH OF WARRANTY, BREACH OF CONTRACT, NEGLIGENCE, STRICT LIABILITY IN TORT OR ANY OTHER LEGAL THEORY EVEN IF Neogen HAS BEEN ADVISED OF THE POSSIBILITY OF SUCH DAMAGES, INCLUDING, BUT NOT LIMITED TO, LOSS OF PROFITS, REVENUE, EQUIPMENT, USE, DATA, OR INFORMATION OF ANY KIND.

1. As a final check before power is applied to the system, ensure all wires are safely covered, and the area is free from any spillage. It is now safe to plug the power cable into the socket.
2. Turn on the power to the Neogen MLS II by using the switch at the rear of the machine, then switch on the computer.
3. Download the Neogen MLS II Software onto the computer following the guidelines on the CD provided. When prompted by the Neogen MLS II Software, connect the PC to the newly set-up Neogen MLS II Instrument via the USB cable.
4. Once the Neogen MLS II Software has been successfully loaded, double-click on the 'MLSII' icon on the Desktop to open the program. For the initial use of the program, the user name will be 'Administrator' and the password section left blank. This will remain unchanged until altered by the user. Press the 'Enter' key.

4.21 Change Password

1. It is recommended that the password be changed as soon as possible to ensure maximum security. To carry out this step, click on the 'Tools' tab and select the Password button (Figure 2).


Figure 2. The Password button

2. An option will appear to allow the user to alter the password as desired. Enter the old password (leave blank when entering a password for the first time), enter the new password, repeat the same new password and click 'OK.' The Neogen MLS II system is now secured.

4.22 Access Control

1. The user can define a user's rights to individuals or groups of users. To carry out this step, click on the 'Tools' tab and select the Access Control button (Figure 3).


Figure 3. The Access Control button

2. The user can manage different user groups by selecting which Neogen MLS II Software functions can be performed.
3. To populate the different user groups, select the 'Users' tab where new user name can be added and maintained.

5. Neogen MLS II Software Functions

5.1. Log-on window

Figure 4. Log-on Window

This allows secured access to the software. Enter the user name and password.

The user can run the Neogen MLS II Software from a PC without connection to the Neogen MLS II Instrument simply by checking the '**Demo mode**' box. This may be useful when training personnel.

5.2 Actual Tab

5.2.1 Live Window

The '**Actual**' tab appears when the user opens the software. The '**Assay Protocol**' displays saved protocols. Choose the assay to run from the drop-down button.


Figure 5. Actual Window

Note: *Reagent Control Assay may or may not be labelled with "Cogent"

To select an assay, click on the name of the assay. It will show a Neogen MLS Microwell Plate graph representing the selected assay.

5.22 Actual Tab with Opened Assay Protocol

The user can select the number of wells, between 1 and 96, to test with the selected assay. Where wells are surplus to requirements, they can be removed from the assay by clicking on the unwanted well(s). Similarly, wells can be reselected using the same method. For selection or de-selection of entire columns or rows, click on the respective column or row identifier.


Figure 6. Actual tab with selected assay

- **Start:** Is used to start the selected assay
- When an assay is running, the user will see an estimated time remaining on the lower left-hand side of the 'Actual' tab window.
- **Pause:** Is used to pause or interrupt a running assay
- **Abort:** Is used to end or exit a running assay

5.23 Actual Tab with Finished Assay

The result status appears in each well as it is generated. The result appears in real-time and is color-coded according to the Pass/Caution/Fail values set for each assay: green (Pass), yellow (Caution), and red (Fail). The quantitative results can be viewed while an assay is running. The result of a specific test is displayed by placing the cursor over the well in question. All results shall be displayed when the assay is completed.



Figure 7. Actual tab with selected assay

5.3 Lot Specific Values For Tab

The window 'Lot Specific Values For' allows the user to include the batch name, batch number and expiration date of user-defined reagent(s) information.

- **Add:** Is used to add additional user-defined reagent(s) information.
- **Remove:** Is used to remove additional user-defined reagent(s) information.

5.6 Load Plate Tab

This dialog box allows the user to change the Plate ID and assure that the correct assay has been selected.

Load Plate

Plate ID: UHT Assay V1.1-2

Assay(s): UHT Assay V1.1 (Strips 1-12)

Barcode: Scan

Plate Layout:

	1	2	3	4	5	6	7	8	9	10	11	12
A	T1	T9	T17	T25	T33	T41	T49	T57	T65	T73	T81	T89
B	T2	T10	T18	T26	T34	T42	T50	T58	T66	T74	T82	T90
C	T3	T11	T19	T27	T35	T43	T51	T59	T67	T75	T83	T91
D	T4	T12	T20	T28	T36	T44	T52	T60	T68	T76	T84	T92
E	T5	T13	T21	T29	T37	T45	T53	T61	T69	T77	T85	T93
F	T6	T14	T22	T30	T38	T46	T54	T62	T70	T78	T86	T94
G	T7	T15	T23	T31	T39	T47	T55	T63	T71	T79	T87	T95
H	T8	T16	T24	T32	T40	T48	T56	T64	T72	T80	T88	T96

OK Cancel

Figure 9. Load Plate Information

- **Plate ID:** Is used as the name of the result file and the export file. The user can manually enter the plate name.
- **Assay(s):** Shows the assay(s) to be run on the displayed plate.
- **Plate Layout:** Shows the layout of the Neogen MLS Microwell Plate
- **Barcode:** Not Enabled
- **OK:** If the plate has been loaded, the run will start. If a scanner is used to load the plate, the bar code is read first and the run is started by clicking 'OK.'
- **Cancel:** Aborts the loading process. The user will be re-directed to the 'Actual' tab.

1.1. Report Tab

The user can access saved test results from the 'Report' tab.


Figure 10. Report tab

- **Available Results:** Displays all available test results defined by the assay
- **Open:** Opens saved reports
- **Print:** Prints saved reports
- **Sign:** Adds electronic signature of current user
- **Recalculate:** Recalculates the results
- **Save Model:** Saves open report

1.2 Protocol Editor Tab


Figure 11. Protocol Editor Tab

The functions in the 'Protocol Editor' tab allow modification of the protocol assays loaded in the Neogen MLS Software. It also allows the creation of new assays if it is specifically required by the user. Please note that it is not recommended that settings be altered without first consulting a Neogen Service Representative.

1.3 Tools Tab

The main functions located in the 'Tools' tab include:

- Plate Command
- Liquid handling
- Plate Editor
- Change Password Option
- Access Control
- Preferences


Figure 12. Tools Tab

1.3.1 Plate Command

Using the plate 'In' or plate 'Out' buttons, the user has the option of positioning the plate carrier to the desired in or out position while an assay is not running.

1.3.2 Liquid Handling

Most requirements will be met by using the 'Wash' and 'Prime' assays, but the three injectors, waste system and bottle monitoring option can be selected and controlled by clicking on the 'Liquid Handling' tab. When selected, the user is presented with a choice of setting preferences that can be applied to the injectors, waste pump or bottle monitoring system.


Figure 13. Liquid Handling Options

The injection system can be controlled by using the options displayed in the 'Prime Injector' field.

- **Injector:** Choose injector you need to prime
- **Volume:** Enter volume in microliters
- **Speed:** Enter speed
- **Cycles:** Enter cycles
- **Prime:** Starts the injection cycle as indicated by the options above

The user can define the 'Waste Pump' parameters.

- **Duration:** Defines the waste pump activity time in seconds
- **Activate:** Starts the waste pump

The user can empty the injector(s) using the 'Injector' field. This function is particularly useful when preparing the Neogen MLS II Instrument for transport or a period of dormancy.

- **Injector:** Choose the injector to be emptied
- **Empty:** Pumps fluid in the line back into the storage bottle
- **Empty all:** Pumps all back in one step
- **Up:** Moves the selected injector (syringe) up
- **All up:** Moves all injectors (syringes) up
- **Down:** Moves the selected injector (syringe) down
- **Down all:** Moves all injectors (syringes) down

The 'Reagent Bottle Monitor' option can be enabled, thereby allowing the user to accurately monitor or count down the volume of each reagent available throughout testing. This will visually alert the user if too little reagent is available to allow an assay to be completed. Refer to Section 8.3.

- **Enabled:** If bottle is empty the Neogen MLS II Instrument will sound a warning beep.

1.3.3 Plate Editor

The 'Plate Editor' tab allows the user to describe the dimension of a plate and links it with a supplier name. Please note that it is not recommended that settings be altered without consulting a member of the Neogen Service Representative. Also note that some options are not changeable.

Figure 14. Plate Editor

Here the user can enter the specifications of the plate. The saved data will show up in the Assay Header ID of the Assay Protocol.

- **Focus Point:** Ideal focus point Z-axis, if known
 ✓ The user can enter the A1 offset and the last well offset of the plate.
- **A1 Offset:** Edit the position of A1 in X-axis and Y-axis
- **Last Well Offset:** Edit the position of the last well in X-axis and Y-axis

1.3.4 Change Password

Refer to Section 4.2.1. The user can set up password-protected access to the Neogen MLS II Software.

1.3.5 Access Control

Refer to Section 4.2.2. The user can define different access rights to the user groups who want to access the Neogen MLS II Software.

5.7.6 Preferences

Upon installing the Neogen MLS II Software, the files are saved in C:\Program Files. The user can modify or browse the 'Directories' if needed.


Figure 15. Preferences

- **Language:** The user can select what language to use
- **Results:** To print test results automatically click on box

2. Cleaning the System for Initial Use

WARNING

To reduce the risks associated with exposure to chemicals:

- Always follow standard lab safety practices, including wearing appropriate protective apparel and eye protection while using the products.

When using any reagent with the device:


- Always secure the waste bottle to the device.
 - Avoid contact with the contents of the waste container.
1. Place 3 amber vials filled with Neogen MLS Injector Cleaning Solution into the short neck bottle holder, and insert each reagent-dispensing tube into each vial.
 2. Each reagent dispensing tube displays a label presenting a letter A, B or C. Please note that the liquid handling system – including the reagent dispensing tubes and syringes – will also be labelled as 1 (= injector A), 2 (= injector B) and 3 (= injector C) throughout this document.
 3. Double-click on the **MLSII** icon on the desktop, click on the **'Actual'** tab and select the **'Wash'** assay from the drop-down list. Place a Neogen MLS Microwell Plate in the plate carrier. During the **'Wash'/'Prime'** assay, the fluid is passed through the reagent tubes, into the waste canal (to the left of the plate carrier) and is then removed via the transparent waste tube at the rear of the instrument into the waste container. It is important that, despite the transfer of the reagents through the waste pipe into the waste container, thereby by-passing the plate area, a Neogen MLS Microwell Plate must be placed onto the plate carrier for each assay run. The instrument will not allow an assay to continue unless a plate has been detected on the plate carrier. As this plate will not be used for any testing, the same plate can be re-used for this purpose.
 4. It is necessary for the user to check that the waste pump is working efficiently and that the system does not have air entering the syringes. To ensure the waste pump is functioning as required, select the **'Tools'** tab, then the **'Liquid Handling'** button. Proceed to initialize the waste pump for 10 seconds. The user should be able to hear the pump running.
 5. If no sound is generated at this time, remove the tower cover carefully, and check the connection of the waste pump.
 6. The presence/absence of air is confirmed by watching the movement of fluid within the syringes during operation. The presence of bubbles indicate a loose connection – in which case each junction associated with the syringe in question will need to be checked and tightened. **Please note that bubbles will inevitably be present for the first set of aspirations, but they should quickly diminish.**

- Once the assay is complete, replace the 3 amber vials containing injector cleaning fluid with 3 amber vials filled with ATP-free water. Refer to step 3 to rinse out the injector cleaning fluid.

Please note that following a ATP Reagent Control assay (which is to be performed before a UHT assay is run; see Neogen MLS Reagent Control Kit (ATP50). Instructions for Use for more information), any contamination within the tubes and/or reagents highlighted by the assay may indicate a requirement for a more thorough cleaning than is usually recommended during routine use – in which case the injector cleaning fluid should be run through the reagent dispensing tubes several times using the ‘Wash’ assay.

- The set-up is now complete.

3. Use of the Neogen MLS II Washing/Priming the Injectors

- Injectors
- Dispensing Vials
- Pump Tower
- Waste Bottle
- Door to plate Carrier

Figure 16. A view of the Neogen MLS II Instrument

- Switch on the Neogen MLS II Instrument first and then the PC.
- Double-click on the **MLSII** icon on the Desktop
- Unless the password was already changed, you will see a password screen with a window showing ‘Administrator’ as the User and a blank dialog box for Password. Press the **ENTER** key. (The only reason to use an individual or lab password is to prevent unauthorized users.) You are now ready to prepare the Neogen MLS II Instrument for use.

3.1. Washing the Injectors

WARNING

To reduce the risks associated with exposure to chemicals:


- Always follow standard lab safety practices, including wearing appropriate protective apparel and eye protection while using the products.

When using any reagent with the device:

- Always secure the waste bottle to the device.
 - Avoid contact with the contents of the waste container.
- Place each reagent dispensing tube into an amber vial filled with Injector Cleaning Fluid.
 - Select the ‘Actual’ tab then the ‘Wash’ assay. Place a Neogen MLS Microwell plate onto the plate carrier. Ensure the waste tube has been securely inserted into the waste bottle.
 - Once the ‘Wash’ assay has been completed, replace the three amber vials containing Injector cleaning fluid with ATP-free water and repeat Step 2 to rinse out the Injector Cleaning Fluid. It is important that gloves are worn during this procedure to avoid contamination should hands make contact with the washed section of tubing.

3.2 Priming the Injectors

1. Replace vials containing ATP-free water with reconstituted ATPase (on injector A), Extractant (on injector B) and reconstituted LL1 (on injector C). Refer to kit instructions for details.
2. In the 'Actual' tab, select the 'Prime' assay. Place a Neogen MLS Microwell Plate onto the plate carrier.
3. If desired, enter lot specific data for reagent kit and click 'OK' – this will bring back the 'Actual' tab. Click 'Start' to begin the assay.
4. An ATP Reagent Control assay should be performed before starting any testing. We recommend this be done routinely before testing is commenced. Refer to Neogen MLS Reagent Control Kit (ATP50) product instructions for details of the Reagent Control assay.

1. Injector A (ATPase reagent)
2. Injector B (Extractant)
3. Injector C (LL1 reagent)

Figure 17. Positions of the reagent bottles

3.3 Reagent Use Calculation for the Neogen® Microbial Luminescence System Ultra-High Temperature (UHT) Screen Kits Reagents

Please refer to the Neogen MLS UHT Screen Kits Instructions for Use provided with each reagent kit.

Marks from 50 to 300 are shown on the left side of the ATPase and LL1 reagent bottles label and represent an estimation of the number of tests remaining in each reagent bottle.

Reagent volumes per kit:

Neogen MLS UHT Dairy Screen Kits #3000DPQCOG and #3060


Extractant	- 35 mL per bottle
ATPase	- 23.5 mL per bottle after rehydrated with ATPase Buffer
LL1	- 35 mL per bottle after rehydrated with LL1 Buffer

Reconstituted vials of ATPase and LL1 with the respective buffer can be used to run up to 300 tests.

Neogen MLS UHT Beverage Screen Kit #BEV600 (using UHT Assay)

Extractant	- 35 mL per bottle
ATPase	- 35 mL per bottle after rehydrated with ATPase Buffer
Beverage LL1	- 35 mL per bottle after rehydrated with LL1 Buffer

Reconstituted vials of ATPase and LL1 with the respective buffer can be used to run up to 300 tests.

 Alternatively, the user can select the '**Reagent Bottle Monitoring**' option to be activated by selecting the '**Enabled**' option within '**Tools**' tab and then the '**Liquid Handling**' dialog (Figure 18).

Figure 18. The Reagent Bottle Monitoring option

Example of Bottle Volume of Neogen MLS UHT Dairy Screen Kit

- Bottle 1 (ATPase): Set the maximum and current volume* to 23.5 mL
- Bottle 2 (Extractant): Set the maximum and current volume to 35 mL
- Bottle 3 (LL1): Set the maximum and current volume to 35 mL

* Volumes pertain to Neogen MLS UHT Dairy Screen Kits (#3000DPQCOG and #3060). Reagent volumes differ for other kits; refer to kit Instruction for Use for bottle volumes.

It is the responsibility of the user to ensure that all settings displayed in the bottle monitoring option are accurate. Edit reagent bottle information each time a new bottle is added to the system.

4. Ultra-High Temperature (UHT) Sterility Testing

WARNING

To reduce the risks associated with exposure to chemicals:

- Always follow standard lab safety practices, including wearing appropriate protective apparel and eye protection while using the products.

When using any reagent with the device:

- Always secure the waste bottle to the device.
- Avoid contact with the contents of the waste container.

Neogen MLS UHT Screen Kits are available for use with this instrument. These kits utilize ATP bioluminescence to detect microbial (contaminant) ATP in sterile or ultra-pasteurized beverages. The kit reagents are designed to exclude ATP from non-microbial sources and measure only microbial ATP. Detection of microbial ATP (outside appropriate Pass/Fail criteria) indicates microbial contamination in the product.

NOTE: This section applies to instruments that have been primed and passed the Reagent Control Assay. Please refer to the Neogen MLS Reagent Control Kit (ATP50) Product Instructions to perform a Reagent Control Assay.

1. You will need a clean 96-well Neogen MLS Microwell Plate (or Neogen MLS Microwell Strips for fewer numbers of tests).
2. Using a fresh pipette tip for each sample/ batch, pipette 50 μ L of UHT sample into the bottom of the wells. **It is important to avoid depositing any sample onto the well wall as this may cause the results to be inaccurate.** Always dispense the samples in the order of columns A1 – H1, A2 – H2, etc.

Figure 19. Pipette the sample into the bottom of the wells

3. In the 'Actual' tab, select the appropriate assay depending on the kit and product being tested in the scroll down menu and press the 'Start' button.
4. Carefully place the Neogen MLS Microwell Plate containing the samples onto the plate carrier (ensuring well 'A1' is positioned at the left corner closest to the user).

Figure 20. Bottom left-hand corner of plate carrier illustrating the 'A1' mark

5. If desired, complete the Lot Specific information and click 'OK'.
6. Complete the 'Load Plate' information section.
7. Click 'OK'. The assay will commence.
8. Refer to kit instructions for interpretation of assay results.

5. Cleaning the Instrument After Testing

With normal use, microbial contamination and biofilm formation can occur in the Neogen MLS II Instrument. The Neogen MLS Injector Cleaning Kit (3005) helps prevent microbial contamination and biofilm formation in the injectors and reagent lines, thus avoiding high background readings, false positive results, and inaccurate dispensing of reagents. For best results, use the Neogen MLS Injector Cleaning Kit at the start and end of each testing day.

WARNING

To reduce the risks associated with exposure to chemicals:

- Always follow standard lab safety practices, including wearing appropriate protective apparel and eye protection while using the products.

When using any reagent with the device:

- Always secure the waste bottle to the device.
 - Avoid contact with the contents of the waste container.
1. When cleaning the instrument, the reagent dispensing tubes and injectors must be rinsed by flushing ATP-free water through the system.
 2. On the bottle rack, replace the reagents, ATPase, Extractant and LL1, with an amber vial containing ATP-free water. Cover and store the LL1 and ATPase in the fridge, or discard appropriately. The Extractant may be stored at ambient temperature.
 3. Reagents may also be flushed by going into the 'Tools' tab, select 'Liquid handling' icon, then select 'Empty all.'
 4. Select the 'Actual' tab and click on the 'Wash' assay in the drop down list. Place a Neogen MLS Microwell Plate or strip holder onto the plate carrier.
 5. It is recommended that the waste container be emptied on a daily basis so as to avoid overflowing and the need to dispose of a large amount of waste fluid at any one time.
 6. Once the 'Wash' assay is complete, place each reagent tube into a vial containing injector cleaning fluid and repeat the 'Wash' assay. Allow the reagent tubes to stand in the container of injector cleaning fluid while not in use.
 7. At the end of the day, the user can switch off the Neogen MLS II Instrument via the button at the rear of the machine. It is now safe to continue to turn off your computer.

6. Weekly Cleaning

With normal use, microbial contamination and biofilm formation can occur in the Neogen MLS II Instrument. The Neogen MLS Injector Maintenance Solution/Protein Removal Fluid (3006 or BMLSCK) removes protein residues and microbial contamination in the injectors and reagent lines, thus avoiding high background readings, false positive results, and inaccurate dispensation of reagents. To properly maintain your Neogen MLS II Instrument, use the Neogen MLS Maintenance Solution/Protein Removal Fluid to clean the injectors a minimum of once a week.

WARNING

To reduce the risks associated with exposure to chemicals:

- Always follow standard lab safety practices, including wearing appropriate protective apparel and eye protection while using the products.

When using any reagent with the device:

- Always secure the waste bottle to the device.
- Avoid contact with the contents of the waste container.

Cleaning the Injectors

1. Remove the tubing from the reagent bottles and replace with vials of ATP-Free Water.
2. Click on the 'Actual' tab select the 'Wash assay' from the menu. Place a Neogen MLS Microwell Plate onto the plate carrier and click 'OK'.

3. Once the 'Wash' assay is complete, replace the ATP-free water with amber vials filled with Neogen MLS Maintenance Solution or Neogen MLS Protein Removal Fluid.
4. Repeat Step 2 above.
5. Allow the Neogen MLS Maintenance Solution or Neogen MLS Protein Removal Fluid to remain in the line(s) for a minimum of 20 minutes, preferably 2 hours. **It is recommended not to exceed 4 hours.**
6. Replace the Neogen MLS Maintenance Solution or Neogen MLS Protein Removal Fluid vials with vials of ATP-Free water and repeat Step 2 twice to ensure the lines are free from the weekly cleaning solution and other contaminants.
7. Once the 'Wash' assay is complete, place each reagent tube into a vial containing injector cleaning fluid and repeat the 'Wash' assay. Allow the reagent tubes to stand in the container of injector cleaning fluid while not in use.

7. Routine and Periodic Maintenance Summary

Maintenance	Before and After Testing Daily	Weekly	Monthly	Annually (or every 60,000 tests if sooner)	When required
Routine Maintenance					
Clean the lines with ATP-free water	●				
Clean the lines using Injector Cleaning Fluid and ATP-free water	●				
Clean Injectors with Weekly Cleaning Solution		●			
Periodic Maintenance – by Neogen Service Technician					
Verify injector volumes				●	
Check tubing lines – possibly replace				●	
Check Syringes – possibly replace				●	
Check valves – possibly replace				●	
Check ribbon cable – possibly replace				●	
Exchange mirror				●	
Clean fan filter – possibly replace				●	

8. Moving the Neogen MLS II Instrument to Another Workbench

Contact your Neogen Food Safety Service Representative to move the Neogen MLS II to another workbench.

9. Instructions for Shipping or Prolonged Storage of the Neogen MLS II Instrument

WARNING

To reduce the risks associated with exposure to chemicals:

- Always follow standard lab safety practices, including wearing appropriate protective apparel and eye protection while using the products.

To reduce the risks associated with hazardous voltage:

- Always turn the instrument off and disconnect from outlet before cleaning, disinfecting or decontaminating.

To reduce the risks of exposure to biohazards:

- Always follow standard lab safety practices, including wearing appropriate protective apparel and eye protection while handling instrument, contaminated samples and reagents.

Before shipping or prolonged storage of the Neogen MLS II Instrument contact your local Neogen Food Safety Representative.

10. Instructions for Installing the Neogen MLS II Instrument After Service Repair

When the Neogen MLS II Instrument returns from the Service Center contact your Neogen Food Safety representative.

11. Instructions for Cleaning External Surfaces of the Neogen MLS II Instrument

WARNING

To reduce the risks associated with exposure to chemicals:

- Always follow standard lab safety practices, including wearing appropriate protective apparel and eye protection while using the products.

To reduce the risks associated with hazardous voltage:

- Always turn the instrument off and disconnect from outlet before cleaning, disinfecting or decontaminating.

To reduce the risks of exposure to biohazards:

- Always follow standard lab safety practices, including wearing appropriate protective apparel and eye protection while handling instrument, contaminated samples and reagents.

△ CAUTION

To reduce the risks associated with environmental contamination, which if not avoided could result in minor or moderate injury:

- Follow current industry standards for disposal of contaminated waste.

Liquid detergents, disinfectants or decontamination liquids **may not be poured or sprayed into the instrument or sprayed inside the system.** Always use saturated disposable paper towels.

- For the cleaning, disinfection or decontamination, a cloth moistened with detergent, disinfectant or decontamination liquid should be used.
- Only approved detergents, disinfectants or decontamination liquids may be used.
- Only approved cleaning, disinfection or decontamination methods may be used.

Persons performing the decontamination process must be familiar with the operation of the instrument.

Cleaning the external surfaces

1. Prepare decontamination solution (70% ethanol). Use prophylactic gloves, lab coat and safety glasses when preparing and using these solutions.
2. Wipe all external shrouds with a disposable paper towel saturated with soapy water.
3. Wipe off the external shrouds with a disposable paper towel saturated with water.
4. Wipe all external shrouds with a disposable paper towel saturated with decontamination solution. Allow the decontamination solution to remain at least 4 minutes on the surfaces to dry.
5. Wipe off the external shrouds with a disposable paper towel saturated with water.

Approved decontamination solutions

Solution	Preparation
70% ethanol	Measure 700 mL of 96% ethanol and add enough water to have a total volume of 1000 mL of solution.

Important: Decontamination solution must be thoroughly wiped off and out of the system!

12. Contact Information

Contact your local Neogen Service personnel. Have a description of the problem with as much detail as possible. Save relevant log files.

For Neogen Food Safety Department Product Information: neogen.com

Appendix I – Troubleshooting Guide:

Failures	Possible Root Cause(s)	Corrective Actions
Air bubbles visible in injector syringe tubes.	Reagent lines are not submerged in the liquid. Reagents not primed before use.	Submerge lines in solution and prime the system
Front door to plate carrier opened during run	User's fault	Re-run assay
Wrong plate position	Guide bar X-axis contaminated or in need of lubrication	Call your Neogen Service Personnel
Wrong plate carrier position	Guide bar Y-axis contaminated or in need of lubrication	Call your Neogen Service Personnel
Wrong plate carrier position	X-Axle or Y-Axle drive belt broken	Call your Neogen Service Personnel
Too low or too high RLU	Injector/tubing contaminated	Clean system with cleaning kit (sections 9 and 10) and run a Reagent control assay
Too low or too high RLU	Clogged injection probe	Call your Neogen Service Personnel
Inaccurate injection volume	Injection probe snapped off	Call your Neogen Service Personnel
Washing station overflows	Blocked waste tube	Call your Neogen Service Personnel
Syringe is lose, system is leaking	Syringe was not completely screwed in	Call your Neogen Service Personnel
Plunger of syringe is lose, inaccurate injection volume	Plunger was not completely locked	Call your Neogen Service Personnel

Neogen Food Safety

Neogen Corporation

620 Lesher Place
Lansing, MI 48912 USA
Neogen.com

Neogen Europe Ltd.

The Dairy School
Auchincruive
Ayr, KA6 5HU
Scotland, UK

Neogen Ireland, Ltd.

Bray Business Park, Bray
Co. Wicklow
A98YV29, Ireland

Neogen Corporation

620 Lesher Place Lansing, MI 48912 USA
neogen.com

© Neogen Corporation 2024. All rights reserved.
Neogen and Petrifilm are registered trademarks of Neogen Corporation.
FS01346A

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The screenshot shows the 'Plate Editor' dialog box with the following fields and values:

- Plate Name: []
- Supplier: []
- Width: 0 mm
- Height: 0 mm
- Depth: 0 mm
- Columns: 1
- Rows: 1
- Volume: 0 ul
- Bottom: []
- Colour: []
- Material: Other
- Focus Point: 0 mm
- A1 Offset X: 0 mm
- A1 Offset Y: 0 mm
- Last Well Offset X: 0 mm
- Last Well Offset Y: 0 mm
- Created by: []
- Last edited: []
- Checksum: []
- Electronic signature: []

Buttons at the bottom: Close, New Plate, Open..., Edit..., Save, Sign, Audit Trail.

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