*Read instructions carefully before starting test*

# Veratox<sup>®</sup> for Gliadin

*Refrigerate at 2–8°C (35–46°F) • Do not freeze*

## **GLIADIN/GLUTEN**

Gliadin is an alcohol-soluble protein found in wheat that belongs to a group of proteins called prolamins. Other prolamins include secalin, found in rye, and hordein, found in barley. Gluten consists of two groups of proteins (prolamins and glutelins) that are found in differing amounts in wheat, barley, rye and oats.

Gliadin and other prolamins have been identified as major causal agents in a number of disorders, including wheat allergy and gluten intolerance (celiac disease). Wheat allergy is a specific immune response to a number of wheat proteins, including gliadin, albumin, globulin, and glutenin. Celiac disease is a chronic reaction to gluten proteins that results in the poor absorption of nutrients in the small intestine.

Those with celiac disease must avoid gluten, and rely upon the correct labeling of food to make appropriate, safe food choices. Testing for the presence of gluten components ensures food manufacturers that an unlabeled—and potentially dangerous—ingredient did not make its way into a food product.

## **INTENDED USE**

Veratox<sup>®</sup> for Gliadin is intended for the quantitative analysis of ingredients, clean-in-place solutions, and finished food products intended to be gluten-free for the presence of gliadins and prolamins found in wheat, barley and rye.

## **INTENDED USER**

This test kit is designed for use by quality control personnel and others familiar with foods possibly contaminated by gluten. Since technique is very important, operators should be trained by a Neogen representative or someone who has completed the Neogen training.

## **ASSAY PRINCIPLES**

Veratox for Gliadin is a sandwich enzyme-linked immunosorbent assay (S-ELISA). Gliadin is extracted from samples with a 40% ethanol solution by shaking in a shaker or rotator. Extract is diluted in Phosphate Buffered Saline (PBS) and diluted samples are added to antibody-coated wells (capture antibody) where gliadin will bind to the antibody during an incubation period. Any unbound gliadin is washed away and a second antibody, which is enzyme labeled (detector antibody) is added. The detector antibody binds to the gliadin during another incubation period. Unbound enzyme-labeled antibody is washed away and a one step substrate is added. Color develops as a result of the presence

of bound-labeled antibody. A stopping reagent is added and the color of the solution is observed. Blue color indicates samples containing high levels of gliadin while purple or red samples contain little or no gliadin. The optical densities of the controls form a standard curve, and the sample optical densities are plotted against the curve to calculate the exact concentration of gliadin in parts per million (ppm).

### STORAGE REQUIREMENTS

The kit can be used until the expiration date on the label when stored refrigerated at 2–8°C (35–46°F).

### MATERIALS PROVIDED

1. 48 antibody-coated microwells
2. 48 red-marked transfer wells
3. 5 yellow-labeled bottles of 0, 5, 10, 20, and 50 ppm gliadin controls (0, 12.5, 25, 50 and 125 ng/mL gliadin)
4. 2 blue-labeled bottles of enzyme-labeled antibody conjugate
5. 1 green-labeled bottle of K-Blue® Substrate
6. 1 red-labeled bottle of Red Stop solution
7. 40 mL of 10 mM PBS-Tween washing reagent in a wide mouth bottle. Each bottle can prepare 1 L in distilled or deionized water (pH 7.4)
8. 1 foil pouch of 10 mM PBS dry powder containing enough powder to prepare 1 L of dilution solvent
9. 50 g of extraction additive in 2 specimen cups
10. Plastic scoop to measure extraction additive

### MATERIALS RECOMMENDED BUT NOT PROVIDED

1. Orbital rotator or shaker to hold 50 cc centrifuge tubes for a 1 g sample, or shaker water bath with clamps adjusted to hold 125 mL (4 oz) extraction bottles for a 2 g sample
2. Oven or water bath adjustable to 50°C if analyzing heat-processed samples
3. Scale capable of weighing  $1 \pm 0.1$  g (Neogen item 9427), or scale capable of weighing  $0.25 \pm 0.01$  g if using gliadin cocktail solution (Neogen item 9395)
4. Pipettor, 50–200 µL adjustable (Neogen item 9276)
5. Pipettor, 12-channel (Neogen item 9273)
6. Pipette tips (Neogen item 9410, 9417, 9407)
7. Three 1 L bottles to prepare washing solution, extract solution, and sample extract dilution solution (Neogen item 9472)
8. Test tubes to perform sample extract dilution
9. Timer (Neogen item 8426, 9452)
10. Microwell reader with a 650 nm filter (Neogen item 9303)
11. Three reagent boats for 12-channel pipettor (Neogen item 9435)
12. Microwell strip holder (Neogen item 9402)
13. Wash bottle (Neogen item 9400)
14. Paper towels or equivalent absorbent material
15. Waterproof marker
16. Distilled or deionized water
17. Laboratory grade ethanol (190 proof)
18. Special extraction additive for dark chocolate, cocoa and tannin (Neogen item 8482)
19. Cocktail solution for analysis of heat-processed samples (Neogen item 8483)

### PRECAUTIONS

1. Ethanol solution is highly flammable. Keep container tightly closed, and keep away from heat, sparks, open flame and those smoking. It is toxic if swallowed, or if vapor is inhaled. Avoid contact with skin.
2. Use chemical hood when using cocktail solution to analyze heat-processed samples.
3. Components of Veratox for Gliadin, such as controls and extraction additive, may contain one or more of the following potentially allergic materials: gluten, casein, almond protein and soy protein. If allergic to any of these compounds, use caution when using this product.
4. Store test kit between 2–8°C (35–46°F) when not in use. Do not freeze test kits, and avoid prolonged storage at ambient temperatures.

5. Bring kits to room temperature (18–30°C, 64–86°F) prior to use.
6. Do not use kit components beyond expiration date.
7. Do not mix reagents from one kit with reagents from a kit with a different serial number.
8. Do not run more than 24 wells per test.
9. Follow proper pipetting techniques (e.g., prime tips and use clean tips).
10. Use only incubation times specified; others may give inaccurate results.
11. Use clean pipette tips and glassware for each sample to avoid cross-contamination. Thoroughly wash all glassware between samples.

## PROCEDURAL NOTES

1. **Dark chocolate, cocoa and tannin.** If testing commodities where dark chocolate, cocoa, and/or tannin are present, such as dark chocolate bars and cocoa powder, contact Neogen for a special extraction additive (Neogen item 8482). Add 1 scoop of this special extraction additive in addition to 1 scoop of the additive supplied with this test kit.
2. **Substrate.** K-Blue Substrate is ready for use. The substrate should be clear to light blue — discard if it has turned dark blue. Only pour the needed volume of substrate into a reagent boat. **Do not return unused substrate to the bottle.** Cover the reagent boat to keep the substrate protected from light until needed.
3. **Sample extract dilution solution (PBS).** Prepare extract dilution solution by adding a foil pouch of dilution solvent, 10 mM PBS, to 1 L distilled or deionized water. Swirl to mix thoroughly.
4. **Wash buffer.** Prepare the wash buffer solution by diluting the wash buffer concentrate in 960 mL of distilled or deionized water in an empty 1 L container. Ensure the transfer of all concentrate by rinsing the bottle several times with water, and swirl to assure thorough mixing.  
**NOTE:** Discard unused portions of extraction solution and wash buffer when the test kit has been used completely.
5. **Antibody wells.** Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after samples are extracted, and the test procedure is set to begin.

## SAMPLE PREPARATION AND EXTRACTION

For analyzing **heat-processed** samples, follow extraction **procedure C**. For all commodities that were **not heat-processed**, follow either extraction procedure **A or B**. Samples of an unknown origin should be extracted using extraction procedure C.

### A. Extraction of non-heat processed samples with orbital shaker or rotator

1. Prepare 40% ethanol extraction solution by combining 4 parts ethanol with 6 parts distilled water. Prepare sample extract dilution solution (PBS) as detailed in procedural note 3.
2. Add 1 g ground sample, or 1 mL liquid sample, to a clean 50 cc tube.
3. Add 1 scoop of extraction additive to the tube (see procedural note 1).
4. Add 10 mL (9 mL for liquid samples) of 40% ethanol to the tube, cap tightly, then shake the tube vigorously by hand for about **20 seconds**, or vortex for **10 seconds**, to ensure complete mixing.
5. Extract by shaking (150 rpm) in an orbital shaker or rotator by laying down the tube on its side over the flat pad of the instrument, and holding it tightly using a rubber band or tape. Rotate or shake for **15 minutes** at room temperature.
6. Remove the tube and let it stand in a rack for about **10 minutes** to enable the sample extract to settle before withdrawing the clear extract.
7. Dilute each sample 1:40 by withdrawing 100 µL of the upper layer of the extract and transferring it to a small tube or vial containing 3.9 mL of PBS, or mix 0.5 mL of extract with 19.5 mL of PBS in a test tube.
8. To mix, vortex the tube for **5 seconds**, or invert several times by hand.
9. Test diluted samples within **2–3 hours** of extraction.

**B. Extraction of non-heat processed samples with shaker or shaker water bath**

1. Prepare 40% ethanol extraction solution by combining 4 parts ethanol with 6 parts distilled water. Prepare sample extract dilution solution (PBS) as detailed in procedural note #3.
2. Add 2 g ground sample, or 2 mL liquid sample, to a 125 mL clean extraction bottle.
3. Add 1 scoop of extraction additive to the bottle (see procedural note 1).
4. Add 20 mL (18 mL for liquid samples) of 40% ethanol, cap the bottle tightly, then shake vigorously by hand for about **20 seconds** to ensure complete mixing.
5. Extract by shaking (150 rpm) in a shaker for **15 minutes** at room temperature (a shaker water bath can work, but do not turn the heat on). Remove the bottle from shaker or bath.
6. Let the bottle stand for about **10 minutes** to enable some of the sample to settle before withdrawing the clear extract.
7. Dilute each sample 1:40 by withdrawing 100 µL of the upper layer of the extract and transferring it to a small tube or vial containing 3.9 mL of PBS, or mix 0.5 mL of extract with 19.5 mL of PBS.
8. To mix, vortex the tube for **5 seconds**, or invert several times by hand.
9. Test diluted samples within **2–3 hours** of extraction.

**C. Extraction of heat processed commodities**

Heat-processed commodities require an extraction cocktail solution (Neogen item 8483) that re-natures the heated gliadin and allows the accurate detection of any possible gliadin in a sample. To extract gliadin from heat-processed samples:

1. Prepare 55% ethanol extraction solution by combining 55 parts ethanol with 45 parts distilled water.
2. Prepare sample extract dilution solution (PBS) as detailed in procedural note 3.  
**NOTE:** Extraction with extraction cocktail solution should be performed under a chemical hood.
3. Weigh out 0.25 g sample to 50 cc screw cap centrifuge tube.
4. Add 2.5 mL of extraction cocktail solution (dilution factor 1:10).
5. Cap and vortex **10–20 seconds** to homogenize cocktail and sample.
6. Incubate **40 minutes** at 50°C (water bath or oven).
7. Remove samples and let cool for **5–10 minutes**.
8. Add 1 scoop of special extraction additive (see procedural note 1).
9. Add 7.5 mL of 55% ethanol and vortex again for **10–20 seconds** (the final concentration of ethanol will be 41% and the sample dilution to this point is 1:40).
10. Shake (150–200 rpm) for **1 hour** at room temperature on a rotator (tube on its side).
11. Centrifuge sample (if necessary) for **5 minutes** at 2500 rpm.
12. Dilute the sample 1:10 into PBS (200 µL sample into 1.8 mL PBS).
13. Samples are ready to run (1:400 final dilution).

**TEST PROCEDURE**

Allow the test kit and all reagents to warm to room temperature (18–30°C, 64–86°F) before using.

1. Remove 1 red-marked mixing well for each sample to be tested plus 5 red-marked wells for controls, and place in the well holder.
2. Remove an equal number of antibody-coated wells. Return antibody wells which will not be used immediately to the foil pack with desiccant. Reseal the foil pack to protect the antibody. Mark one end of the strip with a "1", and place strip in the well holder with the marked end on the left.
3. Mix each reagent by swirling the reagent bottle prior to use.
4. Using a new pipette tip for each, transfer 150 µL of controls and sample extracts to the red-marked transfer wells as shown in the template below. Only run up to two 12-well strips at a time.

0	5	10	20	50	S1	S2	S3	S4	S5	S6	S7
S8	S9	S10	S11	S12	S13	S14	S15	S16	S17	S18	S19

- Place tips on the 12-channel pipettor and transfer 100  $\mu$ L of the controls and sample extracts to the antibody-coated wells. Mix for **20 seconds** by sliding the well holder back and forth on a flat surface.
- Incubate microwells **10 minutes** at room temperature (18–30°C, 64–86°F). Discard the red-marked transfer wells.
- Empty the contents of the wells into a sink. With a wash bottle fill each antibody well with the wash buffer solution and dump out. Repeat the washing 5 times, then turn the wells upside down and tap out on a paper towel until all washing solution is removed.
- Pour the needed volume of conjugate from the blue-labeled bottle into a clean reagent boat.
- Using the 12-channel pipettor and new tips, transfer 100  $\mu$ L of the conjugate into all the wells and mix for **20 seconds** by sliding the well holder back and forth on a flat surface.
- Incubate for **10 minutes** at room temperature (18–30°C, 64–86°F).
- Wash all wells with the wash buffer solution as described in step 7.
- Pour the needed volume of substrate solution from the green-labeled bottle into a clean reagent boat.
- Place new tips on the 12-channel pipettor and transfer 100  $\mu$ L of substrate into each well and mix for **20 seconds**. Do not eject tips.
- Incubate for **10 minutes** at room temperature (18–30°C, 64–86°F).
- Pour the needed volume of Red Stop solution from the red-labeled bottle into a clean reagent boat.
- With the same tips used to dispense the substrate, transfer 100  $\mu$ L of Red Stop into each well and mix for **20 seconds**.
- Wipe the bottom of the microwells and read in a microwell reader with a 650 nm filter.
- Interpret the test's results using the Neogen 4700 microwell reader, or an equivalent strip reader. If using a strip reader, calculate the results using Neogen's Veratox software for Windows.

### INTERPRETATION OF RESULTS

Standard controls were made from wheat gliadin and calculated as gliadin. Approximately 50% of the gluten is available as gliadin. Therefore, to calculate the gluten value of the samples, multiply the ppm gliadin results by 2.

### PERFORMANCE CHARACTERISTICS

**Limit of quantitation:** 5 ppm (Described as the lowest concentration point on the calibration curve that this test can reliably detect gliadin.)

**Range of quantitation:** 5–50 ppm (For quantitating samples above 50 ppm, contact a Neogen representative for dilution instructions).

### CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

### SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at [foodsafety.neogen.com](http://foodsafety.neogen.com), or by calling Neogen at 800.234.5333 or 517.372.9200.

### TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit [www.neogen.com/en/terms-and-conditions](http://www.neogen.com/en/terms-and-conditions).

### WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

## TESTING KITS AVAILABLE FROM NEOGEN

### Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

### Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

### Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

### Food allergens

- Almonds, coconut, crustaceans, eggs, gliadin, hazelnut, milk, multi-treenut, mustard, peanuts, sesame, soy, walnuts

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