
Please note this kit insert has been revised as of June 2020.  
Read instructions carefully before starting test.

# **Reveal<sup>®</sup> 3-D**

## for Hazelnut

### **INTRODUCTION AND INTENDED USE**

Reveal<sup>®</sup> 3-D for Hazelnut is uniquely designed with three lines of detection and can be used virtually anywhere to screen environmental swabs, clean-in-place (CIP) rinses, and food products and ingredients for the presence of significant levels of hazelnut. The test's 3-D technology ensures greater reliability with screening than ever before.

Reveal 3-D for Hazelnut is intended for use only in an industrial food manufacturing/preparation context. Because of the problems of adequate sampling and extraction of hazelnut, it is not suitable for the testing of foods to be consumed in the home or in a restaurant by allergic individuals. The test detects significant (low parts per million, or ppm) levels of hazelnut in CIP rinses, environmental swabs, and foods.

### **LIMIT OF DETECTION**

The test utilizes highly specific antibodies to detect hazelnut protein. Hazelnut can be detected from various surfaces using the provided environmental swabs or from CIP rinses and foods.

Utilizing the environmental swabs supplied, levels of 10  $\mu\text{g}/100 \text{ cm}^2$  of hazelnut extract can be detected.

When analyzing CIP rinses and foods, hazelnut residues are detectable as low as 0.75-1.5 ppm hazelnut protein (5-10 ppm total hazelnut). The presence of cleaners and sanitizers in rinses, as well the type of food matrix being tested, can affect limit of detection.

## **CROSS REACTIVITY**

The test was used to analyze a panel of potentially cross reacting commodities including nuts, seeds, legumes and grains. Of these, whole walnut gave a clear cross reactivity with the test. Further information is available in the Reveal 3-D Hazelnut validation report, which is available on request.

## **TEST PERFORMANCE**

Reveal 3-D for Hazelnut has undergone rigorous validation to evaluate the specificity, sensitivity, robustness and intra- and inter-batch variability of the test method on CIP rinses, environmental swabs, and foods.

## **SAMPLE COMPATIBILITY**

Reveal 3-D for Hazelnut is designed to detect hazelnut on environmental swabs, CIP rinses and foods. Although every effort has been made to validate as many variables as possible, there may be some sample types that are not suitable for testing.

Users should perform in-house matrix-specific spike recovery validation work in conjunction with a validated laboratory assay (e.g., Veratox®) to help confirm or validate results. This procedure will highlight any problematic matrices encountered.

For certain food ingredients or complex food matrices, it is recommended that the use of a more broadly inclusive test method, such as Veratox for Hazelnut Allergen, be used to evaluate. Existing commodity validations should be revalidated when suppliers or the manufacturing process has been changed. Please contact a NEOGEN representative or distributor for additional details.

## **MATERIALS PROVIDED**

The Reveal 3-D for Hazelnut test pack contains the following:

1. An instruction leaflet
2. 1 foil pouch, containing 10 blue Reveal 3-D for Hazelnut devices
3. 10 sachets containing Type 11 Extraction Buffer
4. 10 sample tubes and caps
5. 10 individually packaged, sterile swabs with break-off tips
6. 1 bottle of swab wetting solution

## **MATERIALS REQUIRED FOR FOOD EXTRACTION BUT NOT PROVIDED**

1. 3-D Food Buffer (NEOGEN item 8504)
2. 50 mL graduated tubes, 25 ea (NEOGEN item 9381)

## **SAMPLING TECHNIQUE**

### **1. Clean-in-place rinse water**

As only a small amount of material is required for the Reveal 3-D for Hazelnut test, it is important to test a representative portion of liquid.

### **2. Environmental swabbing**

The swabs supplied are intended to be used for the collection of environmental samples from which the presence of hazelnut can be tested. This method can be used to validate the adequacy of cleaning and/or to identify problem areas (e.g., unwanted build-up of hazelnut in processing equipment).

### **3. Food samples**

As only a small sample of material is required for Reveal 3-D for Hazelnut, it is important to test a representative portion of the sample.

## PRECAUTIONS

1. For testing in an industrial food manufacturing/preparation, or labeling enforcement context only.
2. Do not use any part of the test beyond the expiry date.
3. Do not open the foil pouch until just before use.
4. Ensure the foil bag is tightly sealed after removal of a device.
5. Always store the kit between 2–8°C (35–46°F). Avoid freezing.
6. Bring kit to room temperature 18–30°C (64–86°F) prior to use.
7. The test devices use a pale green loading dye at T (test), O (overload) and C (control) positions for quality purposes. The pale green lines should be visible on the devices prior to use. If the lines are not visible, please do not use the device and contact a NEOGEN representative for assistance.

## TEST STORAGE / SAMPLE STABILITY

Store the Reveal 3-D for Hazelnut kit between 2–8°C (35–46°F), and use within the expiry date stated on the outer label. Extracted samples should be used within 3 hours of extraction.

## LIMITATIONS


A negative test from a surface, CIP rinse, or food cannot exclude the possibility that the sample or swabbed surface has hazelnut present since it may be distributed unevenly on the surface, in the sample or rinse and may be below the limit of detection of the test.

Reveal 3-D for Hazelnut is qualitative and should only be used as a preliminary screen for the presence of hazelnut. The validity of results obtained with the test should preferably be viewed in conjunction with data from a quantitative assay, such as Veratox for Hazelnut Allergen.

## RINSE SAMPLING

Remove the following and allow to equilibrate at room temperature before use (**20–30 minutes** out of refrigerator):



1. 1 sample tube
  2. 1 Type 11 Extraction Buffer sachet
  3. 1 Reveal 3-D for Hazelnut device (in foil pouch)
1. Carefully tear/cut off the top of the Type 11 Extraction Buffer sachet and add the entire contents to the sample tube.
  2. Add the 1 mL sample to the sample tube. If a pipette is unavailable, 1 mL may be crudely estimated by half-filling a white sample tube cap. **NOTE:** Test results will be less accurate.
  3. Secure the white cap and shake for **1 minute**.

## SWAB SAMPLING


Remove the following and allow to equilibrate at room temperature before use (**20–30 minutes** out of refrigerator):

1. 1 bottle of swab wetting solution
  2. 1 sample tube
  3. 1 Type 11 Extraction Buffer sachet
  4. 1 Reveal 3-D Hazelnut device (in foil pouch)
  5. 1 sterile swab
1. Carefully tear/cut off the top of the Type 11 Extraction Buffer sachet and add the entire contents to the sample tube.
  2. Estimate a swabbing area of approximately 10 x 10 cm. Alternatively, use the swab to collect samples of contamination from problem areas (e.g., of processing equipment).
  3. Gather the sample with the swab, using one of the following methods:


**For dry surfaces:** Remove a sterile swab from the packaging and wet with two drops of swab wetting solution. Swab a 10 x 10 cm area by using a crosshatch technique revolving the swab on the surface. Repeat this swabbing procedure using movements at right angles to those used in the first swabbing.

**For wet surfaces:** Remove a sterile swab from the packaging and swab a 10 x 10 cm area by using a crosshatch technique, revolving the swab on the surface. Repeat this swabbing procedure using movements at right angles to those used in the first swabbing.


4. Return the swab to the Extraction Buffer in the sample tube and carefully break off the moistened end at the pre-scored mark so that it remains in the tube.
5. Secure the cap of the sample tube, taking care to ensure that the stem does not prevent the tube from being properly sealed. Shake for **1 minute**.


## SWAB AND RINSE TESTING

1. Remove the lid and fill it with the liquid from the tube. Any froth should remain in the tube.
2. Dip the head of the Reveal 3-D device into the liquid in the lid. Ensure that the cavity is saturated with the liquid.
3. Leave the cavity saturated until liquid is seen running in the test window.
4. Place device on a flat surface and allow test to develop for **5 minutes**.

**NOTE:** It is essential to place the device flat on a level surface as soon as the liquid has entered the test window to stimulate flow through the device. Additionally, the devices are pre-stripped with pale green loading dye in positions T (test), O (overload) and C (control). The loading dye assists with quality and manufacturing checks and does not impact test performance. The loading dye is removed from the test window as the sample flows through the device.



## FOOD EXTRACTION AND TESTING

The extraction of food samples requires the use of NEOGEN's 3-D Food Buffer (NEOGEN item 8504).


From the test kit, remove the following and allow to equilibrate at room temperature before use (**20–30 minutes** out of refrigerator):

1. 3-D Food Buffer (NEOGEN item 8504)
2. Reveal 3-D for Hazelnut device (in foil pouch)

1. Add 20 mL of 3-D Food Buffer into a new 50 mL tube (18 mL for liquid samples) or a similar container.
2. Weigh out 2 g (mL) of sample and add to tube with buffer.
3. Secure the caps and invert to suspend sample. Shake vigorously or vortex for **1 minute**.
4. Tilt the tube at approximately a 45 degree angle. Dip the head of the Reveal 3-D device into the liquid. Ensure the cavity is saturated with the liquid, and that no liquid gets on the test window.
5. Leave the cavity saturated until liquid is seen running in the test window.
6. Place device on a flat surface and allow test to develop for **5 minutes**.

**NOTE:** It is essential to place the device flat on a level surface as soon as the liquid has entered the test window to stimulate flow through the device. Additionally, the devices are pre-stripped with pale green loading dye in positions T (test), O (overload) and C (control). The loading dye assists with quality and manufacturing checks and does not impact test performance. The loading dye is removed from the test window as the sample flows through the device.


## READING RESULTS

Liquid will flow into the test window. Read the result **5 minutes** after dipping.

1. **Negative result**


No line at position T (test): Level of hazelnut is below the limit of detection. (See Limitations section.)

Negative

2. **Positive result**


Any intensity of the red line at position T (test): Level of hazelnut above detection limit.


Positive

3. **High positive results**

No line is visible at position O (overload) and a red line may be faintly visible or absent at position T: Sample is overloaded with hazelnut.

High positive

High positive

4. **Invalid results**

If a red line does not appear at position C (control), then the test may be invalid

## CUSTOMER SERVICE

NEOGEN Customer Assistance and Technical Services can be reached by using the contact information on below. Training on this product, and all NEOGEN test kits, is available.

## SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of NEOGEN's test kits, on NEOGEN's website at [neogen.com](http://neogen.com), or by calling NEOGEN at 800.234.5333 or 517.372.9200.

## TERMS AND CONDITIONS

For NEOGEN's full terms and conditions, please visit [neogen.com/terms-and-conditions/](http://neogen.com/terms-and-conditions/)

## WARRANTY

NEOGEN Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, NEOGEN will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. NEOGEN shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

## TEST KITS AVAILABLE FROM NEOGEN

### Natural toxins

- Aflatoxin, DON, ergot alkaloids, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

### Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*

### Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

### Food allergens

- Almonds, coconut, crustaceans, eggs, gliadin, hazelnut, milk, mustard, peanuts, sesame, soy, walnuts

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