Read instructions carefully before starting test

Veratox[®] **for Histamine Tuna Pack**

REFRIGERATE at 2–8°C (35–46°F). DO NOT FREEZE

HISTAMINE

High levels of histamine may develop in a variety of fish species as they decompose. These species include tuna, mahi-mahi, marlin, bluefish, sardines, anchovy, bonito, herring and mackerel. Ingestion of histamine may cause scombroid poisoning in humans, which may lead to a variety of symptoms, including rash, nausea, vomiting, diarrhea, hypotension, palpitations and muscle weakness. Paralysis and death have also been reported in cases of scombroid poisoning.

Because of its potential for human illness, the U.S. Food and Drug Administration (FDA) has ruled that extensive refrigeration records and/or histamine testing must be included in a hazard analysis and critical control point (HACCP) program for relevant fish species. The FDA has set an action level of 50 parts per million (ppm) for histamine in domestic and imported fish.

INTENDED USE

Veratox for Histamine Tuna Pack is intended for the quantitative analysis of histamine in scombroid species of fish, such as tuna, bluefish and mahi-mahi, fish meal, and dry animal protein.

INTENDED USER

The test kit is designed for use by quality control personnel and others familiar with histamine analysis in fish. Since technique is very important, operators should be trained by a Neogen representative or someone who has completed Neogen training.

STORAGE REQUIREMENTS

The kit can be used until the expiration date on the label when stored refrigerated at 2–8°C (35–46°F). Do not freeze.

ASSAY PRINCIPLES

Veratox for Histamine is a competitive direct enzyme-linked immunosorbent assay (CD-ELISA) that allows the user to obtain exact concentrations of histamine in parts per million (ppm). Free histamine in the samples and controls is allowed to compete with enzyme-labeled histamine (conjugate) for the antibody binding sites. After a wash step, substrate is added, which reacts with the bound conjugate to produce blue color. More blue color means less histamine. The test is read in a microwell reader to yield optical densities. The optical densities of the controls form the standard curve, and the sample optical densities are plotted against the curve to calculate the exact concentration of histamine.

MATERIALS PROVIDED

1. 48 12-well antibody-coated microwell strips
2. 48 12-well red-marked mixing well strips
3. 24 yellow-labeled bottles (4 bottles each of 0, 2.5, 5, 10, 20 and 50 ppm histamine controls)
4. 4 blue-labeled bottles of histamine-HRP conjugate solution
5. 8 foil pouches of sample extract diluent buffer concentrate of 10 mM PBS dry powder
6. 4 foil pouches of wash buffer concentrate of 10 mM PBS-Tween
7. 4 green-labeled bottles of K-Blue[®] Substrate solution
8. 4 red-labeled bottles of Red Stop solution

MATERIALS RECOMMENDED BUT NOT PROVIDED

1. Veratox for Histamine Extraction Kit (Neogen item 9510)—Contains 38 bottles, 38 filter syringes and 38 tubes
2. Disposable bottles with 125 mL capacity (Neogen item 9398)
3. Filter syringes for histamine (Neogen item 9420H), Whatman #1 filter paper (Neogen item 9430), centrifuge or equivalent
4. Sample collection tubes (Neogen item 9421, 9421B)
5. 100 mL graduated cylinder
6. Blender
7. Test tube rack (Neogen item 9440)
8. Tubes, 15 mL, graduated (Neogen item 9474)
9. Scale capable of weighing 10–50 grams (Neogen item 9427)
10. Microwell strip reader with a 650 nm filter (Neogen item 9303)
11. Pipettor, 12-channel (Neogen item 9273)
12. Pipettor, 100 μ L (Neogen item 9272)
13. Tips for 12-channel and 100 μ L pipettors (Neogen item 9410, 9407, 9417)
14. Paper towels or equivalent absorbent material
15. Microwell holder (Neogen item 9402)
16. Timer (Neogen item 9426)
17. Waterproof marker
18. Wash bottle (Neogen item 9400)
19. 1 L bottle with lid (Neogen item 9472)
20. 2 reagent boats for 12-channel pipettor (Neogen item 9435)
21. Distilled or deionized water

PRECAUTIONS

1. Store test kit between 2–8°C (35–46°F) when not in use.
2. Do not use kit components beyond expiration date.
3. Glassware should not be used for extraction purposes. Histamine may adhere to glass, therefore using glassware may affect test results.
4. Do not mix reagents from a kit with one serial number with reagents from a kit with a different serial number.
5. Do not run more than 24 wells per test.
6. Follow proper pipetting techniques, including priming tips by filling and dispensing solution once before use.
7. Use of incubation times other than those specified may give inaccurate results.
8. Kits should be brought to 18–30°C (64–86°F) prior to use.
9. Avoid prolonged storage of kits at ambient temperatures.
10. Do not freeze test kits.
11. To avoid cross-contamination, use clean pipette tips for each sample.
12. Commodity extracts should have a pH of 6–8 before testing. Excessively acidic or alkaline samples should be adjusted. For instructions on adjusting pH contact Neogen.

PROCEDURAL NOTES

1. **Substrate.** K-Blue Substrate is ready for use. The substrate should be clear to light blue — discard if it has turned dark blue. Only pour the needed volume of substrate into a reagent boat. **Do not return unused substrate to the bottle.** Cover the reagent boat to keep the substrate protected from light until it is needed.
2. **Controls.** Six controls are provided with this kit. Neogen recommends using a combination of at least 5 controls with each assay. This combination may vary. One possible combination is to eliminate the 5 ppm control, to yield full quantitative test results in the range of 2.5 to 40 ppm. Another possibility is to eliminate the 2.5 ppm control to yield results in the range of 5 to 40 ppm. All 6 controls may be used with a microwell reader for full quantitation of histamine in the range of 2.5 to 40 ppm. If needed, contact Neogen for more information on the use of controls.
3. **Sample extract diluent buffer.** Prepare by adding foil pouch of extract buffer (10 mM PBS) to 1 L of distilled or deionized water. Swirl to mix. Label bottle appropriately. Store remaining buffer covered at room temperature.
4. **Wash buffer.** Prepare by adding one foil pouch of wash buffer concentrate (10 mM PBS-Tween) to 1 L of distilled or deionized water. Swirl to mix. Do not shake. Label bottle appropriately. Store remaining wash buffer at room temperature.
5. **Antibody wells.** Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after the samples are extracted, and the test procedure is set to begin.

SAMPLE PREPARATION

- A. **Canned and pouched tuna: AOAC 937.07b**
Place entire contents of can or pouch (meat and liquid) into a blender. Blend until homogenous. Store samples at 2–8°C (35–46°F) until analyzed.
- B. **Fresh or thawed frozen raw fish: AOAC 937.07a**
Clean and eviscerate three fish. Cut three cross-sectional pieces 2.5 cm (1 inch) thick, from back of the pectoral fin, halfway to vent and one posterior to the vent. Debone slices and blend or grind combined samples until homogenous. Store samples at 2–8°C (35–46°F) until analyzed.
- C. **Dry Samples (example: Fishmeal)**
The sample to be tested should be collected according to accepted sampling techniques and be representative of the lot. After homogenization, grind part of the sample (minimum of 200 g) so at least 95% of the ground material passes through a 20 mesh sieve (about the particle size of fine espresso).

SAMPLE EXTRACTION

If using Neogen's Veratox for Histamine Extraction Kit (Neogen item #9510), follow the instructions in that kit. Otherwise, continue with the instructions that follow. **NOTE:** Glassware should not be used for extraction purposes as histamine may adhere to the glass.

Liquid/Wet Samples

1. Add **10 g** of the homogenous mixture to a clean disposable extraction bottle containing **90 mL** of distilled or deionized water.
2. Tightly cap and vigorously shake the bottle for **15-20 seconds** in order to suspend the fish tissue in the water.
3. Wait approximately **5 minutes**, then shake the bottle for **15-20 seconds** in order to re-suspend the fish tissue.
4. Wait an additional **5 minutes**, and again shake the bottle for **15-20 seconds** in order to re-suspend the fish tissue. Allow the tissue to settle to the bottom of the bottle for about **30 seconds**.
5. If necessary, filter the contents through folded filter paper or a histamine filter syringe (Neogen #9420H) into a clean container. The sample is now ready for extract dilution.

ALTERNATIVE: Centrifuge the sample and use the clear supernatant as the sample for extract dilution.

Dry Samples

1. Add **10 g** of the homogenous mixture to a clean disposable extraction bottle containing **100 mL** of distilled or deionized water.
2. Tightly cap and vigorously shake the bottle for **15-20 seconds** in order to suspend the sample in the water.
3. Wait approximately **5 minutes**, then shake the bottle for **15-20 seconds** in order to re-suspend the sample.
4. Wait an additional **5 minutes**, and again shake the bottle for **15-20 seconds** in order to re-suspend the sample. Allow the sample to settle to the bottom of the bottle for about **30 seconds**.
5. If necessary, filter the contents through folded filter paper or a histamine filter syringe (Neogen #9420H) into a clean container. The sample is now ready for extract dilution.

ALTERNATIVE: Centrifuge the sample and use the clear supernatant as the sample for extract dilution.

SAMPLE EXTRACT DILUTION

Because of the sensitivity of the ELISA format, the sample extract must be diluted before performing the test (see **Procedural Note 3** for preparing the sample extract diluent buffer).

1. Add 10 mL of sample extract diluent buffer to a clean test tube or bottle.
2. Using a clean pipette tip, add 100 μ L of the extract to the sample extract diluent buffer. Gently swirl to mix.
3. The sample is now ready to test. Repeat for all samples.

NOTE: For optimal performance with expected results greater than 40 ppm, a further dilution of the sample extract using extract diluent buffer is recommended. Begin with a 1:1 dilution, repeating dilution and retesting if necessary until the new result (prior to application of the dilution factor) is less than 40 ppm. Corrected values require multiplying by 2x for total X repetitions of these additional 1:1 dilutions. Samples should be tested within 4 hours of extraction.

TEST PROCEDURE

Allow all reagents to warm to 18–30°C (64–86°F) prior to use.

1. Remove 1 red-marked mixing well for each sample to be tested plus 5 red-marked wells for controls, and place in the well holder.
2. Remove an equal number of antibody-coated wells. Return antibody wells that will not be used immediately to the foil pack with desiccant. Reseal the foil pack to protect the antibody. Mark one end of strip with a “1” and place strip in the well holder with the marked end on the left. Do not mark on the inside or bottom of the wells.
3. Mix each reagent by swirling the reagent bottle prior to use.
4. Add 100 μ L of conjugate from the blue-labeled bottle to each red-marked mixing well.
5. The controls (see **Procedural Note 2**) are supplied ready to use—do not dilute. Using a new pipette tip for each, transfer 100 μ L of controls and diluted samples to the red-marked mixing wells as shown:

If concerned about histamine levels between 2.5 and 40 ppm:

0	2.5	10	20	50	S1	S2	S3	S4	S5	S6	S7	Strip 1
S8	S9	S10	S11	S12	S13	S14	S15	S16	S17	S18	S19	Strip 2

If concerned about histamine levels between 5 and 40 ppm:

0	5	10	20	50	S1	S2	S3	S4	S5	S6	S7	Strip 1
S8	S9	S10	S11	S12	S13	S14	S15	S16	S17	S18	S19	Strip 2

6. Using a 12-channel pipettor, mix the liquid in the wells by pipetting up and down 3 times. Transfer 100 μ L to the antibody-coated wells.
7. Incubate for **10 minutes** at room temperature 18–30°C (64–86°F) mixing for the first **10–20 seconds** by sliding the microwell holder back and forth on a flat surface without splashing

reagents from the wells. Discard red-marked mixing wells.

8. Shake out the contents of the antibody wells. Fill each antibody well with diluted wash buffer and dump them out. Repeat this step 3 times, then turn the wells upside down and tap out the remaining liquid on an absorbent towel.
9. Pour the needed volume of substrate from the green-labeled bottle into the green-labeled reagent boat, and with new tips on the 12-channel pipettor, pipette 100 μ L of substrate into the wells.
10. Incubate for **10 minutes** at room temperature 18–30°C (64–86°F) mixing for the first **10–20 seconds** by sliding back and forth on a flat surface. Discard remaining substrate and rinse the reagent boat with water.
11. Pour Red Stop solution from the red-labeled bottle (same volume as substrate) into the red-labeled reagent boat. Using the same pipette tips on the 12-channel pipettor as were used to dispense substrate, add 100 μ L Red Stop to each well and mix by sliding back and forth on a flat surface. Discard the tips.
12. Wipe the bottom of the microwells with a dry cloth or towel and read in a microwell reader using a 650 nm filter and Neogen Veratox software comparing the standard curve to calculate results. Air bubbles should be eliminated, as they could affect analytical results. Results should be read within 20 minutes of completion of the test.
13. All materials are safe to be disposed of in the trash.

RETESTING

If positives occur in commodities not previously tested, confirm with an additional approved method prior to taking action. Results between 2.5 and 40 ppm are quantitative. For optimal performance, results of greater than 40 ppm should be diluted 1:1 with extract diluent buffer as described in the **Sample Extract Dilution** section.

PERFORMANCE CHARACTERISTICS

Limit of detection: 2 ppm (determined by the mean average of 10 histamine free samples plus 3 standard deviations)

Limit of quantitation: 2.5 ppm (described as the lowest concentration point on the calibration curve that this test can reliably detect histamine)

Range of quantitation: 2.5–40 ppm (for optimal performance, results greater than 40 ppm require a further dilution; see **Sample Extract Dilution** for further instructions.)

Validated matrices: Fresh, canned and pouched tuna in oil or water*, mahi-mahi, bluefish and fish meal.

*AOAC-R1 validated matrix

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at www.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/Corporate/termsconditions.html.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TEST KITS AVAILABLE FROM NEOGEN

Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, crustaceans, eggs, gliadin, hazelnut, lupine, milk, mustard, peanuts, sesame, soy, walnuts, multi-treenut

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