U.S. English

SKU #: 700002592

REF #: 8510

Veratox[®]

for Gliadin R5

Refrigerate at 2°C–8°C (35°F–46°F). Do not freeze.

Veratox[®] for Gliadin R5

SKU #: 700002592 | REF #: 8510

Gliadin/Gluten

Gliadins are a class of alcohol-soluble proteins found in wheat that belong to a group of proteins called prolamins. Other prolamins include secalin, found in rye, and hordein, found in barley. Gluten consists of two groups of proteins (prolamins and glutelins) that are found in differing amounts in wheat, barley, rye, and oats.

Gliadin and other prolamins have been identified as major causal agents in a number of disorders, including wheat allergy and gluten intolerance (celiac disease). Wheat allergy is a specific immune response to a number of wheat proteins, including gliadin, albumin, globulin, and glutenin. Celiac disease is a chronic reaction to gluten proteins that results in the poor absorption of nutrients in the small intestine.

Those who must avoid gluten rely on the correct labeling of food to make appropriate, safe food choices. Testing for the presence of gluten components is a vital part of a food manufacturer's allergen control plan.

Intended Use

Veratox[®] for Gliadin R5 is intended for the quantitative analysis of ingredients, clean-in-place solutions, and finished food products intended to be gluten-free. The presence of gliadin and prolamins of wheat, barley, and rye can be used as a measure of gluten contamination. The detection of wheat gluten in rice flour, cereal, bread, and cooked hamburger was validated by the AOAC Research Institute (AOAC-RI #061201). Clean-in-place solutions and other foods were validated internally by Neogen[®].

Veratox for Gliadin R5 has not been evaluated with all possible food products, food processes, or testing protocols, or with all possible sources of gliadin. Neogen has not documented this product for use in industries other than food and beverage. For example, Neogen has not documented Veratox for Gliadin R5 for testing in water, pharmaceuticals, or cosmetics. Do not use this product for the diagnosis of conditions in humans or animals.

Intended User

This test kit is designed for use by quality control personnel and others familiar with foods possibly contaminated by gliadin. Since technique is very important, operators should be trained by a Neogen representative or someone who has completed the Neogen training.

User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at neogen.com or contact your local Neogen representative or distributor for more information. As with all test methods used for food analysis, the test matrix can influence the results. When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results. The food sample itself may influence results. It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples to satisfy the user that the chosen test method meets the user's criteria. It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements. As with any test method, results obtained from use of any Neogen product do not constitute a guarantee of the quality of the matrices or processes tested.

Assay Principles

Veratox for Gliadin R5 is a sandwich enzyme-linked immunosorbent assay (S-ELISA). Gliadin is extracted from samples with a 60% ethanol solution by shaking in a shaker or rotator. The extract is diluted in Phosphate Buffered Saline (PBS) and diluted samples are added to R5 antibody-coated wells (capture antibody), where gliadin will bind to the antibody during an incubation period. The detector antibody, which is a conjugated R5 antibody, is then added. The detector antibody binds to captured gliadin during another incubation period. Unbound enzyme-labeled antibody is washed away and a one-step substrate is added. The substrate is chemically altered by the enzyme and a blue color develops as a result. A stopping reagent is added and the color of the solution is observed. Blue indicates samples containing high levels of gliadin while purple or red samples contain little or no gliadin. The optical densities of the controls form a standard curve, and the sample optical densities are plotted against the curve to calculate the concentration of gliadin in parts per million (ppm).

Storage Requirements

The kit can be used until the expiration date on the label when stored refrigerated at 2°C–8°C (35°F–46°F).

Materials Provided

1. 48 antibody-coated microwells
2. 48 red-marked transfer wells
3. 6 yellow-labeled bottles of 0, 2.5, 5, 10, 20, and 40 ppm gliadin controls (5, 10, 20, 40, and 80 ng/mL gliadin)
4. 2 blue-labeled bottles of enzyme-labeled antibody conjugate
5. 1 green-labeled bottle of K-Blue® Substrate
6. 1 red-labeled bottle of Red Stop solution
7. 40 mL of 10 mM PBS-Tween washing reagent in a widemouth bottle; each bottle can prepare 1 L in distilled or deionized water (pH 7.4)
8. 1 foil pouch sample diluent concentrate of 10 mM PBS dry powder containing enough powder to prepare 1 L of dilution buffer
9. 50 g of extraction additive in 2 specimen cups
10. Plastic scoop to measure extraction additive

Materials Recommended, Not Provided

1. Allergen extraction kit (Neogen item 700002568)
 - a. 20 disposable plastic extraction bottles
 - b. 20 sample collection tubes (12 x 75 mm) with caps
2. Orbital rotator or shaker to hold 50 cc centrifuge tubes for a 1 g sample, or shaker water bath with clamps adjusted to hold 125 mL (4 oz) extraction bottles for a 2 g sample
3. Oven or water bath adjustable to 50°C if analyzing heat-processed samples
4. Scale capable of weighing 0.25 ± 0.01 g if using gliadin renaturing cocktail solution to analyze heat-processed samples
5. Gliadin renaturing cocktail solution for heat processed samples (Neogen items 700002594, 700002595)
6. Laboratory-grade ethanol (190 proof)
7. Microwell reader with a 650 nm filter (Neogen item 700003947)
8. Pipettor, 50–200 µL adjustable (Neogen item 100001331)
9. Pipettor, 12-channel (Neogen item 700006487)
10. Pipette tips (Neogen items 100001353, 700004009, 100001352)
11. Two 1 L bottles to prepare washing solution and sample extract dilution solution (Neogen item 700006489)
12. Test tubes to perform sample extract dilution
13. Timer (Neogen item 100001358)
14. 3 reagent boats for 12-channel pipettor (Neogen item 700002730)
15. Microwell strip holder (Neogen item 100001351)
16. Wash bottle (Neogen item 700006491)
17. Paper towels or equivalent absorbent material
18. Waterproof marker
19. Distilled or deionized water
20. Centrifuge (optional)
21. Vortex (Neogen item 100001381)

Precautions

1. This method uses the R5 antibody, which is aligned with the Type I method for gluten quantification as defined by the Codex Alimentarius. As is the case with all methods that use this antibody, Veratox Gliadin R5 demonstrates an overestimation of gluten from rye and barley sources compared to the wheat standard.
2. Ethanol solution is highly flammable. Keep container tightly closed, and keep away from heat, sparks, open flame, and those smoking. It is toxic if swallowed, or if vapor is inhaled. Avoid contact with skin.
3. Components of Veratox for Gliadin R5, such as controls and extraction additive, may contain one or more of the following potentially allergic materials: gluten, casein, and soy protein. If allergic to any of these compounds, use caution when using this product.
4. Concentrated food additives, colors, and flavors may cause interferences on ELISA test methods. Contact Neogen Technical Services for validation information.
5. Sponges should not be used for sample collection and gliadin testing. Sample collection swabs other than Neogen swabs should be validated prior to use. General sponges and swabs may contain solutions or materials that may interfere with the test kits.
6. Store test kit between 2°C–8°C (35°F–46°F) when not in use. Do not freeze test kits, and avoid prolonged storage at ambient temperatures.
7. Bring kits to room temperature (18°C–30°C, 64°F–86°F) prior to use.
8. Do not use kit components beyond expiration date.
9. Do not mix reagents from one kit lot with reagents from a different lot.
10. Do not run more than 24 wells per test.
11. Follow proper pipetting techniques (e.g., prime tips and use clean tips).
12. Use only incubation times specified; others may give inaccurate results.
13. Use clean pipette tips and glassware for each sample to avoid cross-contamination. Thoroughly wash all glassware between samples.

Procedural Notes

1. **Sample Extract Dilution Solution (PBS)**
Prepare extract dilution solution by adding a foil pouch of sample dilution concentrate, 10 mM PBS, to 1 L distilled or deionized water. Swirl to mix thoroughly.
2. **Controls**
Six controls containing a known amount of gliadin are provided with this kit to create a standard curve which is utilized when measuring gliadin concentration in a sample. Neogen recommends using all 6 or a combination of at least 5 controls with each assay. This combination may vary. One possible combination is to eliminate the 2.5 ppm control to yield results in the range of 5–40 ppm. If needed, contact your Neogen representative for more information on the use of these controls.
3. **Wash Buffer**
Prepare the wash buffer solution by pouring the wash buffer concentrate into a 1 L container. Add 960 mL of distilled or deionized water. Swirl to ensure thorough mixing.
Note: Discard unused portions of extract dilution solution and wash buffer when the test kit has been used completely.
4. **Extraction Additive**
To be used with the provided scoop with all the samples utilizing extraction procedure A or B (non-heat-processed samples). Heat-processed samples (procedure C or D) should only use the extraction additive if samples contain buckwheat, chestnut flour, or tannins/phenolic compounds such as chocolate, coffee, cocoa, wine, herbs, or fruits.
5. **Substrate**
Substrate is ready for use. The substrate should be clear to light blue — discard if it has turned dark blue. Only pour the needed volume of substrate into a reagent boat. Do not return unused substrate to the bottle. Cover the reagent boat to keep the substrate protected from light until needed.
6. **Antibody Wells**
Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after samples are extracted and the test procedure is set to begin.

Sample Preparation and Extraction

When testing samples that are of unknown origin or are heat processed, use either extraction method C or D. Extraction method D is the methodology used with AOAC Performance Tested Method RI #061201. When testing samples that are not heat processed, use extraction method A or B.

- A. Extraction of non-heat-processed samples with orbital shaker or rotator
1. Prepare 60% ethanol extraction solution by combining 6 parts ethanol with 4 parts distilled water. Prepare sample extract dilution solution (PBS) as detailed in procedural note 1.
 2. Add 1 g ground sample, or 1 mL liquid sample, to a clean 50 cc tube.
 3. Add 1 level scoop of extraction additive to the tube.
 4. Add 10 mL (9 mL for liquid samples) of 60% ethanol to the tube, cap tightly, then shake the tube vigorously by hand for about 1 minute, or vortex for 30 seconds, to ensure complete mixing.
 5. Extract by shaking (150 rpm) in an orbital shaker or rotator by laying the tube down on its side over the flat pad of the instrument and holding it tightly using a rubber band or tape. Rotate or shake for 10 minutes at room temperature.
 6. Centrifuge sample (if necessary such as highly viscous or high particulate powders) for 10 minutes at > 2500 g at room temperature.
 7. Dilute each sample 1:50 by withdrawing 100 μ L of the upper layer of the extract and transferring it to a small tube or vial containing 4.9 mL of sample extract dilution solution (PBS).
 8. To mix, vortex the tube for 5 seconds.
 9. Test diluted samples within 2–3 hours of extraction.
- B. Extraction of non-heat-processed samples with shaker or shaker water bath
1. Prepare 60% ethanol extraction solution by combining 6 parts ethanol with 4 parts distilled water. Prepare sample extract dilution solution (PBS) as detailed in procedural note 1.
 2. Add 2 g ground sample, or 2 mL liquid sample, to a 125 mL clean extraction bottle.
 3. Add 1 level scoop of extraction additive to the bottle.
 4. Add 20 mL (18 mL for liquid samples) of 60% ethanol, cap the bottle tightly, then shake vigorously by hand for about 20 seconds to ensure complete mixing.
 5. Extract by shaking (150 rpm) in a shaker for 10 minutes at room temperature (a shaker water bath can work, but do not turn the heat on). Remove the bottle from shaker or bath.
 6. Centrifuge sample (if necessary such as highly viscous or high particulate powders) for 10 minutes at > 2500 g at room temperature.
 7. Dilute each sample 1:50 by withdrawing 100 μ L of the upper layer of the extract and transferring it to a small tube or vial containing 4.9 mL of sample extract dilution solution (PBS).
 8. To mix, vortex the tube for 5 seconds.
 9. Test diluted samples within 2–3 hours of extraction.
- C. Extraction of heat-processed or unknown origin commodities
- This method requires the gliadin renaturing cocktail as described in the Materials Recommended, Not Provided Section of the instructions.
1. Prepare 80% ethanol extraction solution by combining 8 parts ethanol with 2 parts distilled water.
 2. Prepare sample extract dilution solution (PBS) as detailed in procedural note 1.
 3. Add 0.25 g sample into a 50 cc screw cap centrifuge tube.
 4. Add 2.5 mL of renaturing cocktail solution.
 5. If samples are of unknown composition or are known to contain buckwheat, chestnut flour, or tannins/phenolic compounds such as chocolate, coffee, cocoa, wine, herbs, or fruits, add 1 level scoop of extraction additive. If sample being tested is known not to contain the above ingredients, it is not necessary to use additive.
 6. Cap and vortex 30 seconds to homogenize cocktail and sample.
 7. Incubate 40 minutes at 50°C (water bath or oven).
 8. Remove samples and let cool for 5–10 minutes.
 9. Add 7.5 mL of 80% ethanol and vortex again for 10–20 seconds.
 10. Shake (150–200 rpm) for 1 hour at room temperature on a rotator (tube on its side).

11. Centrifuge sample (if necessary such as highly viscous or high particulate powders) for 10 minutes at > 2500 g at room temperature.
 12. Dilute the sample 1:12.5 into PBS (200 µL sample into 2.3 mL PBS).
 13. To mix, vortex the tube for 5 seconds.
 14. Test diluted samples within 2–3 hours of extraction.
- D. Extraction of heat-processed or unknown origin commodities (AOAC RI PTM™ #061201 method). AOAC method must be run on the StatFax Microwell Reader. This procedure requires the gliadin renaturing cocktail solution described in method C.
1. Prepare 80% ethanol extraction solution by combining 8 parts ethanol with 2 parts distilled water.
 2. Prepare sample extract dilution solution (PBS) as detailed in procedural note 1.
 3. Add 1 g sample into an extraction bottle.
 4. Add 10 mL of renaturing cocktail solution.
 5. If samples are of unknown composition or are known to contain buckwheat, chestnut flour, or tannins/phenolic compounds such as chocolate, coffee, cocoa, wine, herbs, or fruits, add 1 level scoop of extraction additive. If sample being tested is known not to contain the above ingredients, it is not necessary to use additive.
 6. Cap and vortex 30 seconds to homogenize cocktail and sample.
 7. Incubate 40 minutes at 50°C (water bath or oven).
 8. Remove samples and let cool for 5–10 minutes.
 9. Add 30 mL of 80% ethanol and vortex or manually shake for 10-20 seconds.
 10. Shake (150–200 rpm) for 1 hour at room temperature on a rotator (tube on its side).
 11. Centrifuge sample (if necessary such as highly viscous or high particulate powders) for 10 minutes at > 2500 g at room temperature.
 12. Dilute the sample 1:12.5 into PBS (200 µL sample into 2.3 mL PBS).
 13. To mix, vortex the tube for 5 seconds.
 14. Test diluted samples within 2–3 hours of extraction.

Test Procedure

Allow the test kit and all reagents to warm to room temperature (18°C–30°C, 64°F–86°F) before using.

1. Remove 1 red-marked mixing well for each sample to be tested plus 1 red-marked well for each control, and place in the well holder.
2. Remove an equal number of antibody-coated wells. Return antibody wells that will not be used immediately to the foil pack with desiccant. Reseal the foil pack to protect the antibody. Mark one end of the strip with a “1,” and place strip in the well holder with the marked end on the left.
3. Mix each reagent by swirling the reagent bottle prior to use.
4. The controls (see procedural note 2) are supplied ready to use — do not dilute. Using a new pipette tip for each, transfer 150 µL of controls and diluted samples to the red-marked mixing wells as shown in one of the templates below. Only run up to two 12-well strips at a time.

If using all 6 controls for a range of 2.5–40 ppm (AOAC method):

0	2.5	5	10	20	40	S1	S2	S3	S4	S5	S6
S7	S8	S9	S10	S11	S12	S13	S14	S15	S16	S17	S18

If using 5 controls for a range of 5–40 ppm:

0	5	10	20	40	S1	S2	S3	S4	S5	S6	S7
S8	S9	S10	S11	S12	S13	S14	S15	S16	S17	S18	S19

5. Place tips on the 12-channel pipettor and transfer 100 µL of the controls and sample extracts to the antibody-coated wells. Mix for 20 seconds by sliding the well holder back and forth on a flat surface.
6. Incubate microwells 10 minutes at room temperature (18°C–30°C, 64°F–86°F). Discard the red-marked transfer wells.
7. Empty the contents of the wells into a sink. With a wash bottle, fill each antibody well with the wash buffer solution and dump out. Repeat the washing 5 times, then turn the wells upside down and tap out on a paper towel until all washing solution is removed.
8. Pour the needed volume of conjugate from the blue-labeled bottle into a clean reagent boat.

9. Using the 12-channel pipettor and new tips, transfer 100 μ L of the conjugate into all the wells and mix for 20 seconds by sliding the well holder back and forth on a flat surface.
10. Incubate for 10 minutes at room temperature (18°C–30°C, 64°F–86°F).
11. Wash all wells with the wash buffer solution as described in step 7.
12. Pour the needed volume of substrate solution from the green-labeled bottle into a clean reagent boat.
13. Place new tips on the 12-channel pipettor and transfer 100 μ L of substrate into each well and mix for 20 seconds. Do not eject tips.
14. Incubate for 10 minutes at room temperature (18°C–30°C, 64°F–86°F).
15. Pour the needed volume of Red Stop solution from the red-labeled bottle into a clean reagent boat.
16. With the same tips used to dispense the substrate, transfer 100 μ L of Red Stop into each well and mix for 20 seconds.
17. Wipe the bottom of the microwells and read in a microwell reader with a 650 nm filter.
18. Interpret the test results using the StatFax 4700 microwell reader (AOAC RI PTM #061201) or an equivalent strip reader. If using an equivalent reader, it is recommended to calculate the results using Neogen's Veratox Software for Windows.

Interpretation of Results

Standard controls were made from wheat gliadin and calculated as gliadin. Approximately 50% of the gluten is available as gliadin. Therefore, to calculate the gluten value of the samples, multiply the ppm gliadin results by 2.

Performance Characteristics

Limit of quantitation (LOQ) in ppm gliadin: From AOAC validation: Rice flour, 0.84 ppm; soy flour, 0.69 ppm; cooked hamburger, 1.22 ppm; oat cereal, 0.35 ppm; cornbread (bread), 0.24 ppm. Neogen internally validated commodity LOQs were all under 2.5 ppm (described as the lowest concentration point on the calibration curve that this test can reliably measure gliadin).

Limit of detection (LOD) in ppm wheat gliadin: From AOAC validation: Rice flour, 0.62 ppm; soy flour, 0.64 ppm; cooked hamburger, 0.79 ppm; oat cereal, 0.22 ppm; cornbread, 0.17 ppm. Neogen internally validated commodity LODs were all under 2.5 ppm (defined as the mean value of 10 negative samples + 3.3 standard deviations).

Range of quantitation: 2.5–40 ppm gliadin. (For quantitating samples above 40 ppm, contact a Neogen representative for dilution instructions).

Customer Service

Neogen Customer and Technical Services can be contacted through neogen.com and product training is available by request.

Safety Data Sheets (SDS) Information Available

SDS are available for all test kits at neogen.com or by calling 800.234.5333 or 517.372.9200.

Terms and Conditions

Neogen's full terms and conditions are available [online](#).

Warranty

Neogen makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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