*Read instructions carefully before starting test*

# for Total Milk Allergen

## Quantitative Test

*Refrigerate at 2–8°C (35–46°F) • Do not freeze*

### **MILK ALLERGEN**

Food allergens are proteins in food that can create an immune response in sensitive individuals. Once ingested, food allergens can cause a number of reactions, ranging in severity from hives and itching to anaphylaxis. Anaphylaxis is a severe allergic reaction, involving vomiting, diarrhea, difficulty breathing, swelling of the mouth and tongue, and a rapid drop in blood pressure.

An estimated 3.5 to 4 percent of adults, and 6 to 8 percent of children, are sensitive in some degree to food allergens. More than 12 million people in the United States alone are known to have a food allergy.

Food manufacturers protect those with food allergies by clearly labeling their products with a list of ingredients. Testing for the presence of milk components assures food manufacturers that an unlabeled — and potentially dangerous — ingredient did not make its way into a food product.

### **INTENDED USE**

Veratox for Total Milk Allergen is intended for the quantitative analysis of food products, such as juices, cake mixes, cookies, sauces and sorbets, and clean-in-place rinses for the presence of milk proteins casein and whey.

### **INTENDED USER**

The Veratox for Total Milk Allergen test kit is designed for use by quality control personnel and others familiar with foods possibly contaminated by casein or whey. Since technique is very important, operators should be trained by a Neogen representative or someone who has successfully completed the Neogen training.

### **ASSAY PRINCIPLES**

Veratox for Total Milk Allergen is a sandwich enzyme-linked immunosorbent assay (S-ELISA). Milk protein residue is extracted from samples with a buffered salt solution (Phosphate Buffered Saline, or PBS) by shaking in a heated water bath. Milk protein residue is extracted from samples and added to antibody-coated wells (capture antibody) where it binds to the antibody during an incubation. Any unbound protein residue is washed away and a second antibody (detector antibody), which is enzyme labeled, is added. The detector antibody binds to the already bound milk protein residue. After a second wash, the substrate is added. Color develops as a result of the presence of bound detector antibody. Red Stop reagent is added and the color of the resulting solution is observed. The test is read in a microwell reader to yield optical densities. The optical densities of the controls form a standard curve, and the sample optical densities are plotted against the curve to calculate the exact concentration of allergen.

### **STORAGE REQUIREMENTS**

The kit can be used until the expiration date on the label when stored refrigerated at 2–8°C (35–46°F). Do not freeze.

## MATERIALS PROVIDED

1. 48 antibody-coated microwells (24 per pouch)
2. 48 red-marked transfer wells (24 per pouch)
3. 5 yellow-labeled bottles of 0, 2.5, 5, 10 and 25 ppm nonfat dry milk (NFDM) controls
4. 4 blue-labeled bottles of enzyme-labeled antibody conjugate
5. 1 green-labeled bottle of K-Blue® Substrate
6. 1 red-labeled bottle of Red Stop solution
7. 5 foil pouches of 10 mM PBS dry powder extraction solvent; each pouch contains enough powder to prepare 1 L of extraction solvent
8. 2 wide-mouth bottles of 40 mL PBS-Tween wash buffer concentrate; each bottle contains enough concentrate to prepare 1 L of wash buffer
9. 50 g of extraction additive in a specimen cup
10. Plastic scoop to measure extraction additive

## MATERIALS RECOMMENDED BUT NOT PROVIDED

1. Allergen Extraction Kit (Neogen item 8429)
  - a. 20 disposable plastic extraction bottles
  - b. 20 disposable transfer pipettes
2. Shaker water bath adjusted to 60°C ± 1° (140°F) with holder for extraction bottles (Neogen items 9298, 9299)
3. Pipettor, adjustable 50–200 µL (Neogen item 9276)
4. Pipette tips (Neogen item 9410, 9407, 9417)
5. Pipettor, 12-channel (Neogen item 9273)
6. Timer (Neogen item 9452, 9426)
7. Microwell strip reader with a 650 nm filter (Neogen item 9303)
8. 1 L bottle to prepare washing solution (Neogen item 9472)
9. 1 L heat safe bottle to prepare extract solution (Neogen item 9472)
10. Paper towels or equivalent absorbent material
11. Microwell holder (Neogen item 9402)
12. Waterproof marker
13. Wash bottle (Neogen item 9400)
14. Graduated cylinder capable of measuring 125 mL (Neogen item 9368)
15. Scale capable of weighing 5 g ± 0.1 g (Neogen item 9427)
16. 3 reagent multichannel pipettor boats (Neogen item 9435)
17. Distilled or deionized water

## PRECAUTIONS

1. Samples intended to be tested for **milk must be extracted separately** from samples intended to be tested for other food allergens, such as peanut and egg residues. The extraction additives for each type of test are designed specifically for the target food allergen.
2. The controls and extraction reagents of the Veratox for Total Milk Allergen test kit may contain one or more of the following potentially allergic materials: milk, egg protein, peanut protein, soy protein, or tree nut protein. If allergic to any of these compounds, please use caution when handling this product.
3. Concentrated food additives, colors and flavors may cause interferences on ELISA test methods. Contact Neogen Technical Services for validation information.
4. Hydrolyzed and fermented proteins may not be detected using ELISA methods for allergen testing. Although the proteins may be undetectable in the assay because of their properties, there still could be active allergenic protein residue present.
5. Infant formula: Contact Neogen for additional information concerning the testing of infant formulas.
6. Sponges should not be used for sample collection and allergen testing. Sample collection swabs other than Neogen swabs should be validated prior to use. General sponges and swabs may contain solutions or materials that may interfere with the test kit.
7. Disposable extraction bottles and tubes must be used to extract milk samples to avoid cross-contamination.
8. The testing area must be totally free of milk products. A minute amount of milk protein in the environment can affect test results.
9. Store test kit between 2–8°C (35–46°F) when not in use. Do not freeze test kits.
10. Do not use kit components beyond expiration date.
11. Sample extracts must be cooled to room temperature 18–30°C (64–86°F) prior to use.
12. Bring kits to room temperature between 18–30°C (64–86°F) prior to use.
13. Do not mix reagents from kit serial with reagents from a different kit serial.
14. Do not run more than 24 wells per test.
15. Follow proper pipetting techniques (e.g., prime tips and use clean tips).

16. Use only the incubation times specified. Others may give inaccurate results.
17. Avoid prolonged storage of kits at ambient temperatures.

## PROCEDURAL NOTES

1. **Substrate:** K-Blue Substrate is ready for use. The substrate should be clear to light blue — discard if it has turned dark blue. Only pour the needed volume of substrate into a reagent boat. **Do not return unused substrate to the bottle.** Cover the reagent boat to keep the substrate protected from light until it is needed.
2. **Antibody wells:** Keep wells sealed in the foil pouch until needed. Remove wells from the foil pouch only after samples are extracted, and the test procedure is set to begin.
3. **Extraction solution:** Prepare extraction solution by adding a foil pouch of extraction solvent, 10 mM PBS, to 1 L distilled or deionized water at pH 7.4. Cover and store unused portion, refrigerated at 2–8°C (35–46°F).
4. **Wash buffer:** Prepare wash buffer solution by mixing 40 mL of wash buffer concentrate in wide mouth bottle into 960 mL of distilled or deionized water at pH 7.4. Cover and store unused portion refrigerated at 2–8°C (35–46°F).  
**NOTE:** Discard unused portions of extraction solution and wash buffer when test kit has been used completely.

## SAMPLE PREPARATION AND EXTRACTION

The sample to be tested should be collected according to accepted sampling techniques (see Neogen's Food Allergen Handbook). The sample should be ground and thoroughly mixed prior to proceeding with the extraction procedure.

**NOTE:** Glassware and material used for other allergen testing cannot be used for milk residue testing due to the potential of cross-contamination. For this reason, it is also highly recommended that any labware, such as wash bottles, graduated cylinders, and 1 L bottles, be solely dedicated for use with the milk allergen kit.

1. Prepare the extraction solution (PBS, pH 7.4) as described in the procedural notes.
2. Preheat extraction solution to 60°C (140°F) by immersing the bottle containing the solution into the water bath and allowing it to reach 60°C.
3. Using your sampling and collection procedure, obtain a representative sample and grind it to a very fine particle size.
4. Transfer 5 g of sample, or 5 mL of liquid sample, into a 250 mL disposable plastic extraction bottle.
5. Add 1 level scoop of the extraction additive into the sample bottle. (Do not use the extraction additive from another allergen test kit.)
6. Pour 125 mL of the 60°C (140°F) extraction solution into the sample bottle.
7. Cap the sample bottle to prevent contents from splashing during the extraction.
8. Extract by shaking (150 rpm) in a 60°C water bath for **15 minutes**. Remove the bottle from the bath.
9. Let the material settle for **5 minutes** before proceeding to the next step.
10. Use the supernatant (the top liquid portion of the extract) as your sample. **Do not filter.** Begin the test procedure once the sample has cooled to room temperature (at least **15 minutes**).

## TEST PROCEDURE FOR QUANTITATION

Allow the test kit and all reagents to warm to room temperature 18–30°C (64–86°F) before using.

1. Mix each reagent by swirling the reagent bottle prior to use.
2. Using a new pipette tip for each, transfer 150 µL of controls and sample extracts to the red-marked transfer wells as shown in the template below. Run only two 12-well strips at a time.

0	2.5	5	10	25	S1	S2	S3	S4	S5	S6	S7	Strip 1
S8	S9	S10	S11	S12	S13	S14	S15	S16	S17	S18	S19	Strip 2

3. Place tips on the 12-channel pipettor and transfer 100 µL of the controls and sample extracts to the antibody-coated wells. Mix for **10 seconds** by sliding the microwell holder back and forth on a flat surface without splashing reagents from the wells.
4. Incubate microwells **10 minutes** at room temperature 18–30°C (64–86°F). Discard the red-marked transfer wells.
5. Empty the contents of the wells into a sink. With a wash bottle, fill each antibody well with the wash buffer solution and dump out. Repeat the washing 10 times, then turn the wells upside down and tap out on a paper towel until the remaining washing solution is removed.
6. Pour the needed volume of conjugate from the blue-labeled bottle into a clean reagent boat.
7. Using the 12-channel pipettor transfer 100 µL of the conjugate into all the wells and mix for **10 seconds** by sliding the microwell holder back and forth on a flat surface.
8. Incubate for **10 minutes** at room temperature 18–30°C (64–86°F).
9. Wash all wells with the wash buffer solution as described in step 5.
10. Pour the needed volume of substrate solution from the green-labeled bottle into a clean reagent boat.

- Place new tips on the 12-channel pipettor and transfer 100 µL of substrate into each well and mix for **10 seconds** by sliding the microwell holder back and forth on a flat surface.
- Incubate for **10 minutes** at room temperature 18–30°C (64–86°F).
- Pour the needed volume of Red Stop solution from the red-labeled bottle into a clean reagent boat.
- Place new tips on the 12-channel pipettor and transfer 100 µL of Red Stop solution into each well and mix for **10 seconds** by sliding the microwell holder back and forth on a flat surface.
- Wipe the bottom of the microwells with a dry cloth or towel and read in a microwell reader with a 650 nm filter. Air bubbles should be eliminated, as they could affect analytical results. Results should be read within **20 minutes** after the addition of Red Stop Solution.
- Read and calculate the test's results using Neogen's Stat-Fax microwell reader, or an equivalent strip reader. If using a strip reader, calculate the results using Neogen's Veratox software for windows.

### TEST PROCEDURE FOR SCREENING

Allow the test kit and all reagents to warm to room temperature 18–30°C (64–86°F) before using.

- Remove 1 well for each sample to be tested plus 1 well for the control, and place into the well holder.
- Choose one yellow-labeled control bottle to serve as the screening level for the test.
- Mix each reagent by swirling its bottle prior to use.
- Add 100 µL from the yellow-labeled control bottle to the first well. Add 100 µL from each sample extract to a respective well as indicated in the template below. For environmental swabs, add 3 drops from swab tube with dropper tip. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface.

Control	S1	S2	S3	S4	S5
---------	----	----	----	----	----

- Incubate microwells **10 minutes** at room temperature 18–30°C (64–86°F).
- Shake out the contents of the wells. Using a wash bottle filled with wash buffer solution, fill each well and shake out. Repeat 10 times. Remove excess wash buffer by turning wells upside down and vigorously tapping wells on absorbent towel.
- Add 100 µL from the blue-labeled conjugate bottle to each well. Mix for 20 seconds by sliding the microwell holder back and forth on a flat surface.
- Incubate for **10 minutes** at room temperature 18–30°C (64–86°F).
- Shake out the contents of the wells. Using a wash bottle filled with wash buffer solution, fill each well and shake out. **Repeat 10 times.** Remove excess wash buffer by turning wells upside down and vigorously tapping wells on absorbent towel.
- Add 100 µL from the green-labeled substrate bottle to each well. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface.
- Incubate for **10 minutes** at room temperature 18–30°C (64–86°F).
- Add 100 µL from the red-labeled Red Stop bottle to each well. Mix for **20 seconds** by sliding the microwell holder back and forth on a flat surface. The results are now ready to be interpreted.
- Visually compare the color of a sample well to the color of the control well. If the sample well has more blue color than the control well, the sample tests positive for milk contamination of more than the control used. If the sample well has less blue color, or more red color, than the control well, the sample contains less than the control used of milk contamination.

**Alternative:** Read wells (wipe the bottom of wells with a dry cloth or towel first) in a microwell reader with a 650 nm filter. If the sample well has an optical density (OD) higher than the control well, the sample is positive for milk contamination of more than the control used. If the sample well has an OD lower than the control well, the sample contains less than the control used of milk contamination.

### PERFORMANCE CHARACTERISTICS

**Limit of detection:** 1 ppm nonfat dried milk (NFDM) (described as the lowest concentration point on the calibration curve that this test can reliably detect milk protein.)

**Range of quantitation:** 2.5–25 ppm NFDM (for quantitating samples below 2.5 ppm or above 25 ppm, contact a Neogen representative for dilution instructions.)

**Allergen detection:** Test detects casein and whey proteins from cow, goat and sheep and the results are expressed as ppm of NFDM.

**Protein Conversion:** Test results are expressed as ppm NFDM. To express results as ppm protein, multiply the NFDM result by 0.36 (e.g., 2.5 ppm NFDM x 0.36 = 0.9 ppm NFDM protein). \*NFDM contains 36% protein.

\*Source: USDA National Nutrient Database for Standard Reference Release 28, #01091 - Milk, dry, nonfat, regular, without added vitamin A and vitamin B.

## CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

## SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's Food Safety test kits, at [www.neogen.com](http://www.neogen.com), or by calling Neogen at 800/234-5333 or 517/372-9200.

## TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit [www.neogen.com/Corporate/termsconditions.html](http://www.neogen.com/Corporate/termsconditions.html).

## WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

### TEST KITS AVAILABLE FROM NEOGEN

#### Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

#### Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

#### Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

#### Food allergens

- Almonds, crustaceans, eggs, gliadin, hazelnut, milk, mustard, peanuts, sesame, soy, walnuts, multi-treenut

#### Genetic modification

- CP4 (Roundup Ready®)

#### Ruminant by-products

- Meat and bone meal, feed

#### Species Identification

- Raw and cooked meat samples, feed

#### North America Neogen Headquarters

800/234-5333 (USA/Canada)  
[foodsafety@neogen.com](mailto:foodsafety@neogen.com)  
[foodsafety.neogen.com](http://foodsafety.neogen.com)

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©Neogen Corporation, 2018. Neogen, Veratox and K-Blue are registered trademarks of Neogen Corporation. All other brand and product names are trademarks or registered trademarks of their respective companies.

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