-  **(EN) Microbial Luminescence System  
Ultra High Temperature Beverage Screen Kit**
-  **(FR) Microbial Luminescence System  
Kit pour les boissons ultra haute température**
-  **(DE) Mikrobielles Lumineszenz System  
Ultra Hoch Temperatur Getränke Screening Kit**
-  **(IT) Kit di controllo per bevande trattate a temperatura  
elevatissima per sistemi di luminescenza microbica**
-  **(ES) Sistema de Luminiscencia para Detección Microbiana Kit  
de Análisis para Bebidas Ultra Alta Temperatura**
-  **(NL) Kit voor het screenen van Ultra Hoge Temperatuur dranken  
op basis van een microbieel luminescentiesysteem**
-  **(PT) Kit de Triagem de Bebidas em Temperatura  
Ultraelevada para Sistema de Luminescência Microbiana**
-  **(RU) Комплект для скрининга напитков,  
прошедших ультравысокотемпературную  
обработку, в составе микробиологической  
люминисцентной системы**
-  **(JA) 微生物発光システム超高温の飲料スクリーンキット**
-  **(ZH) 微生物发光系统超高温饮料筛选试剂盒**
-  **(TH) ชุดน้ำยาทดสอบจุลินทรีย์ในผลิตภัณฑ์เครื่องดื่มยูเอชที**
-  **(AR) معدة فحص المشروبات بالحرارة العالية جدًا بجهاز اللعنان  
الميكروبي**

## Product Instructions

# Microbial Luminescence System Ultra High Temperature Beverage Screen Kit

### PRODUCT DESCRIPTION AND INTENDED USE

The Neogen® Microbial Luminescence System (MLS) Ultra High Temperature (UHT) Beverage Screen Kit is intended for use with the Neogen® Microbial Luminescence System (MLS) II Instrument. These kits offer a rapid detection system that utilizes the adenosine triphosphate (ATP) bioluminescence technology to detect the presence of microbial ATP in beverages such as; dairy, dairy substitutes, juices, caffeinated and electrolytic drinks that have been processed for commercial sterilization including Ultra-High Temperature (UHT) and Extended Shelf life (ESL) products. After the beverage is subjected to an enrichment step in its original unopened container, the Neogen MLS UHT Beverage Screen Kit has the capability to exclude ATP from non-microbial sources and then measure only the ATP released by microorganisms. In less than 30 minutes, the microbial ATP is measured in a light output format of Relative Light Units (RLU) proportional to the amount of microbial ATP present in the enriched sample. The Neogen MLS UHT Beverage Screen Kit significantly reduces the time for product release compared to traditional microbiology testing methods.

Neogen Food Safety is certified to ISO (International Organization for Standardization) 9001 for design and manufacturing.

### SAFETY

The user should read, understand, and follow all safety information in the instructions for the Neogen MLS UHT Beverage Screen Kit. Retain the safety instructions for future reference.

**⚠ WARNING:** Indicates a hazardous situation, which, if not avoided, could result in death or serious injury and/or property damage.

**⚠ CAUTION:** Indicates a hazardous situation, which, if not avoided, could result in minor or moderate injury and/or property damage.

**NOTICE:** Indicates a potentially hazardous situation which, if not avoided, could result in property damage.

### ⚠ CAUTION

#### To reduce the risk of inaccurate results:

- The dehydrated ATPase and LL1 reagents and corresponding buffers should be stored at 2-8°C.
- Avoid the exposure of the kit's reagents to the direct sunlight.
- Do not use reconstituted reagents beyond 5 days, stored at 2-8°C.
- Do not expose reconstituted reagents to ambient temperatures (20-25°C) for more than 12 hours during the 5 day reconstituted reagents shelf life.
- Do not shake the reconstituted reagents. Follow instructions for use described for reagent preparation.

#### To reduce the risk associated with exposure to chemicals and biohazards:

- Handle enriched samples following standard laboratory safety practices, including wearing appropriate protective apparel and eye protection.

### NOTICE

#### To reduce the risks associated with cross-contamination while preparing the assay

- The use of gloves is recommended.
- ATP is a common substance. It is present on skin, hair and on many surfaces. To avoid ATP contamination, do not touch pipette tips or reagent vial stoppers with bare hands, or any part of the Neogen MLS II Instrument that comes in direct contact with the reagents.

Consult the Safety Data Sheet for additional information.

For information on documentation of product performance, visit our website at [www.neogen.com](http://www.neogen.com) or contact your local Neogen representative or distributor.

### USER RESPONSIBILITY

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at [www.neogen.com](http://www.neogen.com), or contact your local Neogen representative or distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results. The food sample itself may influence results.

It is the user’s responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user’s criteria.

It is also the user’s responsibility to determine that any test methods and results meet its customers’ and suppliers’ requirements.

As with any test method, results obtained from use of any Neogen Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

**LIMITATION OF WARRANTIES / LIMITED REMEDY**

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, NEOGEN DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any Neogen Food Safety Product is defective, Neogen or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify Neogen within sixty days of discovery of any suspected defects in a product and return it to Neogen. Please contact your Neogen representative or authorized Neogen distributor for any further questions.

**LIMITATION OF NEOGEN LIABILITY**

NEOGEN WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall Neogen’s liability under any legal theory exceed the purchase price of the product alleged to be defective.

**STORAGE AND DISPOSAL**

See Table 1 for specific storage conditions. DO NOT FREEZE.

**⚠ DISPOSAL**

After use, Neogen MLS UHT Beverage Screen Kit reagents must be disposed of per local, regional, national, and industry standards for disposal of chemical reagents.

**Table 1.** Kit contents for BEV600.

Item	Description	Quantity	Storage	Comments
ATPase** (blue cap)	Purified lyophilized ATPase enzyme that degrades somatic cells and free ATP to allow detection of only microbial ATP	2 vials with lyophilized enzyme	2-8°C DO NOT FREEZE	Reconstitute with ATPase buffer (read instructions below)
ATPase Buffer (blue cap)	Buffer to reconstitute ATPase enzyme	2 vials	2-8°C DO NOT FREEZE	Used to reconstitute the ATPase lyophilized enzyme (read instructions below)
Extractant (black cap)	Rapidly lyses microorganisms to release microbial ATP	2 bottles	20-25°C DO NOT FREEZE	Ready to use <b>NOTE:</b> Refrigerate until first use. Once in use maintain at 20-25°C. Store in the dark.
LL1** (white cap)	Purified Luciferin/ Luciferase complex that interacts with microbial ATP to produce light captured by the Neogen MLS II Instrument	2 vials with lyophilized enzyme	2-8°C DO NOT FREEZE	Reconstitute with LL1 buffer (read instructions below)
LL1 Buffer (white cap)	Buffer to reconstitute LL1 enzyme	2 vials	2-8°C DO NOT FREEZE	Used to reconstitute the LL1 lyophilized enzyme (read instructions below)
Quick Start Guide		1		

\*NOTE: The BEV600 can be utilized to perform up to 600 tests. Reconstitute one vial of LL1 and ATPase with the respective buffer to run up to 300 tests.

\*\* Marks from 50 to 300 on the left side of the label represent an estimation of the number of samples remaining in the vial. The white space on the right side of the label underneath the “Date Opened:” text can be utilized to write in the date that the enzyme was reconstituted.

The following table describes accessory products related to Neogen MLS II that are not included with the Neogen MLS UHT Beverage Screen Kit.

**Table 2.** Neogen MLS II Products.

Catalog Number	Description
ATP50	Neogen® MLS Reagent Control Kit
3005	Neogen® Microbial Luminescence System (MLS) Injector Cleaning Kit
3006	Neogen® Microbial Luminescence System (MLS) Maintenance Solution
3007	Neogen® Microbial Luminescence System (MLS) Microwell Plate
3008	Neogen® Microbial Luminescence System (MLS) Microwell Strip
3009	Neogen® Microbial Luminescence System (MLS) Microwell Strip Holder
BMLSCK	Neogen® Microbial Luminescence System (MLS) Weekly Cleaning Kit

## INSTRUCTIONS FOR USE

Follow all instructions carefully. Failure to do so may lead to inaccurate results.

The overall work flow of all procedures to run the Neogen MLS UHT Beverage Screen Kit is described in Figure 1 at the end of this document.

### SAMPLE ENRICHMENT

Prior to the screening of UHT or ESL beverages with the Neogen MLS UHT Beverage Screen Kit, samples need to be enriched to allow growth of low populations of microorganisms.

1. Place the beverage in the incubator. To achieve uniform temperature within the beverage container, ensure that containers do not touch each other or the walls of the incubator. Table 3 provides guidance for enrichment of dairy, dairy substitutes, juices and other beverages. **It is the user's responsibility to validate the enrichment protocol to ensure the results meet the user's criteria.**

**Table 3.** Guidance for enrichment conditions of UHT and ESL beverages.

Product type (UHT/ESL)	Examples (This is not a complete a list)	Temperature of enrichment (°C)	Time of enrichment (hours)
Dairy	Milk Flavored milk Dairy based creamers Ice cream mix Caffeinated beverages containing dairy or dairy substitute ingredients	30 ± 1	48 ± 2
Dairy substitutes	Soy milk Rice milk Almond milk Coconut milk	30 ± 1	48 ± 2
Juices	Orange juice Apple juice Peach juice Grape juice Tomato juice Dairy or dairy substitute beverages containing fruit juice	30 ± 1	72 ± 2
Beverages	Caffeinated beverages without sugar, dairy or dairy substitutes. Coconut water	30 ± 1	72 ± 2

2. Following enrichment, remove the beverage container from the incubator.
3. Proceed with sample screening with the Neogen MLS UHT Beverage Screen Kit.

## REAGENT PREPARATION

### ATPase Reagent\*\*

1. Open the ATPase vial (blue cap) containing the lyophilized ATPase enzyme. Carefully remove the rubber stopper; contents are under vacuum.
2. Open the ATPase buffer vial (blue cap) and pour or pipette the entire contents of the buffer into the ATPase vial.
3. Replace the rubber stopper and invert 5 - 10 times to dissolve the lyophilized ATPase, then gently swirl. DO NOT SHAKE.

### LL1 Reagent\*\*

1. Open the LL1 vial (white cap) containing the lyophilized Luciferin/Luciferase complex. Carefully remove the rubber stopper; contents are under vacuum.
2. Open the LL1 buffer vial (white cap) and pour or pipette the entire contents of the buffer into the LL1 vial.
3. Replace the rubber stopper and invert 5 times to dissolve the lyophilized LL1, then gently swirl. DO NOT SHAKE.

**\*\*NOTE:** The reconstituted ATPase and LL1 Reagents have a shelf life of 5 days when stored at 2-8°C including up to 12 hours at 20-25°C. Return the reconstituted ATPase and LL1 Reagents to storage at 2-8°C when not in use. **DO NOT FREEZE THE RECONSTITUTED ATPase OR LL1 REAGENTS.**

### Extractant (Black cap)

The Extractant solution is ready for immediate use. Once in use, it can be maintained at 20-25°C. It is not recommended to return to the refrigerator. When not in use store in a dark place.

## PREPARING THE NEOGEN MLS II INSTRUMENT

### Cleaning the system for initial use

1. Switch on the Neogen MLS II Instrument first and then the computer.
2. Launch the Neogen MLS II Software.
3. Place three vials containing Neogen MLS Cleaning Solution on the Neogen MLS II Instrument.
4. In the Neogen MLS II Software, open the “**Actual**” tab and run a “**Wash Assay**”.
5. Replace the three vials containing Neogen MLS Cleaning Solution with three vials containing ATP-free water and run a second “**Wash Assay**”.

Refer to the Neogen MLS II User’s Manual for detailed information related to cleaning the Neogen MLS II Instrument injectors.

### Priming the Reagents in the Neogen MLS II Instrument Injector lines

**Important:** Ensure all reagents (reconstituted ATPase, Extractant and reconstituted LL1) have reached ambient temperature (20-25°C) before starting the Reagent Control Assay and Beverage Screening Assay.

1. Replace vials containing ATP-free water with reconstituted ATPase (on injector A), Extractant (on injector B) and reconstituted LL1 (on injector C). Retain rubber stoppers and caps for resealing vials for later storage.
2. In the Neogen MLS II Software, open the “**Actual**” tab and run a “**Prime Assay**”.

For more detailed information related to priming of the Neogen MLS II Instrument refer to the Neogen MLS II User’s Manual.

**NOTE:** Priming injectors before running a Reagent Control Assay and a Beverage Screening Assay is needed to ensure the instrument’s injectors and reagent tubes are pre-filled with the Neogen MLS UHT Beverage Screen Kit Reagents. Failure to do so can result in False-Negatives or False-Positive results.

## RUNNING A REAGENT CONTROL ASSAY

Performing a Reagent Control Assay is always recommended prior to product testing, to ensure the Neogen MLS II Instrument and Reagents are performing appropriately. Please refer to the Neogen MLS ATP Reagent Control Kit Instructions for Use for detailed information about the Reagent Control Assay and reconstitution of the ATP Reagent.

1. Place a Neogen MLS Microwell Strip (a minimum of 6 wells is required) in a Neogen Microwell Strip Holder. A Neogen MLS Microwell Plate can be used instead of the Neogen MLS Microwell strips.
2. Pipette 50 µL of the reconstituted ATP Reagent (See Table 2) into the **bottom** of the last 4 wells (C1, D1, E1, and F1). Wells A1 and B1 should be empty.
3. In the Neogen MLS II Software, open the “**Actual**” tab and run a “**Reagent Control Assay**”.

For more detailed information about the Reagent Control Assay and interpretation of results, please refer to the Neogen MLS ATP Reagent Control Kit Instructions for Use.

## RUNNING A BEVERAGE SCREENING ASSAY

1. Homogenize the enriched beverage in its original unopened container by shaking.
2. Aseptically, remove 50 µL of the enriched beverage with a micropipette and deposit the sample into the **bottom** of a well in a Neogen MLS Microwell Plate (a Neogen MLS Microwell Strip can be utilized for fewer number of assays). **Read Special Instructions for orange juice with high content of pulp\*\*\*.**

**NOTE.** Careful pipetting of samples is required to ensure accurate results. Depositing the sample onto the well wall can cause results to be inaccurate. Users should pipette samples directly into the **bottom** of the well.

3. Repeat this step utilizing a fresh pipette tip for each sample. It is recommended to dispense the sample in order of columns A1 through H1, A2 through H2 etc.
4. In the Neogen MLS II Software, open the **“Actual”** tab and select the assay according to Table 4.

**Table 4.** Neogen MLS II Software Assays.

Matrix	Assay Name
<b>Beverages containing fruit juice</b> Examples: coconut water, tomato juice, fruit juices and dairy or dairy substitute mixes with fruit.	<b>“UHT Juice Assay”</b>
<b>All other beverages</b> including dairy and dairy substitute products.	<b>“UHT Assay”</b>

5. If desired, complete the Neogen MLS UHT Beverage Screen Kit Lot information and click **“OK”**
6. Select the wells that contain samples to be analyzed in the 96 well grid displayed on the **“Actual”** tab (unselected wells will not be analyzed).
7. Click **“Start”**
8. Complete the **“Load Plate”** information (name file of the Assay).
9. Carefully place the Neogen MLS Microwell Plate containing the samples onto the plate carrier of the Neogen MLS II Instrument ensuring well **“A1”** is positioned at the left corner of the plate carrier, closest to the user.
10. Click **“OK”** and the assay will start.
11. After the Assay is completed, remove the Neogen MLS Microwell Plate from the Neogen MLS II Instrument and dispose the plate according to current industry standards.

### \*\*\*Special instructions for orange juice with high content of pulp.

High RLU readings from blanks (not enriched samples) may be obtained as a result of pulp present in the well. Proper pipetting technique is important to ensure correct results. The following steps may help in alleviating inconsistent results:

Before loading the sample in the Neogen MLS Microwell Plate.

- Perform a 1:1 ratio dilution (for example: 10 mL of juice with 10 mL of ATP-free water).
- Utilize wide bore pipette tips to collect the sample and deposit it into the well.
- When analyzing juices with high pulp content, include at least 5 replicates of each sample. RLU readings for the blank should be lower than 150 + 15%. Larger RLU values within the replicates can be considered outliers. Calculate the background RLU of the sample with at least three consistent readings.

Example:

Well	RLU
A1	165
B1	172
C1	176
<b>D1</b>	<b>1650****</b>
E1	170

\*\*\*\*The RLU obtained in well D1 is likely to be an outlier and should not be considered to determine the background RLU of the sample.

## INTERPRETATION OF RESULTS

**NOTE. It is the user’s responsibility to validate the Pass/Fail limits to ensure this test method meets the user’s criteria.**

To establish the threshold/background RLU, it is advised to obtain repeated measurements (multiple replicates) of RLU readings from several lots of non-contaminated product. Utilize this data to determine the standard deviation of the background RLU of the analyzed product. A minimum Pass/Fail limit can be defined as the average background RLU plus three times the standard deviation of the background RLU for each specific matrix/product. For specific instructions on changing the Pass/Fail limits in the software, contact your Neogen Technical Service representative.

1. Results will appear after approximately 27 minutes (for a full plate of 96 samples) in the “Report” Tab. The results will be displayed by means of color code as the assay progresses (Green/Pass and Red/Fail). Placing the mouse over the well will show the RLU value in real time as the assay progresses.
2. A Pass or Fail result is determined by the establishment of a Pass/Fail limit based on a RLU threshold value. As a default in the software, a result of <150 RLU will be marked as a Pass; however Pass/Fail limits may vary depending on customer or product requirements.


#### **CLEANING THE NEOGEN MLS II INSTRUMENT AFTER USE**

Following sample testing, the reagents dispensing lines and injectors must be rinsed by flushing ATP-free water through the system.

1. Replace the ATPase vial, Extractant and LL1 vial with vials of ATP-free water. Place rubber stopper/caps back onto their respective vials and store the ATPase and LL1 Reagents at 2-8°C, or discard appropriately. Store Extractant at 20-25°C after first use (Table 1).
2. Select the “**Actual**” tab and click on the “**Wash Assay**” with ATP-free water first and then repeat with Neogen MLS Cleaning solution. Leave the Neogen MLS Cleaning solution on the instrument until next use.

Refer to the Neogen MLS User’s Manual for detailed information related to cleaning the Neogen MLS II Instrument injectors.

**NOTE.** The waste container for the Neogen MLS II Instrument should be emptied on a daily basis.

**Figure 1.** Work flow to perform screening of UHT and ESL beverages.

If you have questions about specific applications or procedures, please visit our website at [www.neogen.com](http://www.neogen.com) or contact your local Neogen representative or distributor.

## EXPLANATION OF SYMBOLS

[info.neogen.com/symbols](http://info.neogen.com/symbols)

# Neogen Food Safety

## Neogen Corporation

620 Leshar Place  
Lansing, MI 48912 USA  
Neogen.com

## Neogen Europe Ltd.

The Dairy School  
Auchincruive  
Ayr, KA6 5HU  
Scotland, UK

## Neogen Ireland, Ltd.

Bray Business Park, Bray  
Co. Wicklow  
A98YV29, Ireland

## Neogen Corporation

620 Leshar Place Lansing, MI 48912 USA  
www.neogen.com

Neogen is a trademark of Neogen Corporation.

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All rights reserved.
© Neogen Corporation 2024.