*Read instructions carefully before using*

# ALERT<sup>®</sup>

## for Histamine Extraction Kit

*Designed for use with  
Neogen's Alert for Histamine test kit*

### INTENDED USE

This Alert<sup>®</sup> for Histamine Extraction Kit contains materials necessary for the extraction of histamine from scombroid species of fish, such as tuna, bluefish and mahi-mahi, and in fish meal

### MATERIALS PROVIDED

1. 20 disposable extraction bottles (see **Precaution 1**)
2. 20 disposable filter syringes
3. 20 sample collection tubes, plastic (12 x 75 mm)
4. Directions for use

### MATERIALS REQUIRED BUT NOT PROVIDED

1. Blender
2. Graduated cylinder capable of measuring 190 mL (Neogen item 9368)
3. Timer (Neogen item 9426)
4. Test tube rack (Neogen item 9440)
5. Distilled or deionized water

### STORAGE

This kit may be stored at room temperature, 18–30°C (64–86°F).

### PRECAUTIONS

1. If using extraction bottles other than those included in the Alert for Histamine Extraction Kit, make sure the bottles are plastic, **NOT GLASS**. Histamine adheres to glass, thus the use of glassware can affect test results.
2. Good laboratory practices should always be followed.
3. Always use a new extraction bottle, sample collection tube and filtering syringe for each sample to avoid cross-contamination.
4. Destroy and dispose of filtering syringe in a manner to prohibit reuse.
5. Do not take filtering material out of syringe. Material can cause skin and eye irritation.

### SAMPLE PREPARATION FOR EXTRACTION

1. **Canned tuna: AOAC 937.07b**  
Place entire contents of can (meat and liquid) into a blender. Blend until homogenous.
2. **Fresh or thawed frozen raw fish: AOAC 937.07a**  
Clean and eviscerate three fish. Cut three cross-sectional pieces 2.5 cm (1 inch) thick, from back of the pectoral fin, halfway to vent and one posterior to the vent. Debone slices and blend or grind combined samples until homogenous.

## SAMPLE EXTRACTION PROCEDURE

1. Add 190 mL of distilled or deionized water to an extraction bottle.
2. Add 10 g of the homogenous mixture to the extraction bottle.
3. Tightly cap and vigorously shake the bottle for **15–20 seconds** to suspend the fish tissue in the water.
4. Wait approximately **5 minutes**, then shake the bottle for **15–20 seconds** to resuspend the fish tissue.
5. Wait an additional **5 minutes**, and again shake the bottle for **15–20 seconds**. After shaking, allow the tissue to settle for about **30 seconds**.
6. Ensure the filtering material is tightly packed in the syringe by depressing the plunger. Remove the plunger from the syringe and hold syringe tip over a sample collection tube. Fill the syringe by pouring the liquid from the extraction bottle into the syringe until it is at least 1/2 full.
7. Replace the plunger in the syringe and apply gentle pressure until the sample collection tube is at least 1/3 full. Discard the syringe. Test tubes and extraction bottles should be labeled for proper sample identification.
8. The sample is now ready for testing with the Alert for Histamine test kit. Proceed to the test kit's **SAMPLE EXTRACT DILUTION** step.

## CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

## MSDS INFORMATION AVAILABLE

Material safety data sheets (MSDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at [www.neogen.com](http://www.neogen.com), or by calling Neogen at 800/234-5333 or 517/372-9200.

## TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit [www.neogen.com/Corporate/termsconditions.html](http://www.neogen.com/Corporate/termsconditions.html).

## WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

### North America Neogen Headquarters

620 Leshler Place, Lansing, MI 48912 USA  
800/234-5333 (USA/Canada) or 517/372-9200  
Fax: 517/372-2006 • [foodsafety@neogen.com](mailto:foodsafety@neogen.com)  
[www.neogen.com](http://www.neogen.com)

### Europe, Middle East and Africa Neogen Europe

The Dairy School, Auchincruive, Ayr  
KA6 5HW Scotland, UK  
+ 44 (0) 1292 525 600  
Fax: + 44 (0) 1292 525 601  
[info\\_uk@neogeneurope.com](mailto:info_uk@neogeneurope.com)  
[www.neogeneurope.com](http://www.neogeneurope.com)

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