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

## A. Environmental Scrub Sampler Stick with 10 mL Wide Spectrum Neutralizer (ESS10WSN)

## B. Environmental Scrub Sampler with 10 mL Wide Spectrum Neutralizer and Gloves (HES10WSN2G)

### Product Description and Intended Use

Neogen® Environmental Scrub Sampler Stick with 10 mL Wide Spectrum Neutralizer and Neogen® Environmental Scrub Sampler with 10 mL Wide Spectrum Neutralizer and Gloves are a ready-to-use environmental scrub sampler devices intended to be used in the food and beverage industry for environmental microbial surface sampling.

Neogen Environmental Scrub Sampler combines an innovative design and high-tech materials into a sampling device. Its function integrates a scouring feature that enables disruption of organic build up including biofilm and effective sample pick-up.

The Wide Spectrum Neutralizer offers a neutralizing capability for a variety of sanitizers used in the food industry. It has been formulated using components accepted for use in the food industry, without any known allergenic components, animal derived materials (ADM), polymerase chain reaction (PCR) inhibitors, or components derived from genetically modified organisms (GMOs). The Wide Spectrum Neutralizer is also compatible with the following: Neogen® Petrifim® Plates: Aerobic Count, Rapid Aerobic Count, Rapid *E. coli* / Coliform Count, Select *E. coli* Count, Coliform Count, *E. coli* / Coliform Count, High Sensitivity Coliform Count, Rapid Yeast and Mold Count, Yeast and Mold Count, and with the Neogen® Molecular Detection System.

Neogen Environmental Scrub Sampler Stick contains a metal additive to transform it into a detectable product that can be identified with a metal detector system\*.

Product Description	Catalog number	Volume of buffer (mL)	Gloves
Neogen Environmental Scrub Sampler Stick with 10 mL Wide Spectrum Neutralizer	ESS10WSN	10	none
Neogen Environmental Scrub Sampler with 10 mL Wide Spectrum Neutralizer and Gloves	HES10WSN2G	10	100**

\*\* Aseptic gloves, 100 individually wrapped pairs per case.

\* The Neogen Environmental Scrub Sampler Stick is designed with a metal additive that enables metal detection. Detectability in metal detection systems is dependent on many factors, including food matrix, food package size, detection system settings, and orientation and size of the metal particle, among others. Users must validate with their own product, production systems and detection system to determine the level of metal detectability of the Neogen Environmental Scrub Sampler Stick.

For Laboratory Use Only. Material Safety Data Sheet available.

## Safety

The user should read, understand, and follow all safety information in the instructions. Retain the safety instructions for future reference.

 **WARNING** Indicates a hazardous situation, which, if not avoided, could result in death or serious injury and/or property damage

 **WARNING**

**The user must train its personnel in current proper methods for testing and surface sampling techniques: for example, Good Laboratory Practices<sup>1</sup>, ISO/IEC 17025<sup>2</sup>, or ISO 18583:2018<sup>3</sup>**

### To reduce the risks associated with environmental contamination:

- The Neogen Environmental Scrub Sampler products are intended to be used for testing for microorganisms on surfaces. Surfaces may potentially contain pathogenic organisms, such as *Listeria monocytogenes* or *Salmonella*.
- Individuals should be trained in accordance with applicable regulatory and company/institution requirements before working with potentially infectious materials.
- All enrichment broths should be sterilized following any culture based confirmatory steps.
- Strict compliance with BSL-2 (Biosafety Level 2) practices and current industry standards/local and federal regulations for disposal of contaminated waste should be followed.

### To reduce the risks associated with exposure to chemicals and biohazards:

- Dispose of samples according to all applicable government regulations, industry standards and applicable laboratory procedures for disposal of biohazardous waste. Strict compliance with BSL-2 practices should be followed.
- Always follow standard laboratory safety practices (GLP<sup>1</sup> or ISO 17025<sup>2</sup>), including proper containment procedures, wearing appropriate protective apparel, disposable gloves and eye protection while handling reagents and contaminated samples.

### To reduce the risk associated with false negatives resulting in the use of contaminated environmental surfaces for food or beverage products:

- Always reference package label for storage instruction and expiration date.
- Always reference product instruction for usage.

### To reduce the risk of false-positive results due to cross contaminated environmental surfaces for food or beverage products that may result in re-testing or the rejection of food or beverage product:

- Do not touch the Neogen Environmental Scrub Sampler device to any unintended surface.
- Do not break the Neogen Environmental Scrub Sampler device while sampling.
- Do not reach into the Neogen Environmental Scrub Sampler device bag.

### To reduce the risk of cross contamination from reuse of Neogen Environmental Scrub Sampler device:

- Do not use the same Neogen Environmental Scrub Sampler device more than once.
- Do not use the same Neogen Environmental Scrub Sampler device for sampling more than one surface area.
- Review that the bag does not have any defect that can compromise the aseptic conditions of the Scrub Sampler device.

### The colors of the Neogen Environmental Scrub Sampler and Neogen Environmental Scrub Sampler Stick are designed to be visible in food production environments.

Consult the Safety Data Sheet for additional information.

For information on documentation of product performance, visit our website at [www.neogen.com](http://www.neogen.com) or contact your local Neogen representative or distributor.

## User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at [www.neogen.com](http://www.neogen.com), or contact your local Neogen representative or distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results. The food sample itself may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any Neogen Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

### Limitation of Warranties / Limited Remedy

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, NEOGEN DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any Neogen Food Safety Product is defective, Neogen or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify Neogen within sixty days of discovery of any suspected defects in a product and return it to Neogen. Please contact your Neogen representative or authorized Neogen distributor for any further questions.

### Limitation of Neogen Liability

NEOGEN WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall Neogen's liability under any legal theory exceed the purchase price of the product alleged to be defective.

### Storage

Refer to package label for storage information and shelf life. Store the Neogen Environmental Scrub Sampler products in foil bag at 2-8°C for longest shelf life. Neogen Environmental Scrub Sampler products can also be stored in foil bag at 15-25°C for up to three months. Neogen Environmental Scrub Sampler products should not be used past their expiration date.

### DISPOSAL

After use, Neogen Environmental Scrub Sampler products may contain microorganisms that may be a potential biohazard. Follow current industry standards and local regulations for disposal of contaminated waste. Consult the Material Safety Data Sheet for additional information.

### Validated Methods

**AOAC® Performance Tested Method<sup>SM</sup>  
#022104**


The Neogen Environmental Scrub Sampler with 10 mL Wide Spectrum Neutralizer recovered *Listeria* or *Salmonella spp.* from several different environmental surfaces including: stainless steel, plastic (polystyrene), and sealed concrete. Using POD analysis, no statistical differences were observed between the number of positive samples detected by the Neogen Environmental Scrub Sampler with 10 mL Wide Spectrum Neutralizer sampling method and the reference sampling method for all samples tested (Dey-Engley Neutralizing Broth). The Neogen Wide Spectrum Neutralizer successfully neutralized a range of sanitizers including quaternary ammonium, high acid, hydrogen peroxide/peroxyacetic acid, and chlorine/bleach and was found to be non-toxic to the target organisms. The Neogen Environmental Scrub Sampler with 10 mL Wide Spectrum Neutralizer was found to be a robust and stable sampling device through robustness and product consistency testing.


## Instructions for Use

Follow all instructions carefully. Failure to do so may lead to inaccurate results.


### A. ESS10WSN


Wearing gloves, tear off the top of the bag along the perforation.


Aseptically open the bag using the red tabs on either side of the bag. Be sure not to touch the inside or edges of the bag.


Squeeze out excess Neutralizer solution so the Neogen Environmental Scrub Sampler device is moist but not dripping.


Working from the outside of the bag, move the device up allowing the stick to protrude from the bag.


Aseptically, using one hand, grasp the stick above the thumb stop (see blue arrow in the image above) and remove the device from the bag, being sure the device does not touch the outside part of the bag.


Practicing aseptic technique, press the Neogen Environmental Scrub Sampler device down firmly and flex the stick to ensure full contact with the sampling surface. Sample in a zigzag motion in one direction across the entire sampling surface. Optionally, scour vigorously in a zigzag motion in one direction across the entire sampling surface to disrupt organic matter if present. Sample an area from 10x10 cm (4x4 in) to 30x30 cm (12x12 in), following appropriate standards or regulatory guidance<sup>3,4,5,6,7</sup>.


Turn the device over to the other side, change sampling direction by 90° and repeat the swabbing procedure described in #6 in the same sampling site.


Return the sampling device back into the bag, without going beyond the thumb stop, and hold the device with one hand from the outside of the bag.

Using the other hand hold the Neogen Environmental Scrub Sampler Stick and twist it to separate from the Neogen Environmental Scrub Sampler device. Allow the Neogen Environmental Scrub Sampler to drop to the bottom of the bag to re-submerge into the Neutralizer. Discard the stick.


Close the bag by rolling the blue wires down and folding in the ends of the wires.


It is recommended to clean/sanitize the area after sampling to remove any remaining neutralizing solution residue from the sampled surface<sup>3</sup>


## B. HES10WSN2G


1 Tear off the top of the bag along the perforation.


2 Aseptically open the bag using the red tabs on either side of the bag. Be sure not to touch the inside or edges of the bag.


3 Work from the outside of the bag and move the Neogen Environmental Scrub Sampler to the top of the bag.


4 Squeeze out excess Neutralizer solution so the Neogen Environmental Scrub Sampler is moist but not dripping.


5 Aseptically put on the provided sterile glove. Be aware that the glove should not touch the outside part of the bag.


6 Using the gloved hand, aseptically remove the Neogen Environmental Scrub Sampler from the bag, being sure it does not touch the outside of the bag.


7 Practicing aseptic technique, press the Neogen Environmental Scrub Sampler device firmly using the side with scrub dots down to ensure full contact with the sampling surface. Sample in a zigzag motion in one direction across the entire sampling surface. Optionally, scour vigorously the sampling surface with the device in a zigzag motion to disrupt organic matter if present. Sample an area of about 30x30 cm (12x12 in), following appropriate standards or regulatory guidance<sup>3,4,5,6,7</sup>.


8 Turn the Neogen Environmental Scrub Sampler device over to the other side, change sampling direction by 90° and sample in a zigzag motion across the entire same sampling surface. Sample an area of about 30x30 cm (12x12 in), following appropriate standards or regulatory guidance<sup>3,4,5,6,7</sup>.

9 Return the Neogen Environmental Scrub Sampler device into the bag without touching the edges of the bag. Move the device to the bottom of the bag to re-submerge into the Neutralizer.

10 Close the bag by rolling the blue wires down and folding in the ends of the wires.

11 It is recommended to clean/sanitize the area after sampling to remove any remaining neutralizing solution residue from the sampled surface<sup>3</sup>.

If you have questions about specific applications or procedures, please visit our website at [www.neogen.com](http://www.neogen.com) or contact your local Neogen representative or distributor.

## References

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
3. ISO 18593:2018. Microbiology of the food chain – Horizontal methods for surface sampling.
4. American Public Health Association - Compendium of Methods for the Microbiological Examination of Foods - Chapter 3: Microbiological Monitoring of the Food Processing Environment, 4th edition.
5. US Food and Drug Administration - Bacteriological Analytical Method (available online at <https://www.fda.gov/food/laboratory-methods-food/bacteriological-analytical-manual-bam>).
6. United States Department of Agriculture - Microbiological Lab Guidebook - Chapters 4.04, 5.04, and 8.07.
7. American Public Health Association - Standard Methods for the Examination of Dairy Products - Chapter 13: Microbiological Tests for Equipment, Containers, Water, and Air, 17th edition.

Refer to the current versions of the standard methods listed above.

## Explanation of Symbols

[info.neogen.com/symbols](http://info.neogen.com/symbols)

## Neogen Food Safety

### Neogen Corporation

620 Leshar Place  
Lansing, MI 48912 USA  
Neogen.com

### Neogen Europe Ltd.

The Dairy School  
Auchincruive  
Ayr, KA6 5HU  
Scotland, UK

### Neogen Ireland, Ltd.

Bray Business Park, Bray  
Co. Wicklow  
A98YV29, Ireland

### Neogen Corporation

620 Leshar Place Lansing, MI 48912 USA  
www.neogen.com

Neogen and Petrifilm are trademarks of Neogen Corporation.

© Neogen Corporation 2024.

Unauthorized use prohibited. All rights reserved.

All other trademarks are the property of their  
respective companies.

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