700002516 | 8188

## Product Instructions

-  **EN** Neogen® Reveal® Q+ MAX for Zearalenone
-  **ES** Neogen® Reveal® Q+ MAX para zearalenona
-  **PT** Neogen® Reveal® Q+ MAX para Zearalenona

## Product Instructions

### Reveal Q+ MAX for Zearalenone

*Quantitative Test with Aqueous Extraction*

#### The Toxin

Zearalenone is primarily produced by the mold *Fusarium graminearum*, which also commonly produces deoxynivalenol (DON). Hence, there is evidence that if zearalenone is detected, there is a high probability that other fusarial mycotoxins may be present. Zearalenone is classified as an estrogenic mycotoxin because it frequently causes estrogenic responses in animals.

When zearalenone-contaminated feed or grain is eaten by livestock, it can cause a wide variety of reproductive problems. In swine, it causes vulvovaginitis, low birth weights, fetal reabsorption, aborted pregnancies, reduced litter sizes, abnormal estrus and feminization of immature males. Zearalenone can delay the breeding process and cost the producer significant economic and physical losses.

Livestock producers are becoming increasingly aware of zearalenone problems and have looked for ways to reduce risks related to contaminated feed.

The best protection against mycotoxins is monitoring for their presence in feeds and foods. That means testing all along the pathway from initial harvest of grains to the finished product.

#### Intended Use/User

Reveal Q+ MAX for Zearalenone is intended for the quantitative analysis of corn, wheat and other commodities for zearalenone. This test utilizes an aqueous extraction, thus eliminating the need for the disposal of hazardous waste. The test kit is designed for use by quality control personnel and others familiar with commodities possibly contaminated by zearalenone.

#### User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at [neogen.com](http://neogen.com), or contact your local Neogen representative or authorized distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, laboratory technique and the sample itself may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any Neogen Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

## Assay Principles

Reveal Q+ MAX for Zearalenone is a single-step lateral flow immunoassay. In the test, the extract is wicked through the reagent zone, which contains antibodies specific for zearalenone conjugated to colloidal gold particles (gold complex). If zearalenone is present, it will be captured by the gold complex. The gold complex, along with any free gold-complex, is then wicked onto a membrane, which contains a zone of zearalenone conjugated to a protein carrier. This zone captures any unbound gold complex, allowing the particles to concentrate and form a visible line.

As the level of zearalenone in a sample increases, free zearalenone will bind with the gold complex allowing less gold complex to be captured in the test zones. Therefore, as the concentration of zearalenone in the sample increases, the test line density decreases. Algorithms programmed into the readers convert the line intensities into a quantitative result displayed in parts per billion (ppb). The membrane also contains a control zone where an immune complex present in the reagent zone is captured by an antibody, forming a visible line. The control line will always form regardless of the presence of zearalenone, ensuring the strip is functioning properly.

## Storage Requirements

Store kit components at room temperature (18–30°C, 64–86°F). Test strips should remain capped in their original tubes until used.

## Materials Provided

1. 25 Reveal Q+ MAX for Zearalenone test strips
2. 25 red conical sample dilution cups
3. 25 clear sample cups
4. 2 bottles of red-labeled sample diluent
5. 1 bottle of yellow-labeled dilution buffer (for samples > 800 ppb)
6. 25 MAX 1 Aqueous Extraction packets

## Materials Recommended, Not Provided

1. Sample collection cups with lids (700004011 | 9428, 700004012 | 9428B)
2. Agri-Grind grinder or equivalent (100001350 | 9401, 700004021 | 9453)
3. Microcentrifuge tubes (700003932 | 9172)
4. Mini centrifuge (700003963 | 9330)
5. Filter syringes (700002724 | 9420)
6. Sample collection tubes with caps (700002726 | 9421, 700002727 | 9421B)
7. Whatman™ #4 filter paper or equivalent (700004035 | 9519, 700006493 | 9429)
8. Pipettor, 250 µL (700004072 | 9692)
9. Pipettor, 100 µL (100001422 | 9860, 100001330 | 9272)
10. Pipette tips, 100 µL to 1 mL (100001373 | 9464, 100001378 | 9487)
11. Scale capable of weighing 5–50 g ±0.1 g (700004010 | 9427)
12. Timer (100001358 | 9426)
13. Reveal sample cup rack (700002734 | 9475)
14. Dispensing pump or graduated cylinder (700004020 | 9448, 100001367 | 9447)
15. Distilled or deionized water
16. Raptor® Integrated Analysis Platform (700002783 | 9680)
17. Raptor cartridges (700002786 | 9681)
18. Raptor 400 µL Exact Volume pipettes (700002787 | 9682)

## Precautions

1. The test strips must remain inside the stay-dry tube before use.
2. Store test kit at room temperature (18–30°C, 64–86°F) when not in use. Do not freeze.
3. Do not use kit components beyond expiration date.
4. Treat all used liquids, including sample extract, and labware as if contaminated with zearalenone. Gloves and other protective apparel should be worn at all times.
5. Commodity extracts should have a pH of 6–8 before testing. Excessively acidic or alkaline samples should be adjusted. For instructions on adjusting pH, contact Neogen Technical Services.
6. Please refer to the SDS information available at [neogen.com](http://neogen.com) for complete safety information on the assay and components.

## Sample Preparation

The sample to be tested should be collected according to accepted sampling techniques (see FGIS sampling protocol or contact your Neogen representative). Obtain a representative sample (minimum 100 g). Grind the sample so at least 95% of the ground material passes through a 20 mesh sieve (about the particle size of fine espresso).

## Test Procedure — Raptor Integrated Analysis Platform

### Sample Extraction

1. Weigh out  $10 \pm 0.1$  g of sample into extraction cup.
2. Add contents of one (1) MAX 1 aqueous extraction packet to the extraction cup.
3. Add 50 mL distilled or deionized water to the extraction cup.
4. Vigorously shake, using hand or mechanical means, for 3 minutes, or blend for 1 minute.
5. Allow the sample to settle, then filter with a filter syringe or Whatman No. 4 filter paper, collecting a minimum of 3 mL filtrate into a sample collection tube. Or, you may also pipette 1 mL of sample into a 2.0 mL microcentrifuge tube, and centrifuge for 1 minute using a microcentrifuge (approx. 2,000 x g).

### Testing Procedure — Corn or Wheat

1. Place the appropriate number of red sample dilution cups into a sample cup rack. Label cups if necessary.
2. Add 250  $\mu$ L of sample diluent (red-labeled bottle) to each red sample dilution cup.
3. Add 250  $\mu$ L of sample extract to the red sample dilution cup containing the sample diluent. Mix by pipetting up and down 5 times. Continue to step 4 below.
4. Fully insert a Reveal Q+ MAX for Zearalenone test strip into a Raptor cartridge.
5. Insert the Raptor cartridge containing the test strip into any of the three ports within the Raptor.
  - a. The bar code on the test strip will be read. The system identifies the type of test strip and the lot number. If the lot number is not found in the system, the bar code reader on the front of the Raptor will turn on automatically.
  - b. Scan the QR code found on the tube containing the test strips. The information will be stored on the reader. (The QR code can also be found on the kit Certificate of Analysis.)
6. Enter the Sample ID if desired.
7. Select Curve 1 for wheat or Curve 2 for corn.
8. Add 400  $\mu$ L of sample from the red sample dilution cup to the Raptor cartridge. The Raptor system will start automatically. Additional samples can be started in the other ports while the first sample is processing.
9. Results will be displayed on the Raptor screen after the 5-minute incubation is complete.

## Dilution Procedure — Raptor Integrated Analysis Platform

Samples greater than 500 ppb will need to be diluted and retested.

### Corn or Wheat

1. Add 100 µL of sample filtrate to a sample collection tube.
2. Add 500 µL of dilution buffer (yellow-labeled bottle) to the sample collection tube. Mix by pipetting up and down 5 times.
3. Add 250 µL of sample diluent (red-labeled bottle) to a red sample dilution cup.
4. Add 250 µL of sample extract to the red sample dilution cup containing the sample diluent. Mix by pipetting up and down 5 times. Continue to step 5 in the following section.
5. Fully insert a Reveal Q+ MAX for Zearalenone test strip into a Raptor cartridge.
6. Insert the Raptor cartridge containing the test strip into any of the three ports within the Raptor.
  - a. The bar code on the test strip will be read. The system identifies the type of test strip and the lot number. If the lot number is not found in the system, the bar code reader on the front of the Raptor will turn on automatically.
  - b. Scan the QR code found on the tube containing the test strips. The information will be stored on the reader. (The QR code can also be found on the kit Certificate of Analysis.)
7. Enter the Sample ID if desired.
8. Select Curve 1 for wheat or Curve 2 for corn.
9. Add 400 µL of sample from the red sample dilution cup to the Raptor cartridge. The Raptor system will start automatically. Additional samples can be started in the other ports while the first sample is processing.
10. Results will be displayed on the Raptor screen after the 5-minute incubation is complete.
11. Multiply your result by 6.

## Performance Characteristics

### Corn

Limit of detection: 13.6 ppb

Range of detection: 50–500 ppb

### Wheat

Limit of detection: 17.0 ppb

Range of detection: 25–500 ppb

## Validated/Verified Matrices

Wheat and corn

**Note:** Neogen continues to validate new commodities. Please contact a representative for the latest validated commodity list.

## Customer Service

Neogen Customer and Technical Services can be contacted through [neogen.com](http://neogen.com) and product training is available by request.

## Safety Data Sheets (SDS) Information Available

SDS are available for all test kits at [neogen.com](http://neogen.com) or by calling 800.234.5333 or 517.372.9200. Dispose of in compliance with all local and national regulations.

## Terms and Conditions

Neogen's full terms and conditions are available [online](#).

## Warranty

Neogen makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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