**U.S. English**

*Product Number:*  
8385

***Reveal*® **Q+****

*for DON*

*Quantitative Test*

*FGIS 2021-151*

*Revised July 2021*

# Reveal<sup>®</sup> Q+ for DON

Product Number: 8385

## The Toxin

Deoxynivalenol (DON) most commonly is produced by the pink mold *Fusarium graminearum*. DON, a member of the trichothecene family, is produced by fungi living on cereal commodities such as wheat, corn, barley and ensilages. The toxicological effects attributed to DON include nausea (vomiting), feed refusal, gastroenteritis, diarrhea, immunosuppression and blood disorders.

Pigs have been shown to be highly sensitive to DON. They will refuse to eat feeds when DON levels of  $\geq 1$  parts per million (ppm) are present. DON and its analogs cause toxic effects in other species as well, with varying degrees of sensitivity. DON has been implicated in causing problems in processed food, including off-flavor in ready-to-eat cereals and adverse effects on dough quality. Accurate determination of the presence of the toxin is of major importance to those monitoring the quality of feed and food in which DON may occur. Testing these commodities for the toxin requires careful sampling, extraction, sanitation, and quantitative analysis.

### The FDA has issued advisory levels for DON as follows:

For	Level	Commodities
Humans	1 ppm	Finished wheat products (flour, bran, and germ)
Ruminating beef, feedlot cattle, and chickens	10 ppm in < 50% of diet (5 ppm total diet)	All grains, grain by-products
Swine	5 ppm in < 20% of diet (1 ppm total diet)	All grains, grain by-products
All other animals	5 ppm in < 40% of diet (2 ppm total diet)	All grains, grain by-products

### European Union regulations for DON as follows:

Foodstuffs	Maximum levels parts per billion (ppb)
Unprocessed cereals other than durum wheat, oats, and maize	1250 ppb
Unprocessed durum wheat and oats	1750 ppb
Cereal flour, including maize flour, maize grits, and maize meal	750 ppb
Bread, pastries, biscuits, cereal snacks, and breakfast cereals	500 ppb
Processed cereal-based foods for infants and young children and baby food	200 ppb

## Intended Use or User

Reveal<sup>®</sup> Q+ for DON is intended for the quantitative analysis of DON levels in grain and grain products, such as wheat, barley, corn, oats, and soybeans. The test kit is designed for use by quality control personnel and others familiar with commodities possibly contaminated by DON.

## Assay Principles

Reveal Q+ for DON is a single-step lateral flow immunochromatographic assay based on a competitive immunoassay format. The extract is wicked through a reagent zone, which contains antibodies specific for DON conjugated to colloidal gold particles. If DON is present, it will be captured by the particle-antibody complex. The DON-antibody-particle complex then is wicked onto a membrane, which contains a zone of DON conjugated to a protein carrier. This zone captures any uncomplexed DON antibody, allowing the particles to concentrate and form a visible line. As the level of DON in a sample increases, free DON will complex with the antibody-gold particles. This, in turn, allows less antibody-gold to be captured in the test zone. Therefore, as the concentration of DON in the sample increases, the test line density decreases. Algorithms programmed into the reader convert these line densities into a quantitative result displayed in parts per million (ppm). The membrane also contains a control zone where an immune complex present in the reagent zone is captured by an antibody, forming a visible line. The control line will always form regardless of the presence of DON, ensuring the strip is functioning properly.

## Storage Requirements

Store kit components at room temperature (18–30°C, 64–86°F) to ensure full shelf life. Test strips should remain capped in their original tubes until used to ensure optimal performance.

## Materials Provided

1. 25 Reveal Q+ for DON test strips
2. 25 red sample dilution cups
3. 25 clear sample cups
4. 2 bottles of sample diluent
5. Instructions for use

## Materials Recommended but not Provided

1. Extraction materials
  - a. Sample collection cups with lids 125 mL (NEOGEN item 9428, 9428B)
  - b. Sample collection tubes with caps (NEOGEN item 9421, 9421B)
  - c. Filter syringes, Whatman #4 filter paper or equivalent (NEOGEN item 9420, 9429, 9519)
  - d. Centrifuge, mini (NEOGEN item 9330)
  - e. Microcentrifuge tubes (NEOGEN item 9172)
  - f. Dispensing pump or graduated cylinder (NEOGEN item 9448, 9447)
  - g. Distilled or deionized water
2. Agri-Grind grinder or equivalent (NEOGEN item 9401, 9453)
3. Scale capable of weighing 5–50 g ± 0.1 g (NEOGEN item 9427)
4. Timer (NEOGEN item 9426)
5. Reveal sample cup rack (NEOGEN item 9475)
6. Pipettor, 100 µL (NEOGEN item 9860)
7. Pipette tips, 1–200 µL (NEOGEN items 9407, 9410, 9417)
8. Pipettor, 100–1,000 µL (NEOGEN items 9463)
9. Pipette tips, 100–1,000 µL (NEOGEN items 9464, 9487)
10. AccuScan<sup>®</sup> Gold reader or AccuScan Pro reader
11. Raptor<sup>®</sup> Integrated Analysis System (NEOGEN item 9680 or 9696)
  - a. Raptor cartridges (NEOGEN item 9681)

NEOGEN.com

## Precautions

1. The test strips must remain inside the stay-dry tube before use.
2. Store test kit at room temperature (18–30°C, 64–86°F) when not in use. Do not freeze.
3. Do not use kit components beyond expiration date listed on kit box label.
4. Treat all used liquids, including sample extract, and labware as if contaminated with DON.
5. Ensure the device lot number and the curve details match the lot ID number selected on the reader. Failure to update the lot-specific QR code within the AccuScan Pro and AccuScan Gold reader will cause inaccurate results.
6. Commodity extracts should have a pH of 6–8 before testing. Excessively acidic or alkaline samples should be adjusted. For instructions on adjusting pH, contact a NEOGEN representative or Technical Services.

## AccuScan Reader Set Up

### AccuScan Pro reader

1. Enter the lot-specific QR code by selecting the QR code icon on the reader. Place the QR code into the cartridge and insert the cartridge into the reader.  
**Note:** For instructions on manually entering sample IDs, see the AccuScan Pro user manual.
2. Return to the home screen and select the test strip icon. Touch the mycotoxin category, then select the DON Q+ test type.

### AccuScan Gold reader

1. Enter the lot specific QR code by selecting scan QR from the main screen. Place the lot specific QR code into the white cartridge adapter labeled Cal/QR and place the cartridge into the reader.
2. The valid code will be scanned by the reader and provide information on the lot number and expiry date. Verify this information is correct and then add the lot ID to the reader by pressing add lot ID.  
**Note:** The lot ID for the current lot will now be stored with the test ID (e.g., aflatoxin, DON) and can be selected when running a test.

## Sample Extraction

1. Extract at a ratio of 1 part sample to 10 parts distilled or deionized water. For example, combine 10 g of ground sample with 100 mL of distilled or deionized water.
2. Vigorously shake, using hand or mechanical means (250 rpm) for 3 minutes, or blend for 1 minute
3. Allow the sample to settle, then filter at least 4 mL with a filter syringe, or Whatman No. 1 filter paper.
4. You may also pipette 1 mL of sample into a 2.0 mL microcentrifuge tube, and centrifuge for 30 seconds.
5. The sample is now ready for testing.

## Sample Extraction (AOAC Method)

1. Combine 50 g of ground sample with 250 mL of distilled or deionized water and shake vigorously by mechanical shaker (250 rpm) for 3 minutes. Allow the sample to settle.
2. Filter the extract by pouring 4 mL through a NEOGEN filter syringe and collecting the filtrate.
3. Dilute the filtrate 1:2 (1 + 1) with distilled or deionized water (e.g., add 1 mL of extract to 1 mL of distilled or deionized water). Vortex for a few seconds before testing.

## Sample Extraction — 50 g (FGIS Method) — Raptor Integrated Analysis Platform

1. Combine 50 g of ground sample with 250 mL of distilled or deionized water into a Whirl-Pak bag. Shake vigorously by mechanical shaker (250 rpm) or by hand with similar shaking action, for 3 minutes. Allow the sample to settle.
2. Filter the extract by pouring 3 mL through a NEOGEN filter syringe and collecting the filtrate.
3. Dilute the filtrate 3:2 by adding 900 µL of filtered extract to 600 µL of distilled or deionized water in a new collection tube. Vortex for 10 seconds. This is the diluted filtered extract and is ready for testing following the instructions in Test Procedure — 50 g (FGIS Method) — Raptor Integrated Analysis Platform.

## Test Procedure — AccuScan Reader


1. Place the appropriate number of red sample dilution cups and clear sample cups into a sample cup rack. Label cups if necessary.
2. Add 100 µL of sample extract to each red sample dilution cup.
3. Add 1,000 µL of sample diluent to red sample dilution cup. Mix by pipetting up and down 5 times.
4. Using a new pipette tip, transfer 100 µL of diluted sample extract into a clear sample cup.
5. Place a new Reveal Q+ for DON test strip with the sample end down into the sample cup and set timer for 3 minutes. Ensure the test strip comes into contact with liquid and begins to wick.
6. Remove the strip from the sample cup after it has developed for exactly 3 minutes. Immediately (within 30 seconds) read test strip in the AccuScan reader.
7. For the AccuScan PRO reader:                 Select category: Mycotoxin  
  Test name: DON Q+  
  
For the AccuScan GOLD reader:                 Select category: Mycotoxin Q Plus  
  Test name: DON Q+

## Dilution Procedure — AccuScan Reader

Samples greater than 6 ppm will need to be diluted and retested.

1. Add 100 µL sample filtrate to a sample collection tube (Use sample filtrate obtained in Sample Extraction Procedure from page 4, step 5, or from Sample Extraction Procedure (FGIS and AOAC Method) on page 5, step 3).
2. Add 500 µL distilled or deionized water to the sample collection tube. Mix well by pipetting up and down 5 times, or cap the tube and shake.
3. Add 100 µL of diluted sample extract (from step 2) to each red sample dilution cup
4. Add 1,000 µL of sample diluent (blue labeled bottle) to the red dilution cup with sample extract. Mix by pipetting up and down 5 times.
5. Using a new pipette tip, transfer 100 µL of diluted sample extract into a clear sample cup
6. Place a new Reveal Q+ for DON test strip with the sample end down into the sample cup and set timer for 3 minutes. Ensure the test strip comes into contact with liquid and begins to wick.
7. Remove the strip from the sample cup after it has developed for exactly 3 minutes. Immediately (within 30 seconds) read test strip in the AccuScan reader.
8. For the AccuScan PRO reader:                 Select category: Mycotoxin  
  Test name: DON Q+  
  
For the AccuScan GOLD reader:                 Select category: Mycotoxin Q Plus  
  Test name: DON Q+

**Note:** The AccuScan reader will not calculate your dilution. Final result displayed will need to be multiplied by 6.


2

3

4

4

5

## Reading Test Results — AccuScan Reader

**Note:** Test strips must be read within 30 seconds of completion of the 3-minute incubation. Refer to AccuScan Reader Set Up section for test selection and setup information.

1. Select the assay type from the menu and ensure the device lot number matches the lot ID number selected on the reader. **Note:** Failure to update the lot-specific QR code will cause inaccurate results.
2. Fully insert the Reveal Q+ test strip into the R-labeled cartridge adapter with the sample end first and results facing out.
3. Insert the cartridge with test strip upside-down into the AccuScan GOLD reader (the test lines will face downward into the reader). If using the AccuScan Pro reader, insert the cartridge right side up. The reader will automatically begin analyzing the cartridge. **Caution:** Removing cartridge prior to completion can result in invalid readings.
4. The AccuScan reader will analyze the test strip and results will be displayed and stored in the reader.

## Notes

1. Ensure device is fully inserted into cartridge.
2. Readings should be made within 30 seconds of the 3-minute incubation time. Readings after 3.5 minutes may be inaccurate due to over-development of the device.

## Test Procedure — Raptor Integrated Analysis Platform

1. Place the appropriate number of red sample dilution cups and clear sample cups into a sample cup rack. Label cups if necessary.
2. Add 100  $\mu$ L of sample extract to each red sample dilution cup.
3. Add 1,000  $\mu$ L of sample diluent to red sample dilution cup. Mix by pipetting up and down 5 times.
4. Fully insert a Reveal DON Q+ test strip into a Raptor cartridge.
5. Insert the Raptor cartridge containing the test strip into any of the three ports within the Raptor.
  - a. The barcode on the test strip will be read. The system identifies the type of test strip and the lot number. If the lot number is not found in the system, the barcode reader in the front of the Raptor will turn on automatically.
  - b. Scan the QR code found on the tube containing the test strips. The information will be stored on the reader.
6. Enter sample ID if desired.
7. Add 400  $\mu$ L of sample extract from the red sample dilution cup to the Raptor cartridge.
  - a. The Raptor system will start automatically.
  - b. Additional samples can be started in the other ports while the first sample is processing.
8. Results will be displayed on the Raptor screen after the 3-minute testing period is complete.

## Dilution Procedure — Raptor Integrated Analysis Platform

Samples greater than 6 ppm will need to be diluted and retested.

1. Add 100 µL sample filtrate to a sample collection tube (Use sample filtrate obtained in Sample Extraction Procedure from page 4, step 5, or from Sample Extraction Procedure(AOAC Method) on page 4, step 3).
2. Add 500 µL distilled or deionized water to the sample collection tube. Mix well by pipetting up and down 5 times, or cap the tube and shake.
3. Add 100 µL of diluted sample extract (from step 2) to each red sample dilution cup.
4. Add 1,000 µL of sample diluent (blue labeled bottle) to the red dilution cup with sample extract. Mix by pipetting up and down 5 times.
5. Insert the Raptor cartridge containing the test strip into any of the three ports within the Raptor.
  - a. The barcode on the test strip will be read. The system identifies the type of test strip and the lot number. If the lot number is not found in the system, the barcode reader in the front of the Raptor will turn on automatically.
  - b. Scan the QR code found on the tube containing the test strips. The information will be stored on the reader.
6. Enter sample ID if desired.
7. Add 400 µL of sample extract from the red sample dilution cup to the Raptor cartridge.
  - a. The Raptor system will start automatically.
  - b. Additional samples can be started in the other ports while the first sample is processing.
8. Results will be displayed on the Raptor screen after the 3-minute testing period is complete.

**Note:** The Raptor reader will not calculate your dilution. Final result displayed will need to be multiplied by 6.

## Test Procedure — 50 g (FGIS Method) — Raptor Integrated Analysis Platform

1. Place the appropriate number of red sample dilution cups and clear sample cups into a sample cup rack. Label cups if necessary.
2. Add 1000 µL of sample diluent to a red sample dilution cup.
3. Add 100 µL of the diluted filtered extract (from sample extraction procedure FGIS Method- step 3) to the red sample dilution cup. Mix by pipetting up and down 5 times.
4. Fully insert a Reveal DON Q+ test strip into a Raptor cartridge.
5. Insert the Raptor cartridge containing the test strip into any of the three ports within the Raptor.
  - a. The bar code on the test strip will be read. The system identifies the type of test strip and the lot number. If the lot number is not found in the system, the bar code reader in the front of the Raptor will turn on automatically.
  - b. Scan the QR code found on the QR Code Card in the kit (FGIS 50g). Do not use the QR code on the device tube (used for the 10 g sample size).
6. Enter sample ID if desired.
7. Add 400 µL of sample extract from the red sample dilution cup to the Raptor cartridge.
  - a. The Raptor system will start automatically.
  - b. Additional samples can be started in the other ports while the first sample is processing.
8. Results will be displayed on the Raptor screen after the 4-minute testing period is complete.

## Dilution Procedure — 50 g (FGIS Method) — Raptor Integrated Analysis Platform

Samples greater than 5 ppm will need to be diluted and retested.

1. Add 100 µL diluted filtered extract (from step 3 of Sample Extraction — 50 g (FGIS Method) Raptor Integrated Analysis Platform procedure) to a sample collection tube.
2. Add 1400 µL distilled or deionized water to the sample collection tube. Mix well. Capped samples should be analyzed within 4 hours. Vortex for a few seconds prior to analysis.
3. Follow the same test procedure detailed in Test Procedure — 50 g (FGIS method) — Raptor Integrated Analysis Platform
4. Multiply the results by 15.

## Performance Characteristics

- 1 Limit of detection: 0.1 ppm
- 2 Range of quantitation: 0.3–6 ppm

**Note:** Samples greater than 6 ppm must be diluted and retested. See dilution procedure on page 5 (AccuScan GOLD) or page 7 (Raptor).

Samples less than the range of quantitation should be reported as less than 0.3 ppm.

## Validated/Varified Matrices

NEOGEN continues to validate new commodities. Please contact a representative for the latest validated commodity list.

## Customer Service

NEOGEN Customer and Technical Services can be contacted through [NEOGEN.com](http://NEOGEN.com) and product training is available by request.

## SDS Information Available

Safety data sheets are available for all test kits at [NEOGEN.com](http://NEOGEN.com) or by calling 800.234.5333 or 517.372.9200.

## Terms and Conditions

NEOGEN's full terms and conditions are available [online](#).

## Warranty

NEOGEN makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, NEOGEN will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product, or of the fitness of the product for any purpose. NEOGEN shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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