Read instructions carefully before starting test

Reveal[®] Q+ MAX

for T-2/HT-2

Quantitative Test with Aqueous Extraction

THE TOXIN

T-2/HT-2 toxins are trichothecene mycotoxins produced by several species of *Fusarium* molds. Since T-2 toxin is readily metabolized to HT-2 toxin and the toxins have been shown to produce numerous adverse effects on many animals, these two mycotoxins are frequently evaluated together.

Animals affected by the toxins include swine, dairy cattle, poultry, dogs, cats and horses. Effects of the toxins include digestive disorders, hemorrhage, edema, oral lesions, dermatitis, and blood disorders. Damage caused by the toxins to the digestive track is irreversible. In the most severe cases, these toxins will cause death. T-2 toxin is the principal causal toxin in the human disease alimentary toxic aleukia.

Poultry studies have shown T-2 intoxication has led to a reduction in weight gain and other problems such as beak lesions, poor feathering, motor function impairment and increased susceptibility to *Salmonella* spp.

The best protection against these mycotoxins is monitoring for their presence in feeds and foods. That means testing all along the pathway from initial harvest of grains to the finished product.

INTENDED USE/USER

Reveal Q+ MAX for T-2/HT-2 is intended for the quantitative analysis of corn and corn products for T-2/HT-2 toxins. This test utilizes a common aqueous extraction used with all Q+ MAX kits allowing the user to test for multiple mycotoxins with one extraction. The test kit is designed for use by quality control personnel and others familiar with commodities possibly contaminated by T-2/HT-2.

ASSAY PRINCIPLES

Reveal Q+MAX for T-2/HT-2 is a single-step lateral flow immunoassay. In the test, the extract is wicked through the reagent zone, which contains antibodies specific for T-2/HT-2 conjugated to colloidal gold particles (gold complex). If T-2/HT-2 is present, it will be captured by the gold complex. The gold complex, along with any free gold-complex, is then wicked onto a membrane, which contains a zone of T-2/HT-2 conjugated to a protein carrier. This zone captures any unbound gold complex, allowing the particles to concentrate and form a visible line. As the level of T-2/HT-2 in a sample increases, free T-2/HT-2 will bind with the gold complex allowing less gold complex to be captured in the test zones. Therefore, as the concentration of T-2/HT-2 in the sample increases, the test line density decreases. Algorithms programmed into the AccuScan® readers convert the line intensities into a quantitative result displayed in parts per billion (ppb). The membrane also contains a control zone where an immune complex present in the reagent zone is captured by an antibody, forming a visible line. The control line will always form regardless of the presence of T-2/HT-2, ensuring the strip is functioning properly.

STORAGE REQUIREMENTS

Store kit components at room temperature (18–30°C / 64–86°F) to ensure full shelf life. Test strips should remain capped in their original tubes until used to ensure optimal performance.

MATERIALS PROVIDED

1. 25 Reveal Q+ MAX for T-2/HT-2 test strips
2. 25 red sample dilution cups
3. 25 clear sample cups
4. 2 bottles of green labeled sample diluent
5. 25 MAX 1 Aqueous Extraction packets
6. Instructions for use

MATERIALS RECOMMENDED BUT NOT PROVIDED

1. Extraction Materials
 - a. Sample collection cups with lids—125 mL (Neogen item 9428, 9428B)
 - b. Sample collection tubes with caps (Neogen item 9421, 9421B)
 - c. Filter syringes, Whatman #4 filter paper or equivalent (Neogen item 9420, 9429, 9519)**OR**
 - d. Centrifuge, mini (Neogen item 9330)
 - e. Microcentrifuge tubes (Neogen item 9172)
 - f. Dispensing pump or graduated cylinder (Neogen item 9448, 9447)
2. Agri-Grind grinder or equivalent (Neogen item 9401, 9453)
3. Scale capable of weighing 5–50 g ± 0.1 g (Neogen item 9427)
4. Timer (Neogen item 9426)
5. Reveal sample cup rack (Neogen item 9475)
6. Pipettor, 100 µL (Neogen item 9272, 9860)
7. Pipette tips, 1–200 µL (Neogen item 9407,9410))
8. Pipettor, 100–1000 µL (Neogen item 9463)
9. Pipette tips, 100 µL–1 mL (Neogen item 9464,9487)
10. Distilled or deionized water
11. AccuScan Pro (Neogen item 9565), or AccuScan Gold (Neogen item 9595) or Raptor Integrated Analysis Platform (Neogen item 9680)
12. Raptor cartridges (Neogen item 9681)
13. Raptor 400 µL exact volume pipettes (Neogen item 9682)
14. MAX 1 Aqueous Extraction packets (Neogen item 8089)

SAMPLE PREPARATION

The sample to be tested should be collected according to accepted sampling techniques (see FGIS sampling protocol or contact your Neogen representative). Obtain a representative sample (minimum 100 g). Grind the sample so at least 95% of the ground material passes through a 20 mesh sieve (about the particle size of fine espresso).

PRECAUTIONS

1. The test strips must remain inside the stay-dry tube before use.
2. Store test kit at room temperature (18–30°C / 64–86°F) when not in use. Do not freeze.
3. Do not use kit components beyond expiration date.
4. Treat all used liquids, including sample extract, and labware as if contaminated with T-2/HT-2 toxins. Gloves and other protective apparel should be worn at all times.
5. Ensure the device lot number and the curve details match the lot ID number selected on the reader. Failure to update the lot-specific QR code within the AccuScan Pro and AccuScan Gold reader will cause inaccurate results.
6. Commodity extracts should have a pH of 6–8 before testing. Excessively acidic or alkaline samples should be adjusted. For instructions on adjusting pH contact a Neogen representative or Technical Services.

ACCUSCAN GOLD READER SETUP

1. Enter the lot specific QR code by selecting **Scan QR** from the main screen. Place the lot specific QR code into the white cartridge adapter labeled Cal/QR and place the cartridge into the reader.
2. The valid code will be scanned by the reader and provide information on the lot number and expiry date. Verify this information is correct and then add the lot ID to the reader by pressing **Add Lot ID**.

NOTE: The lot ID for the current lot will now be stored with the test ID (e.g., T-2/HT-2) and can be selected when running a test.

ACCUSCAN PRO READER SETUP

1. Enter the lot-specific QR code by selecting the QR code icon on the reader. Place the QR code into the cartridge and insert the cartridge into the reader.

NOTE: For instructions on manually entering sample IDs, see the AccuScan Pro user manual.

2. Return to the home screen and select the test strip icon. Touch the Mycotoxin Q+ MAX category, then select the **T2HT2Q+MAX** test type.

SAMPLE EXTRACTION

1. Weigh out 10 g \pm 0.1 g of sample into extraction cup.
2. Add contents of one (1) MAX 1 aqueous extraction packet to the extraction cup.
3. Add 50 mL distilled or deionized water to the extraction cup.
4. Vigorously shake, using hand or mechanical means, for 3 minutes, or blend for 1 minute.
5. Allow the sample to settle, then filter with a filter syringe or Whatman #4 filter paper to collect a minimum of 3 mL filtrate into a sample collection tube. Or, you may also pipette 1 mL of sample into a 2.0 mL micro-centrifuge tube, and centrifuge for 30 seconds using a microcentrifuge (approx. 2,000 x g).

TEST PROCEDURE

1. Place the appropriate number of red sample dilution cups and clear sample cups into a sample cup rack. Label cups if necessary.
2. **Add 100 µL of sample extract to a red sample dilution cup.**
3. Add 1500 µL of green-labeled sample diluent to the red dilution cup with sample extract. Mix by pipetting up and down 5 times.
4. Using a new pipette tip, transfer 100 µL of diluted sample extract into a new **clear** sample cup.
5. Place a new Reveal Q+ MAX for T-2/HT-2 test strip with the sample end down into the clear sample cup and set timer for **5 minutes**. Ensure the test strip comes into contact with liquid and begins to wick.
6. Remove the strip from the sample cup after it has developed for **5 minutes** and read immediately (within 30 seconds).
7. **For the AccuScan Gold reader: Select CATEGORY: Mycotoxin Q+ MAX
TEST NAME: Q+M T2HT2**

**For the AccuScan Pro reader: Select CATEGORY: Mycotoxin Q+ MAX
TEST NAME: T2HT2Q+ MAX**

DILUTION PROCEDURE

Samples **greater than 500 ppb** will need to be diluted and retested.

- a. Add 100 µL sample filtrate to a sample collection tube.
- b. Add 400 µL distilled or deionized water to the sample collection tube. Mix by pipetting up and down 5 times.
- c. **Add 100 µL of diluted sample extract from Step b to a red sample dilution cup.**
- d. Add 1500 µL of green-labeled sample diluent to the red dilution cup with sample extract. Mix by pipetting up and down 5 times.
- e. Using a new pipette tip, transfer 100 µL of diluted sample extract into a new **clear** sample cup.
- f. Place a new Reveal Q+ MAX for T-2/HT-2 test strip with the sample end down into the clear sample cup and set timer for **5 minutes**. Ensure the test strip comes into contact with liquid and begins to wick.
- g. Remove the strip from the sample cup after it has developed for **5 minutes** and read immediately (within 30 seconds).
- h. **For the AccuScan Gold reader: Select CATEGORY: Mycotoxin Q+ MAX
TEST NAME: Q+M T2HT2**

**For the AccuScan Pro reader: Select CATEGORY: Mycotoxin Q+ MAX
TEST NAME: T2HT2Q+ MAX**
- i. Final result displayed will need to be multiplied by 5.

READING TEST RESULTS

NOTE: Test strips should be read within **30 seconds** of completion of the 5 minute incubation. Refer to **AccuScan Reader Set Up** for test selection and set up information.

1. Select the assay type (e.g., T-2/HT-2) from the menu and ensure the device lot number matches the lot ID number selected on the reader.

NOTE: Failure to update the lot-specific QR code will cause inaccurate results.

2. Fully insert the Reveal Q+ MAX test strip into the R-labeled cartridge adapter with the sample end first and results facing out.
3. Insert the cartridge with test strip upsidedown into the AccuScan Gold reader (the test lines will face downward into the reader). If using the AccuScan Pro reader, insert the cartridge right side up. The reader will automatically begin analyzing the cartridge.

CAUTION: Removing cartridge prior to completion can result in invalid readings.

4. The AccuScan reader will analyze the test strip and results will be displayed and stored in the reader.


NOTES:

1. Ensure device is fully inserted into cartridge.
2. Readings should be made between **5 and 5.5 minutes**. Readings after **5.5 minutes** may be inaccurate due to overdevelopment of the device.
3. The strips must be read using Neogen's AccuScan Gold or AccuScan Pro readers.

TEST PROCEDURE — RAPTOR INTEGRATED ANALYSIS PLATFORM

1. Place the appropriate number of red sample dilution cups into a sample cup rack. Label cups if necessary.
2. Add 100 μ L of sample extract to each red sample dilution cup.
3. Add 1,500 μ L of green-labeled sample diluent to red sample dilution cup. Mix by pipetting up and down five times.
4. Fully insert a Reveal Q+ MAX for T-2/HT-2 test strip into a Raptor cartridge.
5. Insert the Raptor cartridge containing the test strip into any of the three ports within the Raptor.
 - a. The bar code on the test strip will be read. The system identifies the type of test strip and the lot number. If the lot number is not found in the system, the bar code reader in the front of the Raptor will turn on automatically.
 - b. Scan the QR code found on the tube containing the test strips. The information will be stored on the reader.
6. Enter Sample ID if desired.
7. Add 400 μ L of sample extract from the red sample dilution cup to the Raptor cartridge.
 - a. The Raptor system will start automatically.
 - b. Additional samples can be started in the other ports while the first sample is processing.
8. Results will be displayed on the Raptor screen after the 5-minute testing period is complete.

**DILUTION PROCEDURE - RAPTOR INTEGRATED ANALYSIS PLATFORM**

Samples **greater than 500 ppb** will need to be diluted and re-tested.

1. Add 100 μ L sample filtrate to a sample collection tube.
2. Add 400 μ L distilled or deionized water to the sample collection tube. Mix well by pipetting up and down five times, or cap the tube and shake.
3. Add 100 μ L of diluted sample extract (from step 2) to each red sample dilution cup.
4. Add 1,500 μ L of sample diluent (green labeled bottle) to the red dilution cup with sample extract. Mix by pipetting up and down 5 times.
5. Insert the Raptor cartridge containing the test strip into any of the three ports within the Raptor.
 - a. The bar code on the test strip will be read. The system identifies the type of test strip and the lot number. If the lot number is not found in the system, the bar code reader in the front of the Raptor will turn on automatically.
 - b. Scan the QR code found on the tube containing the test strips. The information will be stored on the reader.
6. Enter Sample ID if desired.
7. Add 400 μ L of sample extract from the red sample dilution cup to the Raptor cartridge.
 - a. The Raptor system will start automatically.
 - b. Additional samples can be started in the other ports while the first sample is processing.
8. Results will be displayed on the Raptor screen after the 5-minute testing period is complete.

NOTE: The Raptor reader will not calculate your dilution. **Final result displayed will need to be multiplied by 5.**

PERFORMANCE CHARACTERISTICS

1. Limit of quantitation: 50 ppb
2. Range of detection: 50–3000 ppb (with dilutions)

NOTE: For quantification of samples above 500 ppb, follow the dilution protocol. Results below the limit of detection should be reported as less than 50 ppb.

VALIDATED MATRICES

Wheat, oats, and corn

NOTE: Neogen continues to validate new commodities. Please contact a representative for the latest validated commodity list.

CUSTOMER SERVICE

Neogen Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all Neogen test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of Neogen's test kits, on Neogen's website at www.neogen.com, or by calling Neogen at 800/234-5333 or 517/372-9200.

TERMS AND CONDITIONS

For Neogen's full terms and conditions, please visit www.neogen.com/en/terms-and-conditions.

WARRANTY

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

TEST KITS AVAILABLE FROM NEOGEN

Natural toxins

- Aflatoxin, DON, ochratoxin, zearalenone, T-2/HT-2 toxins, fumonisin, histamine

Foodborne bacteria

- *E. coli* O157:H7, *Salmonella*, *Listeria*, *Listeria monocytogenes*, *Campylobacter*, *Staphylococcus aureus*, *Salmonella enteritidis*

Sanitation

- ATP, yeast and mold, total plate count, generic *E. coli* and total coliforms, protein residues

Food allergens

- Almonds, crustaceans, eggs, gliadin, hazelnut, milk, mustard, peanuts, sesame, soy, walnuts, multi-treenut

Genetic modification

- CP4 (Roundup Ready[®])

Ruminant by-products

- Meat and bone meal, feed

Species Identification

- Raw and cooked meat samples, feed

**North America
Neogen Headquarters**
800/234-5333 (USA/Canada)
foodsafety@neogen.com
foodsafety.neogen.com

**Europe, Middle East and Africa
Neogen Europe**
+ 44 (0) 1292 525 600
info_uk@neogeneurope.com
www.neogeneurope.com

**Mexico
Neogen Latinoamerica**
+52 (55) 5254-8235
informacion@neogenlac.com
www.neogenlac.com

**Brazil
Neogen do Brasil**
+55 19 3935.3727
info@neogendobrasil.com.br
www.neogendobrasil.com.br

**China
Neogen Bio-Scientific Technology**
+86 21 6271 7013
info@neogenchina.com.cn
www.neogenchina.com.cn

**India
Neogen Food and Animal Security**
+91 484 2306598, 2301582
info@neogenindia.com
www.neogenindia.com

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©Neogen Corporation, 2018. Neogen, AccuScan, Raptor y Reveal son marcas registradas de Neogen Corporation. Todas las otras marcas y nombres de productos mencionados son marcas registradas o marcas de sus respectivas compañías.