
Read the instructions carefully before starting the test.

NEOGEN[®]
ANSR[®]
for *Campylobacter*

INTENDED USE

ANSR[®] for *Campylobacter* has been found to be an effective procedure for detection of *Campylobacter jejuni*, *C. coli*, and *C. lari* in poultry rinse samples and carcass swabs.

ASSAY PRINCIPLES

ANSR for *Campylobacter* is an isothermal, amplified nucleic acid assay. The ANSR for *Campylobacter* method is based on nicking enzyme amplification reaction (NEAR™) technology preceded by the reverse transcription of 23S ribosomal RNA. Target complementary DNA is amplified through a mechanism of polymerization from the ends of nicks created in double-stranded DNA by the action of a specific endonuclease. Amplified target sequences are detected in real time using fluorescent molecular beacon probes.

A 2-stage lysis reaction is performed, first at $37 \pm 2^\circ\text{C}$ for 10 minutes, then at $80 \pm 2^\circ\text{C}$ for 20 minutes. Next, a portion of the lysed sample is transferred to a strip tube containing lyophilized ANSR reagents. The tubes are sealed and incubated at $56 \pm 2^\circ\text{C}$ on the ANSR reader. Results are generated by the reader and displayed in the ANSR software within 18 minutes. Positive results may be confirmed from the enrichment cultures following standard procedures. Each tube of ANSR reagents contains an internal positive control, ensuring that the reagents are functioning properly.

INTENDED USER

The ANSR® for *Campylobacter* test is designed for use by personnel with appropriate training in microbiology. Training in the use of the ANSR test system is available through NEOGEN®.

MATERIALS PROVIDED

1. 12 strips of 8 cluster tubes, 1.2 mL
2. 12 strips of 8 reaction tubes, 200 µL, containing lyophilized ANSR for *Campylobacter* reagents in 2 sealed foil pouches with desiccant pack
3. 12 strips of 8 permanent caps for the reaction tubes
4. 1 bottle lysis reagent suspension buffer, 60 mL
5. 3 vials containing lyophilized lysis reagents
6. 1 kit insert

EQUIPMENT REQUIRED

1. ANSR reader (NEOGEN item 9828)
2. Computer and software for connection to the ANSR reader (NEOGEN item 9832)
3. 1 dual heater block with aluminum block inserts for 1.2 mL cluster tubes, $80 \pm 2^\circ\text{C}$ and $37 \pm 2^\circ\text{C}$ (NEOGEN items 9836-DUAL48, 9829-48) or 2 single heater blocks with aluminum block inserts for 1.2 mL cluster tubes, $80 \pm 2^\circ\text{C}$ and $37 \pm 2^\circ\text{C}$ (NEOGEN items 9836-48D, 9829-48)
4. Pipettor, 100–1000 µL (NEOGEN item 9463)
5. Pipette tip rack, 100–1000 µL, sterile, 96 tips (NEOGEN item 9487)
6. Pipettor, 10–100 µL, 8-channel (NEOGEN item 9388)
7. Pipette tips, 100 µL, sterile, filtered, 96 tips (NEOGEN item 9389)
8. Vortex mixer, adjustable speed (NEOGEN item 9494)
9. 3 thermometers, traceable (NEOGEN item 9518)
10. Timer, 3-channel (NEOGEN item 9426)
11. Optional-for-use heater block with 0.2 mL reaction tube aluminum block insert, $56 \pm 2^\circ\text{C}$ (NEOGEN items 9836-48D, 9829-48)
12. Stomacher (optional)
13. Webcam (NEOGEN item WEBCAM)
14. ANSR ethernet cable (NEOGEN item 9835)
15. 10 mL pipette pump (NEOGEN item 9277)
16. Pipettes, sterile serological (NEOGEN item 8686)
17. 40-slot, 20 mm test tube rack, autoclavable (NEOGEN item 9553)
18. Pipettor, 20–200 mL (NEOGEN item 9276)

OTHER MATERIALS REQUIRED

1. Stomacher-type bags for sample enrichment. Filtered bags are recommended (NEOGEN item 6827)
2. Flat bottomed 4 oz. Whirl-Pak® sampling bag (NEOGEN Item 9508)
3. Graduated cylinder, 250 mL (NEOGEN item 9368)

MEDIA ENRICHMENT BROTH REQUIRED

1. ANSR for *Campylobacter* Enrichment Broth (NEOGEN Item 9818)
2. Buffered Peptone Water (NEOGEN Item 7418A or equivalent)
3. 1 L purified water
4. Phosphate Buffered Saline (PBS) (NEOGEN item 8425)

STORAGE

Store the ANSR® reagents at 2–8°C. After removing reaction tubes from the foil pouch, promptly re-seal the pouch. Leave the desiccant pack in the pouch at all times.

PRECAUTIONS

1. Use good microbiology laboratory practices.
2. Dispose of used pipette tips in a covered container containing a fresh solution of 10% bleach. The 10% bleach solution should be made fresh each day. Use 9 parts water with 1 part commercially available household-strength bleach to make the 10% bleach solution. Undiluted stock solutions of bleach should be used within 30 days after opening.
3. Discard bleach solution and tips as regular waste at the end of each day.
4. Do not use reagents beyond the expiration date.
5. Enrichment media should be pre-warmed to the incubation temperature (all protocols) before use.
6. Use of enrichment media and incubation times or temperatures other than those specified may lead to erroneous results.
7. Remove the reaction tubes from the foil pouch just before use and keep covered until heating process begins. Reseal the pouch containing the remaining reaction tubes to avoid prolonged exposure to light. More than 15 minutes of total exposure time may lead to erroneous results.
8. Do not, under any circumstance, remove caps from reaction tubes after the assay has been started. This is essential in order to prevent accidental contamination of the environment with amplification products.
9. Exercise care in all pipetting steps to avoid cross-contamination of samples.
10. Complete all assay steps in sequence, avoiding delays between steps.
11. Tap reaction tubes on bench top to make sure the lyophilized reagents are at the bottom of the tube prior to adding lysed sample.

PREPARATION OF ENRICHMENT BROTH (2X STRENGTH)

Dissolve 56 g of the medium into 1 L of purified water prewarmed to $42 \pm 1^\circ\text{C}$. Use the rehydrated media the same day as prepared. Do not autoclave.

SAMPLE ENRICHMENT

Poultry rinse

1. Weigh out 30 mL of sample in a 4 oz. Whirl-Pak® bag.
2. Add 30 mL 2x ANSR for *Campylobacter* Enrichment Broth prewarmed to $42 \pm 1^\circ\text{C}$ to the bag. Swirl briefly to mix.
Note: Attempt to minimize airspace in the bag as well as exposure to oxygen with gentle, brief mixing.
3. Incubate the culture under aerobic or microaerophilic conditions at $42 \pm 1^\circ\text{C}$ for 20–24 hours.

Carcass sponge

1. Place the sponge in a Stomacher-type bag.
2. Add 25 mL of sterile Buffered Peptone Water and 25 mL of 2x ANSR for *Campylobacter* Enrichment Broth to the bag.
3. Grasp the bag tightly at top and shake the bag vigorously using a side-to-side motion or squeeze the sponge several times to ensure complete mixing.
4. Incubate the culture under aerobic or microaerophilic conditions at $42 \pm 1^\circ\text{C}$ for 20–24 hours.

LYSIS REAGENT SOLUTION PREPARATION

Reconstitute 1 vial of the lyophilized lysis reagents with 18 mL of the lysis reagent suspension buffer by adding the buffer to the reagent vial. Swirl gently to mix.

Note: 1 vial of the lysis reagents is enough for approximately 32 samples. Prepared lysis reagent solution can be stored at 2–8°C for up to 30 days.

ANSR TEST PROCEDURE

Before starting the assay:

1. Preheat one lysis heater block to $80 \pm 2^\circ\text{C}$. Preheat the second lysis heater block to $37 \pm 2^\circ\text{C}$. If using the optional single heater, preheat to $56 \pm 2^\circ\text{C}$. Use the thermometer for the temperature reading.
2. Remove the foil pouch containing the reaction tubes from the refrigerator and allow them to warm at room temperature for 15 minutes. To avoid excess light exposure, leave the reaction tubes in the foil pouch until they are needed.
Note: Keep the lysis reagent solution in the refrigerator until ready to use.
3. Connect the ANSR® reader to the computer via USB or ethernet and turn the computer on.
4. Turn on the ANSR reader. The reader will preheat to $56 \pm 2^\circ\text{C}$.
5. Start the ANSR software and click the connect button. Input sample IDs, lot number, and user information.
Note: For instructions on using the reader and software, see the user guide that came with the ANSR reader.

ASSAY PROCEDURE

Note: Mix enriched samples well before beginning assay procedure. Enrichment cultures that are highly turbid or highly colored may require a 1:10 dilution in PBS before lysing the sample. Contact NEOGEN Technical Services for more information. The AOAC approved matrices do not require dilution.

Optional: Pre-pierced plate sealing film (BIO RAD Item 3600040) can be placed over the 1.2 mL cluster tubes if desired.

1. Add 50 μL of enrichment culture (or dilution) to each 1.2 mL cluster tube using a pipettor with 100 μL filtered tips.
Note: Cluster tubes may be pulled apart to provide the number of tubes needed.
2. Add 450 μL of the lysis reagent solution to each cluster tube containing culture. Be sure to switch the pipette tips between samples.
Note: Return the lysis reagent solution to the refrigerator after use (within 1 hour).
3. Incubate each cluster tube(s) at $37 \pm 2^\circ\text{C}$ for 10 minutes.
4. Immediately transfer the cluster tube to the $80 \pm 2^\circ\text{C}$ heater block and incubate for 20 minutes.
Note: The $80 \pm 2^\circ\text{C}$ incubation time may be extended to a total of 60 minutes for the purpose of managing staggered assay start times.
5. For 3–5 minutes before the end of the lysis step, preheat the ANSR reagents to $56 \pm 2^\circ\text{C}$ by placing the reaction tubes in the ANSR reader. Optional: A separate heat block can be used. It should be heated to $56 \pm 2^\circ\text{C}$.
Note: The strip of reaction tubes may be cut to provide the number of tubes needed. Keep all unused tubes in the sealed foil pack. Ensure the pellet in each reaction tube is at the bottom by tapping the tubes gently on the bench top.
6. After the completion of the 20–60 minute lysis incubation, remove and discard the caps from each reaction tube in the ANSR® reader.

Important: Proceed with steps 7–9 without delay. The transfer of the sample from the lysis tubes to the reaction tubes should be completed within 1 minute.

- Using an 8-channel pipette and 100 μ L filter tips, carefully transfer 50 μ L from the top third of each lysed sample in the cluster tube to each reaction tube. Debris may accumulate at the bottom of each lysis tube. This will interfere with the assay performance. Avoid transferring debris by aspirating from the top third of the lysis tube(s). Do not prime the pipette tips and do not mix before aspirating. Place the provided permanent cap(s) on the reaction tube(s).

Note: The lysed sample may be transferred from the same cluster tube a maximum of 3 times.

- Remove the strips(s) of tubes from the reader (or $56 \pm 2^\circ\text{C}$ heat block if one was used) and vortex briefly (about 2 seconds) at medium/high speed (2500 RPM). Quickly visually check each reaction tube to assure that no bubbles are on the bottom or in the middle of the sample. A quick tap of the tubes should release any bubbles from the bottom or middle to the top. Then place into the reader without delay. Close the reader's lid.

Note: The reader will not provide accurate results if the lid is open. Keep the lid closed at all times while the assay is running. Contamination may occur if the permanent caps are not placed on the reaction tubes, or if the permanent caps are removed.

- Click start in the ANSR software to begin the 18 minute assay.
- Results will be displayed as positive, negative or invalid once the assay is finished. If the result displays an invalid, the test must be repeated. A single repeat test can be run starting at step 6. Alternatively, start over from step 1 if lysis from step 6 has been at $80 \pm 2^\circ\text{C}$ for more than 60 minutes.

INTERPRETATION OF RESULTS

The ANSR software will indicate the test results as positive or negative for the presence of *Campylobacter* in the enriched sample. In addition, the real-time fluorescence curve generated from the assay can be viewed.

CONFIRMATION

Samples producing positive ANSR results may be confirmed by streaking the enrichment cultures to Campy-Cefex agar media and continuing with identification of presumptive colonies following standard methods [1].

DISPOSAL

Enrichment cultures and used lysis tubes should be disposed of as biohazard waste. The preferred method of treatment for biohazard waste is autoclaving. Items that cannot be autoclaved, such as used reaction tubes and caps, should be immersed in a fresh 10% bleach solution that is made daily. Consult with the safety advisor of your facility for detailed instructions.

Do not remove the permanent caps, for any reason, from the ANSR reaction tubes once the assay has started, even when disposing of them. Reaction tubes can be disposed of as non-biohazardous waste. It is recommended that they be placed in a sealable container or plastic bag with 10% bleach and immediately disposed of to protect against accidental opening.

CUSTOMER SERVICE

NEOGEN® Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all NEOGEN test kits, is available.

SDS INFORMATION AVAILABLE

Safety data sheets (SDS) are available for this test kit, and all of NEOGEN's test kits, on NEOGEN's website at NEOGEN.com, or by calling NEOGEN at 800.234.5333 or 517.372.9200.

TERMS AND CONDITIONS

For NEOGEN's full terms and conditions, please visit: NEOGEN.com/terms-and-conditions.

WARRANTY

NEOGEN makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, NEOGEN will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. NEOGEN shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

References

[1] USDA-FSIS (2015) Microbiology Laboratory Guidebook, chapter 41.04
<http://www.fsis.usda.gov/wps/wcm/connect/0273bc3d-2363-45b3-befb-1190c25f3c8b/MLG-41.pdf?MOD=AJPERES>

NEOGEN® ANSR® Molecular Diagnostic for Foodborne Pathogen Detection – limited use label license

SYTO® 82

SYTO® 82 contained within this product is provided under an intellectual property license from Life Technologies Corporation, Eugene, OR, and may be used for *in vitro* detection and analysis of (i) food, feeds and beverages, including nutraceuticals, (ii) ingredients for food, feeds and beverages, (iii) process samples from food, feed and beverage preparation, distribution and delivery, and (iv) water from any source for human consumption, all for the purpose of safety and quality assurance. The buyer must not sell or otherwise transfer this product or its components for any other use, including but not limited to: human *in vitro*, veterinary, identity or paternity testing, forensics, or *in vivo* detection of nucleic acid sequences in living beings, or cells. For information on purchasing a license for SYTO 82 for purposes other than food, beverage and water safety, and quality assurance, contact Life Technologies Corporation at outlicensing@lifetech.com.

Molecular beacon probes

One or more molecular beacon probes contained within this product is sold under license from PHRI Properties and may be used under PHRI Properties patent rights only for tests on food products, feeds, beverages and water.

NEAR™ Technology

This product utilizes the patented NEAR isothermal technology and is sold under license from Ionian Technologies, San Diego, CA, and may be used under Ionian Technologies patent rights only for tests on food, beverage and water safety.

PRODUCTS AVAILABLE

- 9465 **Environmental sampling sponges** (with Dey Engley (DE) neutralizing broth) – pack of 100
- 36002 **Sampling sponges, pre-moistened, HiCap** — pack of 100
- 9822 **ANSR® for *E. coli* O157:H7** — up to 96 tests
- 9824 **ANSR for *Listeria monocytogenes*** — up to 96 tests
- 9870 **ANSR for *Salmonella*** — up to 96 tests
- 9871 **ANSR for *Listeria*** — up to 96 tests
- 9872 **ANSR for *Campylobacter*** — up to 96 tests
- 9873 **ANSR for *Listeria* Right Now™** — up to 96 tests
- 9837 **Complete ANSR Equipment Kit** — reader, computer, pipette set, dry bath heaters and inserts, thermometers, vortex, timer, rack

North America

NEOGEN Headquarters
800.234.5333 (U.S./Canada)
foodsafety@neogen.com
NEOGEN.com

Europe, Middle East, and Africa

NEOGEN Europe
+ 44 (0) 1292 525 600
info_uk@neogeneurope.com
NEOGEN.com

Mexico

NEOGEN Latinoamerica
+52 (55) 5254-8235
informacion@neogenlac.com
NEOGEN.com

Brazil

NEOGEN do Brasil
+55 19 3935 3727
info@neogendobrasil.com.br
NEOGEN.com

China

NEOGEN Bio-Scientific Technology
+86 21 6271 7013
info@neogenchina.com.cn
www.neogenchina.com.cn

India

NEOGEN Food and Animal Security
+91 484 2306598, 2301582
info@neogenindia.com
www.neogenindia.com

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