Product Instructions

-  **(EN)** *E. coli* / Coliform Count Plate
-  **(FR)** Test pour la numération des *E.coli* et Coliformes
-  **(DE)** *E. coli* und Coliforme Zählplatte
-  **(IT)** Piastra per il conteggio di *E. coli* e Coliformi
-  **(ES)** Placa para recuento de *E.coli* y Coliformes
-  **(NL)** *E. coli* / Coliform Telplaat
-  **(SV)** Odlingsplatta för *E. coli* och koliforma bakterier
-  **(DA)** *E. coli* / Coliform Tælleplade
-  **(NO)** for *E. coli*- / koliforme bakterier
-  **(FI)** *E. colin* ja koliformien kasvatusalusta
-  **(PT)** Placa para Contagem de *E. coli* / Coliformes
-  **(EL)** Πλακίδιο Καταμέτρησης *E. coli* / Κολοβακτηριδίων
-  **(PL)** Płytko do oznaczania liczby *E. coli* / Coliform
-  **(RU)** Тест-пластина для подсчета *E. coli* и колиформных бактерий
-  **(TR)** *E. coli* / Koliform Sayım Plakası
-  **(JA)** *E.coli*および大腸菌群数測定用プレート（ECプレート）
-  **(ZH)** 大肠杆菌/大肠菌群测试片
-  **(TH)** แผ่นอาหารเลี้ยงเชื้อสำหรับนับจำนวนเชื้ออีโคไล / โคลิฟอร์ม
-  **(KO)** 대장균 / 대장균군 측정용 플레이트

Product Instructions

***E. coli* / Coliform Count Plate**

Product Description and Intended Use

The Neogen® Petrifilm® *E. coli* / Coliform Count (EC) Plate is a sample-ready-culture medium system which contains modified Violet Red Bile (VRB) nutrients, a cold-water-soluble gelling agent, an indicator of glucuronidase activity, 5-bromo-4-chloro-3-indolyl-D-glucuronide (BCIG), and a tetrazolium indicator that facilitates colony enumeration. Neogen Petrifilm EC Plates are used for the enumeration of *Escherichia coli* (*E. coli*) and coliforms in the food and beverage industries. Neogen Petrifilm EC Plate components are decontaminated though not sterilized. Neogen Food Safety is certified to International Standards Organization (ISO) 9001 for design and manufacturing. Neogen Petrifilm EC Plate has not been evaluated with all possible food products, food processes, testing protocols or with all possible microorganism strains.

Safety

The user should read, understand, and follow all safety information in the instructions for Neogen Petrifilm EC Plate. Retain the safety instructions for future reference.

⚠ WARNING: Indicates a hazardous situation, which, if not avoided, could result in death or serious injury and/or property damage.

⚠ WARNING

Do not use this plate for the detection of *E. coli* O157. Because most *E. coli* O157 strains are atypical, for example they are glucuronidase negative, they will not produce a blue color, and will not be detected on Neogen Petrifilm EC Plates.

To reduce the risks associated with exposure to biohazards and environmental contamination:

- Follow current industry standards and local regulations for disposal of biohazardous waste.

To reduce the risks associated with release of contaminated product:

- Follow all product storage instruction contained in the instructions for use.
- Do not use beyond the expiration date.

To reduce the risks associated with bacterial infection and workplace contamination:

- Perform Neogen Petrifilm EC Plate testing in a properly equipped laboratory under the control of a skilled microbiologist.
- The user must train personnel in current proper testing techniques: for example, Good Laboratory Practices¹, ISO 17025² or ISO 7218³.

To reduce the risks associated with misinterpretation of results:

- Neogen has not documented Neogen Petrifilm EC Plates for use in industries other than food and beverage. For example, Neogen has not documented Neogen Petrifilm EC Plates for testing water, pharmaceuticals, or cosmetics.
- Do not use Neogen Petrifilm EC Plates in the diagnosis of conditions in humans or animals.
- The Neogen Petrifilm EC Plates do not differentiate any one microorganism strain from another.

Consult the Safety Data Sheet for additional information.

If you have questions about specific applications or procedures, please visit our website at www.neogen.com or contact your Neogen representative or authorized distributor.

User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at www.neogen.com, or contact your Neogen representative or authorized distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any Neogen Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

Limitation of Warranties / Limited Remedy

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, NEOGEN DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any Neogen Food Safety Product is defective, Neogen or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify Neogen within sixty days of discovery of any suspected defects in a product and return it to Neogen. Please contact your Neogen representative or authorized Neogen distributor for any further questions.

Limitation of Neogen Liability

NEOGEN WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall Neogen's liability under any legal theory exceed the purchase price of the product alleged to be defective.

Storage

Store unopened Neogen Petrifilm EC Plate pouches refrigerated or frozen at temperatures lower than or equal to 8°C (46°F). Just prior to use, allow unopened pouches to come to room temperature before opening. Return unused Neogen Petrifilm EC Plates to pouch. Seal by folding the end of the pouch over and applying adhesive tape. **To prevent exposure to moisture, do not refrigerate opened pouches.** Store resealed pouches in a cool dry place for no longer than four weeks. It is recommended that resealed pouches of Neogen Petrifilm EC Plates be stored in a freezer (see below) if the laboratory temperature exceeds 25°C (77°F) and/or the laboratory is located in a region where the relative humidity exceeds 50% (with the exception of air-conditioned premises).

To store opened pouches in a freezer, place Neogen Petrifilm EC Plates in a sealable container. To remove frozen Neogen Petrifilm EC Plates for use, open the container, remove the plates that are needed and immediately return remaining plates to the freezer in the sealed container for the remainder of the shelf life. Neogen Petrifilm EC Plates should not be used past their expiration date. The freezer that is used for open pouch storage must not have an automatic defrost cycle as this would repeatedly expose the plates to moisture which can damage the plates.

Do not use Neogen Petrifilm EC Plates that show discoloration. Expiration date and lot number are noted on each package of Neogen Petrifilm EC Plates. The lot number is also noted on individual plates Neogen Petrifilm EC Plates.

⚠ Disposal

After use, Neogen Petrifilm EC Plates may contain microorganisms that may be a potential biohazard. Follow current industry standards for disposal.

Instructions for Use

Sample Preparation

1. Use appropriate sterile diluents:

Butterfield's phosphate buffered dilution water⁴, 0.1% peptone water, peptone salt diluent, quarter-strength Ringer's solution, saline solution (0.85-0.90%), Neogen® Wide-Spectrum Neutralizer, bisulfite-free letheen broth or distilled water.

Do not use diluents containing citrate, bisulfite or thiosulfate with Neogen Petrifilm EC Plates; they can inhibit growth. If citrate buffer is indicated in the standard procedure, substitute with one of the buffers listed above, warmed to 40-45°C.

2. Blend or homogenize sample.
3. For optimal growth and recovery of microorganisms, adjust the pH of the sample suspension to 6.6 - 7.2. For acidic products, adjust the pH with 1N NaOH. For alkaline products, adjust the pH with 1N HCl.

Plating

1. Place the Neogen Petrifilm EC Plate on a flat, level surface.
2. Lift the top film and with the pipette perpendicular to the inoculation area dispense 1 mL of sample suspension onto the center of bottom film.
3. Roll the top film down onto the sample to prevent trapping air bubbles.
4. Place the Neogen® Petrifilm® Spreader with the flat side down on the center of the plate. Press gently on the center of the Neogen Petrifilm Spreader to distribute the sample evenly. Spread the inoculum over the entire Neogen Petrifilm EC Plate growth area before the gel is formed. Do not slide the Neogen Petrifilm Spreader across the film.
5. Remove the Neogen Petrifilm Spreader and leave the Neogen Petrifilm EC Plate undisturbed for at least one minute to permit the gel to form.

Incubation

Incubate Neogen Petrifilm EC Plates in a horizontal position with the clear side up in stacks of no more than 20 plates. Several incubation times and temperatures can be used depending on current local reference methods, some of which are listed in the "**Specific Instructions for Validated Methods**" section.

Interpretation

1. Neogen Petrifilm EC Plates can be counted using a standard colony counter or other illuminated magnifier. Do not count colonies on the foam dam since they are removed from the selective influence of the medium. Do not count artifact bubbles that may be present.

The interpretation of *E. coli* colonies on the Neogen Petrifilm EC Plate is as follows:

AOAC Official Methods (998.08 and 991.14) – enumerate blue to red-blue colonies associated with entrapped gas, regardless of size or intensity of color, as confirmed *E. coli*. Blue colonies without gas are not counted as *E. coli*.

Other coliform colonies are red and closely associated (within one colony diameter) with entrapped gas. Colonies not associated with gas (a distance greater than one colony diameter between colony and gas bubble) are not counted as coliforms. The total coliform count consists of both the red and blue colonies associated with gas at 24 hours. Anytime within the validated method incubation period that a blue colony associated with gas appears, it is a confirmed *E. coli*.

▲ WARNING

Do not use this plate for the detection of *E. coli* O157. Because most *E. coli* O157 strains are atypical, for example they are glucuronidase negative, they will not produce a blue color, and will not be interpreted as *E. coli* on Neogen Petrifilm EC Plates.

2. The circular growth area is approximately 20 cm². Estimates can be made on plates containing greater than 150 colonies by counting the number of colonies in one or more representative squares and determining the average number per square. Multiply the average number by 20 to determine the estimated count per plate.
3. When present in large numbers, Neogen Petrifilm EC Plates may have one or more of the following characteristics: a deepening of the gel color with many small, indistinct colonies; and many gas bubbles. High concentrations of *E. coli* will cause the growth area to turn blue, while high concentrations of coliforms (non-*E. coli*) will cause the growth area to turn dark red. When this occurs, record results as too numerous to count (TNTC). When an actual count is required, plate at a higher dilution.
4. Where necessary, colonies may be isolated for further identification. Lift the top film using proper testing technique and pick the colony from the gel. Test using standard procedures.
5. If the plates cannot be counted within 1 hour of removal from the incubator, they may be stored for later enumeration by freezing in a sealable container at temperatures lower than or equal to negative 15°C for no longer than one week.

For further information refer to the “Neogen® Petrifilm® *E. coli*/Coliform Count Plate Interpretation Guide”. If you have questions about specific applications or procedures, please visit our website at www.neogen.com or contact your Neogen representative or authorized distributor.

Specific Instructions for Validated Methods

AOAC® Official MethodsSM (998.08 Confirmed *Escherichia coli* Counts in Poultry, Meats and Seafood, Dry Rehydratable Film Method)

Incubate Neogen Petrifilm EC Plates 24 hours ± 2 hours at 35°C ± 1°C.

AOAC® Official MethodsSM (991.14 Coliform and *Escherichia coli* Counts in Foods, Dry Rehydratable Film Methods)

For coliform results incubate Neogen Petrifilm EC Plates 24 hours ± 2 hours at 35°C ± 1°C.

For *E. coli* results incubate Neogen Petrifilm EC Plates an additional 24 hours ± 2 hours (48 hours ± 4 hours total) at 35°C ± 1°C.

References

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
3. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
4. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
5. ISO 6887. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination.

Explanation of Symbols

info.neogen.com/symbols

AOAC is a registered trademark of AOAC INTERNATIONAL
Official Methods is a service mark of AOAC INTERNATIONAL

Neogen Food Safety

Neogen Corporation

620 Leshar Place
Lansing, MI 48912 USA
Neogen.com

Neogen Europe Ltd.

The Dairy School
Auchincruive
Ayr, KA6 5HU
Scotland, UK

Neogen Ireland, Ltd.

Bray Business Park, Bray
Co. Wicklow
A98YV29, Ireland

Neogen Corporation

620 Leshar Place Lansing, MI 48912 USA
www.neogen.com

© Neogen Corporation 2024. All rights reserved.
Neogen and Petrifilm are registered trademarks of Neogen Corporation.
FS00850A

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REFERANSLAR

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
3. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
4. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.

5. ISO 6887. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination.

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