Product Instructions

- ⓔⓃ Neogen Molecular Detection Assay 2Q - Quantitative *Salmonella* Screening Kit
- ⓕⓇ Kit Neogen de détection moléculaire 2Q – Kit de dépistage de quantification des *Salmonella*
- ⓓⓔ Kit Neogen Molekularer Detektionstest 2Q – Salmonellen, quantitativ
- ⓔⓈ Kit Neogen de Ensayo de Detección Molecular 2Q - Quantificación de *Salmonella*
- ⓔⓟ Kit Neogen Ensaio de Detecção Molecular Neogen 2Q - Quantitativa de *Salmonella*
- ⓙⓐ Neogen 病原菌検出アッセイ 2Q - サルモネラ属菌定量検査スクリーニング キット
- Ⓩⓗ Neogen分子检测分析2Q - 定量 沙门氏菌筛选试剂盒
- ⓓⓗ ชุดทดสอบเชื้อก่อโรคระดับโมเลกุล 2 *ซาลโมเนลลา* แบบนับจำนวน
- ⓕⓐ Neogen 분자 검출 분석 2Q - 정량 *살모넬라* 스크리닝 키트

A large graphic of the letters "2Q" in a bold, green, sans-serif font. The "2" is slightly behind and to the left of the "Q".

Product Instructions

Neogen Molecular Detection Assay 2Q - Quantitative *Salmonella* screening kit

Product Description and Intended Use

The Neogen® Molecular Detection Assay 2Q - Quantitative *Salmonella* (MDA2QSAL96) kit is used with Neogen® Quantitative Rapid Enrichment Dehydrated Media (QRED500) and the Neogen® Molecular Detection System for rapid quantification of *Salmonella* in enriched chicken carcass rinses and raw ground poultry samples. The Neogen® Molecular Detection Assay 2Q – Quantitative *Salmonella* Kit is used in conjunction with a model application that is accessed through the Neogen® Molecular Detection System Software. The resulting output obtained in the Neogen® Molecular Detection System Software is used to assess the quantity of *Salmonella* in chicken carcass rinses and raw ground poultry matrices.

The assay uses loop-mediated isothermal amplification to rapidly amplify nucleic acid sequences with high specificity and sensitivity, combined with bioluminescence to detect the amplification. Presumptive positive results are reported in real-time. Results that are below the limit of detection will be displayed as negative after the assay is completed. Presumptive positive results should be confirmed using your preferred method or as specified by local regulations^(1, 2, 3).

The Neogen Molecular Detection Assay 2Q – Quantitative *Salmonella* is intended for use in a laboratory environment by professionals trained in laboratory techniques. Neogen has not documented the use of this product in industries other than food or beverage. For example, Neogen has not documented this product for testing pharmaceutical, cosmetics, clinical, or veterinary samples. The Neogen Molecular Detection Assay 2Q – Quantitative *Salmonella* has not been evaluated with all possible food products, processes, testing protocols, or with all possible strains of bacteria. As with all test methods, the source, formulation, and quality of enrichment medium can influence the results. Factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may also influence results. Neogen recommends evaluation of the method including enrichment medium, in the user's environment using a sufficient number of samples with particular foods and microbial challenges to ensure that the method meets the user's criteria.

Neogen has evaluated the Neogen Molecular Detection Assay 2Q – Quantitative *Salmonella* with Neogen® Buffered Peptone Water (BPW), Neogen® neutralizing Buffered Peptone Water (nBPW) and Neogen Quantitative Rapid Enrichment Dehydrated Media (QRED500).

The Neogen Molecular Detection Instrument is intended for use with samples that have undergone heat treatment during the assay lysis step, which is designed to destroy organisms present in the sample.

Samples that have not been properly heat treated during the assay lysis step may be considered a potential biohazard and should NOT be inserted into the Neogen Molecular Detection Instrument.

Neogen Food Safety is certified to ISO (International Organization for Standardization) 9001 for design and manufacturing.

The Neogen Molecular Detection Assay 2Q – Quantitative *Salmonella* test kit contains 96 tests, described in Table 1.

Table 1. Kit Components

Item	Identification	Quantity	Contents	Comments
Lysis Solution (LS) tubes	Pink solution in clear tubes	96 (12 strips of 8 tubes)	580 µL of LS per tube	Racked and ready to use
Quantitative <i>Salmonella</i> Reagent tubes	Green tubes	96 (12 strips of 8 tubes)	Lyophilized specific amplification and detection mix	Ready to use
Extra caps	Green caps	96 (12 strips of 8 caps)		Ready to use
Reagent Control (RC)	Clear flip-top tubes	16 (2 pouches of 8 individual tubes)	Lyophilized control DNA, amplification, and detection mix	Ready to use

The Negative Control, not provided in the kit, is a sterile enrichment medium, e.g., BPW, nBPW or QRED500. Do not use water as a Negative Control.

Safety

The user should read, understand, and follow all safety information in the instructions for the Neogen Molecular Detection System and the Neogen Molecular Detection Assay 2Q – Quantitative *Salmonella*. Retain the safety instructions for future reference.

WARNING:	Indicates a hazardous situation, which, if not avoided, could result in death or serious injury and/or property damage.
CAUTION:	Indicates a hazardous situation, which, if not avoided, could result in minor or moderate injury and/or property damage.
NOTICE:	Indicates a potentially hazardous situation which, if not avoided, could result in property damage.

WARNING

The Neogen Molecular Detection Assay 2Q – Quantitative *Salmonella* protocol DOES NOT have interruption points. Samples must be fully analyzed after the sample has been enriched (See protocols in Table 2).

Do not use the Neogen Molecular Detection Assay 2Q – Quantitative *Salmonella* in the diagnosis of conditions in humans or animals.

The user must train its personnel in current proper testing techniques: for example, Good Laboratory Practices, ISO/IEC 17025⁽⁴⁾, or ISO 7218⁽⁵⁾.

To reduce the risks associated with a Below the Limit of Quantitation (LOQ) result leading to the release of contaminated product:

- Do not interpret results below the LOQ of the Neogen Molecular Detection Assay 2Q – Quantitative *Salmonella* assay as a negative *Salmonella* result. Results below the LOQ results do not indicate the complete absence of *Salmonella*, only that the level of *Salmonella*, if present, is below the LOQ for this method.
- Follow the protocol and perform the tests exactly as stated in the product instructions.
- Always use a calibrated micropipette.
- Use the Neogen Molecular Detection Assay 2Q – Quantitative *Salmonella* with foods that have been validated internally or by a third party.

- Store the Neogen Molecular Detection Assay 2Q – Quantitative *Salmonella* as indicated on the package and in the product instructions.
- Always use the Neogen Molecular Detection Assay 2Q – Quantitative *Salmonella* by the expiration date.
- Do not use the sample enrichments outlined in these product instructions in the Neogen® Molecular Detection Assay 2 – *Salmonella* (MDA2SAL96) QUALITATIVE assay. The enrichments included in these product instructions have been developed specifically for this QUANTITATIVE assay.

To reduce the risks associated with exposure to chemicals and biohazards:

- Perform pathogen testing in a properly equipped laboratory under the control of trained personnel. Incubated enrichment media and equipment or surfaces that have come into contact with incubated enrichment media may contain pathogens at levels sufficient to cause risk to human health.
- This procedure uses/detects pathogenic microorganisms and/or their metabolic products above a certain level. Care should be taken to avoid ingestion or inhalation of potentially infectious aerosols or contact with the skin. Always follow standard laboratory safety practices, including wearing appropriate protective apparel and eye protection while handling reagents and contaminated samples.
- Avoid contact with the contents of the enrichment media and reagent tubes after amplification.
- Dispose of enriched samples according to current industry standards.
- Samples that have not been properly heat treated during the assay lysis step may be considered a potential biohazard and should NOT be inserted into the Neogen Molecular Detection Instrument.
- To reduce the risks associated with cross-contamination while preparing the assay:
- Always wear gloves (to protect the user and prevent the introduction of nucleases).

To reduce the risks associated with environmental contamination:

Follow current industry standards for disposal of contaminated waste.

CAUTION

To reduce the risks associated with exposure to hot liquids:

- Do not exceed the recommended temperature setting on heater.
- Do not exceed the recommended heating time.

Use an appropriate, calibrated thermometer to verify the Neogen® Molecular Detection Heat Block Insert temperature (e.g., a partial immersion thermometer or digital thermocouple thermometer, not a total immersion thermometer). The thermometer must be placed in the designated location in the Neogen Molecular Detection Heat Block Insert.

NOTICE

To reduce the risks associated with cross-contamination while preparing the assay:

- Change gloves prior to reagent pellet hydration.
- Use of sterile, aerosol barrier (filtered), molecular biology grade pipette tips is recommended.
- Use a new pipette tip for each sample transfer.
- Use Good Laboratory Practices during centrifugation processing steps and to transfer the sample from the enrichment to the lysis tube. To avoid pipettor contamination, the user may choose to add an intermediate transfer step. For example, the user can transfer each enriched sample into a sterile tube.
- Use a molecular biology workstation containing germicidal lamp where available. Periodically decontaminate laboratory benches and equipment (pipettes, cap/decap tools, etc.) with a 1–5% (v:v in water) household bleach solution or DNA removal solution.

To reduce the risks associated with a false-positive result:

- Never open tubes post amplification.
- Always dispose of the contaminated tubes by soaking in a 1–5% (v:v in water) household bleach solution for 1 hour and away from the assay preparation area.
- Never autoclave reagent tubes post amplification.
- Always use a calibrated micropipette.

Consult the Safety Data Sheet for additional information and local regulations for disposal.

If you have questions about specific applications or procedures, please visit our website at neogen.com or contact your local Neogen representative or authorized distributor.

User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at neogen.com, or contact your local Neogen representative or authorized distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, laboratory technique, and the sample itself may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any Neogen Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

To help customers evaluate the method for various food matrices, Neogen has developed the Neogen® Molecular Detection Matrix Control kit. When needed, use the Matrix Control (MC) to determine if the matrix has the ability to impact the Neogen Molecular Detection Assay 2Q – Quantitative *Salmonella* results. Test several samples, representative of the matrix, i.e. samples obtained from different origin, during any validation period when adopting the Neogen method or when testing new or unknown matrices or matrices that have undergone raw material or process changes.

A matrix can be defined as a type of product with intrinsic properties such as composition and process. Differences between matrices may be as simple as the effects caused by differences in their processing or presentation, for example, raw vs. pasteurized; fresh vs. dried, etc.

Limitation of Warranties / Limited Remedy

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, NEOGEN DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any Neogen Food Safety Product is defective, Neogen or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify Neogen within sixty days of discovery of any suspected defects in a product and return it to Neogen. Please contact your Neogen Food Safety representative or authorized Neogen distributor for any further questions.

Limitation of Neogen Liability

NEOGEN WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall Neogen's liability under any legal theory exceed the purchase price of the product alleged to be defective.

Storage and Disposal

Store the Neogen Molecular Detection Assay 2Q – Quantitative *Salmonella* Kit components at 2–8°C. Do not freeze. Keep the kit away from light during storage. After opening the kit, check that the foil pouch is undamaged. If the pouch is damaged, do not use it. After opening, unused reagent tubes should always be stored in the re-sealable pouch with the desiccant inside to maintain stability of the lyophilized reagents. Store resealed pouches at 2–8°C for no longer than 60 days. Do not use Neogen Molecular Detection Assay 2Q – Quantitative *Salmonella* past the expiration date. Expiration date and lot number are noted on the outside label of the box. After use, the enrichment medium and the Neogen Molecular Detection Assay 2Q – Quantitative *Salmonella* tubes can potentially contain pathogenic materials. When testing is complete, follow local regulations for the disposal of contaminated waste. Consult the Safety Data Sheet for additional information.

Instructions for Use

Follow all instructions carefully. Failure to do so may lead to inaccurate results.

Periodically decontaminate laboratory benches and equipment (pipettes, cap/decap tools, etc.) with a 1–5% (v:v in water) household bleach solution or DNA removal solution.

The user should complete the Neogen Molecular Detection System operator qualification (OQ) training, as described in the “Installation Qualification (IQ) / Operational Qualification (OQ) Protocols and Instructions for Neogen Molecular Detection System” document ⁽⁶⁾.

See Section “Specific Instructions for validated methods” for specific requirements:

Quantification of *Salmonella* in raw ground chicken/turkey

1. Weigh 325 g of raw ground chicken or turkey and place in a sterile bag
2. Add 400 mL of buffered peptone water (BPW).
3. Homogenize thoroughly by stomaching at 180 rpm for 1 minute.
4. Transfer 30 mL of the raw ground chicken or turkey homogenate into a sample bag.
5. Add 30 mL of Neogen Quantitative Rapid Enrichment Dehydrated Media (QRED500) (refer to Quantitative Rapid Enrichment Dehydrated Media Technical Specification Sheet)⁽⁷⁾, that has been pre-warmed to 42°C. This is referred to as the ‘Sample’.
6. Incubate the ‘Sample’ following enrichment conditions in table 2.
7. After enrichment, remove the Sample from the incubator and mix well.
8. Immediately after mixing, aliquot 1.5 mL of the Sample into a 2 mL microcentrifuge tube with conical bottom.
9. Centrifuge in a mini centrifuge at 7,000 x g for 5 minutes to pellet particulates.
10. Carefully discard the supernatant without disturbing the pellet. Gently, invert the tube and tap on paper towel to remove the any remaining excess of supernatant.
11. Resuspend the pellet in 300 µL of QRED500.
12. Vortex until the pellet is completely resuspended in the QRED500.
13. Analyze 50 µL of the resuspended pellet sample (from 12) with the Neogen Molecular Detection Assay 2Q - Quantitative *Salmonella*.

Quantification of *Salmonella* in chicken carcass rinses

1. Place one chicken carcass into a sterile poultry rinse bag.
2. Add 400 mL of Neogen neutralizing buffered peptone water (nBPW) or buffered peptone water (BPW) through the chicken canal (See Table 2).
3. Mix the content for 1 minute to rinse the carcass.
4. Transfer 30 mL of the rinse into a sample bag.
5. Add 30 mL of Neogen Quantitative Rapid Enrichment Dehydrated Media (QRED500) (refer to Quantitative Rapid Enrichment Dehydrated Media Technical Specification Sheet)(6), that has been pre-warmed to 42°C. This is referred to as the ‘Sample’.
6. Incubate the Sample following enrichment conditions in table 2.
7. After enrichment, remove the Sample from the incubator and mix well.
8. Immediately after mixing, aliquot 1.5 mL of the Sample into a 2 mL microcentrifuge tube with conical bottom.
9. Centrifuge in a mini centrifuge at 5,000 x g for 5 minutes to pellet particulates.
10. Carefully discard the supernatant without disturbing the pellet. Gently, invert the tube and tap on paper towel to remove any remaining excess of supernatant.
11. Resuspend the pellet in 100 µL (microliters) of QRED500.
12. Vortex until the pellet is completely resuspended in the QRED500.
13. Analyze 50 µL of the resuspended pellet sample (from 12) with the Neogen Molecular Detection Assay 2Q - Quantitative *Salmonella*.

Table 2 present guidance for general enrichment protocols for chicken carcass rinses and raw ground poultry samples. It is the user’s responsibility to validate alternate sampling protocols or dilution ratios to ensure this test method meets the user’s criteria.

Table 2. General enrichment protocols.

Sample Matrix	Sample Size	Homogenization (first step)	Enrichment	Enrichment Temp	Enrichment Time	Pellet resuspension	Sample Analysis Volume ^(a)
Raw ground chicken	325 g	400 mL BPW	30 mL homogenate + 30 mL QRED	42 ± 1°C	6 hours	300 µL	50 µL
Raw ground turkey	325 g	400 mL BPW	30 mL homogenate + 30 mL QRED	42 ± 1°C	5 hours	300 µL	50 µL
Chicken carcass rinse	30 mL	400 mL nBPW	30 mL nBPW rinse + 30 mL QRED	42 ± 1°C	6 hours	100 µL	50 µL
Chicken carcass rinse	30 mL	400 mL BPW	30 mL BPW rinse + 30 mL QRED	42 ± 1°C	5 hours	100 µL	50 µL

(a) Volume of sample transferred to Lysis Solution tubes. Refer to step 4.5 of Lysis section.

Preparation of the Neogen® Molecular Detection Speed Loader Tray

1. Wet a cloth or disposable towel with a 1–5% (v:v in water) household bleach solution and wipe the Neogen Molecular Detection Speed Loader Tray.
2. Rinse the Neogen Molecular Detection Speed Loader Tray with water.
3. Use a disposable towel to wipe the Neogen Molecular Detection Speed Loader Tray dry.
4. Ensure the Neogen Molecular Detection Speed Loader Tray is dry before use.

Preparation of the Neogen® Molecular Detection Chill Block Insert

Place the Neogen Molecular Detection Chill Block Insert directly on the laboratory bench. Use the block at ambient laboratory temperature (20–25°C).

Preparation of the Neogen® Molecular Detection Heat Block Insert

Place the Neogen Molecular Detection Heat Block Insert in a dry double block heater unit. Turn on the dry block heater unit and set the temperature to allow the Neogen Molecular Detection Heat Block Insert to reach and maintain a temperature of 100 ±1°C.

NOTE: Depending on the heater unit, allow approximately 30 minutes for the Neogen Molecular Detection Heat Block Insert to reach temperature. Using an appropriate, calibrated thermometer (e.g., a partial immersion thermometer, digital thermocouple thermometer, not a total immersion thermometer) placed in the designated location, verify that the Neogen Molecular Detection Heat Block Insert is at 100 ±1°C.



Preparation of the Neogen Molecular Detection Instrument

1. Launch the Neogen Molecular Detection Software and log in. Contact your Neogen Food Safety representative to ensure you have the most updated version of the software that includes Quantitative Assays.
2. Turn on the Neogen Molecular Detection Instrument.
3. Create or edit a run with data for each sample. Refer to the Neogen Molecular Detection System User Manual and the MDA2QSAL96 Software Technical Bulletin for details.

NOTE: The Neogen Molecular Detection Instrument must reach and maintain a temperature of 60°C before inserting the Neogen Molecular Detection Speed Loader Tray with reaction tubes. This heating step takes approximately 20 minutes and is indicated by an ORANGE light on the instrument's status bar. When the instrument is ready to start a run, the status bar will turn GREEN.

Lysis






1. Remove the bottom of Neogen Lysis Solution Rack with a screwdriver or spatula before placing it in the Neogen Molecular Detection Heat Block Insert.



2. Allow the Lysis Solution (LS) tubes to warm up by setting the rack at room temperature (20–25 °C) overnight (16–18 hours). Alternatives to equilibrate the LS tubes to room temperature are to set the LS tubes on the laboratory bench for at least 2 hours, incubate the LS tubes in a 37 ±1°C incubator for 1 hour or place them in a dry double block heater for 30 seconds at 100°C.
3. Invert the capped tubes to mix. Proceed to the next step within 4 hours.
4. One LS tube is required for each sample and the Negative Control (NC) sample (sterile enrichment medium).
 - 4.1 LS tube strips can be cut to desired LS tube number. Select the number of individual LS tubes or 8-tube strips needed. Place the LS tubes in an empty rack.
 - 4.2 Use the Neogen® Molecular Detection Cap/Decap Tool-Lysis to decap one LS tube strip.
 - 4.3 To avoid cross-contamination, decap one LS tubes strip at a time and use a new pipette tip for each transfer step.
 - 4.4 Discard the LS tube cap.
 - 4.5 Transfer 50 µL of the resuspended pellet sample into a LS tube.

Transfer **50 µL** of each resuspended pellet sample into an individual LS tube **first**. Transfer the NC **last**.

5. Repeat steps 2 to 5 until each individual sample has been added to a corresponding LS tube in the strip.

6. Repeat steps 1 to 5 as needed, for the number of samples to be tested.
7. When all samples have been transferred, transfer 50 µL of NC (sterile QRED500) into a LS tube. Do not use water as a NC.
8. Verify that the temperature of the Neogen Molecular Detection Heat Block Insert is at 100 ±1°C.
9. Place the uncovered rack of LS tubes in the Neogen Molecular Detection Heat Block Insert and heat for 15 ±1 minutes. During heating, the LS solution will change from pink (cool) to yellow (hot). Samples that have not been properly heat treated during the assay lysis step may be considered a potential biohazard and should NOT be inserted into the Neogen Molecular Detection Instrument.
10. Remove the uncovered rack of LS tubes from the heating block and allow to cool in the Neogen Molecular Detection Chill Block Insert at least 5 minutes and a maximum of 10 minutes. The Neogen Molecular Chill Block Insert, used at ambient temperature without the Molecular Detection Chill Block Tray, should sit directly on the laboratory bench. When cool, the lysis solution will revert to a pink color.
11. Remove the rack of LS tubes from the Neogen Molecular Detection Chill Block Insert.


Amplification

1. One Neogen Molecular Detection Assay 2Q – Quantitative *Salmonella* is required for each sample and the NC.
 - 1.1 Reagent tube strips can be cut to desired tube number. Select the number of individual Reagent tubes or 8-tube strips needed.
 - 1.2 Place Reagent tubes in an empty rack.
 - 1.3 Avoid disturbing the reagent pellets from the bottom of the tubes.
2. Select one Reagent Control (RC) tube and place it in the rack.
3. To avoid cross-contamination, decap one Reagent Tube strip at a time and use a new pipette tip for each transfer step.
4. Transfer each of the lysates as described below:

Transfer 20 μ L of each sample lysate into individual Reagent tubes first followed by the NC. Hydrate the RC tube last.

5. Use the Neogen Molecular Detection Cap/Decap Tool-Reagent to decap the Neogen Molecular Detection Assay 2Q – Quantitative *Salmonella* Reagent tubes, one Reagent tube strip at a time. Discard cap.
 - 5.1 Transfer 20 μ L of sample lysate from the upper $\frac{1}{2}$ of the liquid (avoid precipitate) in the Neogen Lysis Solution Tube into corresponding Reagent Tube. Dispense at an angle to avoid disturbing the pellets. Mix by gently pipetting up and down 5 times.
 - 5.2 Repeat step 5.1 until all individual Sample lysate has been added to a corresponding Neogen Molecular Detection Assay 2Q – Quantitative *Salmonella* Reagent tube in the strip.
 - 5.3 Cover the Neogen Molecular Detection Assay 2Q – Quantitative *Salmonella* Reagent tubes with the provided extra caps and use the rounded side of the Neogen Molecular Detection Cap/Decap Tool-Reagent to apply pressure in a back and forth motion ensuring that the cap is tightly applied.
 - 5.4 Repeat steps 5.1 to 5.3 as needed, for the number of samples to be tested.
 - 5.5 When all sample lysates have been transferred, vortex the capped reagents tubes for 2 seconds for mixing, and then spin-down for 10 seconds in a mini-centrifuge to make sure all the liquid is at the bottom of the tube.
 - 5.6 Transfer 20 μ L of NC lysate into a Neogen Molecular Detection Assay 2Q – Quantitative *Salmonella* Reagent tube and another 20 μ L of NC lysate into a RC tube. Dispense at an angle to avoid disturbing the pellets. Mix by gently pipetting up and down 5 times.
6. Load capped tubes into a clean and decontaminated Neogen Molecular Detection Speed Loader Tray. Close and latch the Neogen Molecular Detection Speed Loader Tray lid.

7. Review and confirm the configured run in the Neogen Molecular Detection Software Quantitative Assays.
8. Click the Start button in the software and select instrument for use. The selected instrument's lid automatically opens.
9. Place the Neogen Molecular Detection Speed Loader Tray into the Neogen Molecular Detection Instrument and close the lid to start the assay. Results are provided within 40 minutes, although positives may be detected sooner.

10. After the assay is complete, remove the Neogen Molecular Detection Speed Loader Tray from the Neogen Molecular Detection Instrument and dispose of the tubes by soaking in a 1-5% (v:v in water) household bleach solution for 1 hour and away from the assay preparation area.

NOTICE: To minimize the risk of false positives due to cross-contamination, never open reagent tubes containing amplified DNA. This includes Neogen Molecular Detection Assay 2Q – Quantitative *Salmonella* Reagent, Neogen Reagent Control, and Neogen Matrix Control Tubes. Always dispose of sealed reagent tubes by soaking in a 1-5% (v:v in water) household bleach solution for 1 hour and away from the assay preparation area.

Results and Interpretation

The results are determined by the analysis of unique curve parameters by the software.

References:

1. US Food and Drug Administration Bacteriological Analytical Manual. Chapter 5: *Salmonella*, September 2023 Version.
2. US Department of Agriculture (USDA) FSIS Microbiology Laboratory Guidebook 4.14. Isolation and Identification of *Salmonella* from Meat, Poultry, Pasteurized Egg, Carcass, and Environmental Sponges. Effective Date: 05 June 2023.
3. ISO 6579. Microbiology of food and animal feeding stuffs – Horizontal method for the detection of *Salmonella* spp.
4. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
5. ISO 7218. Microbiology of food and animal feeding stuff – General rules for microbiological examination.
6. Neogen Installation Qualification (IQ) / Operational Qualification (OQ) Protocols and Instructions for Neogen Molecular Detection System.
7. Neogen Quantitative Rapid Enrichment Dehydrated Media (QRED500) Technical Specification Sheet.

Explanation of Symbols

info.neogen.com/symbols

Neogen Food Safety

Neogen Corporation
620 Lesher Place
Lansing, MI 48912 USA
[neogen.com](https://www.neogen.com)

Neogen Europe Ltd.

The Dairy School
Auchincruive
Ayr, KA6 5HU
Scotland, UK

Neogen Ireland, Ltd.

Bray Business Park, Bray
Co. Wicklow
A98YV29, Ireland

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Referências:

1. US Food and Drug Administration Bacteriological Analytical Manual. Chapter 5: *Salmonella*, versão de setembro de 2023.
2. US Department of Agriculture (USDA) FSIS Microbiology Laboratory Guidebook 4.14. Isolation and Identification of *Salmonella* from Meat, Poultry, Pasteurized Egg, Carcass, and Environmental Sponges. Data de vigência: 05 de junho de 2023.
3. ISO 6579. Microbiology of food and animal feeding stuffs – Horizontal method for the detection of *Salmonella* spp.
4. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
5. ISO 7218. Microbiology of food and animal feeding stuff – General rules for microbiological examination.
6. Neogen Installation Qualification (IQ) / Operational Qualification (OQ) Protocols and Instructions for Neogen Molecular Detection System.
7. Folha de Especificações Técnicas dos Meios Desidratados de Enriquecimento Rápido Quantitativo (QRED500) da Neogen.

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1. US Food and Drug Administration Bacteriological Analytical Manual. Chapter 5: *Salmonella*, September 2023 Version.
2. US Department of Agriculture (USDA) FSIS Microbiology Laboratory Guidebook 4.14. Isolation and Identification of *Salmonella* from Meat, Poultry, Pasteurized Egg, Carcass, and Environmental Sponges. Effective Date: 05 June 2023.
3. ISO 6579. Microbiology of food and animal feeding stuffs – Horizontal method for the detection of *Salmonella* spp.
4. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
5. ISO 7218. Microbiology of food and animal feeding stuff – General rules for microbiological examination.
6. Neogen 設置適格性 (IQ)/実施適格性 (OQ) プロトコルと Neogen 病原菌自動検出システムの説明書。
7. Neogen 定量迅速増菌粉末培地 (QRED500) 技術仕様書。

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