Product Instructions

-  (EN) High-Sensitivity Coliform Count Plate
-  (FR) Test Haute-Sensibilité pour la numération des Coliformes
-  (DE) Geringe Mengen Coliforme Zählplatte
-  (IT) Piastra ad alta sensibilità per il conteggio dei coliformi
-  (ES) Placa de alta sensibilidad para el recuento de bacterias coliformes
-  (NL) High-Sensitivity Coliform Telplaat
-  (SV) Högkänslig odlingsplatta för koliforma bakterier
-  (DA) High-Sensitivity Coliform Tælleplade
-  (NO) Høysensitiv plate for koliforme bakterier
-  (FI) High-Sensitivity Koliformien kasvatusalusta
-  (PT) Placa Alta Sensibilidade para Contagem de Coliformes
-  (EL) Πλακίδιο Καταμέτρησης Κολοβακτηριδίων Υψηλής Ευαισθησίας
-  (PL) Płytko o wysokiej czułości do zliczania bakterii z grupy coli
-  (RU) Высокочувствительная тест-пластина для подсчета колиформных бактерий
-  (TR) Çok Hassas Koliform Sayım Plakası
-  (JA) 高感度大腸菌群数測定用プレート (HSCCプレート)
-  (ZH) 高灵敏度大肠菌群测试片
-  (TH) แผ่นอาหารเลี้ยงเชื้อสำหรับนับจำนวนโคลิฟอร์มความไวสูง
-  (KO) 고감도 대장균군 플레이트

HSCC

High-Sensitivity Coliform Count

Product Instructions

High-Sensitivity Coliform Count Plate

Product Description and Intended Use

The Neogen® Petrifilm® High-Sensitivity Coliform Count (HSCC) Plate is a sample-ready-culture medium system which contains modified Violet Red Bile (VRB) nutrients, a cold-water-soluble gelling agent, and a tetrazolium indicator that facilitates colony enumeration. Neogen Petrifilm HSCC Plates are used for the enumeration of low levels of coliforms in the food and beverage industries. A 5 mL sample is plated onto one plate. This allows a 1-coliform-per-gram sensitivity if a 1:5 dilution of product made. A 2-coliform-per-gram sensitivity is achieved if a 1:10 dilution is used. Neogen Petrifilm HSCC Plate components are decontaminated though not sterilized. Neogen Food Safety is certified to International Organization for Standardization (ISO) 9001 for design and manufacturing. Neogen Petrifilm HSCC Plates have not been evaluated with all possible food products, food processes, testing protocols or with all possible microorganism strains.

The United States Food and Drug Administration (FDA) Bacteriological Analytical Manual (BAM)¹ define coliforms as Gram-negative rods, which produce acid and gas from lactose fermentation. Coliform colonies growing on the Neogen Petrifilm HSCC Plate produce acid, which causes the pH indicator to deepen the gel color; gas trapped around red colonies indicates coliforms.

ISO defines coliforms by their ability to grow in method-specific, selective media. ISO method 4832⁴, enumerating coliforms by colony count technique, defines coliforms as acid producers on VRB with lactose (VRBL) agar. On Neogen Petrifilm HSCC Plates these acid-producing coliforms are indicated by red colonies with or without gas production. ISO method 4831³, enumerating coliforms by the most probable number (MPN) method, defines coliforms by their ability to grow and produce gas from lactose in a selective broth. On Neogen Petrifilm HSCC Plates these coliforms are indicated by red colonies associated with gas.

AFNOR has validated Neogen Petrifilm HSCC Plate in comparison to ISO method 4831³ for enumeration of total coliforms.

Safety

The user should read, understand, and follow all safety information in the instructions for the Neogen Petrifilm HSCC Plate. Retain the safety instructions for future reference.

⚠ **WARNING:** Indicates a hazardous situation, which, if not avoided, could result in death or serious injury and/or property damage.

⚠ WARNING

To reduce the risks associated with exposure to biohazards and environmental contamination:

- Follow current industry standards and local regulations for disposal of biohazardous waste.

To reduce the risks associated with release of contaminated product:

- Follow all product storage instructions contained in the instructions for use.
- Do not use beyond the expiration date.

To reduce the risks associated with bacterial infection and workplace contamination:

- Perform Neogen Petrifilm HSCC Plate testing in a properly equipped laboratory under the control of a skilled microbiologist.
- The user must train personnel in current proper testing techniques: for example, Good Laboratory Practices⁶, ISO 17025⁷ or ISO 7218⁸.

To reduce the risks associated with misinterpretation of results:

- Neogen has not documented Neogen Petrifilm HSCC Plates for use in industries other than food and beverage. For example, Neogen has not documented Neogen Petrifilm HSCC Plates for testing water, pharmaceuticals or cosmetics.
- Do not use Neogen Petrifilm HSCC Plates in the diagnosis of conditions in humans or animals.
- Neogen Petrifilm HSCC Plates do not differentiate any one coliform strain from another.
- When testing certain unprocessed foods Neogen Petrifilm HSCC Plates may produce a softer gel, which may or may not interfere with the ability to read the plate. Some of the foods that have been found to cause gel softening include some raw shellfish and some raw flour, including mussels, oysters, scallops, shrimp, buckwheat flour, graham flour, and organic whole-wheat flour.
- High sugar content foods may increase growth potential for gas producing non-coliform *Enterobacteriaceae*.
- A high level of background organisms in certain food (for example, some cheeses and some chilled unprocessed foods including raw fish and raw milk) may elevate the number of red colonies **without** gas.

Consult the Safety Data Sheet for additional information.

If you have questions about specific applications or procedures, please visit our website at www.neogen.com or contact your Neogen representative or authorized distributor.

User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at www.neogen.com, or contact your Neogen representative or authorized distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results.

It is the user's responsibility in selecting any test method to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' or suppliers' requirements.

As with any test method, results obtained from use of any Neogen Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

Limitation of Warranties / Limited Remedy

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, NEOGEN DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING, BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any Neogen Food Safety Product is defective, Neogen or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify Neogen within sixty days of discovery of any suspected defects in a product and return it to Neogen. Please contact your Neogen representative or authorized Neogen distributor for any further questions.

Limitation of Neogen Liability

NEOGEN WILL NOT BE LIABLE TO USER OR OTHERS FOR ANY LOSS OR DAMAGE, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING, BUT NOT LIMITED TO, LOST PROFITS. In no event shall Neogen's liability under any legal theory exceed the purchase price of the product alleged to be defective.

Storage

Store unopened Neogen Petrifilm HSCC Plate pouches refrigerated or frozen at temperatures lower than or equal to 8°C (46°F). Just prior to use, allow unopened pouches to equilibrate at room temperature for a minimum of 10 minutes before opening to prevent condensation inside the pouch. Return unused plates to pouch. Seal by folding the end of the pouch over and applying adhesive tape. **To prevent exposure to moisture, do not refrigerate opened pouches.** Store resealed pouches in a cool dry place for no longer than four weeks. It is recommended that resealed pouches of Neogen Petrifilm HSCC Plates be stored in a freezer if the laboratory temperature exceeds 25°C (77°F) and/or the laboratory is located in a region where the relative humidity exceeds 50% (with the exception of air-conditioned premises).

To store opened pouches in a freezer, place Neogen Petrifilm HSCC Plates in a sealable container. To remove frozen Neogen Petrifilm HSCC Plates for use, open the container, remove the plates that are needed and immediately return remaining plates to the freezer in the sealed container for the remainder of the shelf life. Neogen Petrifilm HSCC Plates should not be used past their expiration date. The freezer that is used for open pouch storage must not have an automatic defrost cycle as this would repeatedly expose the Neogen Petrifilm HSCC Plates to moisture which can damage the plates.

Do not use Neogen Petrifilm HSCC Plates that show discoloration. Expiration date and lot number are noted on each package of Neogen Petrifilm HSCC Plates. The lot number is also noted on individual Neogen Petrifilm HSCC Plates.

⚠ Disposal

After use, Neogen Petrifilm HSCC Plates may contain microorganisms that may be a potential biohazard. Follow current industry standards for plate disposal.

Instructions for Use

Sample Preparation

1. Use appropriate sterile diluents:

Butterfield's phosphate buffered dilution water¹, 0.1% peptone water^{1,2}, peptone salt diluent², quarter-strength Ringer's solution, dipotassium hydrogen phosphate, saline solution (0.85-0.90%)², Neogen® Wide-Spectrum Neutralizer, bisulfite-free letheen broth, or distilled water.

Do not use diluents containing citrate, bisulfite or thiosulfate with Neogen Petrifilm HSCC Plates; they can inhibit growth. If citrate buffer is indicated in the standard procedure, substitute with one of the buffers listed above, warmed to 40-45°C.

2. Blend or homogenize sample.

- For optimal growth and recovery of microorganisms, adjust the pH of the sample suspension to 6.5 - 7.5. For acidic products, adjust the pH with 1N NaOH. For alkaline products, adjust the pH with 1N HCl.

Plating

- Place the Neogen Petrifilm HSCC Plate on a flat, level surface.
- Lift the top film and with the pipette perpendicular to the inoculation area dispense 5 mL of sample suspension onto the center of bottom film.
- Roll the top film down onto the sample to prevent trapping air bubbles.
- Place the Neogen® Petrifilm® High Sensitivity Spreader (Catalog #6481) on the center of the plate. Press gently on the center of the Neogen Petrifilm High Sensitivity Spreader to distribute the sample evenly. Spread the inoculum over the entire Neogen Petrifilm HSCC Plate growth area before the gel is formed. Do not slide the Neogen Petrifilm High Sensitivity Spreader across the film.
- Remove the spreader and leave the plate undisturbed for at least two to five minutes to permit the gel to form.

Incubation

Incubate plates in a horizontal position with the clear side up in stacks of no more than 10 plates. Incubate Neogen Petrifilm HSCC Plates 24 hours ± 2 hours. Several incubation times and temperatures can be used depending on current local reference methods, some of which are listed in the section below titled “**Specific Instructions for Validated Methods**”.

Interpretation

- Neogen Petrifilm HSCC Plates can be counted using a standard colony counter or other illuminated magnifier. Do not count artifact bubbles that may be present.

The interpretation of coliform colonies on the Neogen Petrifilm HSCC Plate varies by method depending on current local reference methods, some of which are listed in the section below titled “**Specific Instructions for Validated Methods**”.

For example:

Enumeration of colonies associated with gas:

Coliform colonies are red and closely associated (within one colony diameter) with entrapped gas.

Enumeration of colonies with or without gas for thermotolerant enumeration:

After incubating the Neogen Petrifilm HSCC Plates for 24 hours ± 2 hours at 44°C ± 1°C, thermotolerant coliforms are defined as red colonies with or without gas.

- The circular growth area is approximately 60 cm². Estimates can be made on plates containing greater than 150 colonies by counting the number of colonies in two or more representative squares and determining the average number per square. Multiply the average number by 60 to determine the estimated count per plate.
- When colonies are present in large numbers, Neogen Petrifilm HSCC Plates will have a deepening of the gel color and either or both of the following characteristics: many small, indistinct colonies and/or many gas bubbles. When this occurs, record results as too numerous to count (TNTC). When an actual count is required, plate at a higher dilution.
- Where necessary, colonies may be isolated for further identification. Lift the top film using proper testing technique and pick the colony from the gel. Test using standard procedures.
- If the plates cannot be counted within 1 hour of removal from the incubator, they may be stored for later enumeration by freezing in a sealable container at temperatures lower than or equal to negative 15°C for no longer than one week.

For further information refer to the appropriate “Neogen® Petrifilm® High Sensitivity Coliform Count Plate Interpretation Guide.” If you have questions about specific applications or procedures, please contact your official Neogen representative nearest you or visit our website at www.neogen.com.

Specific Instructions for Validated Methods

AOAC® Official MethodsSM (996.02 Coliform Count in Dairy Products, High-Sensitivity Dry Rehydratable Film Method)

Incubate Neogen Petrifilm HSCC Plates 24 hours ± 2 hours at 32°C ± 1°C.

NF Validation by AFNOR Certification

NF Validation certified method in compliance with ISO 16140-2³ in comparison to ISO 4831³ and ISO 4832⁴ (3M-01/07-03/99)

Use the following details when implementing the above Instructions for Use except for Sample Preparation, refer to the information below:

Scope of the validation: human food products.

Sample preparation:

Use only ISO listed diluents².

Do not use diluents containing citrate, bisulfite or thiosulfate with Neogen Petrifilm HSCC Plates; they can inhibit growth. If citrate buffer is indicated in the standard procedure, substitute with one of the buffers listed above, warmed to 40-45°C.

Incubation:

Incubate Neogen Petrifilm HSCC Plates 24 hours ± 2 hours at 30°C ± 1°C or 37°C ± 1°C.

Interpretation:

Coliforms are defined as red colonies with gas on Neogen Petrifilm HSCC Plates. Calculate the number of microorganisms present in the test sample according to ISO 7218⁸ for one plate per dilution. Estimates are outside of the scope of the NF VALIDATION Certification. Refer to EN ISO 7218 standard for inoculation, colony counting and calculation and expression of results.

3M-01/07-03/99

Alternative Analytical methods for Agribusiness

<http://nf-validation.afnor.org/en>

For more information about end of validity, please refer to NF VALIDATION certificate available on the website mentioned above.

References

1. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
2. ISO 6887: Microbiology of the food chain – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination.
3. ISO 4831: Microbiology of food and animal feeding stuffs – Horizontal method for the detection and enumeration of coliforms – Most probable number technique.
4. ISO 4832: Microbiology of food and animal feeding stuffs – Horizontal method for the enumeration of coliforms – colony-count technique.
5. ISO 16140-2: Microbiology of the food chain – Method Validation – Protocol for the validation of alternative (proprietary) methods against a reference method
6. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
7. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
8. ISO 7218: Microbiology of food and animal feeding stuffs – General rules for microbiological examination.

Explanation of Symbols

info.neogen.com/symbols

AOAC is a registered trademark of AOAC INTERNATIONAL
Official Methods is a service mark of AOAC INTERNATIONAL

Neogen Food Safety

Neogen Corporation

620 Leshar Place
Lansing, MI 48912 USA
Neogen.com

Neogen Europe Ltd.

The Dairy School
Auchincruive
Ayr, KA6 5HU
Scotland, UK

Neogen Ireland, Ltd.

Bray Business Park, Bray
Co. Wicklow
A98YV29, Ireland

Neogen Corporation

620 Leshar Place Lansing, MI 48912 USA
www.neogen.com

© Neogen Corporation 2024. All rights reserved.
Neogen and Petrifilm are registered trademarks of Neogen Corporation.
FS00986A

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