700006963 | NCM1019A

700006964 | NCM1019B

700006965 | NCM1019C

700006966 | NCM1019D

Product Instructions

ⓔ One Plate *Enterobacteriaceae* (OP-EBAC)

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ⓓ One Plate *Enterobacteriaceae* (OP-EBAC)

ⓔ One Plate *Enterobacterias* (OP-EBAC)

ⓓ One Plate *Enterobacteriaceae* (OP-EBAC)

ⓙ One Plate 腸内細菌科 (OP-EBAC)

ⓑ One Plate 腸杆菌科 (OP-EBAC)

Ⓚ One Plate 장내세균 (OP-EBAC)

Product Instructions

One Plate *Enterobacteriaceae* (OP-EBAC)

Intended Use

One Plate *Enterobacteriaceae* (OP-EBAC) offers a rapid method for the enumeration of total *Enterobacteriaceae* using traditional culture methodology.

Product Summary and Explanation

Enterobacteriaceae are a family of gram-negative bacteria typically defined by their ability to ferment glucose to produce acid and gas. There are over 200 recognised species, including all coliform bacteria and *E. coli*. Since the group contains some of the major foodborne pathogens such as *Salmonella* spp. or *E. coli* O157, *Enterobacteriaceae* are considered one of the most significant indicators of microbial safety. Many *Enterobacteriaceae* are non-pathogenic but can pose as opportunistic pathogens in the immunocompromised. The family are present in the intestinal tracts of animals and thus are ubiquitous in soil, water, and sewage, meaning they are a key indicator of hygiene and sanitation in food production.

One Plate for *Enterobacteriaceae* consists of a selective base which permits the detection and enumeration of *Enterobacteriaceae* by their ability to ferment glucose and produce acid visualised by an indicator dye. Non-target organisms are inhibited by a mixture of selective compounds incorporated into the base.

The product can permit a quantitative result for *Enterobacteriaceae* in 18–22 hours for all food products, using only one plate versus 2 or more as described in ISO 7218 and ISO 21528-2:2017. Neither overlays nor confirmations are required for the method.

Intended User

The method is designed for use by qualified personnel with appropriate training.

Product Codes

NOTE: Please refer to ISO 6887 parts 1–5 for suitable sample specific diluents.

Product Name	Format	Pack Size	SKU	
One Plate <i>Enterobacteriaceae</i>	DCM	500 g	700006963	NCM1019A
		5 kg	700006964	NCM1019B
		10 kg	700006965	NCM1019C
		25 kg	700006966	NCM1019D
One Plate <i>Enterobacteriaceae</i>	Ready to Re-melt Agar*	1 x 250 mL bottle	On Request	

*Ready to Re-melt Agar should be melted in a steaming water bath (refer to ISO 7218 for further guidance).

Typical Composition

Formulation may be adjusted and/or supplemented as required to meet performance specifications.


Growth Mix	21.3 g/L
Selective Mix	5.95 g/L
Indicator Mix	0.25 g/L
Agar	10.0 g/L

pH of the prepared media at 25 °C: pH 7.4 ± 0.2.

Preparation

Suspend 37.5 g of the medium in 1 litre of purified water. Heat with frequent agitation and boil for one minute to completely dissolve the medium. DO NOT AUTOCLAVE. Cool to 45–50 °C.

OP-EBAC Flow Diagram

Sample Preparation

Follow the specifications of ISO 6887 or the specific international standard appropriate to the product concerned for the initial suspension and dilutions.

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
1. Transfer 1 mL of the sample if liquid, or 1 mL of the suspension in the case of other products, or any serial dilutions, into **ONE** empty, sterile Petri dish.
2. Pour approximately 15–20 mL of molten OP-EBAC into the plate. Homogenize well by swirling and let solidify on a cool surface.
3. **Overlays are not mandatory** but can help prevent spreading growth on the surface of the agar. In this case, add 4–6 mL of molten agar on top of the set agar and allow to solidify before incubation as described in ISO 7218.

NOTE: Decimal dilutions could be required.

Incubation

Invert the prepared plate and incubate at 37 ± 1 °C for 18–22 hours.

Interpretation

Enterobacteriaceae will present as orange to yellow colonies against a purple background, with or without a yellow halo.

Counting Colony-forming Units

Typical *Enterobacteriaceae* will present as orange to yellow colonies, with or without a yellow halo. If less than one quarter of the dish is overgrown by spreading, count the colonies on the unaffected part and extrapolate for the theoretical total count of the dish.

Select the plate containing 10–150 colonies for accurate enumeration. Refer to ISO 7218 when the results are outside the limits. Apply the dilution factors, which eliminates the need to have two successive dilutions or duplicates.

Confirmation is not required for the expression of the *Enterobacteriaceae* count.

Dilute to Specifications Results

Dilutions are made according to the specification, and a calculated amount of the appropriate dilution is added to the plate.

Examples of specification dilutions:

Dilution	Inoculation	Specifications
1/10	1 mL pour plate	Negative result <10 CFU/mL Positive result ≥ 10 CFU/mL
1/100	1 mL pour plate	Negative result <100 CFU/mL Positive result ≥ 100 CFU/mL
1/1000	1 mL pour plate	Negative result <1000 CFU/mL Positive result ≥ 1000 CFU/mL

Quality Control

Appearance of dehydrated media: Beige, homogeneous, free-flowing powder

Appearance of prepared media: Purple, clear gel

Typical cultural response when incubated aerobically at 37 ± 1°C and examined for growth at 18–22 hours:

Microorganism	WDCM	Approx. Inoculum (CFU)	Expected Results
<i>Escherichia coli</i>	00013 & 00012	50–200	≥50% Recovery, Orange to Yellow*
<i>Klebsiella aerogenes</i>	00175	50–200	≥50% Recovery, Orange to Yellow*
<i>Citrobacter freundii</i>	00006	50–200	≥50% Recovery, Orange to Yellow*
<i>Salmonella Enteritidis</i>	00030	50–200	≥50% Recovery, Orange to Yellow*
<i>Salmonella Typhimurium</i>	00031	50–200	≥50% Recovery, Orange to Yellow*
<i>Enterococcus faecalis</i>	00087 & 00009	>10 ⁴	Complete Inhibition
<i>Pseudomonas aeruginosa</i>	00025	>10 ⁴	Growth

*With or without yellow halo

Precautions and Limitations of the Method

1. Use good microbiology laboratory practices as per ISO 7218 except mandatory duplicate plating for enumeration, as the One Plate method requires only one Petri dish.
2. The minimum time for incubation is 18 hours; if at this point any presumptive colonies present a weak or dubious reaction, the plate can be re-incubated up to a further 4 hours to allow for a stronger development of colour.
3. Overlays are not mandatory (not performed during validation) but can be helpful to prevent spreading growth on the surface of the agar. In this case add 4-6 mL of molten agar on top of the set agar and allow to solidify before incubation as described in ISO 7218.
4. Colony counts in excess of 150 CFU per plate can lead to poorer definition of the diagnostic reactions; it is advised that a higher dilution be used to accurately enumerate higher levels of contamination.
5. Both the limit of detection and quantification are dictated by the subculture volume.
6. An inoculation of 1.0 mL (pour plate) of a 1/10 dilution has a limit of detection of 10 CFU/g and permits the user to quantify (enumerate) > 100 CFU/g.
7. If a plate count results in less than 10 colonies counted, an estimated number should be expressed, e.g., report >10 CFU/g or >100 CFU/g (depending on the inoculation volume).

Verification

Method verification should be done following the protocols described in ISO 16140:3:2021. The user laboratory should follow experimental design and acceptance criteria as described in chapter 6, 'Quantitative methods — Technical protocol for verification'.

Validations

MicroVal (ISO 16140-2:2016)

One Plate *Enterobacteriaceae* has been certified by MicroVal as an alternative to the reference ISO 21528-2:2017 Horizontal method for the detection and enumeration of *Enterobacteriaceae* — Part 2: Colony count technique, according to the reference standard ISO 16140-2:2016, with the scope of enumeration of *Enterobacteriaceae* in a broad range of foods. Refer to the MicroVal certificate for more information.

Safety

Refer to the relevant product safety data sheet (SDS).

Storage

Store dehydrated culture media at 2–30 °C away from direct sunlight. Once it is opened and recapped, place the container in a low-humidity environment at the same storage temperature. Protect from moisture and light by keeping container tightly closed.

Prepared media and/or supplements should be stored at 2–8 °C. Do not freeze. Supplement may be stored at room temperature for up to 5 days.

Expiration

Refer to expiration date stamped on the container. The dehydrated medium should be discarded if not free-flowing, or if appearance has changed from the original colour. Expiry applies to medium in its intact container when stored as directed.

Disposal

Cultures should be disposed of appropriately as biohazard waste. The preferred method of treatment for biohazard waste is autoclaving. Items that cannot be autoclaved may be disinfected with bleach solution. Consult with the safety advisor for

your facility for detailed instructions.

Customer Service

Neogen Customer Services and Technical Services can be reached by using the following contact information:

The Dairy School, Auchincruive, Ayr, KA6 5HU, Scotland UK

+44 (0) 1292 525600

infoUK@neogen.com

Training on this product, and all Neogen test kits, is available.

Terms and Conditions

For full terms and conditions, please visit <https://www.neogen.com/terms-and-conditions>

Warranty

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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