Product Instructions

-  (EN) Staph Express System
-  (FR) Système Staph Express
-  (DE) Staph Express System
-  (IT) Sistema Staph Express
-  (ES) Sistema Staph Express
-  (NL) Staph Express Systeem
-  (SV) Staph Express System
-  (DA) Staph Express system
-  (NO) Staph Express-system
-  (FI) Staph Express-järjestelmä
-  (PT) Sistema Staph Express
-  (EL) Σύστημα Staph Express
-  (PL) System do oznaczania liczby *Staphylococcus*
-  (RU) Система для экспресс-подсчета стафилококков
-  (TR) Staph Express Sistemi
-  (JA) 黄色ブドウ球菌エクスプレス測定用システム
-  (ZH) 金黄色葡萄球菌确认系统
-  (TH) Staph Express System
-  (KO) 황색포도상구균용시스템

Product Instructions

Staph Express System

Product Description and Intended Use

The Neogen® Petrifilm® Staph Express (STX) System consists of a Neogen® Petrifilm® Staph Express Count (STX) Plate and a Neogen® Petrifilm® Staph Express (STX) Disk, which are packaged separately. The Neogen Petrifilm STX Plate is a sample-ready-culture medium system which contains a cold-water-soluble gelling agent. The chromogenic, modified Baird-Parker medium in the plate is selective and differential for *Staphylococcus aureus* but may also indicate *Staphylococcus hyicus* (*S. hyicus*) or *Staphylococcus intermedius* (*S. intermedius*). The Neogen Petrifilm STX Disk contains toluidine blue-O that facilitates the visualization of deoxyribonuclease (DNase) reactions. DNase-positive organisms detected on the Neogen Petrifilm STX Plate include *Staphylococcus aureus* (*S. aureus*). The Neogen Petrifilm STX Plates and Neogen Petrifilm STX Disks are used for the enumeration of DNase positive *Staphylococcus* species in the food and beverage industries. Neogen Petrifilm STX Plate and Neogen Petrifilm STX Disk components are decontaminated though not sterilized.

Neogen Food Safety is certified to International Organization for Standardization (ISO) 9001 for design and manufacturing. Neogen Petrifilm STX System has not been evaluated with all possible food products, food processes, testing protocols or with all possible microorganism strains.

Safety

The user should read, understand, and follow all safety information in the instructions for the Neogen Petrifilm STX Plate and Neogen Petrifilm STX Disk. Retain the safety instructions for future reference.

⚠ **WARNING:** Indicates a hazardous situation, which, if not avoided, could result in death or serious injury and/or property damage.

⚠ WARNING

To reduce the risks associated with exposure to biohazards and environmental contamination:

- Follow current industry standards and local regulations for disposal of biohazardous waste.

To reduce the risks associated with release of contaminated product:

- Follow all product storage instruction contained in the instructions for use.
- Do not use beyond the expiration date.

To reduce the risks associated with bacterial infection and workplace contamination:

- Perform Neogen Petrifilm STX Plate and Neogen Petrifilm STX Disk testing in a properly equipped laboratory under the control of a skilled microbiologist.
- The user must train personnel in current proper testing techniques: for example, Good Laboratory Practices¹, ISO 17025³ or ISO 7218².

To reduce the risks associated with misinterpretation of results:

- Neogen has not documented Neogen Petrifilm STX Plates and Neogen Petrifilm STX Disk for use in industries other than food and beverage. For example, Neogen has not documented Neogen Petrifilm STX Plates and Neogen Petrifilm STX Disk for testing water, pharmaceuticals, or cosmetics.
- Do not use Neogen Petrifilm STX Plates and Neogen Petrifilm STX Disk in the diagnosis of conditions in humans or animals.
- Do not use Neogen Petrifilm STX Plates and Neogen Petrifilm STX Disk for U.S.-recognized laboratory pasteurized counts.
- Acceptance of the Neogen Petrifilm STX Plate and Neogen Petrifilm STX Disk method for the testing of water per an accepted local government regulation is at the sole discretion and responsibility of the end-user.
- Neogen Petrifilm STX Plate and Neogen Petrifilm STX Disk do not differentiate any one microorganism strain from another.

Consult the Safety Data Sheet for additional information.

If you have questions about specific applications or procedures, please visit our website at www.neogen.com or contact your Neogen representative or authorized distributor.

User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at www.neogen.com, or contact your Neogen representative or authorized distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any Neogen Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

Limitation of Warranties / Limited Remedy

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, NEOGEN DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any Neogen Food Safety Product is defective, Neogen or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify Neogen within sixty days of discovery of any suspected defects in a product and return it to Neogen. Please contact your Neogen representative or authorized Neogen distributor for any further questions.

Limitation of Neogen Liability

NEOGEN WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall Neogen's liability under any legal theory exceed the purchase price of the product alleged to be defective.

Storage

Plate Storage

Store unopened Neogen Petrifilm STX Plate pouches refrigerated or frozen at temperatures lower than or equal to 8°C (46°F). Just prior to use, allow unopened pouches to equilibrate at room temperature for a minimum of 10 minutes before opening to prevent condensation inside the pouch. Return unused Neogen Petrifilm STX Plates to pouch. Seal by folding the end of the pouch over and applying adhesive tape. **To prevent exposure to moisture, do not refrigerate opened pouches.** Store resealed pouches in a cool dry place for no longer than four weeks. It is recommended that resealed pouches of Neogen Petrifilm STX Plates be stored in a freezer (see below) if the laboratory temperature exceeds 25°C (77°F) and/or the laboratory is located in a region where the relative humidity exceeds 50% (with the exception of air-conditioned premises).

Disk Storage

Neogen Petrifilm STX Disks are individually packaged within a foil pouch. They are sensitive to both moisture and light. Upon receipt, store unopened pouches of Neogen Petrifilm STX Disks refrigerated or frozen at temperatures lower than or equal to 8°C (46°F). Remove only those individually packaged Neogen Petrifilm STX Disks that will be used immediately; store the remaining Neogen Petrifilm STX Disks in the foil pouch by folding the end of the pouch and applying adhesive tape. **To prevent exposure to moisture, do not refrigerate opened pouches.** Store resealed pouches in a cool, dry place for no longer than six months. It is recommended that resealed pouches of Neogen Petrifilm STX Disks be stored in a freezer (see below) if the laboratory temperature exceeds 25°C (77°F) and/or the laboratory is located in a region where the relative humidity exceeds 50% (with the exception of air-conditioned premises).

Freezer Storage

To store opened pouches of Neogen Petrifilm STX Plates and Neogen Petrifilm STX Disks in a freezer, place Neogen Petrifilm STX Plates and Neogen Petrifilm STX Disks in a sealable container. To remove frozen Neogen Petrifilm STX Plates and Neogen Petrifilm STX Disks for use, open the container, remove the Neogen Petrifilm STX Plates and Neogen Petrifilm STX Disks that are needed and immediately return remaining Neogen Petrifilm STX Plates and Neogen Petrifilm STX Disks to the freezer in the sealed container for the remainder of the shelf life. Neogen Petrifilm STX Plates and Neogen Petrifilm STX Disks should not be used past their expiration date. Do not store open pouches in a freezer with an automatic defrost cycle, as this could damage the Neogen Petrifilm STX Plates due to repeated exposure to moisture.

Do not use Neogen Petrifilm STX Plates and Neogen Petrifilm STX Disks that show discoloration. Expiration date and lot number are noted on each package of Neogen Petrifilm STX Plates and Neogen Petrifilm STX Disks. The lot number is also noted on individual Neogen Petrifilm STX Plates and on the individual Neogen Petrifilm STX Disk packages.

△ Disposal

After use, Neogen Petrifilm STX Plates may contain microorganisms that may be a potential biohazard. Follow current local, regional, national and industry standards for disposal.

Instructions for Use

Follow all instructions carefully. Failure to do so may lead to inaccurate results.

Sample Preparation

1. Use appropriate sterile diluents:

Butterfield's phosphate buffered dilution water⁴, peptone salt diluent⁵, 0.1% peptone water, buffered peptone water⁵, quarter-strength Ringer's solution, saline solution (0.85-0.90%), bisulfite-free letheen broth, or distilled water. See section "**Specific Instruction for Validated Methods**" for specific requirements.

Do not use diluents containing citrate, bisulfite, or thiosulfate with Neogen Petrifilm STX Plates; they can inhibit growth. If citrate buffer is indicated in the standard procedure, substitute warmed to 40-45°C (104-113°F) Butterfield's phosphate buffered dilution water or peptone salt diluent.

Do not use dipotassium hydrogen phosphate as the DNase reaction may be inhibited.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, diluent manufacture, sample preparation, handling, and laboratory technique may influence results. It is important to verify the performance of the Neogen Petrifilm STX Disk with the diluent chosen for sample preparation. Failure to do so, may result in false negatives.

2. Blend or homogenize sample.
3. For optimal growth and recovery of microorganisms, adjust the pH of the sample suspension to 6 - 8. For acidic products, adjust the pH with 1N NaOH. For alkaline products, adjust the pH with 1N HCl.

Plating

1. Place the Neogen Petrifilm STX Plate on a flat, level surface.
2. Lift the top film and with the pipette perpendicular dispense 1 mL of sample suspension onto the center of bottom film.
3. Roll the top film down onto the sample to prevent trapping air bubbles.
4. Place the Neogen Petrifilm Flat Spreader (Catalog #6425) with the flat side down on the center of the plate. Press gently on the center of the spreader to distribute the sample evenly. Spread the inoculum over the entire Neogen Petrifilm STX Plate growth area before the gel is formed. Do not slide the spreader across the film.
5. Remove the spreader and leave the Neogen Petrifilm STX Plate undisturbed for at least one minute to permit the gel to form.

Incubation

Incubate Neogen Petrifilm STX Plates in a horizontal position with the clear side up in stacks of no more than 20 Neogen Petrifilm STX Plates.

Incubate Neogen Petrifilm STX Plates for 24 ± 2 hours at 35°C ± 1°C or 37°C ± 1°C (temperature based on validated references). See section "**Specific Instruction for Validated Methods**" for specific requirements.

Interpretation

1. Count Neogen Petrifilm STX Plates with a standard colony counter or other illuminated magnifier. Do not count colonies on the foam dam since they are removed from the selective influence of the medium.
2. Observe colony colors.
 - a. If no colonies or only red-violet colonies are present after 24 ± 2 hours, count red-violet colonies as *S. aureus*, *S. hyicus*, or *S. intermedius*; the test is complete. Use of the Neogen Petrifilm STX Disk is not necessary.
 - b. If a confirmation test is required or if any colony colors besides red-violet – for example, black or blue-green – are present, use a Neogen Petrifilm STX Disk (see steps 3-11). Black colonies may be stressed microorganisms.

NOTE: High levels of phosphatase, an enzyme that occurs naturally in some foods, such as raw milk products, may cause the entire Neogen Petrifilm STX Plate to turn pink and can obscure the red-violet colonies. If this occurs, dilute your sample further.

3. Remove an individually-packaged Neogen Petrifilm STX Disk from its pouch and allow it to come to room temperature. Then remove the Neogen Petrifilm STX Disk from its individual package by peeling the package to expose the Neogen Petrifilm STX Disk's tab, grasping the tab, and pulling the Neogen Petrifilm STX Disk out.

- Lift the top film of the Neogen Petrifilm STX Plate and place the Neogen Petrifilm STX Disk in the well of the Neogen Petrifilm STX Plate so that the tab remains outside the well.

NOTE: Occasionally the gel may split when the top film is lifted. Performance of the Neogen Petrifilm STX Plate is not affected by gel splitting because the Neogen Petrifilm STX Disk is coated on both sides.

- Lower the top film.
- Apply pressure by sliding a finger firmly across the entire Neogen Petrifilm STX Disk area (including edges) to ensure uniform contact of the Neogen Petrifilm STX Disk with the gel and to eliminate any air bubbles.

NOTE: If too much pressure is applied, the gel may be disturbed, and zones may not be distinct. NOTE: If too little pressure is applied, the gel may dry, and zones may not form.

- Place the disked Neogen Petrifilm STX Plates in stacks of no more than 20 Neogen Petrifilm STX Plates and incubate at 35°C ± 1°C or at 37°C ± 1°C for at least 60 minutes and no longer than 3 hours. Check the Neogen Petrifilm STX Plates periodically—the DNase reaction may give final results in less than 3 hours. If all of the colonies have pink zones, the testing is complete and the colonies can be counted. If some colonies do not have zones then the Neogen Petrifilm STX Plate must be incubated for the full 3 hours.

NOTE: If the Neogen Petrifilm STX Plates cannot be counted within 1 hour of removal from the incubator they may be stored for later enumeration by freezing in a sealable container at temperatures less than or equal to minus 15°C (5°F) for no longer than one week. Allow Neogen Petrifilm STX Plates to reach room temperature, and then use the Neogen Petrifilm STX Disk as described above.

NOTE: Never freeze disked Neogen Petrifilm STX Plates for later interpretation.

- Count all pink zones whether or not colonies are visible. Pink zones are usually associated with *S. aureus* but may indicate *S. hyicus* or *S. intermedius*. Colonies not associated with a pink zone are not DNase producing Staphylococci, and should not be counted. If the entire disked area is pink with no distinct zones, large numbers of DNase producing colonies are present. Record the result as too numerous to count (TNTC) and dilute the sample further to obtain a more accurate count.
- When necessary, colonies may be isolated for further identification. Lift the top film and pick the colony from the gel. If a Neogen Petrifilm STX Disk is covering the gel, peel the Neogen Petrifilm STX Disk away using a forceps and then pick the colony from the gel. Dispose of the Neogen Petrifilm STX Disk by following current industry standards and sterilize the forceps.

For further information refer to the appropriate Neogen Petrifilm Staph Express Count Plate and Disk “Interpretation Guide.” If you have questions about specific applications or procedures, please visit our website at www.neogen.com or contact your Neogen representative or authorized distributor.

Specific Instructions for Validated Methods

AOAC® Official MethodsSM (2003.07 Neogen Petrifilm STX Plate Method for the Enumeration of *Staphylococcus aureus* in Selected Types of Processed and Prepared Foods)

AOAC® Official MethodsSM (2003.08 Neogen Petrifilm STX Plate Method for the Enumeration of *Staphylococcus aureus* in Selected Dairy Foods)

AOAC® Official MethodsSM (2003.11 Neogen Petrifilm STX Plate Method for the Enumeration of *Staphylococcus aureus* in Selected Meat, Seafood, and Poultry)

Incubate 24 hours ± 2 hours at 35°C ± 1°C or 37°C ± 1°C

- If no colonies or only red-violet colonies appear, test is complete, no need to use disk. Count red- violet colonies as *S. aureus*.
- If colony colors besides red-violet appear, insert disk and re-incubate 1 to 3 hours at 35°C ± 1°C or 37°C ± 1°C. Count pink zones as *S. aureus*.

NF Validation by AFNOR Certification

NF Validation certified method in compliance with ISO 16140-2⁶ in comparison to ISO 6888-1⁷ (3M-01/09-04/03A) and ISO 6888-2⁸ (3M-01/09-04/03B)

Scope of the validations: All human foods and pet food

Use the following details when implementing the above Instructions for Use:

Sample preparation

Use only ISO listed diluents.

Incubation

Incubate Neogen Petrifilm STX Plates 24 hours ± 2 hours at 37°C ± 1°C. Incubation of Neogen Petrifilm STX Disks at 37°C ± 1°C for 3 hours.

Interpretation

Calculate the number of microorganisms present in the test sample according to ISO 72182 for one Neogen Petrifilm STX Plate per dilution. Read the Neogen Petrifilm STX Plates after 3 hours of incubation time is complete. Counting range is:

- Less than or equal to 150 red-violet colonies and/or less than or equal to 300 total colonies.
- Less than or equal to 150 pink zones.

Estimates are outside of the scope of the NF VALIDATION Certification. Refer to EN ISO 7218 standard for inoculation, colony counting and calculation and expression of results.

3M 01/09-04/03A and B
ALTERNATIVE ANALYTICAL METHODS FOR AGRIBUSINESS
<http://nf-validation.afnor.org/en>

For more information about end of validity, please refer to NF VALIDATION certificate available on the web-site mentioned above.

References

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Practice Studies.
2. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
3. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
4. U.S. Food and Drug Administration. 1998. Bacteriological Analytical Manual, 8th ed., (Revision A), Appendix 3.64.
5. ISO 6887. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination,
6. ISO 16140-2. Microbiology of the food chain – Method Validation – Protocol for the validation of alternative (proprietary) methods against a reference method
7. ISO 6888-1. Microbiology of food and animal feeding stuffs –Horizontal method for the enumeration of coagulase-positive staphylococci (*Staphylococcus aureus* and other species), Part 1: Technique using Baird-Parker agar medium.
8. ISO 6888-2. Microbiology of food and animal feeding stuffs –Horizontal method for the enumeration of coagulase-positive staphylococci (*Staphylococcus aureus* and other species), Part 2: Technique using rabbit plasma fibrinogen agar medium.

Refer to the current versions of the standard methods listed above.

Explanation of Symbols

info.neogen.com/symbols

AOAC is a registered trademark of AOAC INTERNATIONAL
 Official Methods is a service mark of AOAC INTERNATIONAL

Neogen Food Safety

Neogen Corporation

620 Leshar Place
Lansing, MI 48912 USA
Neogen.com

Neogen Europe Ltd.

The Dairy School
Auchincruive
Ayr, KA6 5HU
Scotland, UK

Neogen Ireland, Ltd.

Bray Business Park, Bray
Co. Wicklow
A98YV29, Ireland

Neogen Corporation

620 Leshar Place Lansing, MI 48912 USA
www.neogen.com

© Neogen Corporation 2024. All rights reserved.
Neogen and Petrifilm are registered trademarks of Neogen Corporation.
FS00985A

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Referências

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Practice Studies.
2. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
3. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
4. U.S. Food and Drug Administration. 1998. Bacteriological Analytical Manual, 8th ed., (Revision A), Appendix 3.64.
5. ISO 6887. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination,

6. ISO 16140-2. Microbiology of the food chain – Method Validation – Protocol for the validation of alternative (proprietary) methods against a reference method
7. ISO 6888-1. Microbiology of food and animal feeding stuffs –Horizontal method for the enumeration of coagulase-positive staphylococci (*Staphylococcus aureus* and other species), Part 1: Technique using Baird-Parker agar medium.
8. ISO 6888-2. Microbiology of food and animal feeding stuffs –Horizontal method for the enumeration of coagulase-positive staphylococci (*Staphylococcus aureus* and other species), Part 2: Technique using rabbit plasma fibrinogen agar medium.

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3M 01/09-04/03A and B

ALTERNATIVE ANALYTICAL METHODS FOR AGRIBUSINESS

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1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Practice Studies.
2. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
3. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
4. U.S. Food and Drug Administration. 1998. Bacteriological Analytical Manual, 8th ed., (Revision A), Appendix 3.64.
5. ISO 6887. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination,
6. ISO 16140-2. Microbiology of the food chain – Method Validation – Protocol for the validation of alternative (proprietary) methods against a reference method

7. ISO 6888-1. Microbiology of food and animal feeding stuffs –Horizontal method for the enumeration of coagulase-positive staphylococci (*Staphylococcus aureus* and other species), Part 1: Technique using Baird-Parker agar medium.
8. ISO 6888-2. Microbiology of food and animal feeding stuffs –Horizontal method for the enumeration of coagulase-positive staphylococci (*Staphylococcus aureus* and other species), Part 2: Technique using rabbit plasma fibrinogen agar medium.

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