Product Instructions

-  **(EN)** Rapid Aerobic Count Plate
-  **(FR)** Test Rapide pour la numération de la Flore Totale Aérobie
-  **(DE)** Rapid Aerobe Gesamtkeimzahl Auszählplatte
-  **(IT)** Piastra per la conta aerobica rapida
-  **(ES)** Placa Recuento Aeróbico Rápido
-  **(NL)** Snelle Aeroob Kiemgetal Telplaat
-  **(SV)** Rapid Aerobic Count odlingsplatta
-  **(DA)** Rapid Aerobic Count Plate
-  **(NO)** Hurtigfilm for totalkim
-  **(FI)** Rapid Aerobien Kasvatusalusta
-  **(PT)** Placa para Contagem Rápida de Aeróbios
-  **(EL)** Πλακίδιο Ταχείας Αερόβιας Απαρίθμησης
-  **(PL)** Płytko do szybkiego oznaczania ogólnej liczby drobnoustrojów
-  **(RU)** Тест-пластина для экспресс-подсчета ОМЧ
-  **(TR)** Hızlı Toplam Canlı Sayım Plakası
-  **(JA)** 生菌数迅速測定用プレート
-  **(ZH)** 快速菌落总数测试片
-  **(TH)** แผ่นอาหารเลี้ยงเชื้อสำเร็จรูปสำหรับนับจำนวนแอโรบิคแบคทีเรียอย่างรวดเร็ว
-  **(KO)** 속성 일반세균 측정용 플레이트

Rapid Aerobic Count Plate

Product Instructions

Product Description and Intended Use

The Neogen® Petrifilm® Rapid Aerobic Count (RAC) Plate is a sample-ready-culture-medium system which contains nutrients, a cold-water-soluble gelling agent, and a dual-sensing indicator technology that facilitates colony enumeration. Neogen Petrifilm RAC Plates are used for the enumeration of aerobic bacteria in the food and beverage industries. Neogen Petrifilm RAC Plate components are decontaminated though not sterilized. Neogen Food Safety is certified to International Organization for Standardization (ISO) 9001 for design and manufacturing. Neogen Petrifilm RAC Plates have not been evaluated with all possible food products, food processes, testing protocols or with all possible microorganism strains.

Safety

The user should read, understand, and follow all safety information in the instructions for the Neogen Petrifilm RAC Plate. Retain the safety instructions for future reference.

⚠ WARNING: Indicates a hazardous situation, which, if not avoided, could result in death or serious injury and/or property damage.

⚠ WARNING

To reduce the risks associated with exposure to biohazards and environmental contamination:

- Follow current industry standards and local regulations for disposal of biohazardous waste.

To reduce the risks associated with release of contaminated product:

- Follow all product storage instructions contained in the instructions for use.
- Do not use beyond the expiration date.

To reduce the risks associated with bacterial infection and workplace contamination:

- Perform Neogen Petrifilm RAC Plate testing in a properly equipped laboratory under the control of a skilled microbiologist.
- The user must train personnel in current proper testing techniques: for example, Good Laboratory Practices¹, ISO 17025² or ISO 7218³.

To reduce the risks associated with misinterpretation of results:

- Neogen has not documented Neogen Petrifilm RAC Plates for use in industries other than food and beverage. For example, Neogen has not documented Neogen Petrifilm RAC Plates for testing water, pharmaceuticals, or cosmetics.
- Do not use Neogen Petrifilm RAC Plates in the diagnosis of conditions in humans or animals.
- Validation of the Neogen Petrifilm RAC Plate method for U.S. -recognized laboratory pasteurized count is the responsibility of the end-user.
- Acceptance of the Neogen Petrifilm RAC Plate method for the testing of water per an accepted local government regulation is at the sole discretion and responsibility of the end-user.
- Neogen Petrifilm RAC Plates do not differentiate any one microorganism strain from another.

Consult the Safety Data Sheet for additional information.

If you have questions about specific applications or procedures, please visit our website at www.neogen.com or contact your local Neogen representative or distributor.

User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at www.neogen.com, or contact your local Neogen representative or distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any Neogen Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

Limitation of Warranties / Limited Remedy

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, NEOGEN DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any Neogen Food Safety Product is defective, Neogen or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify Neogen within sixty days of discovery of any suspected defects in a product and return it to Neogen. Please contact your Neogen representative or authorized Neogen distributor for any further questions.

Limitation of Neogen Liability

NEOGEN WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall Neogen's liability under any legal theory exceed the purchase price of the product alleged to be defective.

Storage

Store unopened Neogen Petrifilm RAC Plate pouches refrigerated or frozen (-20 to 8°C / -4 to 46°F). Just prior to use, allow unopened pouches to come to room temperature before opening (20-25°C / <60% RH). Return unused Neogen Petrifilm RAC Plates to pouch. Seal by folding the end of the pouch over and applying adhesive tape. **To prevent exposure to moisture, do not refrigerate opened pouches.** Store resealed pouches in a cool dry place for no longer than one month. It is recommended that resealed pouches of Neogen Petrifilm RAC Plates be stored in a freezer (see below) if the laboratory temperature exceeds 25°C (77°F) and/or the laboratory is located in a region where the relative humidity exceeds 50% (with the exception of air-conditioned premises).

To store opened pouches in a freezer, place Neogen Petrifilm RAC Plates in a sealable container. To remove frozen Neogen Petrifilm RAC Plates for use, open the container, remove the plates that are needed and immediately return remaining plates to the freezer in the sealed container for the remainder of the shelf life. Neogen Petrifilm RAC Plates should not be used past their expiration date. The freezer that is used for open pouch storage must not have an automatic defrost cycle as this would repeatedly expose the plates to moisture which can damage the plates.

Do not use Neogen Petrifilm RAC Plates that show discoloration. Expiration date and lot number are noted on each package of Neogen Petrifilm RAC Plates. The lot number is also noted on individual Neogen Petrifilm RAC Plates.

⚠ Disposal

After use, Neogen Petrifilm RAC Plates may contain microorganisms that may be a potential biohazard. Follow current local, regional, national and industry standards for disposal.

Instructions for Use

Follow all instructions carefully. Failure to do so may lead to inaccurate results.

Sample Preparation

1. Use appropriate sterile diluents:

Butterfields phosphate buffer, 0.1% peptone water, peptone salt diluent, buffered peptone water, saline solution (0.85-0.90%), Neogen® Wide-Spectrum Neutralizer, bisulfite-free letheen broth or distilled water. See section "Specific Instructions for Validated Methods" for specific requirements.

Do not use diluents containing citrate, bisulfite or thiosulfate with Neogen Petrifilm RAC Plates; they can inhibit growth. If citrate buffer is indicated in the standard procedure, substitute with one of the buffers listed above, warmed to 40-45°C (104-113°F).

2. Blend or homogenize sample per standard.

3. For optimal growth and recovery of microorganisms in acidic products (<pH 5), adjust the pH of the sample suspension to a pH greater than pH 5. For acidic products, adjust the pH with 1N NaOH.

Plating

1. Place the Neogen Petrifilm RAC Plate on a flat, level surface.

2. Lift the top film and with the pipette perpendicular dispense 1 mL of sample suspension onto the center of bottom film.

3. Roll the top film down onto the sample to prevent trapping bubbles.

4. Place the Neogen® Petrifilm® Flat Spreader (6425) on the center of the Neogen Petrifilm RAC Plate. Press on the center of the spreader to distribute the sample evenly. Spread the inoculum over the entire Neogen Petrifilm RAC Plate growth area before the gel is formed. Do not slide the spreader across the film.

5. Remove the Neogen Petrifilm Flat Spreader and leave the Neogen Petrifilm RAC Plate undisturbed for at least one minute to permit the gel to form.

Incubation

Incubate Neogen Petrifilm RAC Plates in a horizontal position with the clear side up in stacks of no more than 40*. Several incubation times and temperatures can be used depending on current local reference methods, some of which are listed in the Specific Instructions section.

*For Standard Methods for the Examination of Dairy Products⁴ maximum stack height is 20.

Interpretation

1. Neogen Petrifilm RAC Plates can be counted using a standard colony counter or other illuminated magnifier. Count all colonies regardless of color, size or intensity.
2. The circular growth area is approximately 30 cm². Gridlines are visible with the use of a backlight to assist with estimated enumeration. Estimates can be made on Neogen Petrifilm RAC Plates containing greater than 300 colonies by counting the number of colonies in two or more representative squares and determining the average number per square. Multiply the average number by 30 to determine the estimated count per plate.
3. High concentrations of colonies on the Neogen Petrifilm RAC Plates will cause the entire growth area to become blue or red. Occasionally, on overcrowded Neogen Petrifilm RAC Plates, the center may lack visible colonies, but many small colonies can be seen on the edges. When any of these occurs, record results as too numerous to count (TNTC). When an actual count is required, plate at a higher dilution.
4. Food samples may occasionally show interference on the Neogen Petrifilm RAC Plates, for example:
 - a. A uniform light blue background color (often seen from the organisms used in cultured products) should not be counted as TNTC.
 - b. Intense, pinpoint blue specs (often seen with spices or granulated products) should not be counted as colonies.
5. Where necessary, colonies may be isolated for further identification. Lift the top film and pick the colony from the gel. Test using standard procedures.
6. If the Neogen Petrifilm RAC Plates cannot be counted within 1 hour of removal from the incubator, they may be stored for later enumeration by freezing in a sealable container at temperatures less than or equal to minus 15°C (5°F) for no longer than one week.

For further information refer to the “Neogen® Petrifilm® Rapid Aerobic Count Plate Interpretation Guide.” If you have questions about specific applications or procedures, please visit our website at www.neogen.com or contact your local Neogen representative or distributor.

Specific Instructions for Validated Methods

AOAC® Official MethodsSM 2015.13

AOAC® Performance Tested Certificate #121403

In AOAC Official Methods of Analysis (OMA) and Performance Tested Method (PTM) studies, the Neogen Petrifilm Rapid Aerobic Count Plate method was found to be equivalent to or better than the average log counts of the FDA/BAM Chapter 3 reference method and the Standard Methods for the Examination of Dairy Products Chapter 6.

Scope of Validation:

Raw ground beef, raw ground pork, ground turkey, chicken carcass rinsate, fresh swai, fresh tuna, fresh tiger shrimp, easy-peel shrimp, cherry tomato wash, frozen blueberries, Mediterranean apricots, creamy salad dressing, fresh pasta, vanilla ice cream, dry milk powder, pasteurized skim milk and environmental samples from stainless steel, rubber and sealed concrete.

Incubation:

Dairy (not including powders) & Seafood:

Incubate Neogen Petrifilm Rapid Aerobic Count Plates 24 hours ± 2 hours at 32°C ± 1°C.

All other foods and environmental surface samples:

Incubate Neogen Petrifilm Rapid Aerobic Count Plates 24 hours ± 2 hours at 35°C ± 1°C.

**Dairy powders including whey powders:**

Incubate Neogen Petrifilm Rapid Aerobic Count Plates 48 hours \pm 3 hours at 32°C \pm 1°C.

NF Validation by AFNOR Certification

NF Validation certified method in compliance with ISO 16140-2⁷ in comparison to ISO 4833-1⁸

Use the following details when implementing the above Instructions for Use:

Scope of the validation:

For testing dairy products

Sample preparation:

Use only ISO listed diluents⁶.

Incubation:

For non-powdered dairy products: Incubate Neogen Petrifilm RAC Plates 28 hours \pm 2 hours at 30°C \pm 1°C.

For powdered milks: Incubate Neogen Petrifilm RAC Plates 48 hours \pm 3 hours at 30°C \pm 1°C.

Interpretation:

Calculate the number of microorganisms present in the test sample according to ISO 7218³ for one plate per dilution. For calculation, take into account only Neogen Petrifilm RAC Plates that contain up to 300 colonies. Estimates are outside of the scope of the NF Validation Certification (cf. interpretation part paragraph 3). Refer to EN ISO 7218 standard for inoculation, colony counting and calculation and expression of results.

3M 01/17 – 11/16

ALTERNATIVE ANALYTICAL METHODS FOR AGRIBUSINESS

<http://nf-validation.afnor.org/en>

For more information about end of validity, please refer to NF VALIDATION certificate available on the website mentioned above.

References

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
3. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
4. Wehr, Michael H. and Joseph F. Frank, *Standard Methods for the Examination of Dairy Products*. American Public Health Association, Jan. 1, 2004. 17th ed.
5. FDA. Bacteriological Analytical Manual (BAM), Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
6. ISO 6887. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination.
7. ISO 16140-2. Microbiology of the food chain - Method Validation - Protocol for the validation of alternative (proprietary) methods against a reference method.
8. ISO 4833-1. Microbiology of the food chain - Horizontal method for the enumeration of microorganisms - Part 1: colony-count technique at 30°C by the pour plate technique.

Please refer to the current versions of the standard methods listed above

Explanation of Symbols

info.neogen.com/symbols

AOAC is a registered trademark of AOAC INTERNATIONAL

Official Methods is a service mark of AOAC INTERNATIONAL

Performance Tested Methods is a service mark of AOAC INTERNATIONAL

Neogen Food Safety

Neogen Corporation

620 Leshar Place
Lansing, MI 48912 USA
Neogen.com

Neogen Europe Ltd.

The Dairy School
Auchincruive
Ayr, KA6 5HU
Scotland, UK

Neogen Ireland, Ltd.

Bray Business Park, Bray
Co. Wicklow
A98YV29, Ireland

Neogen Corporation

620 Leshar Place Lansing, MI 48912 USA
www.neogen.com

Neogen and Petrifilm are trademarks of Neogen.
© 2024, Neogen.

Unauthorized use prohibited. All rights reserved.
All other trademarks are the property of their
respective companies.

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Références

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
3. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
4. Wehr, Michael H. and Joseph F. Frank, *Standard Methods for the Examination of Dairy Products*. American Public Health Association, Jan. 1, 2004. 17th ed.
5. FDA. Bacteriological Analytical Manual (BAM), Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
6. ISO 6887. Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination.
7. ISO 16140-2 (2016). Microbiology of the food chain - Method Validation - Protocol for the validation of alternative (proprietary) methods against a reference method.
8. ISO 4833-1 (2013). Microbiology of the food chain - Horizontal method for the enumeration of microorganisms - Part 1: colony-count technique at 30°C by the pour plate technique.

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