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2 Performance Tested MethodsSM. Program certificate number 022104. AOAC Research Institute.

3 Laboratory procedures for the microbiological analysis of foods. Vol 3. The compendium of analytical methods. MFPL-41: Environmental sampling for the detection of microorganisms. 2010. Available at: <https://www.canada.ca/en/health-canada/services/food-nutrition/research-programs-analytical-methods/analytical-methods/compendium-methods/laboratory-procedures-microbiological-analysis-foods-compendium-analytical-methods.html>.

4 Lethen broth: A neutralizing solution for iodine, chlorine, quaternary ammonium and acid sanitizers. 3M Food Safety. 1985.

5 Sample handling sponges - new sponge qualification performance summary. TB.119.00. 3M Food Safety. 2012.

6 Environmental monitoring procedures. 3M Food Safety. 2018. Available at: <https://multimedia.3m.com/mws/media/2411110/environmental-monitoring-procedures-article.pdf>

7 Ellerbroek, L., Johne, R., Schrader, C., and Schielke, A. (2012). PCR inhibitors - occurrence, properties, and removal. *Journal for Applied Microbiology*. 1365-2672.05384.x.

8 Jersek, B., Klančnik, A., Kovac, M., Piskernik, S., and Toplak, N. (2012). PCR in food analysis, polymerase chain reaction, Patricia Hernandez-Rodriguez and Arlen Patricia Ramirez Gomez. IntechOpen, DOI: 10.5772/38551. Available at: <http://www.intechopen.com/books/polymerase-chain-reaction/pcr-in-food-analysis>

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10 Ignatovich, I., Leishman, O., Podtburg, T., Steinagel, S. (2017). Comparison of neutralizing buffered peptone water and Dey/Engley Broth in the recovery of *Salmonella enterica* from broiler carcass rinsates. *J Food Protection*. 80 (Supplement A): 163.

11 Chen, J., Park, Y.J. (2011). Mitigating the antimicrobial activities of selected organic acids and commercial sanitizers with various neutralizing agents. *Journal of Food Protection*. 74: 820-825. Available at: <https://doi.org/10.4315/0362-028X.JFP-10-447>.

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