Product Instructions

-  (EN) Select *E. coli* Count Plate
-  (FR) Test pour la numération sélective des *E. coli*
-  (DE) selektive *E. coli* (SEC) Zählplatte
-  (IT) Piastra per il conteggio selettivo di *E. coli*
-  (ES) Placa para recuento selectivo de *E. coli*
-  (NL) Select *E. coli* Telplaat
-  (SV) Selektiv odlingsplatta för *E. coli*
-  (DA) Select *E. coli* Tælleplade
-  (NO) Select *E. coli*
-  (FI) Selektiivinen *E. coli* kasvatusalusta
-  (PT) Seleção a Placa para Contagem de *E. coli*
-  (EL) Πλακίδιο Επιλεκτικής Καταμέτρησης *E. coli*
-  (PL) Płytko do selektywnego oznaczania liczby *E. coli*
-  (RU) Тест-пластина для подсчёта *E. coli*
-  (TR) Hızlı *E. coli* Sayım Plakası
-  (JA) 大腸菌選択用プレート (SECプレート)
-  (ZH) 选择性大肠杆菌测试片
-  (TH) แผ่นอาหารเลี้ยงเชื้อสำหรับนับจำนวนอีโคไลแบบจำเพาะ
-  (KO) 선택 대장균 측정용 플레이트

Product Instructions

Select *E. coli* Count Plate

Product Description And Intended Use

The Neogen® Petrifilm® Select *E. coli* Count (SEC) Plate is a sample-ready-culture medium system which contains selective agents, nutrients, a cold-water-soluble gelling agent, and an indicator of glucuronidase activity, 5-bromo-4-chloro-3-indolyl-D-glucuronide (BCIG), which facilitates colony enumeration. Neogen Petrifilm SEC Plates are used for the enumeration of *Escherichia coli* (*E. coli*) in the food and beverage industries. Most *E. coli* strains are thermotolerant and produce beta-glucuronidase, an enzyme that reacts with BCIG indicator in Neogen Petrifilm SEC Plates to produce dark green to blue-green colonies. Neogen Petrifilm SEC Plate components are decontaminated though not sterilized. Neogen Food Safety is certified to International Organization for Standardization (ISO) 9001 for design and manufacturing. Neogen Petrifilm SEC Plates have not been tested with all possible food products, food processes, testing protocols or with all possible strains of *E. coli* or other bacteria.

Safety

The user should read, understand, and follow all safety information in the instructions for the Neogen Petrifilm SEC Plate. Retain the safety instructions for future reference.

⚠ **WARNING:** Indicates a hazardous situation, which, if not avoided, could result in death or serious injury and/or property damage.

⚠ **CAUTION:** Indicates a hazardous situation, which, if not avoided, could result in minor or moderate injury and/or property damage.

⚠ WARNING

Do not use Neogen Petrifilm SEC Plate for the detection of *E. coli* O157. Like most other *E. coli* media, the Neogen Petrifilm SEC Plate will not specifically indicate whether any O157 strain is present. Because most O157 strains are atypical *E. coli*, for example they are glucuronidase negative and they will not produce a blue-green color, they will not be detected on Neogen Petrifilm SEC Plate.

To reduce the risks associated with release of contaminated product:

- Follow all product storage instruction contained in the instructions for use.
- Do not use beyond expiration date.

To reduce the risk associated with bacterial infection and workplace contamination:

- Perform Neogen Petrifilm SEC Plate testing in a properly equipped laboratory under the control of a skilled microbiologist.
- The user must train its personnel in current proper testing techniques: for example, Good Laboratory Practices¹, ISO 17025², or ISO 7218³.

To reduce the risks associated with misinterpretation of results:

- Neogen has not documented Neogen Petrifilm SEC Plates for use in industries other than food and beverage. For example, Neogen has not documented Neogen Petrifilm SEC Plates for testing water, pharmaceuticals, or cosmetics.
- Do not use Neogen Petrifilm SEC Plates in the diagnosis of conditions in humans or animals.
- Neogen Petrifilm SEC Plates do not differentiate any one *E. coli* strain from another.
- High sugar content foods may decrease the color reaction of the glucuronidase and produce pale green to colorless colonies.

To reduce the risks associated with environmental contamination:

- Follow current industry standards and local regulations for disposal of contaminated waste.

⚠ CAUTION

A few strains of bacteria other than *E. coli*, such as some *Shigella*, *Salmonella*, *Enterobacter* and *Klebsiella*, can produce beta-glucuronidase and will produce dark green to blue-green colonies on the Neogen Petrifilm SEC Plate, especially when incubated at lower than recommended temperatures⁴.

If you have questions about specific applications or procedures, please visit our website at www.neogen.com or contact your Neogen representative or authorized distributor.

User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at www.neogen.com, or contact your Neogen representative or authorized distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any Neogen Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

Limitation of Warranties / Limited Remedy

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, NEOGEN DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any Neogen Food Safety Product is defective, Neogen or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify Neogen within sixty days of discovery of any suspected defects in a product and return it to Neogen. Please contact your Neogen representative or authorized Neogen distributor for any further questions.

Limitation of Neogen Liability

NEOGEN WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall Neogen's liability under any legal theory exceed the purchase price of the product alleged to be defective.

Storage

Store unopened Neogen Petrifilm SEC Plate pouches refrigerated or frozen at temperatures lower than or equal to 8°C (46°F). Just prior to use, allow unopened pouches to come to room temperature before opening. Return unused Neogen Petrifilm SEC Plates to pouch. Seal by folding the end of the pouch over and applying adhesive tape. **To prevent exposure to moisture, do not refrigerate opened pouches.** Store resealed pouches in a cool dry place for no longer than four weeks. It is recommended that resealed pouches of Neogen Petrifilm SEC Plates be stored in a freezer (see below) if the laboratory temperature exceeds 25°C (77°F) and/or the laboratory is located in a region where the relative humidity exceeds 50% (with the exception of air-conditioned premises).

To store opened pouches in a freezer, place Neogen Petrifilm SEC Plates in a sealable container. To remove frozen Neogen Petrifilm SEC Plates for use, open the container, remove the plates that are needed and immediately return remaining plates to the freezer in the sealed container. Neogen Petrifilm SEC Plates should not be used past their expiration date. The freezer that is used for open pouch storage must not have an automatic defrost cycle as this would repeatedly expose the Neogen Petrifilm SEC Plates to moisture which can damage the plates.

Do not use Neogen Petrifilm SEC Plates that show discoloration. Expiration date and lot number are noted on each package of Neogen Petrifilm SEC Plates. The lot number is also noted on individual Neogen Petrifilm SEC Plates.

⚠ Disposal

After use, Neogen Petrifilm SEC Plates may contain microorganisms that may be a potential biohazard. Follow local, regional, and national regulations for disposal requirements.

Instructions for Use

Follow all instructions carefully. Failure to do so may lead to inaccurate results.

Sample Preparation

1. Use appropriate sterile diluents: Butterfield's phosphate buffered dilution water⁵, 0.1% peptone water⁵, peptone salt diluent⁶, buffered peptone water⁶, quarter-strength Ringer's solution, dipotassium hydrogen phosphate⁶, saline solution (0.85-0.90%), Neogen® Wide-Spectrum Neutralizer⁶, bisulfite-free letheen broth⁶, or distilled water.

See section "Specific Instructions for Validated Methods" for specific requirements.

Do not use diluents containing citrate, bisulfite or thiosulfate with Neogen Petrifilm SEC Plates; they can inhibit growth. If citrate buffer is indicated in the standard procedure, substitute with one of the buffers listed above, warmed to 40–45°C (104–113°F).

2. Blend or homogenize sample.
3. For optimal growth and recovery of microorganisms, adjust the pH of the sample suspension to 6.5 - 7.5. For acidic products, adjust the pH with 1N NaOH. For alkaline products, adjust the pH with 1N HCl.

Plating

1. Place the Neogen Petrifilm SEC Plate on a flat, level surface.
2. Lift the top film and with the pipette perpendicular to the inoculation area dispense 1 mL of sample suspension onto the center of bottom film.
3. Roll the top film down onto the sample to prevent trapping air bubbles.
4. Place the Neogen® Petrifilm® Spreader with the flat side down on the center of the Neogen Petrifilm SEC Plate. Press gently on the center of the Neogen Petrifilm Spreader to distribute the sample evenly. Spread the inoculum over the entire Neogen Petrifilm SEC Plate growth area before the gel is formed. Do not slide the Neogen Petrifilm Spreader across the film.
5. Remove the Neogen Petrifilm Spreader and leave the Neogen Petrifilm SEC Plate undisturbed for at least one minute to permit the gel to form.

Incubation

Incubate Neogen Petrifilm SEC Plates in a horizontal position with the clear side up in stacks of no more than 20 plates. Incubate Neogen Petrifilm SEC Plates 24 hours ± 2 hours at 42°C ± 1°C. Several incubation times and temperatures can be used depending on current local reference methods, some of which are listed in the section below titled “**Specific Instructions for Validated Methods**”.

Interpretation

1. Neogen Petrifilm SEC Plates can be counted using a standard colony counter or other illuminated magnifier. Do not count colonies on the foam dam since they are removed from the selective influence of the medium. Do not count artifact bubbles that may be present. *E. coli* produce dark green to light-green colonies. Colonies may have gas bubbles associated with them. Count all green to blue-green colonies as *E. coli* whether there are gas bubbles present or not.

⚠ WARNING

Do not use the Neogen Petrifilm SEC Plate for the detection of *E. coli* O157. Because most *E. coli* O157 strains are atypical, for example, they are glucuronidase negative, they will not produce a blue-green color, and will not be interpreted as *E. coli* on the Neogen Petrifilm SEC Plates.

2. The circular growth area is approximately 20 cm². Estimates can be made on Neogen Petrifilm SEC Plates containing greater than 150 colonies by counting the number of colonies in two or more representative squares and determining the average number per square. Multiply the average number by 20 to determine the estimated count per Neogen Petrifilm SEC Plate.
3. When colonies are present in large numbers, Neogen Petrifilm SEC Plates will cause the entire growth area to become gray or blue and either or both of the following characteristics: many small, indistinct colonies and/or many gas bubbles. When this occurs, record results as too numerous to count (TNTC). When an actual count is required, plate at a higher dilution.
4. Where necessary, colonies may be isolated for further identification. Lift the top film using proper testing technique and pick the colony from the gel. Test using standard procedures.
5. If the Neogen Petrifilm SEC Plates cannot be counted within 1 hour of removal from the incubator, they may be stored for later enumeration by freezing in a sealable container at temperatures lower than or equal to negative 15°C (5°F) for no longer than one week*.

For further information refer to the “Neogen® Petrifilm® Select *E. coli* Count Plate Interpretation Guide.” If you have questions about specific applications or procedures, please visit our website at www.neogen.com or contact your Neogen representative or authorized distributor.

Specific Instructions for Validated Methods

NF Validation by AFNOR Certification:

NF Validation certified method 3M 01/08–06/01 is in compliance with ISO 16140-2⁷ in comparison to ISO 16649-2⁸

Use the following details when implementing the above Instructions for Use:

Scope of the validation:

For testing all human food products, pet food and industrial environmental samples.

Sample preparation:

Use only ISO listed diluents.

Incubation:

Incubate Neogen Petrifilm SEC Plates 24 hours $\pm$ 2 hours at 42°C $\pm$ 1°C.

Interpretation:

The Petrifilm SEC Plate has been validated for enumeration using the Petrifilm Plate Reader Advanced (PPRA) and visual interpretation per the general method instructions. Refer to the NF Validation Certificate below for the PPRA software version used in the context of the NF Validation. For instructions on using the PPRA and software, please refer to the user manual. For the latest version visit neogen.com.

*Freezing Petrifilm SEC Plates prior to enumeration is not included in the scope of the NF Validation.

Precautions:

Use good microbiology laboratory practices, such as ISO 7218. Enumeration using a single plate and a single dilution is an option certified in the NF Validation by AFNOR Certification. In this context, the requirements of ISO 7218 concerning the use of two successive dilutions or two plates of the same dilution cannot be applied (see paragraphs “Inoculation” and “Calculation and expression of results”). For sample preparation, refer to ISO 6887 series as specified in the standards. Matrices with a high level of particulate may influence Petrifilm Plate Reader Advanced (PPRA) enumeration and should be reviewed by the user.

For more information about end of validity, please refer to NF Validation certificate available on the website mentioned below.

3M 01/08 –06/01

ALTERNATIVE ANALYTICAL METHODS FOR AGRIBUSINESS

<http://nf-validation.afnor.org/en>

References

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025:2017 General requirements for the competence of testing and calibration laboratories.
3. ISO 7218. Microbiology of the food chain - General requirements and guidance for microbiological examinations.
4. Frampton, E.W. and L. Restaino, 1993. Methods for *Escherichia coli* identification in food, water and clinical samples based on beta-glucuronidase detection. *Journal of Applied Bacteriology* 74, 223-233.
5. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
6. ISO 6887 (all parts) Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination.
7. ISO 16140-2:2016 Microbiology of the food chain – Method validation – Protocol for the validation of alternative (proprietary) methods against a reference method.
8. ISO 16649-2:2001 Microbiology of food and animal feeding stuffs- Horizontal method for the enumeration of β -glucuronidase-positive *Escherichia coli* - Part 2: Colony- count technique at 44°C using 5-bromo-4-chloro-3-indolyl β D-glucuronide.

Refer to the current versions of the standard methods listed above.

Explanation Of Symbols

info.neogen.com/symbols

Neogen Food Safety

Neogen Corporation

620 Leshar Place
Lansing, MI 48912 USA
Neogen.com

Neogen Europe Ltd.

The Dairy School
Auchincruive
Ayr, KA6 5HU
Scotland, UK

Neogen Ireland, Ltd.

Bray Business Park, Bray
Co. Wicklow
A98YV29, Ireland

Neogen Corporation

620 Leshar Place Lansing, MI 48912 USA
www.neogen.com

© Neogen Corporation 2025. All rights reserved.
Neogen and Petrifilm are registered trademarks of Neogen Corporation.
FS00840C-1

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Références

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
3. ISO 7218. Microbiology of food and animal feeding stuffs – General rules for microbiological examination.
4. Frampton, E.W. and L. Restaino, 1993. Methods for *Escherichia coli* identification in food, water and clinical samples based on beta-glucuronidase detection. *Journal of Applied Bacteriology* 74, 223-233.
5. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
6. ISO 6887 Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination.
7. ISO 16140-2. Microbiology of the food chain – Method Validation – Protocol for the validation of alternative (proprietary) methods against a reference method.
8. ISO 16649-2: Microbiology of food and animal feeding stuffs- Horizontal method for the enumeration of β -glucuronidase-positive *Escherichia coli* - Part 2: Colony- count technique at 44°C using 5-bromo-4-chloro-3-indolyl β D-glucuronide.

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Literaturangaben

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
3. ISO 7218. Microbiology of food and animal feeding stuffs – General rules for microbiological examination.
4. Frampton, E.W. and L. Restaino, 1993. Methods for *Escherichia coli* identification in food, water and clinical samples based on beta-glucuronidase detection. *Journal of Applied Bacteriology* 74, 223-233.
5. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM gefunden unter: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
6. ISO 6887 Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination.
7. ISO 16140-2. Microbiology of the food chain – Method Validation – Protocol for the validation of alternative (proprietary) methods against a reference method.
8. ISO 16649-2. Microbiology of food and animal feeding stuffs- Horizontal method for the enumeration of β -glucuronidase-positive *Escherichia coli* - Part 2: Colony- count technique at 44°C using 5-bromo-4-chloro-3-indolyl β D-glucuronide.

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1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
3. ISO 7218. Microbiology of food and animal feeding stuffs – General rules for microbiological examination.
4. Frampton, E.W. and L. Restaino, 1993. Methods for *Escherichia coli* identification in food, water and clinical samples based on beta-glucuronidase detection. *Journal of Applied Bacteriology* 74, 223-233.
5. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
6. ISO 6887 Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination.

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1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
3. ISO 7218. Microbiology of food and animal feeding stuffs – General rules for microbiological examination.
4. Frampton, E.W. and L. Restaino, 1993. Methods for *Escherichia coli* identification in food, water and clinical samples based on beta-glucuronidase detection. *Journal of Applied Bacteriology* 74, 223-233.
5. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.

6. ISO 6887 Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination.
7. ISO 16140-2. Microbiology of the food chain – Method Validation – Protocol for the validation of alternative (proprietary) methods against a reference method.
8. ISO 16649-2: Microbiology of food and animal feeding stuffs- Horizontal method for the enumeration of β -glucuronidase-positive *Escherichia coli* - Part 2: Colony- count technique at 44°C using 5-bromo-4-chloro-3-indolyl β D-glucuronide.

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1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
3. ISO 7218. Microbiology of food and animal feeding stuffs – General rules for microbiological examination.
4. Frampton, E.W. and L. Restaino, 1993. Methods for *Escherichia coli* identification in food, water and clinical samples based on beta-glucuronidase detection. *Journal of Applied Bacteriology* 74, 223-233.
5. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
6. ISO 6887 Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination.
7. ISO 16140-2. Microbiology of the food chain – Method Validation – Protocol for the validation of alternative (proprietary) methods against a reference method.
8. ISO 16649-2: Microbiology of food and animal feeding stuffs- Horizontal method for the enumeration of β-glucuronidase-positive *Escherichia coli* - Part 2: Colony- count technique at 44°C using 5-bromo-4-chloro-3-indolyl βD-glucuronide.

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Źródła

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
3. ISO 7218. Microbiology of food and animal feeding stuffs – General rules for microbiological examination.
4. Frampton, E.W. and L. Restaino, 1993. Methods for *Escherichia coli* identification in food, water and clinical samples based on beta-glucuronidase detection. *Journal of Applied Bacteriology* 74, 223-233.
5. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
6. ISO 6887 Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination.
7. ISO 16140-2. Microbiology of the food chain – Method Validation – Protocol for the validation of alternative (proprietary) methods against a reference method.
8. ISO 16649-2: Microbiology of food and animal feeding stuffs- Horizontal method for the enumeration of β -glucuronidase-positive *Escherichia coli* - Part 2: Colony-count technique at 44°C using 5-bromo-4-chloro-3-indolyl β D-glucuronide.

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1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025:2017 General requirements for the competence of testing and calibration laboratories.
3. ISO 7218. Microbiology of the food chain - General requirements and guidance for microbiological examinations.
4. Frampton, E.W. and L. Restaino, 1993. Methods for Escherichia coli identification in food, water and clinical samples based on beta-glucuronidase detection. Journal of Applied Bacteriology 74, 223-233.
5. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
6. ISO 6887 (all parts) Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination.
7. ISO 16140-2:2016 Microbiology of the food chain – Method validation – Protocol for the validation of alternative (proprietary) methods against a reference method.
8. ISO 16649-2:2001 Microbiology of food and animal feeding stuffs- Horizontal method for the enumeration of β -glucuronidase-positive Escherichia coli - Part 2: Colony- count technique at 44°C using 5-bromo-4-chloro-3-indolyl β D-glucuronide.

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Referanslar

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
3. ISO 7218. Microbiology of food and animal feeding stuffs – General rules for microbiological examination.
4. Frampton, E.W. and L. Restaino, 1993. Methods for *Escherichia coli* identification in food, water and clinical samples based on beta-glucuronidase detection. *Journal of Applied Bacteriology* 74, 223-233.
5. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM şu adreste mevcuttur: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
6. ISO 6887 Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination.

7. ISO 16140-2. Microbiology of the food chain – Method Validation – Protocol for the validation of alternative (proprietary) methods against a reference method.
8. ISO 16649-2: Microbiology of food and animal feeding stuffs- Horizontal method for the enumeration of β -glucuronidase-positive *Escherichia coli* - Part 2: Colony- count technique at 44°C using 5-bromo-4-chloro-3-indolyl β D-glucuronide.

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1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.
3. ISO 7218. Microbiology of food and animal feeding stuffs – General rules for microbiological examination.
4. Frampton, E.W. and L. Restaino, 1993. Methods for *Escherichia coli* identification in food, water and clinical samples based on beta-glucuronidase detection. *Journal of Applied Bacteriology* 74, 223-233.
5. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>.
6. ISO 6887 Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination.
7. ISO 16140-2. Microbiology of the food chain – Method Validation – Protocol for the validation of alternative (proprietary) methods against a reference method.
8. ISO 16649-2: Microbiology of food and animal feeding stuffs- Horizontal method for the enumeration of β -glucuronidase-positive *Escherichia coli*- Part 2: Colony- count technique at 44°C using 5-bromo-4-chloro-3-indolyl β D-glucuronide.

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1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Studies.
2. ISO/IEC 17025 General requirements for the competence of testing and calibration laboratories.
3. ISO 7218 Microbiology of food and animal feeding stuffs – General rules for microbiological examination.
4. Frampton, E.W. and L. Restaino, 1993. Methods for *Escherichia coli* identification in food, water and clinical samples based on beta-glucuronidase detection. *Journal of Applied Bacteriology* 74, 223-233.
5. FDA. Bacteriological Analytical Manual (BAM), 8th Edition, Revision A, 1998. Reagents Index for BAM found at: <http://www.fda.gov/Food/FoodScienceResearch/LaboratoryMethods/ucm055791.htm>
6. ISO 6887 Microbiology of food and animal feeding stuffs – Preparation of test samples, initial suspension and decimal dilutions for microbiological examination.
7. ISO 16140-2. Microbiology of the food chain – Method Validation – Protocol for the validation of alternative (proprietary) methods against a reference method.

8. ISO 16649-2: Microbiology of food and animal feeding stuffs- Horizontal method for the enumeration of β -glucuronidase-positive *Escherichia coli* - Part 2: Colony- count technique at 44°C using 5-bromo-4-chloro-3-indolyl β D-glucuronide.

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