700006967 | NCM1018A
700006968 | NCM1018B
700006969 | NCM1018C
700006970 | NCM1018D

Product Instructions

-  **(EN)** One Plate Total Viable Count (OP-TVC)
-  **(FR)** One Plate Numération viable totale (OP-TVC)
-  **(DE)** One Plate Gesamtkeimzahl (OP-TVC)
-  **(ES)** One Plate Recuento viable total (OP-TVC)
-  **(PT)** One Plate Contagem viável total (OP-TVC)
-  **(JA)** One Plate 総生菌数 (OP-TVC)
-  **(ZH)** One Plate 总活菌数 (OP-TVC)
-  **(KO)** One Plate 총 생균수 (OP-TVC)

Product Instructions

One Plate Total Viable Count (OP-TVC)

Intended Use

One Plate Total Viable Count (OP-TVC) offers a rapid method for the enumeration of total aerobic mesophiles using traditional culture methodology.

Product Summary and Explanation

Total Viable Count (TVC) is the test method for the enumeration of microorganisms that are able to grow and form colonies on/in a solid medium after aerobic incubation at 30 °C. The test is used to determine the microbial bioburden of a broad range of foods.

One Plate Total Viable Count consists of a nutrient agar base, which provides all the required carbon, nitrogen, and growth factors necessary to support the growth of all heterotrophic microorganisms. An additional diagnostic supplement is added post-sterilisation to indicate and colour colony-forming units, expressing metabolic activity. This colouration improves the ability to identify and enumerate colonies on or in the agar.

The product can permit a quantitative result for total viable count 36–48 hours for all food products, using only one plate versus two or more as described in ISO 7218, ISO 4833 (both parts).

Intended User

The method is designed for use by qualified personnel with appropriate training.

Product Codes

NOTE: Please refer to ISO 6887 parts 1–5 for suitable sample specific diluents.

Product Name	Format	Pack Size	SKU	
One Plate TVC	DCM	500 g	700006967	NCM1018A
		5 kg	700006968	NCM1018B
		10 kg	700006969	NCM1018C
		25 kg	700006970	NCM1018D
One Plate TVC Diagnostic Supplement	Liquid-Ready Supplement	1 x 100 mL bottle	700006971	NCM4094-200
One Plate TVC	Prepared Agar Plates*	20 x 90 mm Agar Plates	On Request	
		100 x 90 mm Agar Plates	On Request	
One Plate TVC	Ready to Re-melt Agar**	1 x 250 mL bottle	On Request	

*Prepared Agar Plates do not require supplementation.

**Ready to Re-melt Agar should be melted in a steaming water bath (refer to ISO 7218 for further guidance). Media should be supplemented after tempering.

Typical Composition

Formulation may be adjusted and/or supplemented as required to meet performance specifications.

Peptone	8.5 g/L
Growth Mix	4.65 g/L
Agar	10.0 g/L


pH of the prepared media at 25 °C: pH 6.3 ± 0.1.

Preparation

Suspend 23.15 g of the medium in 1 litre of purified water. Agitate frequently to completely dissolve the medium. Autoclave at 121 °C for 15 minutes. Cool to 45–50 °C. Aseptically add 0.5 mL of the One Plate TVC diagnostic supplement (NCM4094-200), swirl to mix, and pour into Petri dishes.

NOTE: The One Plate TVC diagnostic supplement does not contribute any growth requirements. If the supplement is not added, colonies will present as their natural straw/buff colouration.

OP-TVC Flow Diagram

Sample Preparation

Follow the specifications of ISO 6887 or the specific international standard appropriate to the product concerned for the initial suspension and dilutions.

Surface Inoculation

1. Transfer 0.1 mL of the suspension, or any serial dilutions, onto the surface of **ONE** single plate of prepared or pre-poured OP-TVC agar plate.
2. Spread the inoculum on the surface with the aid of a sterile spreader.

NOTES: Automatic spiral plating deposition techniques can be used for surface inoculation. Decimal dilutions could be required if typical contamination range is unknown. To estimate small numbers, spread 1 mL of the primary dilution over the surface of 3 prepared plates as described in ISO 7218.

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
1. Transfer 1 mL of the sample if liquid, or 1 mL of the suspension in the case of other products, or any serial dilutions, into **ONE** empty, sterile Petri dish.
2. Pour approximately 15–20 mL of molten OP-TVC into the plate. Homogenize well by swirling, and let solidify on a cool surface.

NOTE: Overlays are not mandatory but can help prevent spreading growth on the surface of the agar. In this case, add 4-6 mL of molten agar on top of the set agar and allow to solidify before incubation as described in ISO 7218.

Incubation

Invert the prepared plate and incubate at 30 ± 1 °C for 36–48 hours.

Interpretation

Microbial growth will present as red/pink colonies against a clear background. Some microorganisms may present as weak or no reaction, but all colony-forming units should be interpreted as viable organisms regardless of colour.

Counting Colony-forming Units

Count all colonies regardless of colour, and take care to avoid mistaking particles of matrix debris for colonies. Examine doubtful objects carefully, using higher magnification where required.

Select the plate containing 10–300 colonies for accurate enumeration. Refer to the ISO 7218 when the results are outside the limits. Apply the dilution factors, which eliminates the need for two successive dilutions or duplicates.

Confirmation is not required for expression of the total viable count.

Dilute to Specifications Results

Dilutions are made according to the specification, and a calculated amount of the appropriate dilution is added to the plate.

Examples of specification dilutions:

Dilution	Inoculation	Specifications
1/10	1 mL pour plate	Negative result <10 CFU/mL Positive result ≥ 10 CFU/mL
1/10	0.1 mL surface or pour plate	Negative result <100 CFU/mL Positive result ≥ 100 CFU/mL
1/10	0.01 mL surface (loop) streak plate	Negative result <1000 CFU/mL Positive result ≥ 1000 CFU/mL
1/100	1 mL pour plate	Negative result <100 CFU/mL Positive result ≥ 100 CFU/mL
1/100	0.1 mL surface or pour plate	Negative result <1000 CFU/mL Positive result ≥ 1000 CFU/mL
1/1000	1 mL pour plate	Negative result <1000 CFU/mL Positive result ≥ 1000 CFU/mL

Quality Control

Appearance of dehydrated media: beige, homogeneous, free-flowing powder

Appearance of prepared media (after supplement addition): beige, clear to hazy

Typical cultural response when incubated aerobically at 30 ± 1°C and examined for growth at 48 hours:

Microorganism	WDCM	Approx. Inoculum (CFU)	Expected Results
<i>Bacillus subtilis subsp. spizizenii</i>	00003	50-200	>70 % Recovery, Red/Pink
<i>Escherichia coli</i>	00013	50-200	>70 % Recovery, Red/Pink
<i>Escherichia coli</i>	00012	50-200	>70 % Recovery, Red/Pink
<i>Staphylococcus aureus</i>	00034	50-200	>70 % Recovery, Red/Pink

Precautions and Limitations of the Method

1. Use good microbiology laboratory practices as per ISO 7218, except mandatory duplicate plating for enumeration, as the One Plate method requires only one Petri dish.
2. Some strains may present a weak or absent reaction from the diagnostic supplement. As all colonies are considered target and the supplement is intended to improve counting ability, all colonies should be counted regardless of colour.
3. Overlays are not mandatory (not performed during validation) but can be helpful to prevent spreading growth on the surface of the agar. In this case add 4-6 mL of molten agar on top of the set agar and allow to solidify before incubation as described in ISO 7218.
4. Colony counts in excess of 300 CFU per plate can lead to poorer definition of the diagnostic reactions; it is advised that a higher dilution be used to accurately enumerate higher levels of contamination.
5. Both the limit of detection and quantification are dictated by the subculture volume.
6. An inoculation of 1.0 mL (pour plate) of a 1/10 dilution has a limit of detection of 10 CFU/g and permits the user to quantify (enumerate) > 100 CFU/g.
7. If a plate count results in less than 10 colonies counted, an estimated number should be expressed, e.g., report >10 CFU/g or >100 CFU/g (depending on the inoculation volume).

Verification

Method verification should be done following the protocols described in ISO 16140:3:2021. The user laboratory should follow experimental design and acceptance criteria as described in chapter 6, 'Quantitative methods — Technical protocol for verification'.

Validations

MicroVal (ISO 16140-2:2016)

One Plate TVC has been certified by MicroVal as an alternative to the reference ISO 4833-1:2013 Horizontal method for the enumeration of microorganisms — Part 1: Colony count at 30 °C by the pour plate technique, according to the reference standard ISO 16140-2:2016, with the scope of enumeration of microorganism in a broad range of foods. Refer to the MicroVal certificate for more information.

Safety

Refer to the relevant product safety data sheet (SDS).

Storage

Store dehydrated culture media at 2–30 °C away from direct sunlight. Once it is opened and recapped, place the container in a low-humidity environment at the same storage temperature. Protect from moisture and light by keeping container tightly closed.

Prepared media and/or supplements should be stored at 2–8 °C. Do not freeze. Supplement may be stored at room temperature for up to 5 days.

Expiration

Refer to expiration date stamped on the container. The dehydrated medium should be discarded if not free-flowing, or if appearance has changed from the original color. Expiry applies to medium in its intact container.

Disposal

Cultures should be disposed of appropriately as biohazard waste. The preferred method of treatment for biohazard waste is autoclaving. Items that cannot be autoclaved may be disinfected with bleach solution. Consult with the safety advisor for your

facility for detailed instructions.

Customer Service

Neogen Customer Services and Technical Services can be reached by using the following contact information:

The Dairy School, Auchincruive, Ayr, KA6 5HU, Scotland UK

+44 (0) 1292 525600

infoUK@neogen.com

Training on this product, and all Neogen test kits, is available.

Terms and Conditions

For full terms and conditions, please visit <https://www.neogen.com/terms-and-conditions>

Warranty

Neogen Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, Neogen will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. Neogen shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

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