*Read instructions carefully before starting test*

### **INTENDED USE**

The ANSR<sup>®</sup> for *Salmonella* method provides for rapid and accurate detection of *Salmonella* spp. in a wide variety of foods and environmental samples. Serovars belonging to both *Salmonella enterica* and *Salmonella bongori* can be detected with the ANSR method.

In an AOAC Research Institute Performance Tested<sup>SM</sup> Method study, ANSR for *Salmonella* has been found to be an effective procedure for detection of *Salmonella* spp. in the following sample types: chicken carcass rinse, raw ground turkey, raw ground beef (10% fat), raw shrimp, hot dogs, oat cereal, soy flour, peanut butter, raw almonds, cocoa powder, ice cream, black pepper, pasteurized dried whole egg, pasteurized dried egg white, pasteurized liquid whole egg, pasteurized frozen egg yolk, raw spinach, dry dog food and sponge or swab samples from stainless steel, plastic, sealed concrete, ceramic tile and rubber surfaces. Inclusivity was 99.1% in testing of more than 100 serovars of *S. enterica* and *S. bongori*, with only a single strain of *S. Welasco* testing negative.

The test is approved by the National Poultry Improvement Plan (NPIP) for testing of poultry house environment samples including boot swabs, drag swabs, chick papers, and fluff.

### **ASSAY PRINCIPLES**

ANSR for *Salmonella* is an isothermal, amplified nucleic acid assay. The ANSR method is based on nicking enzyme amplification reaction (NEAR<sup>™</sup>) technology. Target nucleic acid is amplified through a mechanism of polymerization from the ends of nicks created in double-stranded DNA by the action of a specific endonuclease. Amplified target sequences are detected in real time using fluorescent molecular beacon probes.

A 2-stage lysis is performed, first at  $37 \pm 2^\circ\text{C}$  for 10 minutes, then at  $80 \pm 2^\circ\text{C}$  for 20 minutes. Next, a portion of the lysed sample is transferred to a strip tube containing lyophilized ANSR reagents. The tubes are sealed and incubated at  $56 \pm 2^\circ\text{C}$  on the ANSR reader. Results are generated by the reader and displayed in the ANSR software within 10 minutes. Positive results may be confirmed from the enrichment cultures following standard procedures. Each tube of ANSR reagents contains an internal positive control, ensuring the reagents are functioning properly.

### **INTENDED USER**

The ANSR for *Salmonella* test is designed for use by personnel with appropriate training in microbiology. Training in the use of the ANSR test system is available through NEOGEN.

### **MATERIALS PROVIDED**

1. 12 strips of 8 cluster tubes in a rack, 1.2 mL
2. 12 strips of 8 reaction tubes, 200  $\mu\text{L}$ , containing lyophilized ANSR for *Salmonella* reagents in 2 sealed foil pouches with desiccant pack
3. 12 strips of 8 permanent caps for reaction tubes
4. 1 bottle of lysis reagent suspension buffer, 60 mL
5. 3 vials containing lyophilized lysis reagents
6. 1 kit insert

### **EQUIPMENT REQUIRED**

1. ANSR reader (NEOGEN item 9828)
2. Computer and software for connection to ANSR reader (NEOGEN item 9832)
3. 1 dual heater block with aluminum block inserts for 1.2 mL cluster tubes,  $80 \pm 2^\circ\text{C}$  and  $37 \pm 2^\circ\text{C}$  (NEOGEN items 9386-48D, 9829-48) **OR** two single heater blocks with aluminum block inserts for 1.2 mL cluster tubes,  $80 \pm 2^\circ\text{C}$  and  $37 \pm 2^\circ\text{C}$  (NEOGEN items 9386-48D, 9829-48)
4. Pipettor, 100–1000  $\mu\text{L}$  (NEOGEN item 9463)
5. Pipette tip rack, 100–1000  $\mu\text{L}$ , sterile, 96 tips (NEOGEN item 9487)
6. Pipettor, 10–100  $\mu\text{L}$ , 8-channel (NEOGEN item 9388)
7. Pipette tips, 100  $\mu\text{L}$ , sterile, filtered, 96 tips (NEOGEN item 9389)
8. Vortex mixer, adjustable speed (NEOGEN item 9494)
9. Stomacher (optional)
10. 3 thermometers, traceable (NEOGEN item 9518)
11. Timer, 3-channel (NEOGEN item 9426)
12. Optional-for-use heater block with 0.2 mL reaction tube aluminum block insert,  $56 \pm 2^\circ\text{C}$  (NEOGEN items 9386-48D, 9829-48)
13. Webcam (NEOGEN item WEBCAM)
14. ANSR Ethernet cable (NEOGEN item 9835)
15. 10 mL pipette pump (NEOGEN item 9277)
16. Pipettes, sterile serological (NEOGEN item 8686)
17. 40-slot, 20 mm test tube rack, autoclavable (NEOGEN item 9553)
18. Pipettor, 20–200 mL (NEOGEN item 9276)

### **OTHER MATERIALS REQUIRED**

1. Stomacher-type bags for sample enrichment. Filtered bags are recommended (NEOGEN item 6827 or equivalent)
2. Graduated cylinder, 250 mL (NEOGEN item 9368 or equivalent)

## **MEDIA ENRICHMENT BROTH REQUIRED**

1. Novobiocin, sodium salt (NEOGEN item 7985 or equivalent)
2. Rappaport-Vassiliadis Broth (NEOGEN item 7512 or equivalent). For raw meat sample confirmation only.
3. Phosphate Buffered Saline (PBS) (NEOGEN item 8425 or equivalent)
4. Potassium iodide (sigma no. 221945 or equivalent)
5. Iodine (Sigma no. 207772 or equivalent)

## **ENRICHMENT BROTH REQUIRED (DEPENDENT ON SAMPLE TYPE AS DESCRIBED ON PAGES 4 AND 5)**

1. ANSR for *Salmonella* Enrichment Broth 1 (NEOGEN item 9811)
2. ANSR for *Salmonella* Enrichment Broth 2 (NEOGEN item 9812)
3. ANSR for *Salmonella* Enrichment Broth 3 (NEOGEN item 9813)
4. ANSR for *Salmonella* Enrichment Broth 3B (NEOGEN item 9817)
5. Tetrathionate Broth Base, Hajna (TT Hajna) (NEOGEN item 7740A)

## **STORAGE**

Store ANSR reagents at 2–8°C. After removing reaction tubes from the foil pouch, promptly re-seal the pouch. Leave the desiccant pack in the pouch at all times.

## **PRECAUTIONS**

1. Use good microbiology laboratory practices, including biosafety level 2 procedures when handling pathogens such as *Salmonella*.
2. **Dispose of used pipette tips in a covered container containing a fresh solution of 10% bleach. The 10% bleach solution should be made fresh each day. Use 9 parts water with 1 part commercially available household-strength bleach to make the 10% bleach solution. Undiluted stock solutions of bleach should be used within 30 days after opening.**
3. **Discard bleach solution and tips as regular waste at the end of each day.**
4. Do not use reagents beyond the expiration date.
5. Enrichment media should be prewarmed to incubation temperature (all protocols) before use.
6. Use of enrichment media and incubation times or temperatures other than those specified may lead to erroneous results.
7. Remove reaction tubes from the foil pouch just before use and keep covered until heating process begins. Reseal the pouch containing the remaining reaction tubes to **avoid prolonged exposure to light**. More than 15 minutes of total exposure time may lead to erroneous results.
8. Do **NOT**, under any circumstance, remove caps from reaction tubes **AFTER** the assay has been started. This is essential in order to prevent accidental contamination of the environment with amplification products.
9. Exercise care in all pipetting steps to avoid cross-contamination of samples.
10. Complete all assay steps in sequence, avoiding delays between steps.
11. Tap reaction tubes on bench top to make sure lyophilized reagents are at the bottom of the tube prior to adding lysed sample.

**NOTE: IF TESTING USING COLONY PICKS, PROCEED TO PAGE 7 FOR TEST PROCEDURE.**

## PREPARATION OF ENRICHMENT BROTH

**NOTE:** Use rehydrated media the same day as prepared. **Do not autoclave.**

### **ANSR for *Salmonella* Enrichment Broth 1**

Rehydrate 42.1 g of ANSR for *Salmonella* Enrichment Broth 1 with 1 L sterile water prewarmed to  $36 \pm 1^\circ\text{C}$ .

### **ANSR for *Salmonella* Enrichment Broth 2**

Rehydrate 20 g of ANSR for *Salmonella* Enrichment Broth 2 with 1 L sterile water prewarmed to  $42 \pm 1^\circ\text{C}$ .

### **ANSR for *Salmonella* Enrichment Broth 3**

Rehydrate 42.1 g of ANSR for *Salmonella* Enrichment Broth 3 with 1 L sterile water prewarmed to  $36 \pm 1^\circ\text{C}$ .

### **ANSR for *Salmonella* Enrichment Broth 3B**

Rehydrate 39.2 g of ANSR for *Salmonella* Enrichment Broth 3B with 1 L sterile water prewarmed to  $36 \pm 1^\circ\text{C}$  (or  $42 \pm 1^\circ\text{C}$  for pet food).

### **TT Hajna**

Follow manufacturer's instructions for media preparation.

## SAMPLE ENRICHMENT

**For processed meats and other processed foods with a moderate to high microbial load.** Validated for: Hot dogs, raw almonds, peanut butter, black pepper and cocoa powder.

1. Weigh 25 g sample in a Stomacher-type bag.
2. Add 225 mL **ANSR for *Salmonella* Enrichment Broth 1** (NEOGEN item 9811) prewarmed to  $36 \pm 1^\circ\text{C}$  to the bag. Homogenize (Stomacher or equivalent) the sample as appropriate for the sample type.
3. For all products except black pepper, incubate the culture at  $36 \pm 1^\circ\text{C}$  for **16–24 hours**.
4. For black pepper, incubate the culture at  $36 \pm 1^\circ\text{C}$  for **22–26 hours**.

### **For raw meat, poultry carcass rinses, raw leafy greens and flours**

Validated for: raw ground turkey, raw ground beef, chicken carcass rinse, raw frozen shrimp, raw spinach and soy flour.

1. Weigh 25 g sample (or 30 mL poultry carcass rinse) in a Stomacher-type bag.
2. Add 225 mL **ANSR for *Salmonella* Enrichment Broth 2** (NEOGEN item 9812) prewarmed to  $42 \pm 1^\circ\text{C}$  to the bag. For 325 g samples, weigh 325 g sample in a filter-lined Stomacher-type bag, and add 1.3 L ANSR for *Salmonella* Enrichment Broth 2 prewarmed to  $42 \pm 1^\circ\text{C}$  to the bag. Homogenize (Stomacher or equivalent) the sample as appropriate for the sample type.
3. For 25 g ground beef, incubate the culture at  $42 \pm 1^\circ\text{C}$  for **10–24 hours**. For 325 g ground beef and other raw meat samples and poultry carcass rinses, incubate the culture at  $42 \pm 1^\circ\text{C}$  for **12–24 hours**. For spinach and raw frozen shrimp, incubate the culture at  $42 \pm 1^\circ\text{C}$  for **16–24 hours**. For soy flour, incubate the culture at  $36 \pm 1^\circ\text{C}$  for **16–24 hours**.

### For processed foods with a low microbial load

Validated for: oat cereal and ice cream.

1. Weigh 25 g sample in a Stomacher-type bag.
2. Add 225 mL **ANSR for *Salmonella* Enrichment Broth 3** (NEOGEN item 9813) prewarmed to  $36 \pm 1^\circ\text{C}$  to the bag. Homogenize (Stomacher or equivalent) the sample as appropriate for the sample type.
3. Incubate the culture at  $36 \pm 1^\circ\text{C}$  for **16–24 hours**.

### For environmental samples

Validated for: stainless steel, plastic, sealed concrete, ceramic tile and rubber.

1. Place the sponge or swab sample in a Stomacher-type bag or tube.
2. Add an appropriate amount of **ANSR for *Salmonella* Enrichment Broth 3** (NEOGEN item 9813) prewarmed to  $36 \pm 1^\circ\text{C}$ . For sponge samples, an appropriate amount typically is 100–200 mL. For swab samples, it typically is 10 mL. For samples containing high levels of background flora, add novobiocin to a final concentration of 10 mg/L.
3. Incubate the culture at  $36 \pm 1^\circ\text{C}$  for **16–24 hours**.

**NOTE:** Testing of certain foods may require a different sample size (e.g., 100 g, 325 g or 375 g). Generally, for these larger sample sizes, a 4:1 ratio of enrichment medium to sample may be used, unless contraindicated by problems with color, viscosity, etc. Contact NEOGEN for more information.

### For dry pet food

1. For 375 g composite samples, weigh 375 g sample in a filter-lined Stomacher-type bag, and add 3.375 L ANSR for *Salmonella* Enrichment Broth 3B prewarmed to  $42 \pm 1^\circ\text{C}$  to the bag.
2. Homogenize (Stomacher or equivalent) the sample as appropriate for the sample type.
3. Incubate the culture at  $42 \pm 1^\circ\text{C}$  for **16–24 hours**.

### For pasteurized egg products

Validated for: pasteurized dried whole egg, pasteurized dried egg white, pasteurized liquid whole egg and pasteurized frozen egg yolk.

1. Weigh 100 g sample in a Stomacher-type bag.
2. Add 900 mL **ANSR for *Salmonella* Enrichment Broth 3B** (NEOGEN item 9817) prewarmed to  $36 \pm 1^\circ\text{C}$  to the bag.
3. Homogenize (Stomacher or equivalent) the sample as appropriate for the sample type.
4. Incubate the culture at  $36 \pm 1^\circ\text{C}$  for **24–26 hours**. **NOTE:** For pasteurized dried whole egg, incubate for **26 hours**.

### For poultry house environmental samples

**NOTE:** Not part of AOAC Performance Tested Method

1. Place environmental sample into a stomacher-type bag.
2. Aseptically add 100 mL of TT Hajna (prepared earlier) to a Stomacher-type bag containing the environmental sample. Shake bag vigorously in a side-to-side motion at least 10 times to mix.
3. Incubate the sample(s) at  $37 \pm 1^\circ\text{C}$  for  $24 \pm 2$  hours.
4. Before testing in the ANSR assay, dilute a portion of the enrichment culture 1:10 in PBS or a similar buffer.

## LYSIS REAGENT SOLUTION PREPARATION

1. Reconstitute 1 vial of lyophilized lysis reagents with 18 mL of lysis reagent suspension buffer by adding the buffer to the reagent vial. Swirl gently to mix.

**NOTE:** 1 vial of lysis reagents is enough for approximately 32 samples. Prepared lysis reagent solution can be stored at 2–8°C for up to 30 days.

## ANSR TEST PROCEDURE

Prior to starting the assay:

1. Preheat one lysis heater block to  $80 \pm 2^\circ\text{C}$ . Preheat the second lysis heater block to  $37 \pm 2^\circ\text{C}$ . If using the optional single heater, preheat to  $56 \pm 2^\circ\text{C}$ . Use the thermometer for the temperature reading.
2. Remove the foil pouch containing the reaction tubes from the refrigerator and **allow them to warm at room temperature for 15 minutes**. To avoid excess light exposure, leave reaction tubes in foil pouch until they are needed.

**NOTE:** The lysis reagent solution should be set out for 15 minutes prior to use to get to room temperature.

3. Connect the ANSR reader to the computer via USB or Ethernet and turn the computer on.
4. Turn on the ANSR reader. The reader will preheat to  $56 \pm 2^\circ\text{C}$ .
5. Start the ANSR software and click the Connect button. Input sample IDs, lot number and user information.

**NOTE:** For instructions on using the reader and software, see the user guide that came with the ANSR reader.

## ASSAY PROCEDURE

**NOTE:** Enrichment cultures that are highly turbid or highly colored may require a 1:10 dilution in PBS prior to lysing the sample. Examples include peanut butter, black pepper, flours, dry pet food and dried pasteurized egg products. Cocoa powder requires a 1:50 dilution in PBS prior to sample lysis. Contact NEOGEN Technical Services for more information.

**OPTIONAL:** Pre-pierced plate sealing film (BIO RAD Item 3600040) can be placed over the 1.2 mL cluster tubes if desired.

1. Add 50  $\mu\text{L}$  enrichment culture (or dilution) to a 1.2 mL cluster tube(s) using 100  $\mu\text{L}$  filtered tips. **NOTE:** Cluster tubes may be pulled apart to provide the number of tubes needed.
2. Add 450  $\mu\text{L}$  lysis reagent solution to each cluster tube(s) containing culture. Be sure to switch pipette tips between samples. **NOTE:** Return the lysis reagent solution to the refrigerator after use (within 1 hour).
3. Incubate the cluster tube(s) at  $37 \pm 2^\circ\text{C}$  for **10 minutes**.
4. Immediately transfer the cluster tube(s) to the  $80 \pm 2^\circ\text{C}$  heater block and incubate for **20 minutes**. **NOTE:** The  $80 \pm 2^\circ\text{C}$  incubation time may be extended to a total of **60 minutes** for the purpose of managing staggered assay start times.
5. **For 3 to 5 minutes** before the end of the lysis step, preheat the ANSR reagents to  $56 \pm 2^\circ\text{C}$  by placing the reaction tubes in the ANSR reader. Optional: A separate heat block can be used. It should be heated to  $56 \pm 2^\circ\text{C}$ . **NOTE:** The strip of reaction tubes may be cut to provide the number of tubes needed. Keep all unused tubes in the sealed foil pack. Ensure the pellet in the reaction tube(s) is at the bottom by tapping the tubes gently on the bench top.
6. After the completion of the 20–60 minute lysis incubation, remove and discard the caps from the reaction tube(s) in the ANSR reader. **IMPORTANT:** Proceed with steps 7–9 without delay. **The transfer of the sample from the lysis tubes to the reaction tubes should be completed within 1 minute.**
7. Using an 8-channel pipette and 100  $\mu\text{L}$  filter tips, carefully transfer 50  $\mu\text{L}$  from the top third of the lysed sample(s) in the cluster tube(s) to the reaction tube(s). Debris may accumulate at the bottom of the lysis tube(s) that will interfere with assay performance. Avoid transfer of debris by aspirating from the top third of the lysis tube(s). Do not prime the pipette tips and do not mix before aspirating. **Place the provided permanent cap(s) on the reaction tube(s).** **NOTE:** Lysed sample may be transferred from the same cluster tube a maximum of 3 times.
8. Remove the strips(s) of tubes from the reader (or  $56 \pm 2^\circ\text{C}$  heat block if one was used) and vortex briefly (about 2 seconds) at medium/high speed (2500 RPM). Quickly visually check each reaction tube to assure that no bubbles are on the bottom or in the middle of the sample. A quick tap of the tubes should release any bubbles from the bottom or middle to the top. Then place into the reader without delay. Close the reader's lid. **NOTE:** The reader will not provide accurate results if the lid is open. Keep the lid closed at all times while the assay is running. Contamination may occur if the permanent caps are not placed on the reaction tubes and/or if the permanent caps are removed.
9. Click **START** in the ANSR software to begin the **10 minute** assay.
10. Results will be displayed as positive, negative or invalid once the assay is finished. If the result displays an invalid, the test must be repeated. A single repeat test can be run starting at step 6. Alternatively, start over from step 1 if lysis from step 6 has been at  $80 \pm 2^\circ\text{C}$  for more than 60 minutes.

## PROCEED TO PAGE 8 FOR INTERPRETATION OF RESULTS.

### CONFIRMATION

Samples producing positive ANSR results may be confirmed by streaking the enrichment cultures to *Salmonella* selective/differential agar media and continuing with identification of presumptive colonies following standard methods<sup>[1,2]</sup>.

#### For raw meat and other samples with a high microbial load

1. Transfer 0.1 mL enrichment culture to 10 mL Rappaport-Vassiliadis (RV) Broth (NEOGEN item 7512) and incubate for **22–26 hours** at  $42 \pm 0.5^\circ\text{C}$ .
2. Streak RV culture to selective/differential agar plates.

#### FOR USE AS A CONFIRMATION TEST FROM COLONY PICKS (AOAC Official Method 2013.14)

ANSR for *Salmonella* can be used to determine if a presumptive colony on selective/differential or non-selective agar media is *Salmonella* spp. Colonies picked from the following media types have been tested and found to produce accurate results: Hektoen enteric agar, xylose lysine deoxycholate agar, bismuth sulfite agar, brilliant green sulfa agar, xylose lysine tergitol agar (XLT4), double-modified lysine iron agar and tryptic soy agar.

#### Procedure

Pick a colony with an inoculating loop or needle and resuspend the colony in 0.5 mL sterile PBS. Vortex or otherwise mix to ensure the colony is completely resuspended. Proceed with the **ANSR Test Procedure**, followed by the **Colony Picks Assay Procedure** using 50  $\mu\text{L}$  of the PBS resuspension as sample. A positive test result indicates *Salmonella* spp. A negative test result indicates “not *Salmonella* spp.”

---

#### References

1. US Food and Drug Administration (May 2014). Chapter 5. *Salmonella* in *Bacteriological Analytical Manual*: Retrieved from <http://www.fda.gov/food/foodscienceresearch/laboratorymethods/ucm070149.htm>
2. United States Department of Agriculture Food Safety Inspection Service. (January, 2011). Isolation and identification of *Salmonella* from meat, poultry, pasteurized egg and catfish products in *Microbiology Laboratory Guidebook* (4,05): Retrieved from <http://www.fsis.usda.gov/wps/wcm/connect/7f00c05fe-06a2-492a-a6e1-3357f7701f52/MLG-4.pdf?MOD=AJPERES>

## COLONY PICKS ASSAY PROCEDURE

**OPTIONAL:** Pre-pierced plate sealing film (BIO RAD Item 3600040) can be placed over the 1.2mL cluster tubes if desired.

1. Add 50  $\mu$ L colony resuspension to a 1.2 mL cluster tube(s) using 100  $\mu$ L filtered tips.

**NOTE:** Cluster tubes may be pulled apart to provide the number of tubes needed.

2. Add 450  $\mu$ L lysis reagent buffer to each of the cluster tube(s) containing sample. **NOTE:** The colony assay procedure only has one lysis step at  $80 \pm 2^\circ\text{C}$ . There is no  $37 \pm 2^\circ\text{C}$  step needed.

3. Immediately transfer the cluster tube(s) to the  $80 \pm 2^\circ\text{C}$  heater block and incubate for **20 minutes**.

**NOTE:** The  $80 \pm 2^\circ\text{C}$  incubation time may be extended to a total of **60 minutes** for the purpose of managing staggered assay start times.

4. **For 3 to 5 minutes** before the end of the lysis step, preheat the ANSR reagents to  $56 \pm 2^\circ\text{C}$  by placing the reaction tubes in the ANSR reader. Optional: A separate heat block can be used. It should be heated to  $56 \pm 2^\circ\text{C}$ .

**NOTE:** The strip of reaction tubes may be cut to provide the number of tubes needed. Keep all unused tubes in the sealed foil pack. Ensure the pellet in the reaction tube(s) is at the bottom by tapping the tubes gently on the bench top.

5. After the completion of the 20–60 minute lysis incubation, remove and discard the caps from the reaction tube(s) in the ANSR reader.

**IMPORTANT:** Proceed with steps 7–9 without delay. **The transfer of the sample from the lysis tubes to the reaction tubes should be completed within 1 minute.**

6. Using an 8-channel pipette and 100  $\mu$ L filtered tips, carefully transfer 50  $\mu$ L from the top third of the lysed sample(s) in the cluster tube(s) to the reaction tube(s). Debris may accumulate at the bottom of the lysis tube(s) that will interfere with assay performance. Avoid transfer of debris by aspirating from the top third of the lysis tube(s). Do not prime the pipette tips and do not mix before aspirating.

**Place the provided permanent cap(s) on the reaction tube(s).** **NOTE:** Lysed sample may be transferred from the same cluster tube a maximum of 3 times.

7. Remove the strips(s) of tubes from the reader (or  $56 \pm 2^\circ\text{C}$  heat block if one was used) and vortex briefly (about 2 seconds) at medium/high speed (2500 RPM). Quickly visually check each reaction tube to assure that no bubbles are on the bottom or in the middle of the sample. A quick tap of the tubes should release any bubbles from the bottom or middle to the top. Then place into the reader without delay. Close the reader's lid.

**NOTE:** The reader will not provide accurate results if the lid is open. Keep the lid closed at all times while the assay is running. Contamination may occur if the permanent caps are not placed on the reaction tubes and/or if the permanent caps are removed.

8. Click **START** in the ANSR software to begin the **10 minute** assay.

9. Results will be displayed as positive, negative or invalid once the assay is finished. If the result displays an invalid, the test must be repeated. A single repeat test can be run starting at step 5. Alternatively, start over from step 1 if lysis has been at  $80 \pm 2^\circ\text{C}$  for more than 60 minutes.

## **INTERPRETATION OF RESULTS**

The ANSR software will indicate the test results as positive or negative for the presence of *Salmonella* spp. in the enriched sample. In addition, the real-time fluorescence curve generated from the assay can be viewed.

## **DISPOSAL**

Enrichment cultures and used lysis tubes should be disposed of as biohazard waste. The preferred method of treatment for biohazard waste is autoclaving. Items that **CANNOT** be autoclaved, such as used reaction tubes and caps, should be immersed in a fresh 10% bleach solution that is made daily. Consult with the safety advisor of your facility for detailed instructions.

Do **NOT** remove the permanent caps, for **ANY** reason, from the ANSR reaction tubes once the assay has started, even when disposing of them. Reaction tubes can be disposed of as non-biohazardous waste. It is recommended that they be placed in a sealable container or plastic bag with 10% bleach and immediately disposed of to protect against accidental opening.

## **NeoSeek™**

NEOGEN's NeoSeek program provides pathogen detection and identification technology that provides next day, DNA-specific test results for certain pathogenic bacterial strains. Contact NEOGEN for more information.

## **CUSTOMER SERVICE**

NEOGEN Customer Assistance and Technical Services can be reached by using the contact information on the back of this booklet. Training on this product, and all NEOGEN test kits, is available.

## **SDS INFORMATION AVAILABLE**

Safety data sheets (SDS) are available for this test kit, and all of NEOGEN's test kits, on NEOGEN's website at [neogen.com](http://neogen.com), or by calling NEOGEN at 800.234.5333 or 517.372.9200.

## **TERMS AND CONDITIONS**

For NEOGEN's full terms and conditions, please visit [neogen.com/terms-and-conditions/](http://neogen.com/terms-and-conditions/)

## **WARRANTY**

NEOGEN Corporation makes no warranty of any kind, either expressed or implied, except that the materials from which its products are made are of standard quality. If any materials are defective, NEOGEN will provide a replacement of the product. Buyer assumes all risk and liability resulting from the use of this product. There is no warranty of merchantability of this product or of the fitness of the product for any purpose. NEOGEN shall not be liable for any damages, including special or consequential damage, or expense arising directly or indirectly from the use of this product.

**Neogen ANSR Molecular Diagnostic for Foodborne Pathogen Detection – limited use label license**

**SYTO® 82**

SYTO® 82 contained within this product is provided under an intellectual property license from Life Technologies Corporation, Eugene, OR, and may be used for *in vitro* detection and analysis of (i) food, feeds and beverages, including nutraceuticals, (ii) ingredients for food, feeds and beverages, (iii) process samples from food, feed and beverage preparation, distribution and delivery, and (iv) water from any source for human consumption, all for the purpose of safety and quality assurance. The buyer must not sell or otherwise transfer this product or its components for any other use, including but not limited to: human *in vitro*, veterinary, identity or paternity testing, forensics, or *in vivo* detection of nucleic acid sequences in living beings, or cells. For information on purchasing a license for SYTO 82 for purposes other than food, beverage and water safety, and quality assurance, contact Life Technologies Corporation at [outlicensing@lifetech.com](mailto:outlicensing@lifetech.com).

**Molecular beacon probes**

One or more molecular beacon probes contained within this product is sold under license from PHRI Properties and may be used under PHRI Properties patent rights only for tests on food products, feeds, beverages, and water.

**NEAR™ Technology**

This product utilizes the patented NEAR isothermal technology and is sold under license from Ionian Technologies, San Diego, CA, and may be used under Ionian Technologies patent rights only for tests on food, beverage, and water safety.

## PRODUCTS AVAILABLE

- 9465 **Environmental sampling sponges** (with Dey Engley (DE) neutralizing broth) – pack of 100
- 36002 **Sampling sponges, pre-moistened, HiCap** – pack of 100
- 9822 **ANSR for *E. coli* O157:H7** – up to 96 tests
- 9824 **ANSR for *Listeria monocytogenes*** – up to 96 tests
- 9870 **ANSR for *Salmonella*** – up to 96 tests
- 9871 **ANSR for *Listeria*** – up to 96 tests
- 9872 **ANSR for *Campylobacter*** – up to 96 tests
- 9873 **ANSR for *Listeria* Right Now™** – up to 96 tests
- 9837 **Complete ANSR Equipment Kit** – Reader, computer, pipette set, dry bath heaters and inserts, thermometers, vortex, time, rack

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## **Referencias**

1. US Food and Drug Administration (November, 2011). Chapter 5. *Salmonella* in *Bacteriological Analytical Manual*: Retrieved from <http://www.fda.gov/food/foodscienceresearch/laboratorymethods/ucm070149.htm>
2. United States Department of Agriculture Food Safety Inspection Service. (January, 2011). Isolation and identification of *Salmonella* from meat, poultry, pasteurized egg and catfish products in *Microbiology Laboratory Guidebook* (4.05): Retrieved from [http://www.fsis.usda.gov/PDF/MLG\\_4\\_05.pdf](http://www.fsis.usda.gov/PDF/MLG_4_05.pdf)

#### **Norteamérica**

##### **Oficinas Corporativas de NEOGEN**

+1 800.234.5333 (EEUU/Canadá)  
foodsafety@neogen.com  
NEOGEN.com

#### **Europa, Medio Oriente y África**

##### **NEOGEN Europe**

+ 44 (0) 1292 525 600  
info\_uk@neogeneurope.com  
NEOGEN.com

#### **México**

##### **NEOGEN Latinoamérica**

+52 (55) 5254-8235  
informacion@neogenlac.com  
NEOGEN.com

#### **Brasil**

##### **NEOGEN do Brasil**

+55 19 3935.3727  
info@neogendobrasil.com.br  
NEOGEN.com

#### **China**

##### **NEOGEN Bio-Scientific Technology**

+86 21 6271 7013  
info@neogenchina.com.cn  
www.neogenchina.com.cn

#### **India**

##### **NEOGEN Food and Animal Security**

+91 484 2306598, 2301582  
info@neogenindia.com  
www.neogenindia.com