Product Instructions
BPPFV9BPW/BPPFV10BPW/BPPFV9LB/BPPFV10LB
BPPFV9MRS/BPPFV10NB/BPPFV9MRD/BPPFV10DE/BPPFV9BFD

-  **(EN) Mini Flip-Top Vials**
-  **(FR) Mini Flip-Tops**
-  **(DE) Mini Flip-Top Verdünnungsröhrchen**
-  **(ES) Minifrascos Flip-Top**
-  **(PT) Mini Flip-Top**
-  **(JA) ミニフリップトップ**
-  **(ZH) Mini Flip-Top**
-  **(TH) ขวดบรรจุสารละลายบัฟเฟอร์**
-  **(KO) 미니 Flip-Top 바이얼**

Product Instructions

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PRODUCT DESCRIPTION AND INTENDED USE

Neogen® Mini Flip-Top Vials come pre-filled with a variety of buffers that may be used as diluents in the enumeration of microorganisms from food and food process environmental samples. The Neogen Mini Flip-Top Vials offer simplicity in preparing dilutions. Some of the buffers can be used to hydrate the Neogen® Sponge-Stick and the Neogen® Dry-Sponge.

Table 1: Product Description and Storage Conditions

Catalog Number	Product Description	Volume	Storage Symbol
BPPFV9BPW	Neogen Mini Flip-Top Vial with Buffered Peptone Water Broth	9 mL	2-30°C
BPPFV10BPW	Neogen Mini Flip-Top Vial with Buffered Peptone Water Broth	10 mL	20-25°C
BPPFV9LB	Neogen Mini Flip-Top Vial with Letheen Broth	9 mL	20-25°C
BPPFV10LB	Neogen Mini Flip-Top Vial with Letheen Broth	10 mL	20-25°C
BPPFV9MRD	Neogen Mini Flip-Top Vial with Maximum Recovery Diluent	9 mL	20-25°C
BPPFV10DE	Neogen Mini Flip-Top with D/E Broth	10 mL	2-8 °C
BPPFV9BFD	Neogen Mini Flip-Top Vial with Butterfield's Buffer	9 mL	20-25°C
BPPFV10NB	Neogen Mini Flip-Top Vial with Neutralizing Buffer	10 mL	20-25°C

SAFETY

The user should read, understand, and follow all safety information in the instructions for the Neogen Mini Flip-Top Vials. Retain the safety instructions for future reference.

⚠ CAUTION Indicates a hazardous situation, which, if not avoided, could result in minor or moderate injury and/or property damage.

⚠ CAUTION

- Neogen has not documented Neogen Mini Flip-Top Vials diluents for use in industries other than food. For example, Neogen has not documented Neogen Mini Flip-Top Vials for testing water, pharmaceuticals, or cosmetics.
- Neogen Mini Flip-Top Vials have not been tested with all possible food products, food processes, testing protocols or with all possible strains of bacteria.
- Do not use Neogen Mini Flip-Top Vials in the diagnosis of conditions in humans or animals.
- The user must train its personnel in proper testing techniques: for example, Good Laboratory Practices⁽¹⁾, ISO 7218⁽²⁾, or ISO 17025⁽³⁾.
- For information on documentation of product performance, visit our website at www.neogen.com or contact your local Neogen representative or distributor.
- Consult the Safety Data Sheet for additional information.

USER RESPONSIBILITY

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at www.neogen.com, or contact your local Neogen representative or distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results. The food sample itself may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any Neogen Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

LIMITATION OF WARRANTIES / LIMITED REMEDY

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, NEOGEN DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any Neogen Food Safety Product is defective, Neogen or its authorized distributor will, at its option, replace or refund the purchase price of the product. These are your exclusive remedies. Please contact your Neogen representative or authorized Neogen distributor for any further questions.

LIMITATION OF NEOGEN LIABILITY

NEOGEN WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall Neogen's liability under any legal theory exceed the purchase price of the product alleged to be defective.

STORAGE

Upon receipt, store unopened Neogen Mini Flip-Top Vials at the indicated temperatures (See Table 1).

DISPOSAL

After use, Neogen Mini Flip-Top Vials may contain microorganisms that may be a potential biohazard. Follow current local, state and national regulations for disposal.

INSTRUCTIONS FOR USE

Follow all instructions carefully. Failure to do so may lead to inaccurate results.

Wear appropriate protective apparel and follow standard good laboratory practices GLP¹.

Performing a dilution:

1. Open the Neogen Mini Flip Top Vial cap.
2. Aseptically add sample to the vial.
3. Re-seal hinged cap onto vial to close. Firmly press cap to secure.
4. Shake vial to complete the dilution.

Hydrating a sponge or swab:

1. Open Neogen Mini Flip Top Vial cap.
2. Add solution to hydrate a dry sponge or a dry swab.
3. After hydrating, aseptically remove the sponge or swab and use accordingly.

If used with Neogen® Petrifilm® Plates

Refer to the appropriate Neogen Petrifilm Plate "Product Instructions."

If you have questions about specific applications or procedures, please visit our website at www.neogen.com or contact your local Neogen representative or distributor.

REFERENCES

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Practice Studies.
2. ISO 7218. Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
3. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.

Refer to the current versions of the standard methods listed above.

EXPLANATION OF SYMBOLS

info.neogen.com/symbols

Neogen Food Safety

Neogen Corporation

620 Leshar Place
Lansing, MI 48912 USA
neogen.com

Neogen Europe Ltd.

The Dairy School
Auchincruive
Ayr, KA6 5HU
Scotland, UK

Neogen Ireland, Ltd.

Bray Business Park,
Bray Co. Wicklow
A98YV29, Ireland

Neogen Corporation

620 Leshar Place Lansing, MI 48912 USA
www.neogen.com

© Neogen Corporation 2025. All rights reserved.
Neogen is a trademark of Neogen Corporation.
FS00846C

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RÉFÉRENCES

1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Practice Studies.
2. ISO 7218: Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
3. ISO/CEI 17025. General requirements for the competence of testing and calibration laboratories.

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1. U.S. Food and Drug Administration. Code of Federal Regulations, Title 21, Part 58. Good Laboratory Practice for Nonclinical Laboratory Practice Studies.
2. ISO 7218 Microbiology of food and animal feeding stuffs – General requirements and guidance for microbiological examinations.
3. ISO/IEC 17025. General requirements for the competence of testing and calibration laboratories.

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