- ① EN Quantitative Rapid Enrichment Dehydrated Media
- ① FR Milieux déshydratés à enrichissement rapide pour la quantification
- ① DE Quantitative Schnellanreicherung Dehydrierte Medien
- ① ES Medios deshidratados para enriquecimiento rápido en pruebas cuantitativas
- ① PT Meios desidratados de enriquecimento rápido quantitativo
- ① JA 定量的迅速濃縮乾燥培地
- ① ZH 定量快速富集脱水培养基
- ① TH อาหารเลี้ยงเชื้อแบบดื่มน้ำออกเพื่อการเพิ่มจำนวนเชื้ออย่างรวดเร็วเชิงปริมาณ
- ① KO 정량적 신속 증균 탈수 배지

Product Instructions

Quantitative Rapid Enrichment Dehydrated Media

Intended Use

This product is intended to be used with the Neogen® Molecular Detection Assay 2Q - Quantitative *Salmonella* product for quantitation of *Salmonella* in select matrices in a laboratory environment by professionals trained in laboratory techniques.

Neogen has not documented the use of this product in industries other than food or beverage. For example, Neogen has not documented this product for testing pharmaceutical, cosmetics, clinical, or veterinary samples.

Description

This medium is designed to help sub-lethally damaged *Salmonella* recover rapidly for detection. This well-buffered, nutritive broth can be used with the matrices outlined in the Neogen MDA2Q - Quantitative *Salmonella* product instructions.

User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at <https://www.neogen.com/document-search/> or contact your local Neogen representative or authorized distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, laboratory technique, and the sample itself may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any Neogen Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

WARNING

Warnings

Do not use this product in the diagnosis of conditions in humans or animals.

The user must train its personnel in current proper testing techniques: for example, Good Laboratory Practices, ISO/IEC 17025⁽⁴⁾, or ISO 7218⁽⁵⁾.

To reduce the risks associated with a false-negative result leading to the release of contaminated product:

- Follow the protocol and perform the tests exactly as stated in the product instructions.
- Store this product as indicated on the package and in the product instructions.
- Always use this product by the expiration date.
- Use this product with foods that have been validated internally or by a third party.
- A negative result at 6 hours using the Neogen Molecular Detection Assay 2Q - Quantitative *Salmonella* assay does not indicate the complete absence of *Salmonella*, only that the level of *Salmonella*, if present, is below the LOD for this method.
- Do not use the sample enrichments outlined in these product instructions in the Neogen® Molecular Detection Assay 2 – *Salmonella* (MDA2SAL96) QUALITATIVE assay. The enrichments included in these product instructions have been developed specifically for this QUANTITATIVE assay.

**To reduce the risks associated with exposure to chemicals and biohazards:**

- Perform pathogen testing in a properly equipped laboratory under the control of trained personnel. Incubated enrichment media and equipment or surfaces that have come into contact with incubated enrichment media may contain pathogens at levels sufficient to cause risk to human health.
- This procedure uses/detects pathogenic microorganisms and/or their metabolic products. Care should be taken to avoid ingestion or inhalation of potentially infectious aerosols or contact with the skin. Always follow standard laboratory safety practices, including wearing appropriate protective apparel and eye protection while handling reagents and contaminated samples.
- Avoid contact with the contents of the sample plus enrichment media after incubation.
- Dispose of enriched samples according to current industry standards.

To reduce the risks associated with environmental contamination:

- Follow applicable local regulations for disposal of contaminated waste.

Precaution

Refer to SDS

Dehydrated Culture Medium Preparation

Ready-to-use option:

1. Rehydrate 26 grams of the medium with one liter of sterile water preheated to $50 \pm 1^\circ\text{C}$.

Autoclave option:

1. Dissolve 26 grams of the medium into one liter of purified water. Autoclave at 110°C for 15 minutes.
2. Store prepared broth at $15\text{--}30^\circ\text{C}$ for up to seven days, away from direct sunlight. Prepared broth does not need to be stored in the dark.
3. Prewarm the autoclaved, prepared broth at $50 \pm 1^\circ\text{C}$ for at least 15 hours and no longer than 60 hours before use.

Quality Control Specifications**Dehydrated**

Appearance: Medium should be pale to medium orange in color and free flowing; may contain small clumps that break apart easily.

Prepared

Appearance (110°C autoclave cycle): Red, clear to slightly hazy with trace to no precipitate.

pH: 7.0 ± 0.2 .

Performance Provides for the satisfactory cultural recovery of *Salmonella* after 6 hours of incubation at $42 \pm 1^\circ\text{C}$ using 0.5 – 5 CFU/mL of media: ATCC 14028, ATCC 51741, ATCC 9263.

Test Procedure

Incubate prepared Test Samples at $42 \pm 1^\circ\text{C}$. Refer to applicable assay instructions for incubation times and post enrichment sample preparation and assay steps.

Results

Refer to the applicable Pathogen Test product instructions under the section Interpretation of Results for a discussion of the test results.

Storage

Store dehydrated culture media (700006828) at $2\text{--}30^\circ\text{C}$ away from direct sunlight. Once opened and recapped, place container in a low humidity environment at the same storage temperature. Protect from moisture and light by keeping container tightly closed.

Limitations of the Procedure

1. Refer to the Neogen MDA2Q - Quantitative *Salmonella* product instructions.
2. Identification of *Salmonella* must be confirmed according to applicable regulatory guidelines or requirements.
3. Water for hydration or the autoclaved media must be appropriately prewarmed to 50°C before use and used within 5 minutes of being removed from 50°C (otherwise return the bottle of media to 50°C until it can be used).
4. Use non-autoclaved, re-hydrated medium the same day as prepared.
5. Any unused, autoclaved medium in an opened bottle must be used within the same day as the bottle was opened or be discarded.
6. Sample bags must be closed loosely to allow air exchange during incubation, which is vital for organism growth.
7. Do not use expired medium.

Expiration

Refer to expiration date stamped on the container. The dehydrated medium should be discarded if it is not free flowing, or if medium is darker than medium orange. Expiry applies to medium in its intact container when stored as directed.

Neogen Food Safety

Neogen Corporation

620 Leshar Place
Lansing, MI 48912 USA
Neogen.com

Neogen Europe Ltd.

The Dairy School
Auchincruive
Ayr, KA6 5HU
Scotland, UK

Neogen Ireland, Ltd.

Bray Business Park, Bray
Co. Wicklow
A98YV29, Ireland

Neogen Corporation

620 Leshar Place Lansing, MI 48912 USA
www.neogen.com

© Neogen Corporation 2024. All rights reserved.
Neogen is a trademark of Neogen Corporation.

All other trademarks are the property of their respective companies.

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