Clean-Trace® Test

-  (EN) Water Positive Control (3ng)
-  (FR) Témoin positif pour l'eau (3 ng)
-  (DE) Wasser-Positivkontrolle (3 ng)
-  (IT) Controllo positivo dell'acqua (3ng)
-  (ES) Control positivo para agua (3 ng)
-  (NL) Waterpositieve controle (3 ng)
-  (SV) Vatten positiv-test (3 ng)
-  (DA) Vand Positiv Kontrol (3 ng)
-  (NO) Positiv kontroll til vann (3ng)
-  (FI) Water positiivinen kontrolli (3 ng)
-  (PT) Controle Positivo de Água (3 ng)
-  (EL) Θετικής Επιβεβαίωσης (θετικός μάρτυρας) για Νερό (3 ng)
-  (PL) Kontrola dodatnia wody (3 ng)
-  (HU) Víz pozitív kontroll (3 ng)
-  (CS) Pozitivní kontrola vody (3 ng)
-  (RO) Control pozitiv al apei (3 ng)
-  (RU) Положительный контроль воды (3 нг)
-  (TR) Su Pozitif Kontrolü (3 ng)
-  (JA) 水中: ポジティブコントロール (3ng)
-  (ZH) 水质采样棒阳性对照 (3ng)

Clean-Trace® Test

Product Instructions

Water Positive Control (3ng)

Product Description and Intended Use

The measurement of the standard concentration of Adenosine Triphosphate (ATP) by Neogen® Clean-Trace® Water ATP Test Kits is used as an assurance of reagent and instrument function. In addition, the degree of quenching caused by interfering substances in the water sample can be determined using the Neogen Clean-Trace Water ATP Test Kits. The test kit contains 10 ready to use reagent vials. Each vial contains 3.0 ng or 5×10^{-12} moles of ATP. Also included are 10 vials of ATP free water to re-hydrate the ATP.

Neogen Clean-Trace Water ATP Test Kits – A Neogen® Clean-Trace® Water ATP Test Stick is used to collect a sample of re-hydrated ATP positive control. Upon activation of the test device, reagent in the cuvette of the test device reacts with ATP collected on the sample collection rings to produce light. The intensity of light is proportional to the amount of ATP. Measurements are performed in the Neogen® Clean-Trace® Luminometer and results displayed on the digital read out in Relative Light Units (RLU).

The signal obtained should be above a specified level, (see below) which indicates that both the reagents and instrument are operating correctly. The unique format of this product provides the user with a convenient method of assuring that all components used during testing are operating correctly, and that the operator is using the correct method. In some cases, the water samples under test can contain components that reduce or “quench” the activity of the reagents, thus reducing the bioluminescent signal. The degree of quenching can be determined by repeating the assurance test detailed in the descriptions above by replacing the ATP free water with a sample of the water under test. The test indicates quenching if the result obtained is lower than the result achieved during the assurance test. Alternatively, some substances can enhance light output.

This kit has been designed to be used in conjunction with the Neogen® Clean-Trace® Water Plus - Total ATP and Neogen® Clean-Trace® Water - Free ATP Tests to rapidly assess that both reagents and instrument are operating satisfactorily and have been stored and used correctly.

Warnings

Consult the Safety Data Sheet for additional information and local regulations for disposal.

ATP is a common substance. It is present on skin, hair and on many surfaces. Do not touch the swab or the interior surfaces of the positive control vials. Do not use devices past expiry date.

Notice

To reduce the risk of inaccurate results, Neogen® Clean-Trace® Water Positive Control Kit should be stored at 2°-8°C (36°-47°F). To avoid inaccurate results, keep reconstituted Neogen Clean-Trace Water Positive Control Kit at refrigerated temperatures 2°-8°C (36°-47°F). To reduce the risks associated with cross contamination while preparing the assay the use of gloves is recommended. To reduce the risk of inaccurate results do not shake the reconstituted reagents. Follow instructions for reconstituting reagents.

User Responsibility

Users are responsible for familiarizing themselves with product instructions and information. Visit our website at www.neogen.com, or contact your local Neogen representative or distributor for more information.

When selecting a test method, it is important to recognize that external factors such as sampling methods, testing protocols, sample preparation, handling, and laboratory technique may influence results.

It is the user's responsibility in selecting any test method or product to evaluate a sufficient number of samples with the appropriate matrices and microbial challenges to satisfy the user that the chosen test method meets the user's criteria.

It is also the user's responsibility to determine that any test methods and results meet its customers' and suppliers' requirements.

As with any test method, results obtained from use of any Neogen Food Safety product do not constitute a guarantee of the quality of the matrices or processes tested.

Limitation of Warranties / Limited Remedy

EXCEPT AS EXPRESSLY STATED IN A LIMITED WARRANTY SECTION OF INDIVIDUAL PRODUCT PACKAGING, NEOGEN DISCLAIMS ALL EXPRESS AND IMPLIED WARRANTIES, INCLUDING BUT NOT LIMITED TO, ANY WARRANTIES OF MERCHANTABILITY OR FITNESS FOR A PARTICULAR USE. If any Neogen Food Safety Product is defective, Neogen or its authorized distributor will, at its option, replace, or refund the purchase price of the product. These are your exclusive remedies. You must promptly notify Neogen within sixty days of discovery of any suspected defects in a product and return it to Neogen. Please contact your Neogen representative or authorized Neogen distributor for any further questions.

Limitation of Neogen Liability

NEOGEN WILL NOT BE LIABLE FOR ANY LOSS OR DAMAGES, WHETHER DIRECT, INDIRECT, SPECIAL, INCIDENTAL OR CONSEQUENTIAL DAMAGES, INCLUDING BUT NOT LIMITED TO LOST PROFITS. In no event shall Neogen's liability under any legal theory exceed the purchase price of the product alleged to be defective.

For detailed WARNINGS, CAUTIONS, STORAGE AND DISPOSAL information, and complete INSTRUCTIONS FOR USE see product instructions.

Storage and Disposal

For maximum shelf life, store between 2°-8°C (36°-47°F). Expiry date and lot number are noted on each package. After use, may contain microorganisms that may be a potential biohazard. Follow current industry standards for disposal. Do not freeze. Protect from direct sunlight.

Instructions for Use

Before the Neogen Clean-Trace Water Positive Control Kit is used, it is important to ensure your Neogen Clean-Trace Luminometer is switched on and initialized. Refer to the manual provided with the instrument for full details.

Neogen Clean-Trace Water Plus - Total ATP and Neogen Clean-Trace Water - Free ATP Test Applications

1. Place the Neogen Clean-Trace test devices and positive control vials at room temperature for at least 10 minutes before use.
2. Gently tap the positive control vial to ensure that the reagent pellet is at the bottom.
3. Carefully remove the stopper from the vial, avoid contaminating the opening to the vial.
4. Use a disposable Pasteur pipette or laboratory micropipette to transfer 2 mL of the ATP free water into the ATP vial.
5. Swirl the vial gently for 5 seconds.
6. Use a disposable Pasteur pipette or laboratory micropipette to transfer ALL the solution from the ATP vial back into the bottle containing the ATP free water.

Note: Most 1 mL laboratory micropipettes and some disposable 3 mL Pasteur pipettes will not be able to transfer all the solution because they can not reach the bottom of the ATP vial due to their diameter. Therefore, a smaller volume Micropipette or Pasteur pipette would also be required.

7. Swirl to mix the sample. Remove a Neogen Clean-Trace Water Plus - Total ATP test device from the foil pouch. Remove the sample stick from the test device and immerse the sample rings into the liquid under test, tapping the handle gently if bubbles form. Immediately remove the sample stick from the solution to be tested and carefully return the sample stick to the device such that the handle is at its starting position.
8. Process the sampling stick as per the Neogen Clean-Trace Water Plus - Total ATP or Neogen Clean-Trace Water - Free ATP product instructions. A maximum of five test devices may be tested using a single reconstituted vial of ATP positive control. The reconstituted ATP positive control should be used within 6 hours if stored at ambient temperature (21°C, 70°F) or 24 hours if stored at 2-8°C (36° - 47°F).

Test Procedure: Quench Check

- Take a bottle of ATP free water and discard its contents.
- Fill this bottle with 10 mL of the water sample under test.
- Use this sample water to reconstitute an ATP vial, as above.
- Test with a Neogen Clean-Trace Water ATP test device.

Interpretation of Results

Assay Assurance Test

Testing using the Neogen Clean-Trace Water ATP Test Kits should give readings greater than 1,000 RLU.

If the tests fail to reach these levels repeat the test procedure using a fresh vial of ATP positive control reagent. If low results are again achieved it may indicate a possible problem with either the reagent and/or instrument.

Quench Check

Comparing the results of the Neogen Clean-Trace Water ATP with the result for the respective quench check can derive the degree of quenching. The degree of quenching should be taken into account when interpreting your readings for that particular water sample.

Troubleshooting Guide

The following section is a list of potential problems that may cause a low signal:

- Incomplete mixing - shake the device and read again
- Reagents have passed their expiration date
- Reagents stored incorrectly
- Reconstituted ATP is contaminated
- Reagents not allowed to reach room temperature before testing
- Contaminated luminometer
- Damaged luminometer

For further information, contact Neogen Food Safety Technical Services.

Explanation of Symbols

info.neogen.com/symbols

Neogen Food Safety

Neogen Corporation

620 Leshar Place
Lansing, MI 48912 USA
www.neogen.com

Neogen Europe Ltd.

The Dairy School
Auchincruive
Ayr, KA6 5HU
Scotland, UK

Neogen Ireland, Ltd.

Bray Business Park, Bray
Co. Wicklow
A98YV29, Ireland

Neogen Corporation

620 Leshar Place Lansing, MI 48912 USA
www.neogen.com

© Neogen Corporation 2024. All rights reserved.
Neogen and Clean-Trace are registered trademarks of Neogen Corporation.
FS00928A

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