

Sample collection device and neutralizing solution checklist

Neutralizing Solution Features		Neogen® Wide Spectrum Neutralizer ^{1,2}	Current Neutralizing Solution (check all that apply)	Why it Matters
AOAC® Performance Tested Methods SM Program certified?		✓		Certification mark assures users that an independent assessment found that test method performance meets an appropriate standard for its intended use
HACCP International Food Safety Program Certification?		✓		Certification mark assures users that this product is food-safe and suitable for use in food handling facilities that operates in accordance with a HACCP based Food Safety Program
Effective neutralization of multiple sanitizers?				Effective neutralization of residual sanitizers help to ensure sample integrity for accurate downstream testing
Quaternary ammonium compounds		✓		A neutralizing solution that effectively neutralizes multiple types of sanitizers may reduce the need for food processors to match neutralizing solution to the sanitizer currently in use (i.e., one neutralizing solution is now sufficient)
High acid		✓		
Chlorine-based		✓		
Hydrogen peroxide and peroxyacetic acid-based		✓		
Formulation free from known allergenic components?		✓		Reduces the concern of the neutralizing solution leaving residual allergenic material on surfaces in food processing environments
Free from known PCR-inhibitors?		✓		Designed for compatibility with downstream molecular methods
Free from animal derived materials (ADMs)?		✓		Reduces the concern of the neutralizing solution leaving residual Animal Derived Materials on surfaces in food processing environments
Free from Genetically Modified Organism (GMO)-based components?		✓		Reduces the concern of the neutralizing solution leaving residual Genetically Modified Organism-based components on surfaces in food processing environments
Can hold samples for up to 96 hours (2–8°C)?		✓		Allows for more flexibility when shipping and processing samples
Can be stored at room temperature for up to 6 months?		✓		Allows for more flexibility with storage conditions
Designed for compatibility with Neogen® Petrifilm® Plates?		✓		Provides compatibility with Petrifilm Plates: AC, RAC, REC, SEC, EC/CC, CC, HSCC, EB, LAB, YM, and RYM (i.e., one neutralizing solution is now sufficient)
Designed for compatibility with Neogen® Molecular Detection System?		✓		Provides Molecular Detection System users added flexibility (i.e., one neutralizing solution is now sufficient)

Ordering Information

Product Description	SKU Catalog Number	Quantity
Neogen® Environmental Scrub Sampler Stick with 10 mL Wide Spectrum Neutralizer	700002338 ESS10WSN	100/case
Neogen® Environmental Scrub Sampler with 10 mL Wide Spectrum Neutralizer and Gloves	700002337 HES10WSN2G	100/case
Neogen® Sponge Stick with 10 mL Wide Spectrum Neutralizer	700007953 SSL10WSN	100/case

Get results you can trust at
info.neogen.com/samplecollection



¹ Internal Neogen data.
² AOAC® *Performance Tested Methods*™ Program certificate number 022104.
Neogen Corporation, 620 Lesher Place, Lansing, MI 48912 USA.
© Neogen Corporation 2026. All rights reserved. Neogen and Petrifilm are registered trademarks and Molecular Detection System is a trademark of Neogen Corporation. All other marks mentioned are the property of their respective owner.