

Technical Specification Sheet



XLT4 Supplement (NCM4079)

Intended Use

XLT4 Supplement is an anionic surfactant solution used as an inhibitory agent in XLT4 Agar. It is intended for use in culture media for the microbiological analysis of food, animal feed, water and environmental samples, and is for use with:

NCM0100 XLT4 Agar

Typical Formulation by Vial Size

	Per 100 mL Bottle NCM4079-100
Sodium Tetradecyl Sulfate (Tergitol 4)	100 mL

Precaution

Refer to SDS

Preparation

NCM4079 is ready to use.

Quality Control Specifications

Liquid Appearance: Colorless to liquid amber liquid

Storage

Store sealed vial at 15-30°C away from direct sunlight.

Expiration

Refer to expiration date stamped on the container. Expiry applies to medium in its intact container when stored as directed.

Precautions

Discard if supplement demonstrates any precipitate.



620 Leshar Place • Lansing, MI 48912
800-234-5333 (USA/Canada) • 517-372-9200
foodsafety@neogen.com • foodsafety.neogen.com