

Technical Specification Sheet



Potassium Tellurite 1% (NCM4012)

Intended Use

A selective culture media supplement used with multiple media to enhance selectivity. Potassium Tellurite 1% is intended for use in culture media for the microbiological analysis of food, animal feed, water and environmental samples, and is for use with:

NCM0024	Baird-Parker Agar
NCM0184	Modified Giolitti & Cantoni Broth
NCM0282	Vogel Johnson Agar

Typical Formulation by Vial Size

	Per 10 mL Vial NCM4012-10
Potassium Tellurite	0.1 g

Precaution

Refer to SDS

Preparation

NCM4012 is ready to use.

Quality Control Specifications

Liquid Appearance: Clear colorless liquid

Storage

Store sealed vial at 2-8°C away from direct sunlight.

Expiration

Refer to expiration date stamped on the container. Expiry applies to medium in its intact container when stored as directed.

Precautions

Discard if reconstituted supplement demonstrates any precipitate.



620 Leshar Place • Lansing, MI 48912
800-234-5333 (USA/Canada) • 517-372-9200
foodsafety@neogen.com • foodsafety.neogen.com