




Scrub dot technology to
disrupt biofilm and enhance
sample collection



Hydrated with Wide
Spectrum Neutralizer

A new standard of sampling

Do more with the Neogen® Environmental Scrub Sampler

Now you can add a sampling solution to your environmental monitoring program that is designed for compatibility with Neogen Petrifilm® Plates and the Neogen Molecular Detection System™, and is effective at neutralizing sanitizers commonly found in the food industry.

AOAC™ Performance Tested Methods™ Certificate 022104*

A new class of capability

Get to know the Neogen Environmental Scrub Sampler with 10 mL Wide Spectrum Neutralizer, an innovative solution for environmental microbial sampling applications.

- 1 Ergonomic stick design conforms to multiple grip types
- 2 Stick contains an additive for metal detectability
- 3 Bright color for visibility in food processing environments
- 4 1-step twist detachment for ease of use
- 5 Maintains integrity when sampling rough surfaces
- 6 Multilayered non-woven composite product designed for pick-up, release and absorbency
- 7 Hydrated with a proprietary 10 mL Wide Spectrum Neutralizer
- 8 Scrub dot technology to disrupt biofilm and enhance sample collection



Ordering Information

Product Description	Catalog Number	Quantity
Environmental Scrub Sampler Stick with 10 mL Wide Spectrum Neutralizer	ESS10WSN	100/case
Environmental Scrub Sampler with 10 mL Wide Spectrum Neutralizer and Gloves	HES10WSN2G	100/case

Get results you can trust at info.neogen.com/samplecollection



*The Neogen Environmental Scrub Sampler with 10 mL Wide Spectrum Neutralizer is the first and only sample collection device and neutralizing solution to be granted the AOAC *Performance Tested Method* certification.

Neogen Corporation, 620 Leshar Place, Lansing, MI 48912 USA.

© Neogen Corporation 2025. All rights reserved. Neogen and Petrifilm are registered trademarks and Molecular Detection System is a trademark of Neogen Corporation. All other marks mentioned are the property of their respective owner.

Wide Spectrum Neutralizer



- Formulated to effectively neutralize common sanitizers, including quats, high acid cleaners, chlorine-based, hydrogen peroxide and peroxyacetic acid-based sanitizers
- Free from known allergenic components, known PCR-inhibitors, animal-derived materials (ADMs) and components made from genetically modified organisms (GMOs)
- Compatible with Neogen Molecular Detection System
- Compatible with Petrifilm Plates: AC, RAC, REC, SEC, EC/CC, CC, HSCC, EB, LAB, YM, and RYM
- Allows up to a 96-hour sample hold time