



Free up your focus


Neogen® Molecular Detection Assay 2 (MDA 2) - *Campylobacter** and *Campylobacter* Enrichment Broth

There's now a better way to test for bacteria associated with poultry production. Find new efficiencies and simplify *Campylobacter* testing with Neogen MDA 2 - *Campylobacter*.

Featuring all the benefits of the Neogen Molecular Detection System, this assay and one-step enrichment broth are easy to use and accurate.

It's significantly faster than alternatives like PCR, immunoassay, and culture methods.

- Efficiently monitors for *Campylobacter*
- Simplified enrichment prep and protocol
- Easy to use and train
- Streamlined workflow
- Reliable
- Next-day results



Proven performance

Neogen MDA 2 - *Campylobacter*

While *Campylobacter* may not be your most frequent poultry test, it can't be ignored either — this pathogen causes an estimated 1.5 million illnesses each year in the U.S.¹ Now you can test for *Campylobacter* with a solution from Neogen.

The Neogen Molecular Detection System with Neogen MDA 2 - *Campylobacter* and MDA 2 - *Salmonella* is a complete solution that allows you to monitor poultry for both bacteria simultaneously. A single post-enrichment protocol for these pathogens means you can perform up to 96 tests of multiple types in one 60-minute run. And when compared using metrics such as time to result, throughput, unit size, and workflow, Neogen MDA 2 - *Campylobacter* outperforms other methods and manufacturers.

Campylobacter Enrichment Broth

- Easy-to-use
- Simplified enrichment — just add sterile water
- No supplements, microaerophilic incubation, organic solvents, or autoclaving the broth

Certifications and Validations

(view the method for specific details on the validations):

- AOAC® Performance Tested MethodSM (Certificate #111803)



Learn more at info.neogen.com/MDS

¹ US Centers for Disease Control and Protocol Yellow Book, Ch. 4 (2020)




Neogen Corporation, 620 Leshar Place, Lansing, MI 48912 USA.

© Neogen Corporation 2024. All rights reserved. Neogen is a registered trademark of Neogen Corporation.

All other marks mentioned are the property of their respective owners.

Simple Enrichment and Streamlined Workflow



Ordering Information

Description	Cat. No.	Contents
Molecular Detection Assay 2 - <i>Campylobacter</i>	MDA2CAM96	96 tests/kit
<i>Campylobacter</i> Enrichment Broth	CE250	1 bottle
Molecular Detection Assay 2 - <i>Salmonella</i>	MDA2SAL96	96 tests/kit
Bird Rinse with 400 mL Buffered Peptone Water Broth	BRBPW400	35/box
Bird Rinse Plain Sample Bag	BRB3500	500/box