

Sasso Grinder

Sasso Mill is an innovative, commercial grade grain grinder designed to meet industry needs.

Sampling is often cited as the primary contributor to sample variability within mycotoxin testing - Sasso Mill has been designed to advance the efficiency of your workflow and consistency of your test results.

- Grinds ~1 kg of grain per minute
- Hopper size capacity of ~2 kg
- Simple clean out procedure
- Meets industry standard for grind size in mycotoxin testing
- UL certified

Performance (Corn)

Grinder	>20 mesh	< 20 mesh	Total	% passing through 20 mesh	Time min:sec
Sasso	42.25 g	914.94 g	957.19 g	95.59%	1:10
	31.30 g	949.41 g	980.71 g	96.81%	1:09
	41.10 g	936.24 g	977.34 g	95.79%	1:08
Mean	38.22 g	933.53 g	971.75 g	96.06%	1:09

