



## Meet the New *Listeria* Right Now™

**Your environmental testing just got an upgrade.**

Now integrated with the Neogen® Molecular Detection system, *Listeria* Right Now™ features improved performance and a more efficient protocol. Make faster critical decisions with real-time results from user-friendly software.

Learn more at: [info.neogen.com/ListeriaRightNow](https://info.neogen.com/ListeriaRightNow)






### Robust Performance

- Detection where it counts. Collect samples on a range of surfaces, like sealed concrete, stainless steel, and plastic
- Tolerant to a range of sanitizers. chlorine-based, peroxyacetic acid, high acid cleaners, and quaternary ammonium
- Results in ~2 Hours. That includes samples to actionable results
- Unique sampler. Can sample up to a 12" x 12" surface area
- Same trusted sensitivity. Can detect 2-5 CFUs

### Neogen® Molecular Detection System

- User-friendly software with intuitive user interface for quick and easy setup
- Real-time results help make critical decisions faster
- Color-coded assays by organism and matching software screens allow for better tracking and easy identification

### Streamlined Workflow



### Ordering Information

Description	Cat. No.	SKU	Contents
Neogen® Molecular Detection Assay – <i>Listeria</i> Right Now™	MDALRN96	700007660	96 Tests/Kit
Neogen® Molecular Detection Right Now™ Matrix Control	RNMC96	700007661	96 Tests/Kit
Neogen® Molecular Detection Instrument and Accessory Kit	MDS100	700002195	1 Each
Neogen® Molecular Detection Heat/Chill Block Insert	MDSHBIN	700000145	1 Each

Request a free trial at: [info.neogen.com/ListeriaRightNow](https://info.neogen.com/ListeriaRightNow)

