

Technical Specification Sheet



Tryptose No. 2 (NCM0260)

Intended Use

Tryptose No. 2 is a special blend of peptones developed for use in the production of Foot and Mouth vaccine, and is not intended for use in the diagnosis of disease or other conditions in humans.

Description

A special blend of peptones originally developed for use in the production of Foot and Mouth vaccine. Other peptone blends are also available within the catalogue, some of which have been evaluated to demonstrate low cell toxicity and the ability to enhance the growth of cell lines (eg NCM0264).

Precaution

Refer to SDS

Test Procedure

Refer to appropriate references for specific procedures using Tryptose No. 2.

Quality Control Specifications

PHYSICAL/CHEMISTRY:	EXPECTED RESULTS
Appearance of Powder:	Free flowing, fine, homogeneous, beige powder
Clarity:	Clear, pale yellow color
Solubility in water:	Total

Expiration

Refer to expiration date stamped on container. Tryptose No. 2 should be discarded if not free flowing, or if the appearance has changed from the original color. Expiry applies to Tryptose No. 2 in its intact container when stored as directed.

Storage

Store powder at 2-30°C away from direct sunlight. Once opened and recapped, place container in a low humidity environment at the same storage temperature. Protect from moisture and light by keeping container tightly closed.

Record ID: FS-TSS-0105 Revision Number: 1.0 Effective Date: 2019-02-21 12:00 AM EST



620 Leshar Place • Lansing, MI 48912
800-234-5333 (USA/Canada) • 517-372-9200
foodsafety@neogen.com • foodsafety.neogen.com