

Technical Specification Sheet



Vancomycin (10mg/L) (NCM4004)

Intended Use

A selective culture media supplement for the isolation of *Cronobacter* spp. Vancomycin is intended for use in culture media for the microbiological analysis of food, animal feed, water and environmental samples, and is for use with:

NCM0227 Cronobacter Selective Broth

Typical Formulation

Vancomycin 10.0 mg/L

Typical Formulation by Vial Size

	Per 500mL Vial NCM4004-0.5
Vancomycin	5.0 mg

Precaution

Refer to SDS

Preparation

1. Using aseptic technique, remove the crimp cap and rubber stopper. Aseptically retain the stopper.
2. Reconstitute the contents of the vial using the appropriate volume of sterile deionized/RO water as indicated in the table below.
3. Replace the rubber stopper and gently shake or invert the contents to dissolve.
4. Add the vial contents to the cooled medium as per the label directions.

	NCM4004-0.5
Supplement Reconstitution Volume	5mL
Media Prep Volume	500 mL

Quality Control Specifications

Lyophilized Appearance: Homogenous white dried pellet; may break up with agitation.

Reconstituted Appearance: Clear colorless solution without a precipitate.

Storage

Store sealed vial at 2-8°C away from direct sunlight.

Expiration

Refer to expiration date stamped on the container. Expiry applies to supplement in its intact container when stored as directed.

Limitations

Discard if lyophilized supplement demonstrates any moisture uptake or appears to have liquefied; discard if reconstituted supplement demonstrates any precipitate.



620 Leshar Place • Lansing, MI 48912
800-234-5333 (USA/Canada) • 517-372-9200
foodsafety@neogen.com • foodsafety.neogen.com