

Ideal® D3™ Needles

Durable • Detectable • Dependable

Specifically designed for production animal use and processing plant detectability, **D3 Detectable Needles** also address packer HACCP and quality assurance programs.



D3 Needles - packaged 100 per box



Detectability

100% detectable in metal detection systems vs. 15% for standard needles.



Durability

Superior durability; thick wall cannula and plated brass hubs reduces breakage sixfold over conventional needles.



Sharpness

Precision ground tip stays sharper longer, allowing for multiple injections without burring.



Made to be easy

Gauge stamped on hub and on cartridge. Packed in red cartridges and boxes for easy identification. Tamper-evident cartridges are labeled with needle size and lot number.

Item No.	Size	Hub Type	Packaging	Colour Code
9418AL	14 G x 1"	Aluminium	Cartridge	Red
9419AL	14 G x 1 ½"	Aluminium	Cartridge	Red
9415	16 D x 5/8"	Brass	Cartridge	Red
9416	16 D x ¾"	Brass	Cartridge	Red
9420	16 D x 1"	Brass	Cartridge	Red
9417	16 D x 1 ½"	Brass	Cartridge	Red
1001NA	18 G x ¼"	Brass	Cartridge	Red
1001NA-AL	18 G x ¼"	Aluminium	Cartridge	Red
9421	18 D x ¾"	Brass	Cartridge	Red
9422	18 D x 5/8"	Brass	Cartridge	Red
9423	18 D x 1"	Brass	Cartridge	Red
9424	18 D x 1 ½"	Brass	Cartridge	Red
9411AL	20 G x 1"	Aluminium	Cartridge	Red
9412AL	20 G x ½"	Aluminium	Cartridge	Red



neogenaustralasia.com.au | 07 3736 2134 | naa-as@neogen.com | 14 Hume Drive, Bundamba QLD 4304

LK00192

© Neogen Corporation, 2024. Neogen and Ideal are registered trademarks and D3 is a trademark of Neogen Corporation, Lansing, MI. All other trademarks are property of their respective owners.