


# Simplify Your Workflow

## Neogen® Petrifilm® *Bacillus cereus* Count (BC) Plate

The preparation of MYP agar for *B. cereus* can be laborious and requires the handling of hazardous substrates that are perishable and need refrigeration. In addition, ISO 7932 recommends drying steps for this agar before and after inoculation of the sample which adds time to the process. With Petrifilm BC Plates, you can obtain results in just a few easy steps.



**It's easy to use. Easy to read.\* Easy to interpret.\* To request a demo or learn more, visit [info.neogen.com/Petrifilm](http://info.neogen.com/Petrifilm)**

\*Compared to reference method agar.