

Quick Connectors (NCM3201)

Applicable to:

NCM3202	Buffered Peptone Water (BPW) (ISO) (10 x 20L)
NCM3203	Buffered Peptone Water (BPW) (10 x 20L)
NCM3205	Half Fraser (Demi-Fraser) Broth (5 x 20L)
NCM3206	Less Plus Medium (5 x 20L)
NCM3207	Buffered Peptone Water (BPW) (HQ) (10 x 20L)

Intended Use

Used with NEOGEN's ready-to-reconstitute media bags, the quick connectors may be attached to tubing sets for the dispensing of reconstituted media. As with the rest of the ready-to-reconstitute range, the Quick Connectors are not intended for use in the diagnosis of disease or other conditions in humans.

Product Performance and Limitation of the Procedure

Performance is guaranteed as long as manufacturer's guidelines are followed. Aseptic technique must be observed during reconstitution and use. The quick connectors may be rinsed, re-sterilized and re-used.

NCM3201 Quick Connectors **must** be attached to the dispensing tubing and sterilized prior to use.

NCM3201 Quick Connectors are sold in boxes of 10.

Storage

2-30 °C