

Technical Specification Sheet



Letheen Broth Base Modified without Tween SKU: 700003340, 700003341, 700003342, 700003343 NCM0129

Intended Use

Letheen Broth Base Modified without Tween is used with Polysorbate 80 for the recovery of microorganisms from cosmetics in a laboratory setting. Letheen Broth Base Modified without Tween is not intended for use in the diagnosis of disease or other conditions in humans.

Description

In 1948, Weber and Black described the value of a highly nutritional solid medium containing neutralizing agents for quaternary ammonium compounds in sanitizers. The addition of Lecithin and Polysorbate 80 to Tryptone Glucose Extract (TGE) Agar resulted in a medium that effectively neutralizes quaternary ammonium compounds in testing of germicidal activity. Total neutralization of disinfectants is critical. Disinfectant residues can result in a false negative (no-growth) test.

Letheen Broth Base Modified without Tween is based on the formula described in FDA Bacteriological Analytical Manual, and a modification of Letheen Broth Base. Letheen Broth Base Modified is recommended by the FDA for use in the microbiological testing of cosmetics.

Typical Formulation

Letheen Broth Base	20.7 g/L
Enzymatic Digest of Casein	5.0 g/L
Enzymatic Digest of Animal Tissue	10.0 g/L
Yeast Extract	2.0 g/L
Sodium Bisulfite	0.1 g/L

Final pH: 7.2 ± 0.2 at 25°C

Formula is adjusted and/or supplemented as required to meet performance specifications.

Supplement

NCM4081 Tween 80

Precaution

Refer to SDS

Preparation

1. Dissolve 37.8 grams of the medium and 5 grams of NCM4081-100 or 700004920 Tween 80 in one liter of purified water.
2. Heat with frequent agitation to completely dissolve the medium.
3. Autoclave at 121°C for 15 minutes.

Quality Control Specifications

Dehydrated Appearance: Powder is homogeneous, free flowing, and light beige.

Prepared Appearance: Prepared medium is orange to amber and clear to lightly hazy.

Expected Cultural Response: Cultural response in Letheen Broth Base Modified without Tween (supplemented with NCM4081 or 700004920 Tween 80), incubated aerobically at $35 \pm 2^{\circ}\text{C}$ and examined for growth after 18 - 24 hours.



Technical Specification Sheet



Microorganism	Approx. Inoculum (CFU)	Expected Results
<i>Enterococcus faecalis</i> ATCC® 29212	10 - 300	Growth
<i>Escherichia coli</i> ATCC® 25922	10 - 300	Growth
<i>Pseudomonas aeruginosa</i> ATCC® 27853	10 - 300	Growth
<i>Staphylococcus aureus</i> ATCC® 25923	10 - 300	Growth

The organisms listed are the minimum that should be used for quality control testing.

Test Procedure

1. Prepare and dilute samples in Letheen Broth Base Modified (with Tween) in accordance with established guidelines. Proceed with procedure using Letheen Agar Base Modified NCM0130 or 700003344, 700003345, 700003346.
2. Using the spread plate technique, inoculate in duplicate 0.1 mL of the diluted samples onto Letheen Agar Base Modified, Potato Dextrose Agar containing Chlortetracycline, Baird Parker Agar, Anaerobic Agar, and a second set of Letheen Agar Base Modified plates.
3. Incubate one set of Letheen Agar Base Modified (with Tween) at 30°C for 48 hours and the other set at 35°C under anaerobic conditions for 2 - 4 days. Incubate the Potato Dextrose Agar plates at 30°C for 7 days and the Baird Parker Agar plates, if inoculated, at 35°C for 48 hours.
4. Incubate diluted samples from Step 1 at 35°C for 7 days. Subculture enriched samples onto Letheen Agar Base Modified only if there is no growth on the primary Letheen Agar Base, Modified.

Results

Refer to appropriate references for results

Expiration

Refer to expiration date stamped on the container. The dehydrated medium should be discarded if not free flowing, or if appearance has changed from the original color. Expiry applies to medium in its intact container.

Limitation of the Procedure

Due to nutritional variation, some strains may grow poorly or fail to grow on this medium.

Storage

Store dehydrated culture media at 2-30°C away from direct sunlight. Once opened and recapped, place container in a low humidity environment at the same storage temperature. Protect from moisture and light by keeping container tightly closed.

References

1. Weber, G. R., and L. A. Black. 1948. Relative efficiency of quaternary inhibitors. Soap and Sanit. Chem. 24:134-139.
2. www.fda.gov/Food/ScienceResearch/LaboratoryMethods/BacteriologicalAnalyticalManualBAM/default.htm.
3. Quisno, R., I. W. Gibby, and M. J. Foter. 1946. A neutralizing medium for evaluating the germicidal potency of the quaternary ammonium salts. Am. J. Pharm. 118:320-323.
4. Erlandson, A. L., Jr., and C. A. Lawrence. 1953. Inactivating medium for hexachlorophene (G-11) types of compounds and some substituted phenolic disinfectants. Science. 118:274-276.
5. Brummer, B. 1976. Influence of possible disinfectant transfer on *Staphylococcus aureus* plate counts after contact sampling. Appl. Environ. Microbiol. 32:80-84.
6. Favero (chm.). 1967. Microbiological sampling of surfaces-a state of the art report. Biological Contamination Control Committee, American Association for Contamination Control.

