



LIVESTOCK PRODUCTION
VETERINARY SUPPLIES
PROTECTIVE WEAR

ANIMAL SAFETY

CONTENTS



1 Veterinary Consumables

- 1 Wound Care
- 2 Bandages
- 3 General Products
- 4 Surgical Supplies

6 Instruments

- 6 Castration
- 7 Feeding and Watering
- 9 Hoof Care
- 10 Identification and Marking
- 13 Needles
- 19 Syringes
- 27 Syringe Accessories

31 Obstetrics and Reproduction

- 31 OB Equipment
- 35 Semen Collection
- 36 Insemination

37 Protective Wear

- 37 Disposable Boots
- 39 Gloves
- 41 OB Sleeves
- 43 Apparel

44 Index



Global leader in developing and manufacturing animal and food safety solutions

Neogen Australasia, based in Queensland, Australia, is a subsidiary of the Neogen Corporation, and is focussed on advancing human and animal wellbeing through science and technology to fuel a brighter future for global food security.

Neogen Australasia prides itself on being the market leader in the provision of premium animal safety solutions for the Australian and New Zealand livestock industry.

Wound Care

From simple scratches and lacerations to deep abscesses and surgical incisions, wounds are inevitable with animals. Providing timely and proper wound care is essential for quick recovery. Neogen provides a complete line of wound care products for animals for easy and proper wound maintenance. Neogen's products include topical antiseptics, FDA-approved antibacterial ointments, coagulants, and aluminum spray bandages. To keep dressings in place and wounds clean, Neogen has a complete line of cohesive wraps that can be found on the next page.

ALUSHIELD™

Aerosol Bandage

AluShield is a convenient, water-resistant, aerosol bandage that creates a protective barrier that protects wounds from dirt, debris, bacteria, and insects, while allowing the wound to breathe. For use on wounds in small and large animals. Contains aluminum powder, 4%.

Item No.	Description	Case Ct.
79100C	AluShield 75 g Can	12



GAUZE STERILE ROLLS

Individually wrapped sterile rolls of 6 ply gauze. Gauze is 3.4 in wide and 3.6 yds in length.

Item No.	Description	Case Ct.
TA4144	Gauze 6 Ply 3.4 in x 3.6 yd 1/roll	96



HOOF WRAP

Economical cohesive bandage for bovine applications. Sold by 100 count box only.

Item No.	Description	Box Ct.
TA3400HW-HP	SyrFlex Hoof Wrap 4 in Hot Pink	100
TA3400HW-O	SyrFlex Hoof Wrap 4 in Orange	100
TA3400HW-R	SyrFlex Hoof Wrap 4 in Red	100



Bandages

Neogen offers a wide variety of bandage products designed for animal health professionals. From gauze to cohesive bandages, we've got you covered.



SYRFLEX 2 IN COHESIVE BANDAGES

Self-adhesive, maintains original shape and size. Provides consistent lightweight compression; does not tighten with wear. 8 pound tensile strength. Individually wrapped. 2 in x 5 yds (stretched); 2 in x 2.2 yds (unstretched). Sold in eaches (36/box) and 6 Pack.

Item No.	Description	Box Ct.	Master Case Ct.
TA3200BL-E	SyrFlex Bandage 2 in Blue Each	36	216
TA3200HP-E	SyrFlex Bandage 2 in Hot Pink Each	36	216
TA3200NAT-E	SyrFlex Bandage 2 in Natural Each	36	216
TA3200RED-E	SyrFlex Bandage 2 in Red Each	36	216
TA3200LTGRN-E	SyrFlex Bandage 2 in Green Each	36	216
TA3200PUR-E	SyrFlex Bandage 2 in Purple Each	36	216
TA3200YEL-E	SyrFlex Bandage 2 in Yellow Each	36	216
TA3200PET-6pk	Syrflex Elastic Bandage 6 Pack	6	36



SYRFLEX 4 IN COHESIVE BANDAGES

Self-adhesive, maintains original shape and size, offers Starting Line Technology for ease of unrolling. Provides consistent lightweight compression; does not tighten with wear. 8 pound tensile strength. Individually wrapped. 4 in x 5 yds (stretched); 4 in x 2.2 yds (unstretched). Sold in eaches (18/box), 4 Packs, and Bulk (100/Box).

Item No.	Description	Box Ct.	Master Case Ct.
TA3400BL-E	SyrFlex Bandage 4 in Blue	18	108
TA3400BLK-E	SyrFlex Bandage 4 in Black	18	108
TA3400GRN-E	SyrFlex Bandage 4 in Green	18	108
TA3400HP-E	SyrFlex Bandage 4 in Hot Pink	18	108
TA3400NAT-E	SyrFlex Bandage 4 in Natural	18	108
TA3400OR-E	SyrFlex Bandage 4 in Orange	18	108
TA3400PUR-E	SyrFlex Bandage 4 in Purple	18	108
TA3400RED-E	SyrFlex Bandage 4 in Red	18	108
TA3400WH-E	SyrFlex Bandage 4 in White	18	108
TA3400YEL-E	SyrFlex Bandage 4 in Yellow	18	108
TA3400CP*	SyrFlex Bandage 4 in Colour Pack	18	108
TA3400TB-E	SyrFlex Bandage 4 in Teal	18	108




SYRFLEX 4 IN COHESIVE BANDAGES 4 PACK

Item No.	Description	Box Ct.	Master Case Ct.
TA3400BLK-4PK	SyrFlex Cohesive Bandage 4 in Black 4 pk	12	108
TA3400RBW-4PK	SyrFlex Cohesive Bandage 4 in Rainbow 4pk	12	108
TA3400BLUP-4PK	SyrFlex Cohesive Bandage 4 in Blue Crown 4pk	12	108
TA3400GRNP-4PK	SyrFlex Cohesive Bandage 4 in Green Horse Shoes and Ribbons 4pk	12	108
TA3400HPINP-4PK	SyrFlex Cohesive Bandage 4 in Hot Pink Star 4pk	12	108
TA3400PINP-4PK	SyrFlex Cohesive Bandage 4 in Pink Flamingo 4pk	12	108
TA3400TEALP-4PK	SyrFlex Cohesive Bandage 4 in Teal Boot and Cactus 4pk	12	108
TA3400UNIP-4PK	SyrFlex Cohesive Bandage 4 in Light Blue Unicorn 4pk	12	108
TA3400WHIP-4PK	SyrFlex Cohesive Bandage 4 in White Hearts 4pk	12	108
TA3400BLK-BULK	Syrflex Cohesive Bandage 4 in Black 100 count	100	N/A
TA3400WH-BULK	Syrflex Cohesive Bandage 4 in White 100 count	100	N/A



General Products

Along with quality veterinary instruments and animal handling supplies, Neogen offers a complementary line of products to round out your animal care needs. Our offering includes products used in veterinary practices, in livestock operations, and by livestock showmen.



ELASTIC ADHESIVE TAPE

Stretches easily and conforms to any shape without restriction of movement. Porous, allows skin to breath. Tearable. Centre marking line for ease of evenly overlapping tape. 2 in x 5 yds. Sold by box only

Item No.	Description	Case Ct.
TA3102E	Elastic Adhesive Tape 2 in 6/bx	6
TA3104E	Elastic Adhesive Tape 4 in 6/bx	6



CADAVER BAGS

Designed for the disposal of animals and are made from black 3.0 mm polyethylene film. Packaged and sealed in a polyethylene sleeve with twist-tie closures included in each package. Large Bags (30 in x 50 in). Small Bags (24 in x 30 in). Biodegradable and safe for incineration in all states.

Item No.	Description	Case Ct.
C300-R	Cadaver Bags L 50/roll (30 in x 50 in)	1
C240-R	Cadaver Bags S 50/roll (24 in x 30 in)	1



SLEEVE STOPPERS

Natural and anti-pop sleeve stoppers for use on serum bottles.

Item No.	Description	Case Ct.
45-2041	Anti-pop Natural 20 mm 100 ct	10
49-3061	Sleeve Stopper Natural 30 mm 100 ct	10



Surgical Supplies

Neogen's line of medical and surgical supplies consists of a variety of sutures and scalpels. Neogen's scalpels come in a wide range of sizes and blade materials, with sterile and non-sterile options available.



CARBON STEEL STERILE SCALPEL BLADES

Sterile blades are packaged 100/box, individually wrapped.

Item No.	Description
TA608	Carbon Steel Scalpel Blades #10 100/bx
TA611	Carbon Steel Scalpel Blades #12 100/bx
TA604	Carbon Steel Scalpel Blades #20 100/bx
TA606	Carbon Steel Scalpel Blades #22 100/bx

STAINLESS STEEL STERILE SCALPEL BLADES

These Sterile Stainless Steel surgical blades can be used with Ideal Scalpel Handle or Ideal Castration Kit.



Item No.	Description
2034	Scalpels Disposable #10 5pk
2035	Scalpels Disposable #12 5pk
2031	Surgical Blades #10 Curved 10pk
2032	Scalpel Blades #12 Hooked 10pk

SCALPEL HANDLES

Stainless steel. Can be sterilized. For blades up to size 15.

Item No.	Description	Case Ct.
TA627	Scalpel Stainless Steel Handle #4 (Suitable for blades compatible with handles sized #4.) 1/pk	12
2033	Scalpel Handle Stainless Steel #3 (Suitable for blades compatible with handles sized #3.) 1/pk	12



How to properly attach and remove a scalpel blade

Attaching blade to handle



1. Open package and slide handle towards blade.



2. Slide handle up into slot on blade.



3. Once grooves are lined up, push tip of handle down to snap blade into place.




4. Remove blade from package.


Removing blade from handle



5. To remove blade from handle, push up on the base of the blade.



6. Slide blade forward.



7. Remove blade completely and properly dispose.



DISPOSABLE SCALPEL

Disposable stainless steel scalpels are individually wrapped and come with a pre-attached Acrylonitrile Butadiene Styrene (ABS) plastic handle.

Item No.	Description	Case Ct.
TA619	Carbon Steel Disposable Scalpel #10 10/bx	12
TA623	Carbon Steel Disposable Scalpel #12 10/bx	12
TA614	Carbon Steel Disposable Scalpel #22 10/bx	12



IDEAL® MINI SCALPELS

The Ideal Mini Scalpels feature surgical stainless steel blades attached to a uniquely designed, easy-to-grasp handle. Sizes 10, 12, or 22 blade.

Item No.	Description	Case Ct.
PN110	Scalpel Mini Non-Sterile #10 40/bx	6
PN112	Scalpel Mini Non-Sterile #12 40/bx	6
TA642	Scalpel Mini Non-Sterile #12 2/pk	12
PN122	Scalpel Mini Non-Sterile #22 12/bx	6



VICRYL SUTURE

Ethicon size 2-0 coated Vicryl (violet braided) suture, 27 in (70 cm), with a tapered CP-1 needle, sterile, synthetic, and absorbable suture. Packaged on a retail card with a UPC.

Item No.	Description	Case Ct.
5020	Suture Kit	6



CHROMIC SUTURE

Sterile, single-use catgut sutures that are absorbable.

Item No.	Description	Case Ct.
ZA003M	Suture Chromic 2/0 with needle 50 m 3/pk	4
CA4	Suture Catgut Chromic USP 4 30 m	1
CA3	Suture Catgut Chromic USP 3 50 m	1



Castration

Neogen's line of castration and docking instruments are designed to help producers and veterinarians handle animals safely and easily.



CASTRATING BANDS

Strong and dependable, these latex castrating bands can be used with Ideal's economy or premium band castrators. Bands are packaged in bags of 100.



Item No.	Description	Case Ct.
2007	Castrating Bands Green Latex 100 bag	400

ECONOMY BAND CASTRATOR

Ideal's durable Economy Band Castrators were designed for the cost-conscious producer. These are to be used with the above listed bands for castration of lambs, calves, goats, and other animals. They may also be used for tail docking or prolapses. The economy castrator will open to a 1.75 inch diameter.



Item No.	Description	Case Ct.
2008	Plastic Band Castrating Pliers	12
TA131	Metal Band Castrating Pliers	12

PREMIUM BAND CASTRATOR

This easy-to-use instrument was designed for the one-person operation of castrating lambs, calves, goats, and other animals quickly and easily. It may also be used for docking and prolapses. The premium castrator has a 1.75 inch diameter opening.



Item No.	Description	Case Ct.
2005	Premium Metal Band Castrating Pliers	12

Feeding and Watering

Feeding and watering livestock is essential for any successful production operation. Take care of your animals with Neogen's feeding and watering supplies. Our selection includes everything from the popular red screw on teat to feeding buckets with multiple teats, watering bowls, and calf weaners.



CALF BOTTLE

Nursing bottle holds 4 pints (1/2 gallon/2 liters) and comes with pint calibration markings. Bottle is made from strong, durable polyethylene and includes an X-pierced, self-closing, rubber teat. These snap-on type nipples are strong, soft, and pliable, with a tight fit to help prevent leakage.

Item No.	Description	Case Ct.
FW1216-SYR	Calf Bottle Red Teat Snap-On	12
FW1216A-SYR	Red Teat Snap-On	50
FW1222A-SYR	Premium Pink Teat Snap-On	50



FLOATING TEAT

Easy to take apart and clean, the floating teat aids in teaching calves to drink from a pail or bucket. The extra large plate allows the nipple to float and function like a straw, without the risk of tipping over.

Item No.	Description	Case Ct.
FW1207	Floating Teat Complete w/ Plate	12
FW1207A	Replacement Teats 2 pk	12



ORAL FLUID FEEDERS

Safe and practical way to fluid feed calves with a one piece flexible probe to deliver nutrients directly to the stomach. Recommended method for rehydration of stressed cattle.


Item No.	Description	Case Ct.
TA20006	Oral Fluid Feeder 64 oz	1
TA20009	Oral Fluid Feeder 128 oz	1



CALF WEANERS

These anti-sucking devices impede calves' ability to suck. Durable and easy to apply.

Item No.	Description	Case Ct.
TA913	Weaner Calf Plastic Yellow	50
TA916	Weaner Calf Plastic Orange	50
TA918	Weaner Calf Plastic Orange 50 pk	10





LAMB BOTTLE WITH RED SCREW-ON TEAT

Lambs and kids quickly learn to nurse when using Red Screw-On Teats. The teats are made of soft red latex, which closely resemble the ewe's teat, and screws easily onto most standard soda bottles. Threading on cap prevents leakage. Special flutter valve prevents air intake, and the tip of teat can be cut off to increase milk flow. Wash with warm, soapy water. Comes complete with 16 oz bottle.

Item No.	Description	Case Ct.
7004	Lamb Bottle 16 oz with Red Screw-On Teat	24

RED SCREW-ON TEAT

Ideal's Red Screw-On Teat is made of soft latex to closely resemble the ewe's teat. Screws easily onto most standard soda bottles. Teats prevent air intake.

Item No.	Description	Case Ct.
TA1932A	Red Screw-On Teat Each	100
7006	Red Screw-On Teat 2 pk	24



CONTROLLED FLOW LAMB AND KID NIPPLES

Our controlled flow lamb nipples are made of clear isoprene with extremely soft, yet strong walls. The hollow cone tips allow the flow of milk to be regulated from slow to fast. Nipples fit most standard soda bottles by snapping on.

Item No.	Description	Case Ct.
7007	Controlled Flow Nipple 3 pk	24



GUM RUBBER LAMB NIPPLES

Ideal's gum rubber nipples will fit most standard soda bottles by easily snapping on.

Item No.	Description	Case Ct.
7005	Gum Rubber Lamb and Foal Nipples 3 pk	24



Hoof Care

Hoof care is an essential part of animal health and Neogen provides a variety of quality hoof care tools to keep your sheep or goats moving soundly. Whether you are a producer or hobby farmer you can rely on Neogen's range of hoof rot shears, which includes the trusted brand name — Burgon & Ball — known for their quality in every sense of the word.



BURGON AND BALL EZE-TRIM FOOT ROT SHEARS

These super-sharp foot rot shears are made by Burgon and Ball and have non-stick coated, high carbon-steel blades that retain the sharpest edge. Blades have been heat-treated to give the optimum combination of strength and hardness. The blade tension can be adjusted and the blades lock in place. Resin handles are both durable and extremely comfortable to use. Shears come with a 5-year manufacturer's guarantee.



Item No.	Description	Case Ct.
7020	Foot Rot Shear B and B Eze-Trim	10

HOOF ROT SHEARS

Ideal Hoof Rot Shears have sharp, pointed double stainless steel blades. The top blade is serrated to cut hooves with ease. Blue powder coated handle. Use to trim cattle, sheep, and goat hooves.



Item No.	Description	Case Ct.
9801	Hoof Rot Shears	12



Identification and Marking

Being able to quickly and easily identify livestock is essential for livestock producers. Neogen provides a wide array of solutions for both permanent and non-permanent identification. Permanent livestock identification solutions include plastic ear tags, tattoo kits with a wide variety of accessories, and two sizes of ear notchers. For temporary identification to mark animals for treatment, breeding, calving and lambing, etc. we offer reusable plastic and Velcro leg bands, marking sticks, livestock marking paints, and dyes.



PLASTIC LEG BANDS

Reusable, tough leg bands for marking cows in a group pen. The bands have a ratchet lock for easy adjustment on different size ankles. Available in several colours for flexibility in herd identification. Use to indicate a temporary condition such as mastitis, or to number animals individually. Easy to attach and easy to remove and reuse. Sold in eaches.



Item No.	Description	Case Ct.
TA1515	Leg Band Blank White	100
TA1519	Leg Band Blank Orange	100
TA1522	Leg Band Blank Red	100
TA1525	Leg Band Blank Blue	100
TA1528	Leg Band Blank Green	100
TA1531	Leg Band Blank Yellow	100

TATTOO KITS

Ideal Tattoo Kits are versatile, economical instruments with heavy-duty construction used for permanent animal identification. The large pliers hold four 3/8 inch letters or digits. Aluminum alloy pliers.

Item No.	Description	Case Ct.
7502	Tattoo Kit Standard 3/8 in Complete	6
7503	Tattoo Digits (0-9) 3/8 in	6
7504	Tattoo Digits (0-9) 5/16 in	6
7510	Tattoo Kit Small 5/16 in	12
7530	Tattoo Letters (A-Z) 3/8 in	12
7531	Tattoo Letters (A-Z) 5/16 in	12
7577	Tattoo Ink Black 3 oz	12
7578	Tattoo Paste Green 1 oz	12
7580	Tattoo Ink Black 2 oz Roll-on	12
7590	Tattoo Ink Green 2 oz Roll-on	12



PRIMA® RETRACTABLE MARKING STICKS

Designed for easy and efficient animal marking, these markers come in a twistable plastic dispenser with cap to help keep hands and clothing clean. Highly visible, long lasting, and weather and fade resistant. Will mark wet or dry animals. Non-toxic.

Item No.	Description	Case Ct.
333643	Prima Retractable Marking Stick Black	50
333650	Prima Retractable Marking Stick Orange	50
333667	Prima Retractable Marking Stick Red	50
333674	Prima Retractable Marking Stick Green	50
333681	Prima Retractable Marking Stick Blue	50
333704	Prima Retractable Marking Stick Purple	50
333711	Prima Retractable Marking Stick Fluorescent Pink	50



PRIMA MARC™ METERED TIP PAINTS AND DYES

Prima Marc paints and dyes come in 1.75 oz cans that are easy to use and carry. A semi-permanent spray that produces easily visible marks on livestock, and for other uses. Designed for use with Prima Marc Vaccinators. They produce a metered spray that marks 1000+ animals.

Item No.	Description	Case Ct.
369017	Prima Marc Metered Spray Orange	24
369024	Prima Marc Metered Spray Green	24
369031	Prima Marc Metered Spray Red	24
369048	Prima Marc Metered Spray Blue	24
369055	Prima Marc Metered Spray Pink	24



PRIMA MARC PAINTS

This product is designed to work with the Prima Marc Vaccinator. Once mounted to the Prima Marc Syringe, a squeeze of the syringe handle activates the mounted paint can and marks each animal simultaneously with injection. This encourages safe administration practices, improves quality of control efforts, and saves costs of missed animals. The Prima Marc paint can be used as a stand-alone product. It is pocket sized and can be used for treating and sorting livestock, marking fence lines, landscaping, trail marking, and other activities where a quick, easy-to-carry marking solution is needed. Available in 1.75 oz cans.

Item No.	Description	Case Ct.
334077	Prima Marc Orange	24
334091	Prima Marc Green	24
338686	Prima Marc Pink	24
341280	Prima Marc Red	24
341297	Prima Marc Blue	24
342812	Prima Marc Black	24



PRIMA CONCENTRATE PAINT

- Economical and easy to use
- For use in sprayer systems
- Adjust concentration depending on duration needed

Item No.	Description	Case Ct.
333865	Concentrate Pink 16 oz	12
333872	Concentrate Green 16 oz	12
333889	Concentrate Blue 16 oz	12
342669	Concentrate Pink 128 oz	4
342676	Concentrate Blue 128 oz	4
342683	Concentrate Green 128 oz	4



PRIMA SPRAY-ON

- Colours are durable and stable
- Reliable in all temperatures
- Strong colour
- Sprays with can inverted
- Coloured Spray-On available in 12 oz cans
- Black Spray-On available in 13 oz cans

Item No.	Description	Case Ct.
333513	Spray-On Red	12
333520	Spray-On Blue	12
333537	Spray-On Green	12
338747	Spray-On Black	12



PRIMA SPRAY-ON GLO

- Colours are durable and stable
- Reliable in all temperatures
- Sprays with can inverted
- Strong colour with high luminosity
- 13 oz cans

Item No.	Description	Case Ct.
338693	Spray-On Glo Fluorescent Green	12
338389	Spray-On Glo Fluorescent Orange	12
339690	Spray-On Glo Fluorescent Purple	12
339942	Spray-On Glo Fluorescent Pink	12
339959	Spray-On Glo Fluorescent Yellow	12



PRIMA SPRAY-ON II

- Longer lasting than regular Spray-On
- Reliable in all temperatures
- Strong colour
- Sprays with can inverted
- 12 oz cans

Item No.	Description	Case Ct.
354228	Spray-On II Red	12
354235	Spray-On II Blue	12
354242	Spray-On II Green	12



PISTOL GRIP


For use with Spray-On, Spray-On Glo, and Spray-On II.

Item No.	Description	Case Ct.
344458	Pistol Grip for Spray-On Can	1



Needles

Neogen provides the highest quality and broadest range of needles for livestock producers and veterinarians. Our needles are ultra-sharp and designed specifically for animal use with metal inserts for strength and computer ground tri-bevelled point and anti-coring tips for sharpness and durability. Needles are available with polypropylene and aluminium hubs and are colour-coded following ISO standards for easy identification of the gauge size.



IDEAL D3™ DETECTABLE NEEDLES

Durable. Detectable. Dependable.

Neogen's D3 Detectable needle is the leader when it comes to meat industry safety. With a thick wall cannula and brass hub, needle breakage is greatly reduced over conventional and stainless needles.

Durability

- Superior durability; thick wall cannula and plated brass hubs reduce breakage 6-fold over conventional needles

Detectability

- 100% detectable in metal detection systems vs. only 15% for standard needles


Sharpness

- Precision-ground tip stays sharper longer, allowing for multiple injections without burring

Packaging Features

- Gauge stamped on hub and packed in red cartridges and boxes for easy identification
- Cartridges labels are tamper evident and ISO colour-coded to easily identify gauge and length

Item No.	Description	Colour	Case Ct.	Master Case Ct.
Brass Hub D3 Needles				
9415	D3 16 G x 5/8 in Detectable HP 100/bx	Red	10	40
9416	D3 16 G x 3/4 in Detectable HP 100/bx	Red	10	40
9420	D3 16 G x 1 in Detectable HP 100/bx	Red	10	40
9417	D3 16 G x 1 1/2 in Detectable HP 100/bx	Red	10	40
1001NA	D3 18 G x 1/4 in Detectable HP 100/bx	Red	10	40
9428	D3 18 G x 1/2 in Detectable HP 100/bx	Red	10	40
9422	D3 18 G x 5/8 in Detectable HP 100/bx	Red	10	40
9421	D3 18 G x 3/4 in Detectable HP 100/bx	Red	10	40
9423	D3 18 G x 1 in Detectable HP 100/bx	Red	10	40
9424	D3 18 G x 1 1/2 in Detectable HP 100/bx	Red	10	40
Aluminum Hub D3 Needles				
9418AL	D3 AH 14 G x 1 in Detectable HP 100/bx	Red	10	40
9419AL	D3 AH 14 G x 1 1/2 in Detectable HP 100/bx	Red	10	40
9412AL	D3 AH 20 G x 1/2 in Detectable HP 100/bx	Red	10	40
9411AL	D3 AH 20 G x 1 in Detectable HP 100/bx	Red	10	40
9430AL	D3 AH 20 G x 1 1/2 in Detectable HP 100/bx	Red	10	40



IDEAL D3X™ DETECTABLE NEEDLES

Dual Detection

D3X needles are Neogen's patented dual-detection needles featuring our trusted detectable metal alloy and highly visible extraction collar. Made from the same proprietary metal alloy, D3X needles are up to 6 times more detectable when compared to 304 stainless steel needles and allow needle fragments to be detected at line speed using magnetic technology. The extraction collar is designed to prevent a broken needle from submerging into the skin of the animal. With the easy-to-grip extraction collar adhered to the needle, this allows needle removal to happen immediately following a broken needle. The extraction collar is highly visible and recognizable so it can be seen in all lighting conditions.

Durability

D3X Needles have thick-walled cannulas, which are crimped and pressed into sturdy, aluminum hubs. These stronger sidewalls increase load strength and decrease bending and breaking more than 6-fold over some conventional needles. The needle tips stay sharper longer and makes D3X needles the superior option for production needle selection.

Dependability


Neogen has been the trusted manufacturer for production farmers for the past 20 years. D3X needles have ultra-sharp, tri-beveled tips that allows for easier penetration into the animal. The sharpness of the needle increases the possible number of injections per needle and reduces the risk of injection site abscesses while decreasing stress on the animal while working.

Packaging

D3X Needles are packed in trademarked red, colour-coded cartridges for safety and convenience while transporting. Cartridges are tamper evident and ISO colour-coded for easy gauge and length recognition. D3X needles have thick walled detectable cannulas and hubs. D3X Needles are packaged 100 per box.



Product Number	Product Description
9461	NDL D3X 14 G x 5/8 in Detectable HP 100/bx
9456	NDL D3X 14 G x 1 in Detectable HP 100/bx
9448	NDL D3X 14 G x 1 1/2 in Detectable HP 100/bx
9459	NDL D3X 16 G x 1/2 in Detectable HP 100/bx
9446	NDL D3X 16 G x 3/4 in Detectable HP 100/bx
9445	NDL D3X 16 G x 5/8 in Detectable HP 100/bx
9449	NDL D3X 16 G x 1 in Detectable HP 100/bx
9447	NDL D3X 16 G x 1 1/2 in Detectable HP 100/bx
9452	NDL D3X 18 G x 1/2 in Detectable HP 100/bx
9450	NDL D3X 18 G x 3/4 in Detectable HP 100/bx
9451	NDL D3X 18 G x 5/8 in Detectable HP 100/bx
9453	NDL D3X 18 G x 1 in Detectable HP 100/bx
9454	NDL D3X 18G x 1 1/2 in Detectable HP 100/bx
9444	NDL D3X 20 G x 1/2 in Detectable HP 100/bx
9457	NDL D3X 20 G x 1 in Detectable HP 100/bx
9458	NDL D3X 20 G x 1 1/2 in Detectable HP 100/bx



1. Needle intact
2. Bent needle with 1 or 2 of the red arms flared out
3. Broken needle with all 4 arms flared



Broken needle with extraction collar in an animal



Disposable

- A full line of disposable needles manufactured specifically for the animal health market
- Metal insert increases needle strength for animal health applications
- Tri-beveled, anti-coring tips are computer ground for consistent sharpness
- Anti-friction coating for minimal penetration resistance and animal discomfort
- Regular wall stainless steel cannula for durability and breakage prevention
- Tamper-evident cartridges are colour-coded to ISO standards
- Gauge, length, and lot number printed on each cartridge for easy tracking
- Packaged 100 needles/box, 10 boxes/case

ALUMINUM HUB NEEDLES

Aluminum hub for large animal applications that require a durable needle.

Item No.	Description	Colour	Case Ct.	Master Case Ct.
9400AL	AH 14 G x 1 in HP 100/bx	Green	10	40
9401AL	AH 14 G x 1 1/2 in HP 100/bx	Green	10	40
9402AL	AH 14 G x 2 in HP 100/bx	Green	10	40
9403AL	AH 16 G x 5/8 in HP 100/bx	White	10	40
9404AL	AH 16 G x 3/4 in HP 100/bx	White	10	40
9405AL	AH 16 G x 1 in HP 100/bx	White	10	40
9406AL	AH 16 G x 1 1/2 in HP 100/bx	White	10	40
9407AL	AH 18 G x 5/8 in HP 100/bx	Pink	10	40
9408AL	AH 18 G x 3/4 in HP 100/bx	Pink	10	40
9409AL	AH 18 G x 1 in HP 100/bx	Pink	10	40
9410AL	AH 18 G x 1 1/2 in HP 100/bx	Pink	10	40
06-2043B	AH 20 G x 3/4 in HP 100/bx	Yellow	10	40
9413AL	AH 20 G x 1 in HP 100/bx	Yellow	10	40
9414AL	AH 20 G x 1 1/2 in HP 100/bx	Yellow	10	40



POLYPROPYLENE HUB NEEDLES

Clear polypropylene hubs for easy viewing of flashback.

Item No.	Description	Colour	Case Ct.	Master Case Ct.
9319	PH 16 G x 1/2 in HP 100/bx	White	10	40
9325	PH 16 G x 1 in HP 100/bx	White	10	40
9330	PH 16 G x 1 1/2 in HP 100/bx	White	10	40
9340	PH 18 G x 3/4 in HP 100/bx	Pink	10	40
9345	PH 18 G x 1 in HP 100/bx	Pink	10	40
9345K	PH 18 G x 1 in SP 100/bx	Pink	10	40
9350	PH 18 G x 1 1/2 in HP 100/bx	Pink	10	40
9350K	PH 18 G x 1 1/2 in SP 100/bx	Pink	10	40
9360	PH 20 G x 1/2 in HP 100/bx	Yellow	10	40
9365	PH 20 G x 3/4 in HP 100/bx	Yellow	10	40
9370	PH 20 G x 1 in HP 100/bx	Yellow	10	40
9370K	PH 20 G x 1 in SP 100/bx	Yellow	10	40
9375	PH 20 G x 1 1/2 in HP 100/bx	Yellow	10	40
9375K	PH 20 G x 1 1/2 in SP 100/bx	Yellow	10	40
9380K	PH 21 G x 1 in SP 100/bx	Dark Green	10	40
9390	PH 22 G x 3/4 in HP 100/bx	Gray	10	40
9390K	PH 22 G x 3/4 in SP 100/bx	Gray	10	40
9395	PH 22 G x 1 in HP 100/bx	Gray	10	40
9395K	PH 22 G x 1 in SP 100/bx	Gray	10	40
9399	PH 25 G x 5/8 in HP 100/bx	Orange	10	40
9399K	PH 25 G x 5/8 in SP 100/bx	Orange	10	40
02-1605B	PH 16 G x 1/2 in SP 100/bx	White	10	40
02-1615B	PH 16 G x 1 1/2 in SP 100/bx	White	10	40
02-1805B	PH 18 G x 1/2 in SP 100/bx	Pink	10	40
02-2005B	PH 20 G x 1/2 in SP 100/bx	Yellow	10	40



ISO Needle Gauge Colour Code Chart

Green	White	Pink	Yellow	Dark Green	Gray	Blue	Orange
14 G	16 G	18 G	20 G	21 G	22 G	23 G	25 G

Retail Packs

Ideal needles, syringes, and combos offered in small quantities and packaged for the retail shelf.

ALUMINUM HUB CARTRIDGE PACKED NEEDLES

Item No.	Description	Colour	Case Ct.	Master Case Ct.
9317	AH 14 G x 1 in HP 5 pk	Green	150	600
9323	AH 16 G x 1 in HP 5 pk	White	150	600
9328	AH 16 G x 1 1/2 in HP 5 pk	White	150	600
9329	AH 16 G x 1 in HP 25 pk	White	50	200
9342	AH 18 G x 5/8 in HP 5 pk	Pink	150	600
9343	AH 18 G x 3/4 in HP 5 pk	Pink	150	600
9347	AH 18 G x 1 in HP 5 pk	Pink	150	600
9348	AH 18 G x 1 in HP 25 pk	Pink	50	200
9353	AH 18 G x 1 1/2 in HP 5 pk	Pink	150	600
9371	AH 20 G x 1 in HP 5 pk	Yellow	150	600
9376	AH 20 G x 1 1/2 in HP 5 pk	Yellow	150	600



POLY HUB CARTRIDGE PACKED NEEDLES

Item No.	Description	Colour	Case Ct.	Master Case Ct.
9318	PH 16 G x 1/2 in HP 5 pk	White	150	600
9326	PH 16 G x 1 in HP 5 pk	White	150	600
9327	PH 16 G x 1 in HP 25 pk	White	50	200
9331	PH 16 G x 1 1/2 in HP 5 pk	White	150	600
9341	PH 18 G x 3/4 in HP 5 pk	Pink	150	600
9344	PH 18 G x 3/4 in HP 25 pk	Pink	50	200
9346	PH 18 G x 1 in HP 5 pk	Pink	150	600
9351	PH 18 G x 1 1/2 in HP 5 pk	Pink	150	600
9352	PH 18 G x 1 1/2 in HP 25 pk	Pink	50	200
9355	PH 18 G x 1 in HP 25 pk	Pink	50	200
9361	PH 20 G x 1/2 in HP 5 pk	Yellow	150	600
9362	PH 20 G x 1/2 in HP 25 pk	Yellow	50	200
9366	PH 20 G x 3/4 in HP 25 pk	Yellow	50	200
9367	PH 20 G x 3/4 in HP 5 pk	Yellow	150	600
9372	PH 20 G x 1 in HP 5 pk	Yellow	150	600
9373	PH 20 G x 1 1/2 in HP 5 pk	Yellow	150	600
9374	PH 20 G x 1 in HP 25 pk	Yellow	50	200
9391	PH 22 G x 3/4 in HP 5 pk	Gray	150	600
9392	PH 22 G x 3/4 in HP 25 pk	Gray	50	200
9396	PH 22 G x 1 in HP 25 pk	Gray	50	200
9397	PH 22 G x 1 in HP 5 pk	Gray	150	600

Need Help Selecting the Right Needle?

22G

Length	Type	Application
3/4 in	SQ	Small Animals
1 in	IM	Small Animals



18G

Length	Type	Application
1/4 in	SQ	Lambs, Kids, Sheep, Goats
1/2 in	SQ	Cattle, Piglets, Nursery Pigs, Sheep, Goats
5/8 in	SQ	Cattle, Piglets, Sheep, Goats
3/4 in	SQ	Cattle, Sheep, Goats
1/2 in	IM	Piglets
5/8 in	IM	Piglets, Nursery Pigs
3/4 in	IM	Cattle, Nursery Pigs
1 in	IM	Cattle, Sheep, Goats, Horses
1 1/2 in	IM	Horses



20G

Length	Type	Application
1/2 in	SQ	Piglets, Sheep, Goats
5/8 in	SQ	Piglets, Sheep, Goats
3/4 in	SQ	Sheep, Goats
1/2 in	IM	Piglets
5/8 in	IM	Piglets
3/4 in	IM	Calves < 300 lbs
1 in	IM	Calves < 300 lbs, Horses
1 1/2 in	IM	Horses



16G

Length	Type	Application
1/2 in	SQ	Cattle, Nursery Pigs, Sheep, Goats
5/8 in	SQ	Cattle, Sheep, Goats
3/4 in	SQ	Cattle, Market Pigs, Sheep, Goats
1 in	SQ	Breeding Pigs
5/8 in	IM	Nursery Pigs
3/4 in	IM	Cattle > 300 lbs, Nursery Pigs
1 in	IM	Cattle > 300 lbs, Market Pigs, Breeding Pigs, Sheep, Goats
1 1/2 in	IM	Breeding Pigs



The needle selection information provided is only a guide to possible uses of different-sized needles. It is not intended to supersede the advice of a licensed veterinarian. Always consult with your veterinarian on choosing the correct needle size, injection location, and injection type for the medication or vaccine you are administering.

Durable

IDEAL STAINLESS STEEL NEEDLES

Each Ideal Stainless Steel Needle is tri-beveled and precision ground to maintain a sharper, longer-lasting penetration surface. The cannulas are crimped and pressed into a sturdy plated-brass hub. Non-sterile. Needles are available in clamshell packs of 3 or 12.

Item No.	Description	Case Ct.
1178	14 G x 3/4 in SS 3 pk	24
1179	14 G x 1 in SS 3 pk	24
1180	14 G x 1 1/2 in SS 3 pk	24
1184	16 G x 1/2 in SS 3 pk	24
1185	16 G x 3/4 in SS 3 pk	24
1186	16 G x 1 in SS 3 pk	24
1187	16 G x 1 1/2 in SS 3 pk	24
1190	18 G x 1/2 in SS 3 pk	24
1191	18 G x 3/4 in SS 3 pk	24
1192	18 G x 1 in SS 3 pk	24
1193	18 G x 1 1/2 in SS 3 pk	24
1195	20 G x 1/2 in SS 3 pk	24
1197	20 G x 1 in SS 3 pk	24
1202	12 G x 1 in Premium SS 12 pk	12
1204	12 G x 2 in Premium SS 12 pk	12
1205	12 G x 3 in Premium SS 12 pk	12
1210	14 G x 1/2 in Premium SS 12 pk	12
1211	14 G x 3/4 in Premium SS 12 pk	12
1212	14 G x 1 in Premium SS 12 pk	12
1214	14 G x 2 in Premium SS 12 pk	12
1220	16 G x 1/2 in Premium SS 12 pk	12
1221	16 G x 3/4 in Premium SS 12 pk	12
1222	16 G x 1 in Premium SS 12 pk	12
1223	16 G x 1 1/2 in Premium SS 12 pk	12
1230	18 G x 1/2 in Premium SS 12 pk	12
1231	18 G x 3/4 in Premium SS 12 pk	12
1232	18 G x 1 in Premium SS 12 pk	12
1233	18 G x 1 1/2 in Premium SS 12 pk	12
1240	20 G x 1/2 in Premium SS 12 pk	12
1241	20 G x 3/4 in Premium SS 12 pk	12
Transfer Needle		
1280	Transfer 6 pk	12

FLOATING NEEDLE

This attachment for all syringes enables easier injection of large animals. One end is attached to the syringe, the other has a handle with a luer lock adapter for needles.

Item No.	Description	Case Ct.
TA180	Floating Needle	12



IV SETS

The Dexco IV set is one of the easiest to use and most durable intravenous systems ever designed. Its pure gum rubber funnel with air vent fits any standard 500 mL bottle, and comes with stainless steel needle.

The Dexco IV set is also available with latex tubing. Tubing measures 5 feet. The Simplex IV set has a molded finger, clog-stop air vent, and pure gum rubber tubing. It fits any standard 500 mL bottle, and comes complete with a stainless steel needle.

Item No.	Description	Case Ct.
1905	IV Set Dexco w/Gum Rubber Tubing 16 G x 1½ in	12
1906	IV Set Simplex w/Gum Rubber Tubing 16 G x ½ in	12
1907	IV Set Dexco w/Latex Tubing 16 G x 1½ in	12
TA8145	IV Set Simplex w/14 G x 2 in	12



NEOGEN IDEAL BLOOD COLLECTION TUBES

Our disposable blood collection vacuum tubes are made of shatter-resistant plastic and are latex free. The tubes are labeled for easy marking, and the interiors are coated with water soluble silicone. Three varieties available.

Sterile.

Features: Standard Blood Collection Tubes: Red top tubes contain no additives and are used primarily for serum lab tests. EDTA Blood Collection

Tubes: Lavender top tubes contain the additive EDTA/K2 and are used primarily for CBC lab tests. SST Blood Collection Tubes: Yellow/gold top tubes contain the additive Polymer Gel. After centrifugation, the gel forms a barrier between the clot and serum.

Standard Blood Collection Tubes

Item No.	Description	Case Ct.
BCT-3P	Blood Collection Tubes Plain 3 mL Red W/13 mm Cap 100/TRAY	18
BCT-5P	Blood Collection Tubes Plain 5 mL Red W/13 mm Cap 100/TRAY	12
BCT-7P	Blood Collection Tubes Plain 7 mL Red W/16 mm Cap 100/TRAY	12
BCT-10P	Blood Collection Tubes Plain 10 mL Red W/16 mm Cap 100/TRAY	12



Standard Blood Collection Tubes



EDTA Blood Collection Tubes

EDTA Blood Collection Tubes

Item No.	Description	Case Ct.
BCT-4EDTA	Blood Collection Tubes EDTA 4 mL Pur w/13 mm Cap 100/TRAY	18
BCT-6EDTA	Blood Collection Tubes EDTA 6 mL Pur w/13 mm Cap 100/TRAY	12
BCT-10EDTA	Blood Collection Tubes EDTA 10 mL Pur w/16 mm Cap 100/TRAY	12



SST Blood Collection Tubes

SST Blood Collection Tubes

Item No.	Description	Case Ct.
BCT-6SST	Blood Collection Tubes SST 6 mL Yel w/13 mm Cap 100/TRAY	12
BCT-9SST	Blood Collection Tubes SST 9 mL Yel w/16 mm Cap 100/TRAY	12

INFUSION PUMP

Small black rubber bulb-type infusion pump. 12 inch tubing with removable stopper.

Item No.	Description	Case Ct.
1908	Infusion Pump Bulb Type	12



Syringes

Neogen's line of syringes offers a broad spectrum of solutions to the animal health industry. From disposable syringes used in one time settings, to durable syringes built to last in the most rigorous animal production environments, Neogen has a variety of syringes available for any job.



PREMIUM (HARD PACKED) DISPOSABLE SYRINGES

Premium Syringe Features

- Latex-free
- Clear polypropylene barrel for easy viewing of fluid levels
- Large, easy-to-read graduations
- Positive plunger stop allows for maximum filling and eliminates costly spills
- Serrated thumb press for positive grip

Premium Needle Features

- Ultra-sharp, tri-beveled, anti-coring tips
- Metal inserts to increase strength for animal health applications

Premium Packaging Features

- Hard pack design for veterinary office or field use
- ISO colour-coded cartridge caps for easy gauge identification
- Sonic-welded end cap provides tamper-evident seal and ensures sterility
- Product information sticker on box and case box for easy inventory control

Item No.	Needle Gauge / Length	Style
1 cc Tuberculin Syringes (Box of 100; case of 10 boxes)		
8859	25 G x 5/8 in	Luer Slip
3 cc Syringes (Box of 100; case of 10 boxes)		
8877	3 cc Syringe with 20 G x 1 in	Luer Lock
8863	3 cc Syringe with 20 G x 1 in	Luer Slip
8868	3 cc Syringe with 22 G x 1 in	Luer Lock
8865	3 cc Syringe with 22 G x 3/4 in	Luer Lock
8864	3 cc Syringe with 22 G x 3/4 in	Luer Slip
8875	3 cc Syringe Only	Luer Lock
8861	3 cc Syringe Only	Luer Slip
6 cc Syringes (Box of 50; case of 10 boxes)		
8882	6 cc Syringe Only	Luer Lock
8881	6 cc Syringe Only	Luer Slip
12 cc Syringes (Box of 80; case of 6 boxes)		
8886	12 cc Syringe with 18 G x 1 in	Luer Lock
8884	12 cc Syringe with 18 G x 1 in	Luer Slip
8885	12 cc Syringe Only	Luer Lock
8883	12 cc Syringe Only	Luer Slip
20 cc Syringes (Box of 50; case of 6 boxes)		
8890	20 cc Syringe Only	Luer Lock
8889	20 cc Syringe Only	Luer Slip
35 cc Syringes (Box of 30; case of 6 boxes)		
8892	35 cc Syringe Only	Luer Lock
8891	35 cc Syringe Only	Luer Slip
60 cc Syringes (Box of 20; case of 6 boxes)		
8895	60 cc Syringe Only	Catheter Tip
8894	60 cc Syringe Only	Luer Lock
8893	60 cc Syringe Only	Luer Slip

Disposables

Ideal's disposable syringes are made specifically for the animal health industry. The syringes and combos are offered in Premium Hard Packs and Standard Soft Pack packaging. Ideal syringes feature a clear polypropylene barrel and large easy to read graduations to help minimize dosing errors. The combos offer ultra-sharp needles and are ISO colour-coded to indicate gauge for easy identification.



STANDARD (SOFT PACKED) DISPOSABLE SYRINGES

Standard Syringe Features

- Clear polypropylene barrel for easy viewing of fluid levels
- Large, easy-to-read graduations
- Positive plunger stop allows for maximum filling and eliminates costly spills
- Medical grade silicone lubricant inside barrel for smooth, accurate delivery
- Serrated thumb press for positive grip

Standard Needle Features

- Ultra-sharp, tri-beveled, anti-coring tips
- Translucent ISO colour-coded hubs for easy viewing of flashback and quick gauge identification in or out of the packaging
- Anti-friction coating for minimal penetration resistance and animal discomfort

Standard Packaging Features

- Paper-to-plastic packaging is tamper-evident and ensures sterility
- Product information sticker on box and case box for easy inventory control
- Each syringe/combo has UPC, lot number, and manufacture and expiration date printed on package



Item No.	Description	
U-100 Insulin Syringes (Box of 100; case of 10 boxes)		
9103	29 G x 1/2 in	NA
U-40 Insulin Syringes (Box of 100; case of 10 boxes)		
9101	1 cc w/ 29 G x 1/2 in	NA
9102	1/2 cc w/ 29 G x 1/2 in	NA
1 cc Syringes (Box of 100; case of 24 Boxes)		
9178	1 cc Syringe Only	Luer Slip
1 cc Tuberculin Syringes (Box of 100; case of 24 boxes)		
9279	1 cc Syringe with 25 G x 5/8 in	Luer Slip
3 cc Syringes (Box of 100; case of 24 boxes)		
9280	3 cc Syringe with 18 G x 1 1/2 in	Luer Lock
9281	3 cc Syringe with 20 G x 1 in	Luer Lock
9200	3 cc Syringe with 20 G x 1 in	Luer Slip
9282	3 cc Syringe with 20 G x 1 1/2 in	Luer Lock
05-2911B	3 cc Syringe with 21 G x 1 in	Luer Lock
9284	3 cc Syringe with 22 G x 3/4 in	Luer Lock
9210	3 cc Syringe with 22 G x 3/4 in	Luer Slip
9285	3 cc Syringe with 22 G x 1 in	Luer Lock
9215	3 cc Syringe with 22 G x 1 in	Luer Slip
9216	3 cc Syringe with 22 G x 1 1/2 in	Luer Slip
9290	3 cc Syringe with 25 G x 5/8 in	Luer Lock
9291	3 cc Syringe with 25 G x 5/8 in	Luer Slip
9170	3 cc Syringe Only	Luer Lock
9100	3 cc Syringe Only	Luer Slip
6 cc Syringes (Box of 100; case of 24 boxes)		
9171	6 cc Syringe Only	Luer Lock
9110	6 cc Syringe Only	Luer Slip
12 cc Syringes (Box of 100; case of 12 boxes)		
9287	12 cc Syringe with 18 G x 1 in	Luer Lock
9250	12 cc Syringe with 18 G x 1 in	Luer Slip
9288	12 cc Syringe with 20 G x 1 in	Luer Lock
9289	12 cc Syringe with 20 G x 1.5 in	Luer Lock
9172	12 cc Syringe Only	Luer Lock
9120	12 cc Syringe Only	Luer Slip
20 cc Syringes (Box of 50; case of 18 boxes)		
9173	20 cc Syringe Only	Luer Lock
9130	20 cc Syringe Only	Luer Slip
35 cc Syringes (Box of 50; case of 6 boxes)		
9174	35 cc Syringe Only	Luer Lock
9140	35 cc Syringe Only	Luer Slip
60 cc Syringes (Box of 25; case of 8 boxes)		
9175	60 cc Syringe Only	Luer Lock
9150	60 cc Syringe Only	Luer Slip

SOFT PACKED DISPOSABLE SYRINGES

Item No.	Size	Case Ct.
Luer Lock		
9263	3 cc SP Syringe 6 pk	100
9264	6 cc SP Syringe 6 pk	72
9265	12 cc SP Syringe 4 pk	60
9266	20 cc SP Syringe 4 pk	60
9267	35 cc SP Syringe 2 pk	60
9268	60 cc SP Syringe 2 pk	60
Luer Lock with Poly Hub Needles		
9184	3 cc SP w/ PH 22 G x 3/4 in 3 pk	100
9185	3 cc SP w/ PH 22 G x 1 in 3 pk	100
9187	12 cc SP w/ PH 18 G x 1 in 3 pk	60
9188	12 cc SP w/ PH 20 G x 1 in 3 pk	60
Luer Slip		
9276	3 cc SP Syringe 6 pk	100
9271	6 cc SP Syringe 6 pk	72
9272	12 cc SP Syringe 4 pk	60
9273	20 cc SP Syringe 4 pk	60
9274	35 cc SP Syringe 2 pk	60
9275	60 cc SP Syringe 2 pk	60
Catheter Tip (Soft Packed)		
9495	60 cc 2 pk	24




PRIMA-SHOT™

The Prima-Shot 50 mL repeater syringe adjusts from 0.5 mL to 5 mL in 0.5 mL increments, and features an easy-to-grip dose adjuster dial to make dosing changes simple. Features a heavy duty Tritan Barrel that is BPA-free, highly chemical and scratch-resistant, and has chemical compatibility that surpasses both polycarbonate and polypropylene.

- 10 fixed dose settings ranging from 0.5 – 5 mL
- Durable rugged construction for repetitive use
- Ergonomic design is easy on hands
- Tritan barrel for superior chemical compatibility
- Large numbers on barrel that are easy to read
- Amber-coloured barrel for UV Protection

Item No.	Description	Case Ct.
377210	Prima-Shot Syringe 25 mL	12
367723	Prima-Shot Syringe 50 mL	12
Replacement Parts		
374172	Barrel with O-Ring Prima-Shot	
374196	Replacement O-Ring 12 pk	
374158	Pull Knob Inserts Multi-Colour 4 pk	



PRIMA INJECTORS

The versatile Prima Injector dials from 1 mL to 12.5 mL in .5 mL increments. It has an ultra-durable design with an elongated needle nut for extra application protection. Available with polycarbonate or polypropylene barrels for compatibility across all compounds.

Item No.	Description	Case Ct.
214102	12.5 mL Adj Injector Premium Polycarbonate Barrel	12
358844	12.5 mL Adj Injector Premium Polypropylene Barrel	12
336569	12.5 mL Adj Injector Premium 5 in Drench Nozzle	12
343499	12.5 mL Adj Injector Value	32
375124	12.5 mL Adj Injector with Hook Drench Nozzle	12
Replacement Barrel		
206251	Barrel w/ O-ring for 12.5 mL Injector	12

PRIMA SHOOTER SYRINGES

The Prima Shooter Syringe combines a sturdy rugged design with ease of use. It features a durable barrel housing that protects from breakage in unpredictable situations and is designed for quick and easy cleaning. Available in 2 mL and 5 mL adjustable dose.

Item No.	Description	Case Ct.
332240	Shooter Premium 2 mL w/ Unilock	12
332264	Shooter Premium 5 mL w/ Unilock	12
Replacement Barrels		
211903	Barrel w/ O-Ring for 2 mL Shooter	12
212153	Barrel w/ O-Ring for 5 mL Shooter	12



PRIMA BMV SYRINGES

The Prima Bottle Mount Vaccinator is designed to allow for efficient and effective application of all animal health injectables. The self-contained, one-handed unit allows for full mobility and automatic refilling of syringes without returning to the bottle. The line feed convertible barrels allow you to change from a bottle mount to a line feed vaccinator.

Item No.	Description	Case Ct.
351241	0.3 mL Adj Premium	12
351883	0.3 mL Fixed Premium	12
332370	0.5 mL Adj Premium	12
347947	0.5 mL Fixed Premium	12
341112	1 mL Adj Premium	12
341129	1 mL Fixed Premium	12
206558	2 mL Adj Premium	12
339591	2 mL Fixed Premium	12
206657	6 mL Adj Premium	12
364166	0.5 mL Adj Premium w/ Line Feed Convertible Barrel	32
348159	2 mL Adj Premium w/ Line Feed Convertible Barrel	12
348364	5 mL Adj Premium w/ Line Feed Convertible Barrel	12
351906	0.3 mL Adj Value	32
343413	0.5 mL Adj Value	32
351937	0.5 mL Fixed Value	32
356352	1 mL Adj Value	32
344724	1 mL Fixed Value	32
344762	2 mL Adj Value	32
361004	2 mL Fixed Value	32
343444	6 mL Adj Value	32
Replacement Barrels		
352842	Barrel w/ O-Ring Adj 0.3 mL BMV	12
351616	Barrel w/ O-Ring Adj 0.5 mL BMV	12
344779	Barrel w/ O-Ring Adj 1 mL BMV	12
206800	Barrel w/ O-Ring Adj 2 mL BMV	12
206855	Barrel w/ O-Ring Adj 5 mL BMV	12
346506	Barrel w/ O-Ring Fixed 2 mL BMV	12



Premium BMV





BMV w/ line feed convertible barrel



PRIMA VACCINATORS

The Prima Vaccinator has a rugged construction and an advanced design, making it more comfortable to use. The syringe has an ergonomically designed handle to better fit the contour of your hand. The large dosage adjuster on the adjustable syringes makes changing dosages effortless.

Item No.	Description	Case Ct.
351777	0.5 mL Adj Premium	12
351784	0.5 mL Fixed Premium	12
351722	1 mL Fixed Premium	12
205803	2 mL Adj Premium	12
351012	2 mL Fixed Premium	12
206008	6 mL Adj Premium	12
359698	0.5 mL Adj Value w/ Draw-Off Caps	32
351739	1 mL Adj Value	32
359704	1 mL Adj Value w/ Draw-Off Caps	32
351746	1 mL Fixed Value	32
366191	2 mL Fixed Value	32
367754	1.5 mL Fixed Value	32
212757	2 mL Adj Value	32
359711	2 mL Adj Value w/ Draw-Off Caps	32
213655	6 mL Adj Value	32
359728	6 mL Adj Value w/ Draw-Off Caps	32
Replacement Barrels		
353344	Barrel w/ O-Ring Adj 0.5 mL Vaccinator	12
353306	Barrel w/ O-Ring Adj 1 mL Vaccinator	12
206404	Barrel w/ O-Ring Adj 2 mL for Vaccinator	12
217400	Barrel w/ O-Ring Adj 6 mL for Vaccinator	12



Prima Premium vs. Prima Value line

Don't be mistaken by premium or value, the quality of these syringes is identical. The difference between the two lines is packaging, colour, and included accessories. See the difference below.

BMV	Premium Line Includes	Value Line Includes
Packaging	Premium Retail Box	Heat Sealed Bag
Handle Colour	Blue	Gray
Bottle Collars	20 mm, 30 mm, 33 mm	20 mm
Fast Fit	100 mL, 250 mL	100 mL
Needle Nut	Metal	Metal

Vaccinator and Injector	Premium Line Includes	Value Line Includes	Value Line with Draw-Offs
Packaging	Premium Retail Box	Heat Sealed Bag	Heat Sealed Bag
Handle Colour	Blue	Gray	Gray
Draw-Off	Unilock Draw-Off	None	20 mm, 30 mm, 33 mm slip Draw-Offs
Tubing	Soft Medical Grade	Soft Medical Grade	Soft Medical Grade
Needle Nut	Metal	Metal	Metal

PRIMA MARC VACCINATORS

The Prima Marc Vaccinator design is unique to the animal health world allowing the user to simultaneously inject and mark the animal. When the dosage is administered, the squeezing of the syringe activates the spray can, which releases a mark of paint onto the animal. This encourages safe administration practices, improves quality control efforts, and saves the cost of missed animals. The Prima Marc BMV allows for convenience and precision all in the palm of your hand. Prima Marc Paints are located on page 11.

Item No.	Description	Case Ct.
366108	0.5 mL Adj Prima Marc BMV	12
366115	1 mL Adj Prima Marc BMV	12
366122	2.0 mL Adj Prima Marc BMV	12
366139	6.0 mL Adj Prima Marc BMV	12
351678	0.5 mL Fixed Vaccinator	12
351685	1 mL Adj Vaccinator	12
352019	1 mL Fixed Vaccinator	12
334039	2 mL Adj Vaccinator	12
339812	2 mL Fixed Vaccinator	12
334053	6 mL Adj Vaccinator	12
Replacement Barrels		
353375	Barrel w/ O-Ring Adj 0.5 mL Marc	12
353399	Barrel w/ O-Ring Adj 1 mL Marc	12
212153	Barrel w/ O-Ring Adj 2 mL Marc	12
211903	Barrel w /O-Ring Adj 5 mL/6 mL Marc	12



Vaccinator



BMV

DRENCH SYRINGES

Automatic drenchers are built to precise standards and suitable for many drenching applications where accuracy and speed are needed. The applicators have extremely rapid fill time and the select-a-dose system allows changes from large to small doses without any product wastage. Accessories include angled, hooked, and floating hook nozzles.

Item No.	Description	Case Ct.
345103	Prima 70 mL w/ 6.6 in Angled Drench	12
346568	Prima 70 mL w/ 23 in Floating Hook	12
346926	Prima 70 mL w/ Hook Drench	12
365903	Prima 80 mL w/ Angled Drench	12
365934	Prima 80 mL w/ Floating Hook	12
365743	Prima 80 mL w/ Long Angled Tip	12
HA070	70 cc w/ Nozzle	1
HA006	200 cc w/ Nozzle	1
HA300	300 cc w/ Nozzle	1



345103



346568



HA300

PRIMA POUR-ONS

Prima Pour-On guns allow for adjustable dosage increments with a select-a-dose system that is accurate and easy to use. The applicators have a shower head nozzle which avoids wastage from splashing, the nozzles are angled for easy administration of pour on medication, and a T-bar applicator is also available. Applicators are packed with soft medical grade tubing, hose armor springs, and nozzle.

Item No.	Description	Case Ct.
330901	30 mL Adj w/ 8mm Hose Barb	12
351098	30 mL Adj w/ 12 mm Hose Barb	12
352910	30 mL Adj w/ 8 mm Barb/T-BAR/Plastic Nozzle	12
201058	50 mL Adj w/ P Barrel and soft tubing	12
332585	70 mL Adj w/ P barrel and soft tubing	12
346797	70 mL Adj w/ Header Card	12
351630	80 mL Adj w/ Tritan Barrel and soft tubing	12
Replacement Barrels		
353627	Barrel 30 mL Pour-On Applicator	12
346179	Barrel 70 mL Pour-On Applicator	12



330901



201058



352910



351630

SAFE-T-FLEX™ POLE SYRINGES

Ideal pole syringes provide a safe, fast, and simple means of injecting large and small animals in situations where adequate restraint facilities, or space, are not available. The syringe consists of a nylon barrel attached to a flexible plunger. Single-dose syringes are available in 10 cc and 20 cc sizes. Replacement kits are available; each kit includes a barrel, nut, O-ring, and luer lock.

Item No.	Description	Case Ct.
S105	42 in 10 cc	6
S205	42 in 20 cc	6
S105R	Pole Syringe Replacement Kit 10 cc	1
S205R	Pole Syringe Replacement Kit 20 cc	1



Syringe Accessories

Neogen offers a variety of accessories to improve or modify the functionality of our syringes. From accessories that enhance our already dependable syringes to replacement parts, Neogen provides more than just a syringe.



PRIMA MIST SPRAYER

The Prima Mist sprayer is a versatile addition to our Prima line of syringes. The mist produced from the sprayer has implications of use ranging from nasal vaccines, to topical application of medicines and anesthetics. The sprayer fits all Prima syringes, it is not suggested for use with syringes over 2 mL or the BMV Fast fit system.

Item No.	Description	Case Ct.
353597	Prima Mist Sprayer	12
364876	Prima Mist Replacement Foam Tip	12



PRIMA NEEDLE GUARD

The Prima Needle Guard is committed to safety in the agriculture industry. Over 80% of farm workers using syringes have accidentally stuck themselves. The Needle Guard keeps the needle covered at all times, and is easy to lock in times of non-use to prevent accidental sticks. It is designed so that needles can be mounted and dismounted safely without directly handling the needle, and without having to remove the Needle Guard to access the needle.

Item No.	Description	Case Ct.
368089	Prima Needle Guard	1



NEEDLE EXCHANGER

Needle exchanger; for enhanced safety in quickly changing needles.

Item No.	Description	Case Ct.
353177	Needle Exchanger Tool	12



NEEDLE NUT

For use with Prima Syringes.

Item No.	Description	Case Ct.
330840	Needle Nut Luer Lock Metal	12
334671	Needle Nut Replacement Metal 12.5 mL Injector	12



VAC-PACS INSULATED BOTTLE HOLDERS

The Hip-Hugger uses a belt that straps around the waist.

The On-Arm straps to the upper arm with two comfortable Velcro straps that are extremely durable.

The Hanging version has a heavy Velcro strap that can be used to hang the unit on cattle chutes, processing carts, etc.



Hip-Hugger

Hanging

On-Arm

Item No.	Description	Case Ct.
120007	On-Arm L 500-1000 mL	12
120151	Hanging L 500-1000 mL	12
120205	Hanging M 250 mL	12
120304	Hip Hugger L 500-1000 mL	12
120359	Hip Hugger M 250 mL	12

THREADED COLLAR WITH BOTTLE

For use with Prima BMV syringes and any line fed syringe.

Item No.	Description	Case Ct.
347251	Bottle 250 mL w/ Travel Cap /Large Threaded Collar	50



PRIMA BOTTLE COLLARS

Standard bottle collars. For use with Prima Syringes.

Item No.	Description	Case Ct.
218902	Bottle Collar, Small 125 mL / 20 mm Necks	12
332080	Bottle Collar, Medium 250 mL / 30 mm Necks	12
339980	Bottle Collar, Large 250 mL / 33 mm Necks	12
349262	Unilock System w/ Draw-Off Spike / 20, 30, 33 mm Collar and Tube	12



Unilock System

PRIMA DRAW-OFFS

For use with tube-fed syringes. Various sizes to fit all bottle types.

Item No.	Description	Case Ct.
141859	Draw-Off 20 mm Small w/o Tubing	12
141804	Draw-Off 30 mm Large w/o Tubing	12
220356	Draw-Off 33 mm X-Large w/o Tubing	12
343567	Draw-Off Vented for Backpacks for Pour-On Tubing	1



PRIMA FAST-FIT

Bottle holder for 100 and 250 mL bottles. For use with Prima BMV Syringes.

Item No.	Description	Case Ct.
336811	Fast-Fit 100 / 125 mL / 20 mm	12
342713	Fast-Fit 250 mL / 30 mm	12
342720	Fast-Fit 250 mL / 33 mm	12



PRIMA SPIKE ACCESSORIES

For use with Prima BMV Syringes, the filtrating spike is over 99% effective in filtering debris from bottles.


Item No.	Description	Case Ct.
366153	Filtrating Spike	12
340795	Spike Removal Tool	1
340931	Spike Buffer Seal BMV	12



PRIMA DRENCH SYRINGE ACCESSORIES

For use with Prima drench applicators.

Item No.	Description	Case Ct.
345110	Floating Hook and Adapter 23 in x 8 in for Pour-On	12
342737	Hooked Drench 8.6 in x 5.9 in	12
333025	Drench Nozzle 3 in Metal w/ Metal Nut	12
332882	Drench Nozzle 4.75 in Metal w/ Metal Nut	12
334602	Drench Nozzle 5 in Metal w/ Metal Nut	12



PRIMA POUR-ON ACCESSORIES

For use with Prima Pour-on guns to evenly distribute liquids.

Item No.	Description	Case Ct.
352507	Nozzle Modular, T-bar	12
365927	Fan Tip Kit	12



PRIMABACK PACKS

Back packs have a hard plastic construction with vented cap to prevent vacuum from forming, for use with Prima Syringes.

Item No.	Description	Case Ct.
341662	1 L w/ Transit Cap / Stepped Cap / Hang Strap	1
341679	2.5 L w/ Transit Cap / Stepped Cap / Hang Strap	1
341686	5 L w/ Transit Cap / Stepped Cap / Hang Strap	1



PRIMA LUBRICATING OIL

This food grade manufacturing hydraulic oil is an anti-wear lubricant. Lubricates syringe O-rings.

Item No.	Description	Case Ct.
338013	Lubrication Oil Food Grade 2 oz	12



PRIMA O-RINGS

Silicone replacement O-rings for Prima Syringes.

Item No.	Description	Case Ct.
351272	0.3 mL	12
349712	0.5 mL	12
358769	1 mL	12
203304	2 mL	12
203502	5-6 mL	12
206305	12.5 mL	12
377890	25 mL	12
374172	50 mL Prima-Shot	12



OB Equipment

Neogen's various styles of calf pullers, including the well known Dr. Frank's® puller, help assist with difficult birthing situations. Our comprehensive line of calving and lambing equipment is ideal for large livestock operations, small farms, and veterinary practices alike.



DR. FRANK'S CALF PULLER

Dr. Frank's Calf Puller is the original obstetrical instrument developed by B.N. Frank, DVM, a specialist in veterinary obstetrics.

Dr. Frank's is a precision instrument that has proven its effectiveness year after year in thousands of calving cases.

The uniform tension of this unique fetal extractor does not permit slipping or backing-up of the jack.

Traction on the calf and the counter-thrust against the cow induce spreading of the pelvic bones, elevation of the abdomen, and arching of the animal's back, all assisting in the natural birthing process.

The entire unit can easily be assembled in minutes. The smooth rods unscrew for easy storage and cleaning. Dr. Frank's comes packaged with a 45 inch OB chain (item 3102, on page 63). Made in the USA.

Item No.	Description	Pallet Qty.
3003	Dr. Frank's Calf Puller	40
3005	Jack for Dr. Frank's	NA
3006	72 in Rods (pair)	NA
3008	Breechen 23 in	NA
3009	Web Strap	NA
3011	Jack Handle for Dr. Frank's Calf Puller	NA
9784	Dr. Frank's Cams and Springs Kit	NA



CALF EZE™ CALF PULLER

Calf Eze is engineered for the user that prefers a dual action pull, or zig-zag motion, instead of a single action pull instrument.

Calf Eze is superior to ratchet-style pullers that rely on grooved shafts with sharp edges and hard-to-clean grooves. Calf Eze utilizes the same non-slip cam traction system as the Dr. Frank's Calf Puller, but has two pulling points that alternately advance to shift the pull from one leg to the other.

The dual-action jack body is made of special strength aluminum and weighs only 17 lbs. The jack assembly mechanism allows the user to easily position the jack on the rods and lock it in place. Steel cams prevent slipping and provide longer jack life. A third chain hook is also provided to give a single-pull option. The Calf Eze is suggested for use with two 45 inch OB chains (not included; available as item 3102 - see next page). Made in the USA.

Item No.	Description	Pallet Qty.
3022	Calf Eze Calf Puller, complete unit	40
3024	Calf Puller Jack for Calf-Eze	NA
3006	72 in Rods (pair)	NA
3008	Breechen 23 in	NA
3009	Web Strap	NA
9785	Calf Eze Cams and Springs Kit	NA



RATCHET-STYLE CALF PULLER

The Ratchet-Style Calf Puller can be your solution to difficult deliveries. Crafted with the finest quality steel, it will stand up to the toughest deliveries without letting you down. The Ratchet-Style Calf Puller features steel rods that easily disassemble for cleaning and storage. The rods have smooth, clean grooves that lessen wear on the jack and give a much quieter ratchet than competitors' pullers. The jack is extremely durable, and has a large, comfortable handle that allows for single person operation. Large handles also greatly reduce the incidence of breakage, or pinched fingers and hands. The rods screw into the breechen, which provides stability during the birthing process. The Ratchet-Style Calf Puller comes with a standard 66 inch rod and has an additional attachment to extend the rod to a full 78 inch. The jack conveniently offers both a zig-zag pull and a single-pull option. Many prefer the zig-zag pull as it helps avoid locking of the calf's shoulders in the pelvic opening.




Item No.	Description	Pallet Qty.
3020	Calf Puller Ratchet Style	40

OB CHAIN AND HANDLE

Neogen's high quality OB accessories are a comfortable fit for any delivery. The OB handle is 1/2 inch diameter and has a poly coated grip. The larger diameter and comfort grip make those difficult births easier on your hands. The durable, high strength OB chains come in 30, 45, and 60 inch lengths. Both the chains and handles are zinc-plated with an added rust inhibitor.

Item No.	Description	Case Ct.
3100	OB Chain 30 in	88
3101	OB Chain 60 in	88
3102	OB Chain 45 in	88
3104	OB Handle w/ Poly Grip	24



OB STRAP

Heavy-duty 1 inch nylon makes our straps weather and mildew resistant. Nickel-plated, welded steel D-rings on both ends. Economical alternative to OB chains. For use with OB Handles.

Item No.	Description
3111	OB Strap Nylon 1 in x 30 in
3112	OB Strap Nylon 1 in x 60 in



CALF SNARE

Neogen's Save-a-Calf Snare loop is easily made larger or smaller by simply sliding the yoke. The snare is fabricated from high grade steel cable and has a thimble for use with Neogen's OB handle.

Item No.	Description
3108	Calf Snare Save-A-Calf Steel Cable



WIRE SAW

The 27-strand stainless steel wire saw provides easy and quick cutting. Packaged in a reusable plastic dispenser to prevent tangling. Handles come as a pair and are nickel-plated for easy clean-up. Thumb screws lock the wire in place.

Item No.	Description
3115	OB Wire Saw Stainless Steel 32 ft
3116	OB Wire Saw Handles 1 pair



COW LIFT

For cows that cannot stand up on their own, such as right after giving birth. The gummed clamps are adjustable to the cow's size. The lift is placed over the cow's hips and tightened as necessary by cranking the handle, which secures the cow as she is lifted. This enables the cow to regain the use of her legs so she can be mmked or treated.

Item No.	Description
TA1800	Cow Lift - up to 2,000 lbs / 900 kg
TA1805	Cow Lift - up to 3,300 lbs / 1,500 kg



PRO-FIX® TUBES

Designed for the repair of rectal prolapse in swine, ovine, and bovine, Pro-Fix Tubes eliminate the stricture problems sometimes associated with reduction methods of treatment. These precision molded, disposable plastic tubes are made in 5/8 in, 3/4 in, 1 in, 1 1/4 in, 1 1/2 in, and 2 in OD sizes. The tubes are fully tapered on one end for easier insertion and have extra deep centre grooves for better retention. The 1 1/2 in and 2 in OD sizes have holes in the centre groove for suturing in place. Pro-Fix Tubes are available in various sizes.

Item No.	Description
R050	Pro-Fix Tubes 5/8 in 24 ct
R075	Pro-Fix Tubes 3/4 in 16 ct
R100	Pro-Fix Tubes 1 in 12 ct
R125	Pro-Fix Tubes 1.25 in 16 ct
R150	Pro-Fix Tubes 1.5 in 16 ct
R200	Pro-Fix Tubes 2 in 9 ct





Semen Collection

Ideal's ElectroJac® 6 combines the best of more than 30 years experience in automated semen collection instrument design and manufacturing, with state-of-the-art technology, featuring its unique flexibility of automatic and manual settings, ElectroJac 6 is the perfect field unit for even hard-to-collect bulls.

- Unparalleled simplicity of use
- Improved electronics to prevent overcharging
- Powerful gel cell battery
- Bright numeric display
- Trade-in program available

ELECTROJAC 6

The Complete ElectroJac 6 System, with your choice of 1 inch, 2 inch, 2.5 inch, or 3 inch probe.

Item No.	Description
Complete Kits:	
6906-1	ElectroJac 6 Complete w/1 in probe
6906-2	ElectroJac 6 Complete w/2 in probe
6906-25	ElectroJac 6 Complete w/2.5 in probe
6906-3	ElectroJac 6 Complete w/3 in probe
The Complete Kits include:	
6901	ElectroJac 6 Carrying Case
6912	ElectroJac 6 Probe Cable Assembly 20 ft
6915	ElectroJac 6 Wireless Dwell Remote
6916	ElectroJac 6 European/Universal Battery Charger
6983	Semen Collection Director Cone
6985	Semen Collection Handle Assembly
6986	Scrotal Tape Metric
Replacement Probes:	
6989	ElectroJac 6 Ram/Boar Probe 1 in (25 mm)
6990	ElectroJac 6 Bull Probe 2 in (51 mm)
6992	ElectroJac 6 Bull Probe 2.5 in (64 mm)
6993	ElectroJac 6 Bull Probe 3 in (77 mm)
ElectroJac 6 Accessories	
6902	ElectroJac Replacement Battery
C101	Semen Collection Vials 100/bag
C100	Semen Collection Cone 100/pk
6910	ElectroJac 6 Dwell Cable Assembly 20 ft
6903	ElectroJac 3, 4, 5 Battery Charger AC



Insemination

Neogen provides a complete line of artificial insemination and breeding supplies that reproduction professionals trust. With over 30 years of experience in providing the tools for livestock semen collection and artificial insemination, we have built product lines that meet the Neogen quality standards professionals in the field count on to get the job done.



INFUSION TUBES

Drilled infusion tubes have one uniquely drilled end which allows for a tight fit with any standard luer-tip syringe, eliminating the need for connector adapters. Made from high quality polymers with round, smooth, and polished ends for easier insertion. 25 tubes/bag, 40 bags/case.



Item No.	Description	Case Ct.
AI2100	Infusion Tubes Drilled 21 in 25/bag	40

CATHETERS

For artificial insemination in swine.

Item No.	Description	Case Ct.
333117	Catheter Sow Green w/ Handle and Plug 5 pk	500
333124	Catheter Sow Green w/o Handle 5 pk	500
334848	Catheter Foam White w/o Handle 25 pk	500
334855	Catheter Foam White w/ Handle and Plug 5 pk	500
334862	Catheter Gilt Green w/o Handle 5 pk	500
334879	Catheter Gilt Green w/ Handle and Plug 5 pk	500
334886	Catheter Spiral White w/ Handle and Plug N2 Spiral 5 pk	500
334893	Catheter Spiral Blue 25 pk	500
352316	Catheter Uterine deep Soft/Flex w/o Ball Tip/handle	500
369826	Catheter PCAI Gilt Green/White 5 pk	500
369833	Catheter PCAI Sow Pink/White 5 pk	500



SEMEN BOTTLE

For artificial insemination in swine.

Item No.	Description	Case Ct.
334909	Semen Bottle 100 mL w/ Break-Away Tip and Caps 5 pk	500



Disposable Boots

Developed for protective and sanitary use in the farm and industrial applications, Neogen's disposable boots are the solution to keeping your shoes clean and facilities biosecure. Whether used on the farm or in the processing plant, Neogen disposable boots are of the highest quality and meet the needs of the most demanding user.



ELASTIC BAND BOOT

An economical disposable boot designed for protective and sanitary use in farm and industrial operations, the disposable boot is 15 in tall and available in 2.0, 3.0, 4.0, and 6.0 mm thicknesses. The boot has an elastic band sealed at the top for easy application and better protection. Available in two sizes, medium/large (up to shoe size 12) and extra large (up to shoe size 16). Boots are packaged in bags of 50 boots (25 pair).

Item No.	Description	Case Ct.
BC300	Boot w/Elastic Band 3.0 mm M/L 25 pair	10
BC300-XL	Boot w/Elastic Band 3.0 mm XL 25 pair	10
BC400	Boot w/Elastic Band 4.0 mm M/L 25 pair	10
BC400-XL	Boot w/Elastic Band 4.0 mm XL 25 pair	10
BC600-XL	Boot w/Elastic Band 6.0 mm XL 25 pair	8



MAXIBOOT

Developed for protective and sanitary use in farm and industrial applications, the disposable MaxiBoot is made from a dark brown blended polymer film. It is available in 4.0 mm (heavy duty) film thickness. The calf length boot is extra tall (18 in) and is available in two sizes: medium/large (up to shoe size 11) and extra large (up to shoe size 14). The MaxiBoot has a unique press-tab fastener closure. Boots are packaged in resealable plastic bags of 50 boots (25 pair).

Item No.	Description	Case Ct.
MT400-XL	MaxiBoot 4.0 mm XL 25 pair	8



POLYBOOT

Designed for use in farm and industrial applications, the disposable PolyBoot is available in 3.0 mm film thickness and is 18 inches tall. This calf length boot is available in one size, medium/large, is shaped to fit a wide range of shoe sizes, and features press-tab fasteners. Boots are packaged in resealable plastic bags of 50 boots (25 pair).

Item No.	Description	Case Ct.
BT300	PolyBoot 3.0 mm M/L 25 pair	12



POLYETHYLENE OVERBOOT

Clear plastic disposable over-the-shoe boots with tie closure. Easy to slip on and off.

Item No.	Description	Case Ct.
TA724	Polyethylene Overboots 3.0 mm 50 pair	10
TA723	Polyethylene Overboots 4.0 mm 25 pair	10



Gloves

Neogen gloves offer a broad selection of style and materials to provide the user superior protection, durability, and sensitivity at an affordable price. Most gloves have a fully textured palm and fingers, allowing a strong grip for both wet and dry applications.



BLUE NITRILE GLOVES RETAIL PACKS

Our nitrile gloves are 100% latex free. Fully textured palm and fingers allow a strong grip for both wet and dry applications. They provide superior puncture, tear, and chemical resistance with excellent tactile sensitivity and dexterity. Gloves are powder-free and ambidextrous, with a thickness of 4.0 mm and length of 9 inches. Packs of 10.

Item No.	Description	Case Ct.
3132	Blue Nitrile Powder-Free Gloves M 10/pk	24
3120	Blue Nitrile Powder-Free Gloves L 10/pk	24
3122	Blue Nitrile Powder-Free Gloves XL 10/pk	24



LATEX GLOVES

Ideal's powder-free, latex gloves provide superior protection at an affordable price. Gloves are 4.0 mm in thickness and are available in four sizes. 100 per box.

Item No.	Description	Case Ct.
AT300PF-S	Latex Gloves Powder-Free S 100/bx	10
AT300PF-M	Latex Gloves Powder-Free M 100/bx	10
AT300PF-L	Latex Gloves Powder-Free L 100/bx	10
AT300PF-XL	Latex Gloves Powder-Free XL 100/bx	10



PRO-TECT 6 MM NITRILE GLOVES

Ideal's PRO-TECT Nitrile Gloves are designed to provide extra protection at 6 mm in thickness. Gloves are powder-free and latex-free. 100 per box.

Item No.	Description	Case Ct.
TA586-PRO	PRO-TECT Nitrile Gloves Powder-Free 6 mm M 100/bx	10
TA596-PRO	PRO-TECT Nitrile Gloves Powder-Free 6 mm L 100/bx	10



TRUEBLACK PREMIUM POWDER-FREE NITRILE GLOVES

Ideal’s TrueBlack Premium Nitrile Gloves provide strength and durability for extended use and increased hand protection. Gloves are 5.0 mm in thickness in the palm, 5.5 mm in thickness in the fingertip, and have a textured finish. Available in four sizes. 100 per box.

Item No.	Description	Case Ct.
AT455-S	Nitrile Gloves Black S 100/bx**	10
AT450-M	TrueBlack Premium Nitrile Gloves M 100/bx	10
AT450-L	TrueBlack Premium Nitrile Gloves L 100/bx	10
AT450-XL	TrueBlack Premium Nitrile Gloves XL 100/bx	10

** = Not TrueBlack branded



TRUEBLUE PREMIUM POWDER-FREE NITRILE GLOVES

Ideal’s TrueBlue Premium Powder-Free Nitrile Gloves provide superior protection, durability, and sensitivity at an affordable price. Gloves are 4.0 mm in thickness in the palm, 4.5 mm in the fingertip, and have a textured finish. Available in four sizes. 100 per box.

Item No.	Description	Case Ct.
CB400-S	TrueBlue Premium Powder-Free Nitrile Gloves S 100/bx	10
CB400-M	TrueBlue Premium Powder-Free Nitrile Gloves M 100/bx	10
CB400-L	TrueBlue Premium Powder-Free Nitrile Gloves L 100/bx	10
CB400-XL	TrueBlue Premium Powder-Free Nitrile Gloves XL 100/bx	10



OB Sleeves

Neogen is the standard of quality when it comes to breeding sleeves. Neogen offers multiple combinations of lengths, thicknesses, and styles, including sterile options, whether used in artificial insemination, pregnancy checking, or surgery. With Neogen's numerous choices of OB sleeves to choose from, we provide the perfect fit for everyone.



GENIA EXTRA-PROTECTION

Genia OB sleeves provide extra protection by covering the shoulder, chest, and underarm area. The sleeve is held firmly in place by a neck strap. Sleeve measures 48 inches from fingertip to neck strap.

Item No.	Description	Case Ct.
307620	Genia Extra-Protection OB Sleeve w/ Neck Strap 50/bx	10



BREEDER-SLEEVE®

Affordable, quality shoulder length sleeves, packaged in a convenient dispenser box.

Red, 1.0 mm, 36 in Length

Item No.	Description	Case Ct.
BS100	Breeder-Sleeve 100/bx	10



BREEDER-SLEEVE PLUS

Economical extra-long sleeves, packaged in a convenient dispenser box. Orange, 1.0 mm, 39 in Length

Item No.	Description	Case Ct.
BSP100	Breeder-Sleeve Plus 100/bx	10



EQUISLEEVE™

Disposable EquiSleeve OB Sleeves are made of a thin, clear 1.0 mm copolymer film, providing the ultimate in sensitivity and softness. Due to the unique copolymer film, the EquiSleeve has extremely smooth seams, reducing irritation of the sensitive mucosal lining, while providing the strength necessary for the most difficult exams. Packaged in resealable, pressure-sealed bags (PolyPack). 1.0 mm, 36 in Length.

Item No.	Description	Case Ct.
ES100	EquiSleeve 100/bx	12
SES100	EquiSleeve Clear Sterile 12/pk	12



MAXI SLEEVE®

Disposable Maxi Sleeve OB sleeves are made from a blended polymer film for an optimum strength to sensitivity ratio. Packaged in resealable, pressure-sealed bags (PolyPack). 39 in Length.

Item No.	Description	Case Ct.
M100	Maxi Sleeve 1.1 mm Yellow 100/pk	10
M125	Maxi Sleeve 1.25 mm Brown 100/pk	10
M125-F	Maxi Sleeve 1.25 mm Brown Fingerless 100/pk	10
SM100	Maxi Sleeve 1.1 mm Yellow Sterile 12/pk	12
SM125	Maxi Sleeve 1.25 mm Brown Sterile 12/pk	12
SM125-F	Maxi Sleeve 1.25 mm Brown Sterile Fingerless 12/pk	12



POLYSLEEVE®

Disposable, polyethylene PolySleeve OB Sleeves are packaged in resealable, pressure-sealed bags (PolyPack).

Item No.	Description	Case Ct.
P100	PolySleeve 1.1 mm Buff 100/pk	12
P125	PolySleeve 1.25 mm Blue 100/pk	10
P125-S	PolySleeve 1.25 mm Blue Small Hand 100/pk	10
SP100	PolySleeve 1.1 mm Sterile 12/pk	12
SP125	PolySleeve 1.25 mm Sterile 12/pk	12
SP125-S	PolySleeve 1.25 mm Sterile Small Hand 12/pk	12



STANDARD OB SLEEVES

Designed for OB and AI use, Ideal's shoulder-length sleeve is a premium sleeve at an economy price. Made from a high quality polyethylene fiber for maximum strength and optimal sensitivity. 1.25 mm, 36 in Length.

Item No.	Description	Case Ct.
3106	Standard OB Sleeves 36 in Clear 100/bx	20
3110	Standard OB Sleeves 36 in Clear 10/pk	80
TA843	Standard OB Sleeves 36 in Orange 100/bx	10



Apparel

Designed for veterinarian use when performing large animal procedures, Ideal's veterinary apparel is made from a urethane coated nylon oxford material. The garments are navy blue in colour, highly water resistant, breathable, and machine washable. Custom sizes of all garments may be special ordered.



OB GOWN

Full length navy blue veterinary OB Gown with short sleeves for maximum coverage and protection when performing cesarean sections and other large animal procedures. Made from urethane-coated nylon oxford material that is highly water resistant, breathable, and machine washable. Has Velcro fasteners in the back for a secure fit. Available in Large and XLarge, or custom size.

Item No.	Description	Case Ct.
VA500-L	OB Gown L	1
VA500-XL	OB Gown XL	1



OB JACKET

Short sleeve OB jacket with button front closure that covers down to the mid-thigh and is breathable for comfort. Navy blue coloured urethane-coated nylon that is highly water resistant and machine washable. Available in Large and XLarge, or custom size.

Item No.	Description	Case Ct.
VA301-L	OB Jacket L	1
VA301-XL	OB Jacket XL	1



OB SUIT

Get full protection for large animal procedures with this highly water resistant and breathable urethane-coated nylon oxford OB Suit. Short sleeve pull over top and pants bottoms sold separately. Navy blue colour and machine washable. Available in Medium, Large, and XLarge, or custom size.

Item No.	Description	Case Ct.
VA101-M	OB Suit Top M	1
VA101-L	OB Suit Top L	1
VA101-XL	OB Suit Top XL	1
VA102-M	OB Suit Bottom M	1
VA102-L	OB Suit Bottom L	1
VA102-XL	OB Suit Bottom XL	1



SHOULDER PROTECTOR

Urethane coated nylon oxford material that is highly water resistant for protection from fluids and manure during OB and AI procedures. Has an elastic sleeve support to keep OB sleeves in place. May be worn on either arm and has adjustable straps for a secure fit. Navy blue in colour and machine washable. Available in Large.

Item No.	Description	Case Ct.
VA900-L	OB Shoulder Protector L	1



Index

A

Aluminum Hub Needles15, 16
AluShield.....1

B

BandB Eze-Trim Foot Rot Shears...9
Back Packs30
Band Castrator.....6
Bandages2
Blades, Scalpel4,5
Blue Nitrile Gloves39
BMV Syringes23,24,25
Bottle Collars28
Bottle, Calf7
Bottle, Lamb8
Breeder- Sleeve Plus41
Breeder-Sleeve41

C

Cadaver Bags3
Calf Bottle7
Calf-Eze Calf Puller31
Calf Pullers.....31,32
Calf Snare.....32
Calf Weaners7
Carbon Steel Scalpel Blades4
Castrator Bands.....6
Castrators.....6
Catheters36
Chromic Suture5
Cohesive Bandages2
Concentrate Paint.....12
Controlled Flow Lamb Nipples...8
Cow Lift33

D

D3 Detectable Needles14
D3X Detectable Needles15
Disposable Boot37,38
Disposable Needles.....16,17
Disposable Scalpels6
Disposal Cadaver Bags.....4
Dr. Frank's Calf Puller31
Draw-Offs29
Drench Syringe Accessories29
Drench Syringes.....25
Durable Needles17
Dyes.....11

E

Elastic Band Boot37
ElectroJac 635
EquiSleeve41

F

Fan Tip Kit.....29
Fast-Fit29
Floating Needle17
Floating Teat.....7
Fluid Feeders, Oral7
Foot Rot Shears, BandB Eze-Trim.....9

G

Gauze1
Genia Extra- Protection41
Gloves39,40
Gum Rubber Lamb Nipples.....8

H

Hoof Rot Shears9
Hoof Wrap1

I

Infusion Pump18
Infusion Tubes36
Injectors22
IV sets18

L

Lamb Bottle with Red Screw
On Teat8
Lamb Nipples.....8
Latex Gloves.....39
Leg Bands10
Livestock Markers.....11
Lubricating Oil, Prima30

M

Marking Sticks, Retractable11
Maxi Sleeve42
MaxiBoot.....37
Metered Tip Paints and Dyes.....11
Mini Scalpels.....5
Mist Sprayer27

N

Needle Exchanger.....27
Needle Guard.....27
Needle Nut.....27
Needles13-17
Nipples.....8

O

OB Chain and Handles32
OB Gown43
OB Jacket.....43
OB Sleeves41,42
OB Strap.....32
OB Suit.....43
OB Wire Saw33
Oral Fluid Feeders7
O-Rings30

P

Paints11,12
Plastic Leg Bands.....10
Pole Syringes26
PolyBoot37
Polyethylene Overboot38
Polypropylene Hub Needles 15,16
PolySleeve42
Pour-On Accessories.....29
Pour-Ons26
Premium Marking Sticks
and Paints.....11,12
Prima Concentrate Paint.....12
Prima Marc.....11

Prima Marc Metered Tip Paints
and Dyes11
Prima Marc Vaccinators.....25
Prima BMV Syringes23
Prima Bottle Collars28
Prima Injectors22
Prima Draw Off28
Prima Lubricating Oil30
Prima Mist Sprayer27
Prima Needle Guard.....27
Prima O-Rings.....30
Prima Pour-On Accessories.....29
Prima Pour-Ons26
Prima Shooter Syringes.....22
Prima Vaccinators.....24
Prima Back Packs30
Prima-Shot.....21
Pro-Fix Tubes33
PRO-TECT 6 mm Nitrile Gloves....39

R

Ratchet-Style Calf Puller32
Red Screw-On Teat8
Retractable Marking Sticks11

S

Safe-T-Flex Pole Syringes26
Scalpel4,5
Scalpel Blades4
Scalpel Handles.....4
Semen Bottles36
Shooter Syringes22
Shoulder Protector.....43
Sleeve Stoppers.....3
Spike Accessories29
Spray-On12
Spray-On Glo12
Spray-On II.....12
Stainless Steel Needles17
Standard OB Sleeves42
Surgical Blades4
Sutures.....5
SyrFlex Cohesive Bandages.....2
Syringes19-26

T

Tape, Elastic Adhesive3
Tattoo Kits.....10
T-Bar Nozzle.....26
Teat, Floating7
Threaded Collar with Bottle.....28
Treader Disposable Boot.....37
TrueBlack Premium Powder-Free
Nitrile Gloves40
TrueBlue Premium Powder-Free
Nitrile Gloves40
Tubes, Infusion36

V

Vaccinators24,25
Vac-Pacs.....28
Vicryl Suture5

W

Wire Saw33
Wrap, Hoof.....1



©Neogen Corporation 2024. All rights reserved. Neogen, Dr. Frank's, ElectroJac, EqStim, Fura-Zone, GeneSeek, Horse Sense, Ideal, Igenity, JOLT, Maxi Sleeve, PolySleeve, Prima, Pro-Fix, Stress-Dex, SureBond, and UNIPRIM are registered trademarks and AluShield, Calf Eze, D3, D3X, EquiSleeve, MegaShot, PolyPetite, PolyShield, Prima Marc, Prima-Shot, Pro-Shot, and Safe-T-Flex are trademarks of Neogen Corporation. Demotec, Easy Bloc, and FuturaPad are registered trademarks of Demotec. All other marks mentioned are the property of their respective owners.



Questions?

To find out how we can help your business grow,
contact your local Neogen territory manager today



neogenaustralasia.com.au

07 3736 2134 | naa-as@neogen.com

14 Hume Drive, Bundamba QLD 4304