


Neogen One Plate *Enterobacteriaceae*

One plate. One test. No Overlays.

The One Plate Solution is the only single-test method that's ISO-16140-2 validated for all food types.¹


One Plate *Enterobacteriaceae* (OP EBAC) uses less media and requires less labor, with faster time-to-result compared to VRBGA per ISO 21528:2017. In addition, no overlay or confirmation step is required, which adds more efficiency.



¹ Validation studies performed on a broad range of foods.

Less media. Less labor. More savings.

Not only does OP-EBAC offer a single test vs. two tests, but the method also eliminates an additional 5-ml overlay, which means you save even more on media and labor.



³ Plates can be incubated for up to 22 hours.

One Plate EBAC: One and done.



Ordering Information

Description	Format	Cat. No.	Pack Size
One Plate <i>Enterobacteriaceae</i>	DCM	NCM1019A	500 g
		NCM1019B	5 kg
		NCM1019C	10 kg
		NCM1019D	25 Kg

